

HY202

The ideal solution for sit-to-stand hoisting

HY202 sit-to-stand hoist is intended for a user with some weight bearing capacity. The hoist is excellent to use for transfer to and from a toilet. It's low height makes it easy for users to place their feet onto the footplate and the centre of gravity is optimised to provide stability. A large lifting span allows it to be used by both children and adults.

Easy to manoeuvre

HY202 is a compact hoist with a small footprint which facilitates turning in narrow spaces. The big wheels roll easily when moving which provide excellent manoeuvre ability for the caregiver. Altogether it contributes to smooth and more comfortable transfers for the user.



TECHNICAL DATA

1200 mm Length	660 mm Width	880 mm Max fork Distance
880 mm Lifting Distance	60 mm Ground Clearance	24V 8000n Motor
60-80 Times Battery Capacity	60 Kg Net Weight	180 Kg Max Loading Weight