

Document No.	Revision No.	Effective Date	Expiration Date
STLSOP-0008rev02	2	3/24/22	3/24/23

## Sampling Procedure for Cannabis Flower and Cannabis Products

Reference Method(s):

Oklahoma Medical Marijuana Authority

Approved by Lab Director:



Revision No.	Effective Date	Description of Change(s)
0	1/6/21	New Document
1	4/13/21	Updated SOP Numbering
2	3/24/22	Combined sampling procedures for all matrices, added procedure for sampling prerolls, added field sampling log as Appendix A

#### 1. Introduction:

- a. This procedure describes the OMMA-approved sampling requirements for cannabis harvest batches and cannabis production batches for full compliance testing. This includes both the Test sample (TS) and Retention Sample (RS).
- All samples, regardless of matrix, should be submitted in their final (ready-for-sale) form.
   i.e. flower should be dried and cured, edibles should be in their finished form (flavors, frostings, etc) and weight, and concentrates should be in their finished form as they are intended to be sold (example samples being sold as vape carts should be submitted in the cartridge)
- c. A field sampling log is provided in this SOP but growers and producers can opt to use an internally generated form if preferred.

#### 2. Equipment needed:

- a. Stainless Steel Bowl
- b. Table top balance (able to accommodate the weight of your bowl and a maximum of 23g).
- c. Stainless Steel Tongs
- d. Fresh 10% Bleach solution (in a spray bottle)
- e. 70% Ethanol or 70% Isopropyl Alcohol solution (in a spray bottle)
- f. Clean sampling containers
- g. Clean gloves
- h. Clean paper towels
- i. Sampling jars (two per each sample)
- j. Heat gun or hair dryer

#### 3. Cleaning Procedure:

- a. Wearing gloves, prepare 10% bleach solution by mixing 9 parts water to 1 part bleach. Mix well.
- b. Spray sampling tools (mixing bowl and tongs) liberally with 10% bleach solution and allow it to sit for 5 minutes.
- c. Dry with clean paper towels.
- d. Spray sampling tools with 70% Ethanol (or Isopropyl Alcohol) and allow it to sit for 5 minutes.
- e. Wipe dry with clean paper towels.
- f. Set aside on clean surface.

#### 4. Sampling Procedure for Harvest Batches:

- a. Samples should be submitted to the laboratory in TWO sample containers with equal weight of sample in each. One sample is the Test Sample (TS), the other is the Retention Sample (RS).
- b. Selecting sampling sizes.
  - i. Harvest batch size: Less than 6lbs
    - 1. Place clean stainless steel bowl on the table top balance and tare balance.
    - 2. Using the clean tongs randomly select flower from the harvest batch until 10g has been weighed in the bowl.
    - 3. Gently mix flower with tongs.
    - 4. Equally separate the sample into two separate sample containers (5g in each).
    - 5. Label sample containers with the Strain Name, Batch ID, Weight (g), and TS (Test Sample) or RS (Retention Sample).
    - 6. Seal both samples with a custody seal or a piece of tape over the lid.
    - 7. A total of 10g will be submitted.

- ii. Harvest batch size: Greater than 6 lbs
  - 1. Use the following calculation to determine your sampling size for the test and retention samples:

Harvest batch size (lbs) x 2.268 = sampling size (grams) Example: 8.5lb harvest batch size Sampling Size = 8.5 x 2.268 = 19.28g

- 2. Place clean, stainless steel bowl on the table top balance and tare balance.
- 3. Using the clean tongs, randomly select flower from the harvest batch until the weight calculated from the above formula has been weighed into the bowl. (Using the example above, 19.28g would be placed in the bowl).
- 4. Remove bowl from balance and gently mix flower with tongs.
- 5. Place sampling jar on balance and tare.
- 6. From the bowl, randomly select 5g of sample and place into each of the two separate sampling containers.
- 7. Label sample containers with the Strain Name, Batch ID, Weight (g), and TS (Test Sample) or RS (Retention Sample).
- 8. Seal both samples with a custody seal or piece of tape over the lid.
- 9. The remaining flower in the bowl can be returned to the harvest batch.
- 10. A total of 10g will be submitted for testing.

#### 5. Sampling Procedure for Prerolls:

- a. Select the appropriate number of prerolls to be submitted for testing according to Appendix E of the current Title 310 Chapter 681 regulations.
- b. A minimum of 4 one gram prerolls or 8 half gram prerolls must be submitted.
  - i. Label sample containers with the Sample Name, Batch ID, Weight (g), and TS (Test Sample) or RS (Retention Sample).
  - ii. Randomly select the appropriate number of prerolls and place in the sampling container.
  - iii. Label sample containers with the Sample Name, Batch ID, Weight (g), and TS (Test Sample) or RS (Retention Sample).
  - iv. Seal both samples with a custody seal or piece of tape over the lid.
- 6. Sampling Procedure for Concentrate Production Batches:
  - Samples should be submitted to the laboratory in two sample containers with equal weight of sample in each. One sample is the Test Sample (TS) the other is the Retention Sample (RS).
  - b. A total of 6 grams of concentrate must be submitted to the lab for testing. (3g for the TS and 3g for the RS).
  - c. Samples should be submitted in their final sale-ready form.
    - i. Carts submit 6 one gram carts or 12 half gram carts.

- ii. Distillate/Crumble/Shatter/Batter/Etc. submit 6g (3g in each container).
- d. Distillate samples may need to be warmed using the heat gun or hair dryer to more easily remove sample from the jar.
- e. Procedure:
  - i. Place first sampling container on the table top balance and tare balance.
  - ii. Using the clean spatula/tongs, weigh 3g of concentrate into the first sample container.
  - iii. Using the sampling tools, weigh an additional 3g of concentrate into the second sample container.
  - iv. Label sample containers with the Sample Name, Batch ID, Weight (g), and TS (Test Sample) or RS (Retention Sample).
  - v. Seal both samples with a custody seal or piece of tape over the lid.
  - vi. A total of 6g will be submitted for testing.

#### 7. Sampling Procedure for Edible and Topical Production Batches:

- a. Select the appropriate number of units to be submitted for testing according to Appendix D of the current Title 310 Chapter 681 regulations.
- b. Samples should be submitted to the laboratory in two sample containers with equal weight of sample in each. One sample is the Test Sample (TS), the other is the Retention Sample (RS).
- c. A total of 10grams of product must be submitted to the lab for testing. (5g for the TS and 5g for the RS).
- d. If a single product exceeds 10g, then submit a single finished product for testing. The retention sample will be prepared in the lab.
- e. Sampling Procedure
- a. Place first sampling container on the table top balance and tare balance (additionally, edible and topical samples can be submitted in their final product packaging).
- b. Using the clean spatula/tongs, weigh 5g of finished product into the first sample container.
- c. Using the sampling tools, weigh an additional 5g of finished product into the second sample container.
- d. Label sample containers with the Sample Name, Batch ID, Weight (g), and TS (Test Sample) or RS (Retention Sample).
- e. Seal both samples with a custody seal or piece of tape over the lid.
- f. A total of 10g of product will be submitted for testing.

#### 8. Storage of Samples

a. After each sampling event, samples should be placed in a refrigerator or on ice in a cooler until transportation to the laboratory.

#### 9. Transportation of Samples

a. Samples should be transported on ice for preservation.

### Field Sample Log • Transport Manifest • Chain of Custody

www.scissortaillabs.com

405-788-0247

Scissortail Laboratories

Sampling SOP: STLSOP-0008rev02				omma: Laaa-lhy	W-6BUV	OBNDD: 62042	2	24	408 N	√W 10	TH :	St Or	(C, (	)K 73	3107				-
				BUSINE	SS INFORMATION														
Business Name:			OMMA License Number: Business Address:																
Contact Name:			Phone Number:			Email Address:													
				TRANSP	ORT INFORMATION														
Transporter Name:			OMMA License Nu			Transporter Address:													
Make and Model of Transport Vehicle:			License Plate of Tr	ansport Vehicle:		Departure Time:		Arrival Time:											
Driving Directions:			1			1													
			SAMPLING	INFORMATION								TES	STIN	G RE	QUE	ST			
Sampler Name and Title:			Sampled date:		Start Time:	End Tin	ne:		iţe									$\Box$	
				Ambient Temperature and Other Sampling Conditions:  All listed samples are representative of the associated products & batches.				OMMA Compliance Suite		e e	vents	Filth and Foreign Matter	tent						
Sampler Signature:							\ Comp	ا اخ	Terpene Profile	Residual Solvents	nd Fore	Moisture Content	Water Activity	Heavy Metals	ides	jal	Suixc		
Sample Name		Batch ID		Batch Size (weight or units)	Primary Samp Size (weight or units)	Retention Samp Size (weight or units)	Matrix (F•C•E/T)	R&D?	OMM	Potency	Terper	Resid	Filth a	Moist	Water	Heavy	Pesticides	Microbial	Mycotoxins
Relinquished by:		Date:	Time:	Received by:	•	Remark	s:												
Relinquished by:	Received by:																		





# Cannabis and Hemp Sampling and Sample Handling Certificate

This Certifies that is a	trained to
sample harvest batches and production batches for compliance	e testing
and to handle test samples according to the following SO	Ps:

STLSOP-0008rev02

Laboratory Director: Date:	