CIRCAID

reduction kit[™]

A revolutionary product for lymphedema treatment

medi. I feel better.

medi

An introduction to **circaid[®] reduction kit**[™]

The medi[®] company has nearly a century of experience manufacturing compression therapy garments for the venous health and lymphedema markets. Through collaboration with the certified lymphedema therapist community, medi developed a unique, trim-to-fit application with circaid[®] technology products to treat lymphedema.

As the health care environment changes, so does the technology and treatment plan. The circaid reduction kit delivers both patient and therapist benefits to help enhance long-term results through the delivery of more consistent compression and stronger patient engagement.



Phase

Decongestion

Adjustable inelastic com-pression device – for individual, effective and fast reduction • circaid reduction kit medi pcs

circaid reduction kit

Phase

Transition

- medi pcs circaid reduction kit

circaid reduction kit & medi's lymphedema therapy concept

The medi lymphedema therapy concept is a new approach in the treatment of lymphedema with the goal of long-term care for the patient. This approach helps create patient independence early in the treatment process that extends into the maintenance phase, providing more flexibility and freedom in day to day life.

circaid juxtafit

Phase

Maintenance

- circaid graduate
- medi pcs
- circaid reduction kit


mediven plus

mediven 550

circaid reduction kit for the clinician

Measurable therapeutic compression

- The BTS[™] (Built-in-tension system) provides necessary targets for limb reduction throughout the decongestion phase
- Provides patients with guidance and measurable goals for self treatment



Adjustable Compression

- The juxtaposed bands accommodate for limb volume fluctuations to maintain the required compression level
- Easy to perform minor fit adjustments for ambulation and additional comfort
- Can be adjusted for night-time use when less compression may be acceptable





Trimmed to fit patients' specific requirements:

- Each reduction kit piece can be trimmed in clinic to fit the patient's limb and repeated throughout the decongestion phase
- The reduction kit enhances longterm function of the lymphatics through the maintenance phase by limiting deterioriation of swelling and gaining long-term control and selfmanagement by the patient
- Specific kits developed for each body part for treatment for most upper and lower extremity patients
- Kits developed for specific body parts upper and lower extremities, torso and head & neck kits

Maximize clinic time with patients:

- Spend more clinic time on patient education
- Extend time between treatments as patients perform more self treatment
- Save patient out-of-pocket costs due to more production in decongestive phase




circaid reduction kit for the patient

Early results drives motivation

When the reduction kit is properly fit and patient is trained, results can be accelerated. Patients can better comprehend the objective compression settings for better therapy execution.

Progress is easily seen and felt as the reduction kit does its work!



Patient empowerment through engagement

The reduction kit design with juxtaposed bands is easy for patients to learn and perform. As patients become part of the solution there is a greater sense of empowerment.



Applying FLEC foam under reduction kit nighttime vest



Self donning reduction kit nighttime vest



Self donning reduction kit head and neck

Compliance unlike ever before

"Compression needs to be part of your life, not run your life" - Sandra Shaw circaid co-founder.

The reduction kit's inelastic wraps allow patients to take that occasional break, reduce the tension while in bed, or even take time for personal hygiene. Conversely, tension can be added to areas that have been decongested or to prevent product migration.

Setting the stage for long-term management

The reduction kit concept was developed with the end goal of creating a long-term, productive lifestyle for the patient.

The reduction kit moves through all stages - decongestion, transition and maintenence - with the end goal being patient independence.



Repeatable application with the Built-in-Tension system



Self donning on upper extremity



reduction kit glove is a trim-to-fit individualized option

circaid reduction kit for lower extremity

Lower extremity lymphedema can be the most challenging and requires some creativity for successful treatment. The circaid reduction kit is composed of multiple components that allow therapists to implement the pieces they choose for the best overall treatment of the patient.


All lower extremity products are comprised kits that includes the compression garment, the necessary assembly pieces, liner, Built-In-Tension System and instructions.

BOTTOM

50

TOP

circaid reduction kit for lower extremity





reduction kit lobe

circaid reduction kit - upper leg

			•	
Circum in	ference cm	regular	wide	
C	G	< 80	< 125	
cE	2	< 65	< 110	
Length	in cm	short	standard	long
<i>ℓ</i> D1	- G	< 37	37 -41	> 41
ltem	reduction kit description			
CRK2H001	upper leg regular width short length 30cm			
CRK2S001	upper leg regular width standard length 35cm			
CRK2L001	upper leg regular width long length 40cm			
CRK2H011	upper leg wide width short length 30cm			
CRK2S011	upper leg wide width standard length 35cm			
CRK2L011	upper leg wide width long length 40cm			

circaid reduction kit - knee

Circumference in cm		regular
	cE2	< 110
	cE	< 102
	cC	< 80
(one length available: standard 30cm	
ltem	reduction kit description	
CRK3S001	knee regular width standard length 30cm	

circaid reduction kit - lower leg and lobe

	Circumference in cm		wide	
cC		< 65	< 80	
сВ		< 55	< 65	
Length in	ст	short	standard	long
ℓ mall - l	D1	< 32	32-36	> 36
ltem	reduction kit description			
CRK1H001	lower leg regular width short length 30cm			
CRK1S001	lower leg regular width standard length 35cm			
CRK1L001	lower leg regular width long length 40cm			
CRK1H011	lower leg wide width short length 30cm			
CRK1S011	lower leg wide width standard length 35cm			
CRK1L011	lower leg wide width long length 40cm			
CRA7S001	lobe			


circaid reduction kit - whole leg

ltem	circaid reduction kit description
CRWHL011	whole leg kit; upper leg short length and wide width, lower leg long length and wide width
CRWHL031	whole leg kit; upper leg short length and wide width, lower leg long length and regular width
CRWHL001	whole leg kit; upper leg short length and regular width, lower leg long length and regular width
CRWHS011	whole leg kit; upper leg short length and wide width, lower leg standard length and wide width
CRWHS031	whole leg kit; upper leg short length and wide width, lower leg standard length and regular width
CRWHS001	whole leg kit; upper leg short length and regular width, lower leg standard length and regular width
CRWLL011	whole leg kit; upper leg long length and wide width, lower leg long length and wide width
CRWLL031	whole leg kit; upper leg long length and wide width, lower leg long length and regular width
CRWLS011	whole leg kit; upper leg long length and wide width, lower leg standard length and wide width
CRWLS031	whole leg kit; upper leg long length and wide width, lower leg standard length and regular width
CRWSH001	whole leg kit; upper leg standard length and regular width, lower leg short length and regular width
CRWSH031	whole leg kit; upper leg standard length and wide width, lower leg short length and regular width
CRWSL001	whole leg kit; upper leg standard length and regular width, lower leg long length and regular width
CRWSL011	whole leg kit; upper leg standard length and wide width, lower leg long length and wide width
CRWSL031	whole leg kit; upper leg standard length and wide width, lower leg long length and regular width
CRWSS001	whole leg kit; upper leg standard length and regular width, lower leg standard length and regular width
CRWSS011	whole leg kit; upper leg standard length and wide width, lower leg standard length and wide width
CRWSS031	whole leg kit; upper leg standard length and wide width, lower leg standard length and regular width
CRWHH001	whole leg kit; upper leg short length regular width, lower leg short length and regular width
CRWHH031	whole leg kit; upper leg short length and wide width, lower leg short length and regular width



circaid reduction kit whole leg

	circaid reduction kit - lower leg								
	Circumference in cm		regular				wide		
Step 1: Select	сC			< 65				< 80	
lower leg width	сВ			< 55				< 65	
Step 2: Select	Length in cm	shor	t sta	andard	Ion	g	standar	d lo	ng
lower leg length	ℓ mall - D1	< 32	3	2 -36	> 3	6	≤ 36	>	36
	circaid reduct	ion kit -	upper leg	7					
	Circumference in cm	reg	ular		wide			wide	
Step 3: Select	cG	<	80		< 125			< 125	
upper leg width	cE2	<	65		< 110			< 110	
Step 4: Select	Length in cm	short	standard	short	standard	long	short	standard	long
upper leg length	ℓ D1 -Groin	< 37	≥37	< 37	37 - 41	> 41	< 37	37 - 41	> 41




foot options for use with reduction kit



circaid pac band



circaid customizable interlocking ankle foot wrap



circaid reduction kit foot options						
Circumfere	ence			pac band®	single band	
in cm		interlocking afw	standard	large	x-large	
cA1		<44	19 - 36	22 - 41	<44	>44
сВ		<56	19 - 27	25 - 50	<56	>56
ltem	Description					
CPB00013	circa	id pac band xl				
CA001111	circa	id compression a	anklet - pair (star	ndard)		
CA001112	circaid compression anklet - pair (large)					
CFW5S003	S003 circaid customizable interlocking afw					
CFW2S001 circaid single band ankle foot wrap 3" wide						




circaid reduction kit toe cap



circaid reduction kit - toe cap

small	medium	large	X-large	XX-large
20-22	22-24	24-26	26-28	28-30
ltem	circa	id reductio	n kit descrij	ption
CRA1S001	circaid reducti	on kit toe cap,	beige, small	
CRA1S002	circaid reducti	on kit toe cap,	beige, medium	
CRA1S003	circaid reducti	on kit toe cap,	beige, large	
CRA1S004	circaid reduction kit toe cap, beige, x-large			
CRA1S005	circaid reduction kit toe cap, beige, xx-large			
CRA1S501	circaid reduction kit toe cap, black, small			
CRA1S502	circaid reduction kit toe cap, black, medium			
CRA1S503	circaid reduction kit toe cap, black, large			
CRA1S504	circaid reducti	on kit toe cap,	black, x-large	
CRA1S505	circaid reducti	on kit toe cap,	black, xx-large	

circaid reduction kit solutions for upper extremity





circaid reduction kit glove



circaid reduction kit - arm

	ference cm	regular	wide	x-wide
C	F	F < 50 50-65 65-80		
C	В	< 30	< 30	< 30
Length	h in cm	standard	long	
<i>e</i> C	- G < 45		> 45	
ltem	circaid reduction kit description			n
CRK4S001	arm regular width standard length 45cm			
CRK4L001	arm regular width long length 55cm			
CRK4S011	arm wide width standard length 45cm			
CRK4L011	arm wide width long length 55cm			
CRK4S021	arm x-wide width standard length 45cm			
CRK4L021	arm x-wide width long length 55cm			
CRA6S001	adjustable hand wrap			



circaid reduction kitarm and hand circaid

adjustable hand wrap

circaid reduction kit - hand

	circumferences in cm.					
small	medium	large	X-large	XX-large		
18-21	20-23	22-25	24-27	26-29		
ltem	circ	aid reductio	on kit descrip	otion		
CRA2S001	circaid reduc	tion kit glove, b	eige, small			
CRA2S002	circaid reduction kit glove, beige, medium					
CRA2S003	circaid reduction kit glove, beige, large					
CRA2S004	circaid reduction kit glove, beige, x-large					
CRA2S005	circaid reduc	circaid reduction kit glove, beige, xx-large				
CRA2S501	circaid reduc	circaid reduction kit glove, black, small				
CRA2S502	circaid reduction kit glove, black, medium					
CRA2S503	circaid reduction kit glove, black, large					
CRA2S504	circaid reduc	circaid reduction kit glove, black, x-large				
CRA2S505	circaid reduction kit glove, black, xx-large					


circaid reduction vest - nighttime



circaid reduction kit vest -nighttime

	st Width Axilla contour)	Chest cm.	Size	
	≤ 3 9	70 - 120	Small	
	35-39	121 - 180	Large	
	> 39	80 - 180	Large	
Item	circaid reduction kit description			
CRK5S001	circaid reduction kit vest standard small			
CRK5S003	circaid reduction kit vest standard large			
CRK5L001	circaid reduction kit vest long small			
CRK5L003	circaid reduction kit vest long large			

Axilla-Waist Length	Length
20 - 33	Standard
34 -47	Long





circaid reduction head and neck

Contraction of the second

circaid reduction kit - head & neck	
-------------------------------------	--

ltem	description
CRK6S001	reduction kit head & neck, one size

circaid trimmable foam

circaid - foam pads

ltem	description	
CRA8ST02	circaid textured edema control (TEC) foam pad 9"x15", 2 pack	
CRA8SF02	02 circaid flat localized edema control (FLEC) foam pad 16"x16", 2 pac	



FLEC pad 2-pack of FLEC pad, flat localized edema control



NEW

circaid reduction kit · 13



Compression solutions to meet the need

Progress can be tracked with each visit and, as decongestion occurs, the components can easily be adjusted. Patients can continually tweak the tension as the limb reduces.

circaid reduction kit – the technology inside



Breathe-O-Prene® Breathable fabric that improves comfort and compliance



Inelastic

Giraffe skin is inelastic – this is the core technology of all circaid products



Variety

A wide variety of reduction kit components enables trim-to-fit options on a per-patient basis. Choose from leg, arm, hand, foot and lobule solutions.



Easy handling

Easy handling when applying and removing the components fosters patient engagement and independence not possible before



Built-In-Tension System

BTS supports a repeatable level of tension when the patient applies the reduction kit away from the clinic



Juxtaposed band system Ensures an individual fit every time, even after decongestion starts



Juxtalock band Juxtaposed bands for enhanced fit and adjustablilty

CIRCAID circaid reduction kit Case Study¹

Clinicians and patients rate as "Excellent"

A case study conducted through eight facilities in the US utilized the reduction kit on 49 patients to measure the therapeutic and product results.

• The average total volume reduction for all cases in the trial was 10.8%.

- 70% of patients describe the reduction kit as "Excellent"
- 67% of clinicians gave the reduction kit an "Excellent" rating for its ability to control edema.

Parameter	Mean Score - Clinician	Mean Score - Patient
Ease of initial application	1.9 (Excellent)	N/A
Application time	2.1 (Excellent)	N/A
Fit	2.1 (Excellent)	2.1 (Excellent)
Comfort	N/A	2.7 (Good)
Ease of use	N/A	2.3 (Excellent)
Reduction of edema	2.0 (Excellent)	2.4 (Excellent)
Patient hygiene	1.7 (Excellent)	1.7 (Excellent)
Compliance	2.3 (Excellent)	N/A
Overall Quality	1.7 (Excellent)	2.2 (Excellent)

Clinicians and patients rated their experience with the reduction kit on a score from 1 -excellent to 6 -poor.

¹circaid[®] reduction kit[™] Case Study 2014. Use of circaid reduction kit for decongestive lymphatic therapy.

Learn more!





Please join us @mediusaifeelbetter



Please join @mediUSA



Please join us @mediusaifeelbetter

medi USA 6481 Franz Warner Parkway Whitsett, NC 27377 T 800-633-6334 F 888-570-4554 info@mediusa.com mediusa.com