




PROCTOLOGY INSTRUMENTS



FEATURES

- Made from stainless steel which makes them rust proof, smooth and reliable.
- Rounded surface facilitate easy to clean and disinfect with conventional solutions
- Angled tube-end, ensures a good view and easy access with instruments.
- Rounded proximal shield ensures Easy insertion and rotation.
- Opening viewing window with insufflations port.


PROCTOLOGY INSTRUMENTS



FIBER OPTIC ANOSCOPES WITH ANNULAR F.O. ILLUMINATION

FEATURES:


- Annular F.O. illumination enables brighter illumination of the area to be examined.
- Detachable tip enables easier cleaning and disinfection.
- Comes with ACMI and WOLF standard connections, STORZ connector supplied separately.
- Illumination by Fiber Optic cable and projector, Lamp handle or hand held power source.
- Air tight magnifying glass window.
- Insufflations adaptor for pneumatic test.
- 2.0x Swivel lens attachment for performing small procedures.
- Graduated in Centimeter.



Diameter Ø	Working length	Complete with Obturator	Anoscope without Obturator	Obturator Only
20 mm	130 mm	819-088-30	819-088-01	819-088-02

FEATURES:

- Annular F.O. illumination eliminates the bulb in center, thereby increasing the viewing area.
- Bigger Fiber bundle diameter enables brighter illumination of the area to be examined.
- Supplied with ACMI and WOLF connectors as standard STORZ connector supplied separately.
- Illumination by Fiber Optic cable and projector, lamp handle or hand held power source.
- Tubes can be dismantled from the light chamber for cleaning and disinfection purpose.
- Tubes can also be replaced.



Diameter Ø	Working length	Complete with Obturator	Anoscope without Obturator	Obturator Only
15 mm	70mm	819-080-15	819-080-01	819-080-02
20 mm	70mm	819-082-20	819-082-01	819-082-02
25 mm	70mm	819-084-25	819-084-01	819-084-02

SPARES AND ACCESSORIES:

- End Tip (819-088-05)
- 2.0x Swivel Lens (819-099-20)



PROCTOLOGY INSTRUMENTS

FIBER OPTIC SIGMOIDOSCOPE WITH DISTAL ANNULAR ILLUMINATION

FEATURES:

- Stainless steel construction, making it easier to clean and disinfect with conventional solution.
- Distal Annular F.O. illumination enables brighter illumination of the area to be examined.
- Detachable tip.
- Illumination by Fiber Optic cable and projector, Lamp handle or hand held power source.
- Comes with ACMI and WOLF standard Connection, STORZ connector supplied separately.
- Air tight Glass window.
- Insufflations adaptor for pneumatic test.
- 1.5x swivel lens attachment for performing small procedures.
- Graduated in centimeters.



Instruments	Diameter Ø	Working length	Complete	Sigmoidoscope Only	Obturator Only	End Tip
Sigmoidoscope Adult	21 mm	300 mm	819-090-21	819-090-01	819-090-02	819-090-03
Sigmoidoscope Adult	21 mm	250 mm	819-092-21	819-092-01	819-092-02	819-092-03
Sigmoidoscope Adult	21 mm	200 mm	819-094-21	819-094-01	819-094-02	819-094-03
Sigmoidoscope Adolescent	16 mm	250 mm	819-096-16	819-096-01	819-096-02	819-096-03
Sigmoidoscope Child	12 mm	200 mm	819-098-12	819-098-01	819-098-02	819-098-03

SPARES AND ACCESSORIES:

- 1.5x Magnification Swivel Lens (819-099-15)


PROCTOLOGY INSTRUMENTS



DISPOSABLE TUBES

FEATURES:

- Single use
- Made from patient-friendly material
- Light transmission through tube walls.
- Graduated in centimeters.



F.O. INSTRUMENT HEAD

FEATURES:

- Single Instrument Head for all types of Disposable Tubes, Compatible to all manufacturers.
- Stainless steel construction making it easier to clean and disinfect with conventional solutions.
- Annular F.O. illumination enables brighter illumination.
- Illumination by Fiber Optic cable and projector, lamp handle or hand held power source.
- Comes with ACMI and WOLF standard connections, STORZ connector supplied separately.
- Air tight glass window.
- Insufflations adaptor for pneumatic test.
- 1.5x swivel lens for Sigmoidoscope and 2.0x swivel lens for Proctoscope and Anoscope.



Instruments	Diameter Ø	Working length	Complete	Pack of 25	Box of 100
F.O. Head Complete with swivel lens	-	-	819-103-01	-	-
Disposable Anoscope	20 mm	100 mm	-	819-102-25	819-102-50
Disposable Anoscope	20 mm	130 mm	-	819-104-25	819-104-50
Disposable Anoscope	20 mm	250 mm	-	819-106-25	819-106-50

SPARES AND ACCESSORIES:

- 1.5x Magnification swivel lens (819-099-15)
- 2.0x Magnification swivel lens (819-099-20)




PROCTOLOGY INSTRUMENTS

LLOYD DAVIES F.O. SIGMOIDOSCOPIES & PROCTOSCOPES

FEATURES:

- Stainless steel construction, making it easier to clean and disinfect with conventional solutions.
- Single light chamber for all sizes of sigmoidoscopes.
- Comes with ACMI and WOLF standard connections, STORZ connector supplied separately.
- Illumination by Fiber Optic cable and projector, lamp handle or hand held power source.



Instruments	Diameter Ø	Working length	Complete	Proctoscope Only	Obturator Only	Bulbs/Light Guide
Sigmoidoscope Adult	20 mm	300 mm	819-100-20	819-100-01	819-100-02	-
Sigmoidoscope Adult	16 mm	250 mm	819-102-16	819-102-01	819-102-02	-
Proctoscope Child	10 mm	75 mm	819-104-10	819-104-01	819-104-02	-
F.O. Light Chamber	-	-	819-108-00	-	-	819-108-03
Standard Light Chamber	-	-	819-110-00	-	-	833-009-25
Connecting Electrical cable	-	-	819-110-01	-	-	-
Replacement bulb	-	-	833-045-25	-	-	-

SPARES AND ACCESSORIES:

- Glass window (819-108-01)
- Chamber Only (819-108-02)
- Spare light guide (819-108-03)

PROCTOLOGY INSTRUMENTS

STRAUSS RECTOSCOPE



Instruments	Diameter Ø	Working length	Complete with Obturator	Proctoscope Only	Obturator Only	Light Guide
Strauss Recto scope	20 mm	300 mm	819-156-30	819-156-01	819-156-02	819-156-03
Strauss Recto scope	20 mm	250 mm	819-158-25	819-158-01	819-158-02	819-158-03
Strauss Recto scope	20 mm	200 mm	819-160-20	819-160-01	819-160-02	819-160-03
Handle	-	-	819-152-01	-	-	-
Main chamber	-	-	819-150-01	-	-	-
Glass window	-	-	819-154-01	-	-	-
Telescope	-	-	819-162-01	-	-	-


NAUTON MORGAN PROCTOSCOPE



Nauton Morgan Proctoscope	Diameter Ø	Working length	Complete with Obturator	Proctoscope Only	Obturator	Light Attachment	Handle
F.O. Light	17 mm	65 mm	819-112-17	819-112-01	819-112-02	819-140-01	-
F.O. Light	21 mm	65 mm	819-114-21	819-114-01	819-114-02	819-140-01	-
F.O. Light	24 mm	65 mm	819-116-24	819-116-01	819-116-02	819-140-01	-
F.O. Light	17 mm	80 mm	819-118-17	819-118-01	819-118-02	819-140-01	-
F.O. Light	21 mm	80 mm	819-120-21	819-120-01	819-120-02	819-140-01	-
F.O. Light	24 mm	80 mm	819-122-24	819-122-01	819-122-02	819-140-01	-
Standard Light	17 mm	65 mm	819-124-17	819-112-01	819-112-02	819-142-01	808-003-25
Standard Light	21 mm	65 mm	819-126-21	819-114-01	819-114-02	819-142-01	808-003-25
Standard Light	24 mm	65 mm	819-128-24	819-116-01	819-116-02	819-142-01	808-003-25
Standard Light	17 mm	80 mm	819-130-17	819-118-01	819-118-02	819-142-01	808-003-25
Standard Light	21 mm	80 mm	819-132-21	819-120-01	819-120-02	819-142-01	808-003-25
Standard Light	24 mm	80 mm	819-134-24	819-122-01	819-122-02	819-142-01	808-003-25

PROCTOLOGY INSTRUMENTS

PROCTOSCOPE & ANOSCOPE



Instruments	Diameter Ø	Working length	Complete with Obturator
Kelly Proctoscope, Taper	21 mm	50 mm	819-170-21
Kelly Proctoscope, Taper	27 mm	50 mm	819-172-27
Kelly Proctoscope, Straight	22 mm	50 mm	819-180-02
Kelly Proctoscope, Straight	22 mm	100 mm	819-182-04
Kelly Proctoscope, Straight	22 mm	140 mm	819-184-06
Kelly Proctoscope, Straight	22 mm	200 mm	819-186-08
Hirschman Ana scope	16 mm	65 mm	819-190-16
Hirschman Ana scope	19 mm	65 mm	819-192-19
Hirschman Ana scope	22 mm	65 mm	819-194-22

ANOSCOPE & SPECULUM



Instruments	Diameter Ø	Working length	Complete with Obturator
Ives Fansler Ano scope, Small Handle	22 mm	80 mm	819-200-22
Ives Fansler Ano scope, Large Handle	22 mm	80 mm	819-202-22
Chelsea-Eaton Anal Speculum	22 mm	65 mm	819-210-22
Chelsea-Eaton Anal Speculum	30 mm	65 mm	819-212-30
Bensaude Ano scope	19 mm	70 mm	819-220-19
Bensaude Ano scope	22 mm	70 mm	819-222-22

PROCTOLOGY INSTRUMENTS



ACCESSORIES OF PROCTOSCOPY



BIOPSY FORCEPS

FEATURES:

- Stainless steel construction making it durable, easier to clean and disinfect with conventional solutions.
- Precisely finished smooth cutting edges.
- Ergonomically formed handle which is comfortable to hold.
- Single handle for all types of forceps.
- Shaft rotatable through 360°.



SUCTION TUBE

FEATURES:

- Stainless steel construction making it durable, easier to clean and disinfect with conventional solutions.
- Equipped with a clip for attachment to the chamber.



SPONGE HOLDING FORCEPS

FEATURES:

- Stainless steel construction making it durable, easier to clean and disinfect with conventional solutions.



TWIN BELLOW INSUFFLATOR

Instruments	Working length	Angle	Complete	Shaft only	Handle only
Biopsy forceps, 3.5x8mm Cup Jaws	420 mm	Straight	819-600-42	819-600-01	819-600-00
Biopsy forceps, 3.5x8mm Cup Jaws	350 mm	Straight	819-602-35	819-602-01	819-600-00
Biopsy forceps, 3.5x8mm Cup Jaws	420 mm	30°	819-604-42	819-604-01	819-600-00
Biopsy forceps, 3.5x8mm Cup Jaws	350 mm	30°	819-606-35	819-606-01	819-600-00
Biopsy forceps, 4mm cutting Jaws	420 mm	Straight	819-610-42	819-610-01	819-600-00
Biopsy forceps, 4mm cutting Jaws	350 mm	Straight	819-612-35	819-612-01	819-600-00
Biopsy forceps, 4mm cutting Jaws	420 mm	30°	819-614-42	819-614-01	819-600-00
Biopsy forceps, 4mm cutting Jaws	350 mm	30°	819-616-35	819-616-01	819-600-00
Handle only for biopsy forceps	-	-	819-600-00	-	-
Sponge holding forceps	300 mm	-	819-620-30	-	-
	400 mm	-	819-610-40	-	-
Suction tube	400 mm	-	819-620-40	-	-
Twin bellow insufflator	-	-	819-630-00	-	-



AMINO SCOPES

FEATURES:

- Stainless steel construction making it easier to clean and disinfect with conventional solutions.
- Illumination by fiber optic cable and projector, lamp handle or hand held power source.
- Matt inner surface eliminates internal reflection.
- Graduated in centimeters.
- F.O Light clamp can be used with all 3 sizes of tubes.



Instruments	Distal tip Ø	Working length	Complete	Tube only	Obturator only
F.O. Clamp	-	-	822-001-01	-	-
Amino Scope Tube	12 mm	200 mm	822-005-12	822-005-01	822-005-02
Amino Scope Tube	16 mm	200 mm	822-007-16	822-007-01	822-007-02
Amino Scope Tube	20 mm	200 mm	822-009-20	822-009-01	822-009-02

PROCTOLOGY KITS



ALPHA I ILLUMINATION KIT FOR DISPOSABLE ANOSCOPES AND PROCTOSCOPES



Description	Ref
Complete Kit	819-500-01
Mark III Lamp Handle with 6.0V bulb	808-025-60
Plug-in Transformer	808-071-60
Illumination Head for Disposable Proctoscope	819-103-01
Bulb for lamp handle, lens end, 6.0V halogen	833-049-60

ALPHA II ILLUMINATION KIT FOR DISPOSABLE SIGMOIDOSCOPE




Description	Ref
Complete Kit	819-502-01
Mark III Lamp handle with 6.0V bulb	808-025-60
Illumination head for Disposable Proctoscope	819-103-01
Plug in Transformer	808-071-60
Twin bellow insufflators	819-630-00
Sponge holding forcep	819-610-40
Bulb for lamp handle, 6.0V halogen	833-049-60
Rexin case for Sigmoidoscope Set	832-204-01



PROCTOLOGY KITS

ALPHA III F.O ILLUMINATION KIT FOR DISPOSABLE SIGMOIDOSCOPE



Description	Ref
Complete Kit	819-504-01
Mark III add-on F.O Handle	808-027-00
MD F.O Cable instruments End adaptor for WOLF fitting	830-002-02
MD F.O Cable for Halogen Light Source, Bundle dia 4.8mm, length 1.8mm	830-048-18
Illumination Head for Disposable Proctoscope	819-103-01
Twin bellow Insufflators	819-630-00
Sponge Holding Forcep	819-610-40
Rexin Case for Sigmoidoscope Set	832-204-01

ALPHA IV SIGMOIDOSCOPE & PROCTOSCOPE KIT WITH F.O ILLUMINATION



Description	Ref
Complete Kit	819-506-01
Mark III Add-on F.O handle	808-027-00
Biopsy forcep 4.0mm cutting Jaws 30°	819-614-42
MD F.O Cable instruments End adaptor for WOLF fitting	830-002-02
MD F.O Cable for Halogen Light Source Bundle dia 4.8mm, length 1.8mm	830-048-18
Twin bellow Insufflators	819-620-00
Suction tube for Sigmoidoscope	819-620-40
Proctoscope with Distal Annular F.O illumination, 130x21mm Complete with Obturator and Lens	819-088-30
Sigmoido scope with Distal Annular F.O illumination, 21x300mm Complete with Obturator and Lens	819-090-21
Sponge Holding forcep, with spring handle	819-620-40
Rexin Case for Sigmoidoscope Set	832-204-01

PROCTOLOGY KITS



ALPHA V SIGMOIDOSCOPE & PROCTOSCOPE KIT WITH F.O ILLUMINATION



Description	Ref
Complete Kit	819-508-01
Mark III Lamp Handle with 6.0V bulb	808-025-60
Plug-in Transformer	808-071-60
Biopsy Forceps 4.0mm cutting Jaws, 30°	819-614-42
Twin bellow insufflators	819-630-00
Proctoscope with Distal Annular F.O illumination, 130x21mm Complete with Obturator and Lens	819-088-30
Sigmoido scope with Distal Annular F.O illumination 21x300mm Complete with Obturator and Lens	819-090-21
Sponge Holding forcep, with spring handle	819-610-40
Suction tube for Sigmoido Scope	819-620-40
Rexin Case for Sigmoidoscope Set	832-204-01

ALPHA VI ADULT SIGMOIDOSCOPE KIT WITH FIBER OPTIC ILLUMINATION



Description	Ref
Complete Kit	819-516-01
Mark III Add-on F.O Handle	808-027-00
Twin bellow insufflator	819-630-00
MD F.O Cable Instruments End Adaptor for WOLF/HSW/Dyonics/Midicon Fitting	830-002-02
MD F.O Cable for Halogen Light Source, bundle Dia 4.8mm, Length 1.8mm	830-048-18
Biopsy Forceps 4.0mm cutting Jaws, 30°	819-614-42
Sigmoidoscope with Distal Annular F.O illumination, 16x250mm Complete with Obturator and Lens	819-096-16
Sigmoidoscope with Distal Annular F.O illumination, 21x300mm Complete with Obturator and Lens	819-090-21
Sponge Holding forcep, with spring handle	819-610-40
Suction tube for Sigmoido Scope	819-620-40
Rexin Case for Sigmoidoscope Set	832-204-01



PROCTOLOGY KITS

ALPHA VII ADULT SIGMOIDOSCOPE KIT



Description	Ref
Complete Kit	819-518-01
Mark III lamp Handle with 6.0V bulb	808-025-60
Plug-in Transformer	808-071-60
Twin bellow insufflators	819-630-00
Biopsy Forceps 4.0mm cutting Jaws, 30°	819-614-42
Sigmoidoscope with Distal Annular FO illumination, 16x250mm Complete with Obturator and Lens	819-096-16
Sigmoidoscope with Distal Annular FO illumination, 21x300mm Complete with Obturator and Lens	819-090-21
Sponge Holding forcep, with spring handle	819-610-40
Suction tube for SigmoidoScope	819-620-40
Rexin Case for SigmoidoScope Set	832-204-01

ALPHA VIII CHILD SIGMOIDOSCOPE KIT



Description	Ref
Complete Kit	819-520-01
Mark III Add-on FO Handle	808-027-00
Twin bellow insufflators	819-630-00
MD FO Cable Instruments End Adaptor for WOLF Fitting	830-002-02
MD FO Cable for Halogen Light Source, bundle Dia 4.8mm, Length 1.8mm	830-048-18
Sigmoidoscope with Distal Annular FO illumination, 12x200mm Complete with Obturator and Lens	819-098-12
Sponge Holding forcep, with spring handle	819-610-40
Suction tube for SigmoidoScope	819-620-40
Rexin Case for SigmoidoScope Set	832-204-01

PROCTOLOGY KITS



ALPHA XI CHILD SIGMOIDOSCOPE KIT



Description	Ref
Complete Kit	819-522-01
Mark III lamp Handle with 6.0V bulb	808-025-60
Plug-in Transformer	808-071-60
Twin bellow insufflators	819-630-00
Sigmoidoscope with Distal Annular F.O illumination, 12x200mm Complete with Obturator and Lens	819-098-12
Sponge Holding forcep, with spring handle	819-610-40
Suction tube for SigmoidoScope	819-620-40
Rexin Case for SigmoidoScope Set	832-204-01

BETA VI ANOSCOPE KIT



Description	Ref
Complete Kit	819-524-01
Mark III Add-on FO Handle	808-027-00
MD FO Cable Instruments End Adaptor for WOLF/HSW/Dyonics/Midicon Fitting	830-002-02
MD FO Cable for Halogen Light Source, bundle Dia. 4.8mm, Length 1.8mm	830-048-18
FO Anoscope with Annular illumination, 70x15mm, complete with Obturator	819-080-15
FO Anoscope with Annular illumination, 70x20mm, complete with Obturator	819-082-20
FO Anoscope with Annular illumination, 70x25mm, complete with Obturator	819-084-25
Plastic Case XL for Anoscope Set, Black color	832-206-01

BETA VII ANOSCOPE KIT



Description	Ref
Complete Kit	819-526-01
Mark III lamp Handle with 6.0V bulb	808-025-60
Plug-in Transformer	808-071-60
F.O Anoscope with Annular illumination, 70x15mm, complete with Obturator	819-080-15
F.O Anoscope with Annular illumination, 70x20mm, complete with Obturator	819-082-20
F.O Anoscope with Annular illumination, 70x25mm, complete with Obturator	819-084-25
Plastic Case XL for Anoscope Set, Black color	832-206-01

LLOYD DAVIES SIGMOIDOSCOPE & PROTOSCOPE KIT



Description	Ref
Complete Kit	819-100-90
Lloyd Davies Sigmoidoscope Adult W/O Light, 20.0mm Diameter	819-100-20
Lloyd Davies Sigmoidoscope Adolescent W/O light, 16.0mm Diameter	819-102-16
Lloyd Davies Proctoscope, 10.0mm Diameter	819-104-10
Lloyd Davies Sigmoidoscope Chamber with F.O light guide + window	819-108-00
MD F.O Cable Projector End Adaptor for WOLF Fitting	830-002-01
MD F.O Cable Instruments End Adaptor for WOLF Fitting	830-002-02
MD F.O Cable for Halogen Light Source, bundle Dia 4.8mm, Length 1.8mm	830-048-18
Twin bellow insufflators	819-630-00
Nauton Morgan proctoscope with F.O Light, 21x65mm	819-114-21
Rexin Case for SigmoidoScope Set	832-204-01

PROCTOLOGY KITS



LLOYD DAVIES SIGMOIDOSCOPE KIT



Description	Ref
Complete Kit	819-100-91
Lloyd Davies Sigmoidoscope Adult W/O Light, 20.0mm Diameter	819-100-20
Lloyd Davies Sigmoidoscope Adolescent W/O light, 16.0mm Diameter	819-102-16
Lloyd Davies Sigmoidoscope Chamber with FO light guide + window	819-108-00
MD FO Cable Projector End Adaptor for WOLF Fitting	830-002-01
MD FO Cable Instruments End Adaptor for WOLF Fitting	830-002-02
MD FO Cable for Halogen Light Source, bundle Dia 4.8mm, Length 1.8mm	830-048-18
Twin bellow insufflators	819-630-00
Rexin Case for Sigmoidoscope Set	832-204-01

LLOYD DAVIES SIGMOIDOSCOPE & PROCTOSCOPE KIT



Description	Ref
Complete Kit	819-110-90
Lloyd Davies Sigmoidoscope Adult W/O Light, 20.0mm Diameter	819-100-20
Lloyd Davies Sigmoidoscope Adolescent W/O light, 16.0mm Diameter	819-102-16
Lloyd Davies Proctoscope, 10.0mm Diameter	819-104-10
Lloyd Davies Sigmoidoscope Chamber with 2.5V Bulb + window	819-110-00
Twin bellow insufflators	819-630-00
Nauton Morgan proctoscope with Standard Light, 21x65mm	819-126-21
Electrical cord for conventional light Sigmoidoscopes / Proctoscopes	819-110-01
Battery handle Mark I-SV, 2.5V "C" Type	808-001-25
Rexin Case for Sigmoidoscope Set	832-204-01



PROCTOLOGY KITS

LLOYD DAVIES SIGMOIDOSCOPE KIT



Description	Ref
Complete Kit	819-110-91
Lloyd Davies Sigmoidoscope Adult W/O Light, 20.0mm Diameter	819-100-20
Lloyd Davies Sigmoidoscope Adolescent W/O light, 16.0mm Diameter	819-102-16
Lloyd Davies Sigmoidoscope Chamber with 2.5V Bulb + window	819-110-00
Electrical cord for conventional light Sigmoidoscopes / Proctoscopes	819-110-01
Battery handle Mark I-SV, 2.5V "C" Type	808-003-25
Twin bellow insufflators	819-630-00
Rexin Case for Sigmoidoscope Set	832-204-01