

# What are OCEAN Clear Aligners?

OCEAN Clear Aligners appliance involves a series of aligners made from a transparent, thin (typically less than 1mm) plastic material called Essix and formed with CAD-CAM laboratory techniques. These aligners are similar to the splints that cover the clinical crowns and the marginal gingiva. Each aligner is designed to move the teeth a maximum of about 0.25 to 0.3mm over a -2week period and is worn in a specific sequence.



# When and how to use OCEAN Clear Aligners?

OCEAN Clear Aligner appliance is currently recommended for adults and for adolescents with fully erupted permanent teeth who meet an acceptable standard of compliance. It can treat all varieties of cases, including extraction cases, cross bites, rotations, edge to edge, bi max etc. Excellent compliance is mandatory since the appliance has to be worn a minimum of 20 to 22 hours a day and each aligner should be worn 400 hours to be effective



# Why to Choose OCEAN Clear Aligners?

#### 1. CLEAR SHOW OFF THE SMILE, NOT THE ALIGNER

When a patient decides to straighten their teeth, don't let them compromise their daily routine. OCEAN is the most attractive choice for orthodontic treatment because our aligners are crystal clear.

#### 2. PREDICTABLE GET AN ACCURATE AND EFFECTIVE TREATMENT PLAN


OCEAN uses highly advanced software to produce a 3D model of the treatment. This gives you and your patient an excellent example of how their teeth are going to shift. Your treatment plan with OCEAN is predictable and gives you the security which you and your patient would expect and get from the treatment.

#### 3. COMFORTABLE THE MOST PLEASANT WAY TO STRAIGHTEN TEETH


OCEAN aligners are very comfortable to wear. The material is flexible and very thin, designed to stop just before the gums. There is no risk of any emergency, cuts or sores.

#### 4. HIGHLY PRECISE AUTOMATED AND SOPHISTICATED PROCEDURE

The Aligner quality, clarity and finish are unmatchable and it is made by vacuum forming as well as positive pressure technique which makes it snug free.



### 5. **DON'T REQUIRE ATTACHMENT**



# Why to Choose OCEAN Clear Aligners?

#### 6. **REMOVABLE** CLEAN, FRESH, AND FLEXIBLE

Traditional braces may leave spots on the tooth enamel but with OCEAN, this does not happen. The patient can remove their aligners at any time without difficultiesalthough this should not be done too often.

#### 7. TIME SAVING GET THE TIME FOR MORE CASES AND PATIENTS

With OCEAN aligners only 10 minutes of chair side time are required per 6 weeks per patient.

# Requirements:

3D scanning (STL file)

- \*Photos\*
- 1) upper + lower occlusal
- 2) fronted intra oral
- 3) right side intra oral in occlusion
- 4) left side intra oral in occlusion
- 5) face fronted view clips closed + smile
- 6) face profile

OPG x-ray.

www.oceanaligners.com support@oceanaligners.com +962 7 9230 0999

