

*FEHLING...*

*... THE DIFFERENCE*

**PC**

SET FOR PEDIATRIC  
CARDIAC SURGERY

SET FÜR KINDER-  
HERZCHIRURGIE




1690 Stone Village Lane NW • Suite 721 • Kennesaw, GA 30152


Toll Free: 1-800-FEHLING (334-5464) • Fax: 770-794-0122 • Email: [info@fehlingurgical.com](mailto:info@fehlingurgical.com)

CERAMO®


ASD/VSD RETRACTOR  
HERZDEFEKTHAKEN




LANGENBECK



ZENKER



MOBERG



175 mm

BLADE SIZE: DEPTH (a) x WIDTH (b)  
BLATTGRÖSSE: TIEFE (a) x BREITE (b)




1/2

| WIDTH<br>BREITE |       |
|-----------------|-------|
| 4 mm            | EOI-6 |
| 8 mm            | EOI-7 |
| 10 mm           | EOI-8 |
| 14 mm           | EOI-4 |
| 16 mm           | EOI-5 |
| 19 mm           | EOI-9 |
| 24 mm           | EOK-0 |

CERAMO® SUPERPLAST  
SPATULA  
SPATEL  
180 mm

— A —



| SIZE<br>GRÖSSE<br>A x B |       |
|-------------------------|-------|
| 10 x 35 mm              | EOS-1 |
| 14 x 45 mm              | EOS-2 |

CERAMO® SUPERPLAST JAQUISS  
ORGAN RETRACTOR  
ORGANRETRAKTOR  
100 mm

**SUPERPLAST** instruments are made of shape-memory metal. This metal has unique qualities under different thermal conditions. The instruments are malleable in normal OR temperatures. Soft pressure allows "permanent" deformation of shape according to the surgeon's needs. However in higher temperature, such as those found in the autoclave, the instruments regain their initial shape.

**SUPERPLAST** Instrumente sind aus Metall mit Formgedächtnis. Sie bieten verschiedene Eigenschaften je nach thermischen Bedingungen. Sie sind bei normaler OP-Temperatur biegsam und können durch sanften Druck in die gewünschte Form gebracht werden. Unter höheren Temperaturen, wie bei der Sterilisation, nehmen sie ihre ursprüngliche Form wieder an.

blades (included)  
Blätter (inklusive)




10 x 20 mm (a x b)  
(PAIR / PAAR)

MFU-8



10 x 10 mm (a x b)  
(PAIR / PAAR)


MFU-9



**TITANIUM**

STERNAL RETRACTOR  
STERNUMSPERRER  
SPREADING WIDTH = 52 mm  
SPREIZWEITE = 52 mm


MFU-6



**TITANIUM**



MICRO RIB SPREADER  
MIKRORIPPENSPREIZER  
SPREADING WIDTH = 30 mm  
SPREIZWEITE = 30 mm  
STRAIGHT/GERADE

MBW-0




FINOCHIETTO  
MINI-RETRACTOR  
MINI-SPREIZER  
SPREADING WIDTH = 50 mm  
SPREIZWEITE = 50 mm  
CURVED/GEBOGEN

MBT-3



| BLADE SIZE<br>BLATTGRÖSSE<br>a x b | SPREADING WIDTH<br>SPREIZWEITE |       |
|------------------------------------|--------------------------------|-------|
| 14 x 12 mm                         | 40 mm                          | MSB-1 |
| 14 x 12 mm                         | 65 mm                          | MSB-2 |
| 20 x 15 mm                         | 90 mm                          | MSB-6 |
| 23 x 20 mm                         | 100 mm                         | MSB-3 |
| 20 x 15 mm                         | 110 mm                         | MSB-7 |
| 23 x 20 mm                         | 130 mm                         | MSB-4 |
| 23 x 20 mm                         | 160 mm                         | MSB-5 |




| BODY SIZE<br>RAHMENGRÖSSE<br>A / B | BLADE SIZE<br>BLATTGRÖSSE<br>a x b | PATIENT WEIGHT<br>KÖRPERGEWICHT |        |
|------------------------------------|------------------------------------|---------------------------------|--------|
| 60 / 100 mm                        | 12 x 25 mm                         | ≤ 5 kg                          | MSE-4  |
| 75 / 120 mm                        | 18 x 35 mm                         | 5 – 15 kg                       | MSE-5  |
| 100 / 140 mm                       | 21 x 55 mm                         | 15 – 30 kg                      | MSE-6  |
| 150 / 280 mm                       | 30 x 100 mm                        | 30 – 60 kg                      | MSE-9* |

\* Stainless Steel / Rostfreier Stahl




MINI THORACOTOMY RETRACTOR  
MINI-THORAKOTOMIESPREIZER  
115 mm  
MSD-5

V-CUT  
TITANIUM



*exemplary configuration*  
*Konfigurationsbeispiel*



RETRACTOR, BODY  
SPERRER, RAHMEN

| BLADE SIZE<br>BLATTGRÖSSE<br>A x B | PATIENT WEIGHT<br>KÖRPERGEWICHT |       |
|------------------------------------|---------------------------------|-------|
| 80 x 130 mm                        | 3 – 7 kg                        | MSX-1 |
| 90 x 150 mm                        | 7 – 15 kg                       | MSX-2 |
| 120 x 200 mm                       | 15 – 40 kg                      | MSX-3 |



$\frac{1}{2}$

STANDARD  
MSX-5B



$\frac{1}{2}$

RIDER  
SCHLITTEN

ANGLE ADJUSTABLE  
WINKELVERSTELLBAR  
MSX-5C



$\frac{1}{2}$

BLADE GUIDE  
HAKENFÜHRUNG

MSX-5A




$\frac{1}{2}$

HORKE BLADE SWING  
HORKE SCHAUKEL

| HOLE DISTANCE<br>LOCHABSTAND |       |
|------------------------------|-------|
| 20 mm                        | MTH-2 |
| 25 mm                        | MTH-0 |
| 35 mm                        | MTH-1 |




V-CUT  
TITANIUM




| a \ b | 15 mm | 20 mm |
|-------|-------|-------|
| 15 mm | MTH-3 | MTH-4 |
| 21 mm | MTH-5 | MTH-6 |
| 26 mm | MTH-7 | MTH-8 |

STERNAL BLADE  
STERNUMHALTEBLATT



| a \ b | 15 mm | 24 mm |
|-------|-------|-------|
| 40 mm | MSY-6 | MSZ-6 |
| 50 mm | MSY-7 | MSZ-7 |
| 60 mm | MSY-8 | MSZ-8 |
| 70 mm | MSY-9 | MSZ-9 |

DIAPHRAGM RETRACTOR  
ZWERCHFELLHAKEN



| a \ b | 10 mm |
|-------|-------|
| 20 mm | MSX-6 |
| 30 mm | MSX-7 |
| 40 mm | MSX-8 |

LID RETRACTOR  
LIDHAKEN





225 mm  
MCE-9



ANDREWS PYNCHON  
240 mm

| Ø      |       |
|--------|-------|
| 3,0 mm | MNO-4 |
| 4,2 mm | MNO-5 |



245 mm  
MCE-8




210 mm  
MNE-5



SEBENING  
180 mm  
MNE-8

DETACHABLE BASKET  
KORB ABNEHMBAR

GRECO




DETACHABLE BASKET / KORB ABNEHMBAR

| Ø<br>TUBE / ROHR | 180 mm | 160 mm |
|------------------|--------|--------|
| 3,3 mm           | MNM-1  | MNM-2  |
| 5,0 mm           | MNM-3  | MNM-4  |

DETACHABLE BASKET / KORB ABNEHMBAR

| Ø<br>TUBE / ROHR | 95 mm | 80 mm |
|------------------|-------|-------|
| 3,3 mm           | MNM-5 | MNM-6 |
| 5,0 mm           | MNM-7 | MNM-8 |



| Ø<br>TUBE / ROHR | 210 mm | 180 mm | 150 mm |
|------------------|--------|--------|--------|
| 3,3 mm           | MNL-2  | MNL-1  | MNL-3  |
| 5,0 mm           | MNL-5  | MNL-4  | MNL-6  |

GRECO



DETACHABLE BASKET / KORB ABNEHMBAR

|                  |        |
|------------------|--------|
| Ø<br>TUBE / ROHR | 230 mm |
| 6,0 mm           | MNL-7  |




DETACHABLE BASKET / KORB ABNEHMBAR

|                  |        |
|------------------|--------|
| Ø<br>TUBE / ROHR | 205 mm |
| 6,0 mm           | MNL-8  |



VENTRICLE SUCTION CANNULA  
VENTRIKELSAUGER



CARDIOPLEGIA CANNULA  
KARDIOPLEGIEKANÜLE  
30 mm  
MZI-3



| Ø<br>OLIVE |       |
|------------|-------|
| 3,0 mm     | MSF-0 |
| 4,0 mm     | MSF-1 |
| 4,5 mm     | MSF-2 |
| 5,0 mm     | MSF-3 |
| 5,5 mm     | MSF-4 |
| 6,0 mm     | MSF-5 |
| 7,0 mm     | MSF-6 |
| 8,0 mm     | MSF-7 |
| 9,0 mm     | MSF-8 |
| 10,0 mm    | MSF-9 |
| 11,0 mm    | MSG-1 |
| 12,0 mm    | MSG-2 |
| 13,0 mm    | MSG-3 |
| 14,0 mm    | MSG-4 |
| 15,0 mm    | MSG-5 |
| 16,0 mm    | MSG-6 |
| 17,0 mm    | MSG-7 |
| 18,0 mm    | MSG-8 |
| 19,0 mm    | MSG-9 |
| 20,0 mm    | MSH-1 |
| 21,0 mm    | MSH-2 |
| 22,0 mm    | MSH-3 |
| 23,0 mm    | MSH-4 |
| 24,0 mm    | MSH-5 |
| 25,0 mm    | MSH-6 |
| 26,0 mm    | MSH-7 |
| 27,0 mm    | MSH-8 |
| 28,0 mm    | MSH-9 |

**BOULITO SUPERPLAST**  
TUBULAR PROBE  
TUBULARSONDE  
180 mm




| SIZE<br>GRÖSSE |       |
|----------------|-------|
| 2/7 mm         | MSJ-1 |
| 3/8 mm         | MSJ-2 |
| 4/9 mm         | MSJ-3 |

**ZIEMER SUPERLAST**  
PROBE  
SONDE  
180 mm

**SUPERPLAST** instruments are made of shape-memory metal. This metal has unique qualities under different thermal conditions. The instruments are malleable in normal OR temperatures. Soft pressure allows "permanent" deformation of shape according to the surgeon's needs. However in higher temperature, such as those found in the autoclave, the instruments regain their initial shape.

**SUPERPLAST** Instrumente sind aus Metall mit Formgedächtnis. Sie bieten verschiedene Eigenschaften je nach thermischen Bedingungen. Sie sind bei normaler OP-Temperatur biegsam und können durch sanften Druck in die gewünschte Form gebracht werden. Unter höheren Temperaturen, wie bei der Sterilisation, nehmen sie ihre ursprüngliche Form wieder an.



**STORAGE AND STERILIZATION TRAY  
FOR BOULITO PROBES  
STERILISIER- UND LAGERSIEB  
FÜR BOULITO-SONDEN**

|                                         |       |
|-----------------------------------------|-------|
| max. OLIVE Ø 17 mm<br>160 x 210 x 26 mm | MSG-0 |
| max. OLIVE Ø 22 mm<br>160 x 210 x 30 mm | MSH-0 |



CHORDAE GAUGE AND  
KNOT ASSISTANT  
CHORDAE-LEHRE UND  
KNOTENASSISTENT  
200 mm  
MNV-9



CERAMO® BOULITO RONDO  
MICRO VALVE FORCEPS  
MIKRO KLAPPENPINZETTE

|        | 180 mm | 210 mm | 230 mm |
|--------|--------|--------|--------|
| 1,0 mm | COS-0  | —      | —      |
| 1,5 mm | COS-7  | COS-8  | COR-0  |




BOULITO SUPERPLAST  
3D GAUGE FOR AV-RECONSTRUCTION  
3D LEHRE FÜR AK-REKONSTRUKTION  
190 mm

| Ø<br>ANULUS |       |
|-------------|-------|
| 8 mm        | MST-2 |
| 9 mm        | MST-3 |
| 10 mm       | MST-4 |
| 11 mm       | MST-5 |
| 12 mm       | MST-6 |
| 13 mm       | MST-7 |
| 14 mm       | MST-8 |
| 15 mm       | MST-9 |
| 16 mm       | MSU-0 |
| 17 mm       | MSU-1 |
| 18 mm       | MSU-2 |
| 19 mm       | MSU-3 |
| 20 mm       | MSU-4 |
| 21 mm       | MSU-5 |
| 22 mm       | MSU-6 |
| 23 mm       | MSU-7 |
| 24 mm       | MSU-8 |
| 25 mm       | MSU-9 |



STERILIZATION TRAY FOR  
BOULITO SUPERPLAST 3D-GAUGE  
STERILISATIONSTRAY FÜR  
BOULITO SUPERPLAST 3D-LEHRE  
MST-0



WATERSTON

|                                                                                                      |                                                                                                      |                                                                                                      |                                                                                                      |
|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|
| <br>180 mm<br>DAY-7 | <br>210 mm<br>DAY-9 | <br>170 mm<br>DAY-6 | <br>200 mm<br>DAY-8 |
|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------|



BABY-MIXTER  
140 mm  
DBA-1



BABY-MIXTER  
180 mm  
DBA-2



BABY-ADSON  
185 mm  
DAZ-8



BABY-ADSON  
185 mm  
DCG-8



BABY-ADSON  
140 mm  
DAZ-7