## Rapid-Air®Series

### Alternating Air Floatation Mattress & Overlay System with Gentle Low Air Loss

### MODEL 4643 RAPID-AIR® OVERLAY

The **Rapid-Air®** Overlay Alternating Air Floatation Mattress Systems with Gentle Low Air Loss incorporate the features that clinicians demand for their patients in a portable system. The mattress utilizes 18 independent 5.5" deep air cells to reduce interface pressure for the prevention and treatment of Stages I-IV pressure sore ulcers. User friendly controls and ease of setup make these systems ideal for institutional use and the homebound patient. The mat-

tress is designed to fit all standard hospital beds.

CPR Valve quickly deflates mattress system for medical emergencies.

Vyvex-III® multi-stretch cover wicks moisture away from patients skin. The cover is low shear, anti-microbial, fluid resistant, fire retardant and washable. Cover is loose fitting and will conform to the patient.

Ventilated bladders include 20 lasered air holes in sacral area to help reduce skin maceration by reducing excess moisture under patient. Air cells are constructed of TPU, thermal polyurethane fabrics. This combination insures strength, durability and patient comfort.

Optional pump is available to quickly inflate and deflate mattress through CPR valve in under 4 minutes.





Optional 3" PPF® Pressurized Polymer Foam Base with a zippered Vyvex-I Nylon cover. The redistribution foam base protects the patient from bottoming out during power outages.



#### PUMP CONTROL FEATURES

- Alternating mode helps to redistribute pressure well below 32mmHg.
- Static Mode for Patient transfer and care
- Adjustable comfort settings adjusts mattress from soft to firm for patient comfort.
- Variable therapy cycle time 10, 15 & 20 minutes
- Fall Back Mode to prevent patient injury
- Lock Out Mode
- Normal Pressure indicator.
- Low Pressure Alert notifies caregiver of accidental hose disconnection.

# Rapid-Air® Series

### **SPECIFICATIONS**

### MATTRESS ONLY

**Rapid Air:**  $35'' \times 80'' \times 8''$  (including foam base)

Rapid Air Overlay:  $35" \times 80" \times 5.5"$ 

Weight: 20 lbs. Air Cells: 18, 5.5" sacs

15 Alternating cells

10 Ventilated cells with 20 lasered air holes

3 Static cells in head section

Cover: Vyvex-III®-removable zippered, low shear,

moisture/vapor permeable cover

Adjustable anchor straps secures mattress to bed frame. Fire evacuation safety handles standard on all models.



**Dimensions:** 10.75"x 7.5"x 3.5"

Weight: 6.8 lbs. Electrical: 120 60Hz

**Cycle Time:** 10, 15, 20 Minutes **Air Flow:** 7.5 Liters Maximum

Listings: ⊕ (€

Ground Resistance: 0.1 ohms maximum

Current Usage: less than 100 micro amps

**Power Cord:** 12 ft.

**Construction:** ABS fire retardant UL 94V-0

### RAPID-AIR® SERIES ORDERING INFORMATION

Model: 4643 Rapid Air® Mattress 35" x 80" x 8" & Pump Control Unit

Model: 4327 3" PPF® Base with zippered Nylon cover

Model: 4601 Pump Control Unit only

Model: 4320 Mattress only

Model: 4322Carry Bag with StrapModel: 9623Electric Inflator & Deflator

Model: 4324 Replacement 5.5" Ventilated Bladders
Model: 4325 Replacement 5.5" Cell Bladders
Model: 4326 Replacement Zippered Cover

Model: 4326-5R5 Replacement SRS safety rail system cover

\*WARRANTY

Power Control Unit:1 year limitedMattress:1 year limitedCover:6 months limited



Call: 774-236-0322

MediPro Email: info@mediprosupply.com

www.mediprosupply.com