



# ResMed

## AirFit™ F30i

Full face mask

Actual Size: 9 3/4" (Width) x 4" (Height)

Ruler is for scaling purposes.

How to use the fitting template:

1. Please ensure the template is printed to scale
2. Cut out along dotted line
3. Refer to diagram on template to determine appropriate cushion size
4. Refer to diagram on right to determine appropriate frame size

### Frame Sizing

Small frame

Standard frame

Large frame



### Cushion Sizing



**S**



**M**



**SW**



**W**

**1**



**2**



