

Implant Rx

Rx DATE _____

DOCTOR'S NAME _____

DOCTOR'S ADDRESS _____

OFFICE PHONE # _____

DUE DATE (FOR DELIVERY BY 5PM) • 12 WEEKDAY MINIMUM

PATIENT APPT. DATE _____

PATIENT ID/NAME (FIRST/LAST) _____

DENTIST SIGNATURE/LICENSE # (REQUIRED) _____

RUSH (Rush Fee may apply)

Adjustment

Remake

SHADE (Incisal/Body/Cervical) _____

Abutment Type

CEMENT RETAINED:

- Zirconia Custom Abutment
- Titanium Custom Abutment

– OR –

SCREW RETAINED:

- Ti-Base (Stock Abutment)
- Screwmentable – Lab will cement
(One-Piece Custom Titanium Abutment & Crown w/ Access Hole)
- Screwmentable – Doctor will cement
(Two-Piece Custom Titanium Abutment & Crown w/ Access Hole)

Abutment Preference

- Generic
- Dentsply Sirona (Atlantis) - OEM
- Straumann/Neodent (CARES) - OEM
- Nobel Biocare (NobelProcera) - OEM
- BioHorizons (Vulcan) - OEM
- ZimVie (BellaTek) - OEM
- Implant Direct (CustomDirect) - OEM

Jig

- Verification Jig – Mandatory for ≥3 implant screw retained bridge/splinted units
- Seating Jig – Recommended for a Biomet 3i Certain and Astra EV custom abutment

Single Unit: Circle

Bridge/Splinted Units: Use Line to Connect Units (Pontics)



Digital Scans

Date: _____

Type: _____

Current Bite

- Class I
- Class II
- Class III
- Edge to Edge
- Open Bite
- Crossbite

ASC/Dynamic Driver


Does your office have the OEM ASC driver?

- Yes
- No, please provide!

Does your office have the Preat Dynamic (ASC) driver?

- Yes
- No, please provide!


Emergence Profile



- Straight
- Concave
- Convex

Crown Design

PONTIC DESIGN:



OCCLUSAL CLEARANCE:

- Out of occlusion
- Light in centric
- Heavy w/ eccentric

MATERIAL:

- Zirconia ★
- Lithium Disilicate
- PMMA (Temp.)
- Full Cast
- PFM

ESTHETIC CHARACTERIZATIONS:

- Hypo-Calcification
- Craze lines
- High translucency
- Keep diastema
- Close diastema

IF INSUFFICIENT ROOM:

- Trim opposing
- Reduction coping
- Metal occlusion
- Metal island
- Call

Tooth #	Implant Brand/Manufacturer	Implant Interface Platform	Platform Size	Platform Color

Rx SPECIFIC INSTRUCTIONS:

★ SDA Recommended
 ASC: Angled Screw Channel
 OEM: Original Equipment Manufacturer

CONTINUED ON BACK