
Item 25-325 Used MBraun™ Labmaster Pro SP – 4 Port Glovebox with MB 200B Purifier, -40 C Freezer and VOH 150 Vacuum Oven.

Price: Make offer

This used MBraun Labmaster Pro SP glovebox . The glovebox includes large & mini antechamber, MB 200B gas purifier with touch screen controls, freezer, vacuum oven, inline solvent trap, automated evacuation and refill of large antechamber, O2 & H2O sensors and foot switch. Two Vacuum pumps are required but not included. Gloves are not included

90 Day Warranty for Parts.












-40 degree C Freezer



150 degree Vacuum Oven









Purifier with Inline Solvent Trap







Purifier electronic panel



Vacuum oven controller

M BRAUN**MB200B**

Serial No.: 18-410

Purifier Electrical Connection Q1: 115V 20A
Purifier Electrical Connection Q2: 115V 15A

Oven Electrical Connection: 230V 1ph 15A

Manufactured: December 2018

Glove Box

Box Pressure, Max: + 20 mbar

Box Pressure, Min: - 20 mbar

Integral Leak Rate: <0.05 vol % / h

Antechamber

Total Vacuum: 10^{-3} mbarIntegral Leak Rate: < 10^{-5} mbar l/s

Gas Purification System

Processing Gas: N₂ / Ar / HeGas Purification, Max: <1ppm H₂O / O₂Integral Leak Rate: < 10^{-6} mbar l/s

Perm. Ambient Temperature: 15° C - 30° C

Blower Type: Single Stage

Blower Flow Rate, Max: 85 m³/hM.Braun Inc.
14 Station Way
Stoughton MA 03585


UsedGloveBox.com

We Buy, Sell & Service Used Gloveboxes

702-343-0297 info@usedglovebox.com



Glovebox leak check will be conducted with an *Inficon* H2 leak detector using a 5% H2 95% N2 gas



On-Site Installation, Start Up Assistance and Training -We provide the option of on-site installation assistance. Our personnel have combined inert glovebox work experience of over 40 years. We can train your staff on the operation of the glovebox system. We also will provide long term technical support by email and phone at no charge.