

Item 25-325 Used MBraun™ Labmaster Pro SP – 4 Port Glovebox with MB 200B Purifier, -40 C Freezer and VOH 150 Vacuum Oven.

Price: Make offer

This used MBraun Labmaster Pro SP glovebox . The glovebox includes large & mini antechamber, MB 200B gas purifier with touch screen controls, freezer, vacuum oven, inline solvent trap, automated evacuation and refill of large antechamber, O2 & H2O sensors and foot switch. Two Vacuum pumps are required but not included. Gloves are not included

90 Day Warranty for Parts.





UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS









UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS



-40 degree C Freezer

Web: www.usedglovebox.com

Email: info@usedglovebox.com

Tel: 702-343-0297



UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS



150 degree Vacuum Oven

Web: www.usedglovebox.com


Email: info@usedglovebox.com

Tel: 702-343-0297





UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS





UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS



Purifier with Inline Solvent Trap




UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS






UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS





Purifier electronic panel



Vacuum oven controller



UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS

M BRAUN

MB200B

Serial No.:	18-410
Purifier Electrical Connection Q1:	115V 20A
Purifier Electrical Connection Q2:	115V 15A
Oven Electrical Connection:	230V 1ph 15A
Manufactured:	December 2018
Glove Box	
Box Pressure, Max	+ 20 mbar
Box Pressure, Min:	- 20 mbar
Integral Leak Rate:	<0.05 vol % / h
Antechamber	
Total Vacuum:	10^{-3} mbar
Integral Leak Rate:	< 10^{-5} mbar l/s
Gas Purification System	
Processing Gas:	N ₂ / Ar / He
Gas Purification, Max:	<1ppm H ₂ O / O ₂
Integral Leak Rate:	< 10^{-6} mbar l/s
Perm. Ambient Temperature	15° C - 30° C
Blower Type:	Single Stage
Blower Flow Rate, Max:	85 m ³ /h

M. Braun Inc.
14 Martin Way
Statham, NE 68885



UsedGloveBox.com

We Buy, Sell & Service Used Gloveboxes


702-343-0297 info@usedglovebox.com



UsedGloveBox.com
SCIENTIFIC MARKET SOLUTIONS



Glovebox leak check will be conducted with an *Inficon* H2 leak detector using a 5% H2 95% N2 gas



Optional On-Site Installation, Start Up Assistance and Training -We provide the option of on-site installation assistance. Our personnel have combined inert glovebox work experience of over 40 years. We can train your staff on the operation of the glovebox system. We also will provide long term technical support by email and phone at no charge.