

Customer Account Number: _____

Account Manager: _____

Phone Number: _____



DJO SOUTH AFRICA BRACING & SUPPORTS **CATALOGUE 2021**

DJO SOUTH AFRICA BRACING & SUPPORTS CATALOGUE




DJO Orthopaedic South Africa (Pty) Ltd | Unit 6 | 2 on London | 2 London Circle | BrackenGate Business Park | Brackenfell | South Africa | 7560

Tel: +27 (0) 21 2762968 | orders.zaf@DJOglobal.com

www.DJOglobal.eu

© 2020 DJO - 000-000-04527-ZA - Rev D



AIRCAST®

chattanooga®

COMPEX®

DONJOY®

Dr.Comfort®

PROCARE®



MISSION & VISION

OUR MISSION

Life changing innovation and care

CORE VALUES

- Integrity
- Excellence
- Compassion
- Resilience
- Teamwork

DJO SOUTH AFRICA & SUB SAHARA AFRICA DISTRIBUTION



DJO ORTHOPAEDIC SOUTH AFRICA (PTY) LTD

Unit 6 | 2 on London | 2 London Circle
BrackenGate Business Park | Brackenfell
South Africa | 7560

Contents

Upper Limb

UltraSling Quadrant™	6
UltraSling® Pro	6
UltraSling® III	7
UltraSling® III ER	7
UltraSling® III AB	8
Quick-Fit Shoulder Immobilizer	8
Healthcare Bear Arm Sling	9
Clinical Shoulder Immobiliser	9
X-Act ROM Elbow	9
Sully® Shoulder Stabiliser	9
Humeral Stabilising System (HSS)	10
Humeral Fracture Brace/ Over the Shoulder	10
Spider™ Pad Elbow	10
Pneumatic Armband	10

Wrist & Hand

Short Arm Fracture Brace with Boa® - Open Thumb	13
Boxer's Fracture Brace with Boa®	13
Wrist Brace with Boa®	13
Radial Gutter Fracture Brace with Boa®	13
Wrist Brace	14
Thumb Spica Fracture Brace with Boa®	14
Long Thumb Spica II with Boa®	14
Short Thumb Spica II	14
Hand Based Radial Gutter Brace	15
Hand Based Ulnar Gutter Brace	15
Paediatric Short Arm Fracture Brace with Boa® - Open Thumb	15
Extended Short Thumb Spica	16
Humeral Fracture Brace with Boa®	16
Sheet Materials	16
Boa Adjustable Tensioning Kit	16
Straps 12.5cm	17
Top Loading Oven	17
Oven Extender Ring Kit	17
Edging Tape	17
RespiForm® Wrist	18
RespiForm+® Wrist	18
DigiForm® Wrist	19
DigiForm+® Wrist	19
PolyForm® Wrist	20
ErgoForm® Wrist	20
CarpaForm® Wrist	20
RhizoForm® Wrist	20
ComfortFORM™ Wrist	21
ComfortFORM™ Wrist & Thumb	21
Quick Fit Wrist II	21
Quick Fit WTO	21
Frog Splints / Baseball Splints	22
10" Lace up Wrist	22
Finger Splints	22

Back & Neck Supports

Exos FORM™ II	25
DonJoy Back Brace II TLSO	27
Immostrap™	27
Bodystrap	28
Actistrap	29
Abdostrap II	29
MyBabystrap™ Fusion	30
ComfortFORM™ Back Support	30
Premium Panel Elastic Binder	31
Hernia Belt	31
Form Fit Collar	32
California Collar	32

Knee Ligament Bracing

Defiance™	36
Defiance™ Colour Options	37
Armor™ CI	38
FullForce™ with Force Point	38
Playmaker Xpert	39
Quickfit Hinged Knee	39
Reddie Brace	40
Stabilized Knee Support	40

Knee OA

Adjustable OA Defiance™	43
OA Nano™	44
OA Reaction™ Web	44
OA Adjuster™ 3	44

*Products marked "available on request" are non stock items but can be brought in on special order with a 3-4 week delivery period.

Post-Op

XROM	46
X-Act™ ROM Lite	46
Cool IROM®	47
TROM™ Advance	47
Immo AT4 3V	47
VersaROM™ Hip	48
Tibial Fracture Brace	48

Patella Femoral

Reaction™ Web	50
Tru-Pull® Lite	51
Lateral "J"	52
Infrapatellar Band	52

Walkers

AirSelect™ Elite	55
AirSelect™ Standard	56
AirSelect™ Short	57
XP Diabetic Walker™	57
Heel Wedges	58
Nextep™ Contour Walker	58
XcelTrax™ Tall	59
XcelTrax™ Ankle	59
MiniTrax™	60

Ankle

AirSport+	62
Stabilising Speed Pro	63
Stabilising Pro	63
A60™ Ankle Support	64
Air-Stirrup®	64
Air-Stirrup® Universe™	64

Diabetic, Offloading & Foot Supports

Squared Toe Post-Op Shoe	66
Off-Loading Diabetic Shoe	66
Podalux II	67
Remedy Pro™	67
Medical/Surgical Shoe	68
Super-Lite A.F.O.	68
Podous Boot	69
Plantar Fasciitis Night Splint	69
Silicone Heel Cups	69
Podalib	70
Hallux Valgus	70

Post-Op

IceMan® CLEAR3™	72
Universal Cuff	72
McGuire Cuff	72
Cryo/Cuff™ IC Cooler	73
Knee Cuff	73
Ankle Cuff	73
Shoulder Cuff	74
Elbow Cuff	74
Hand/Wrist Cuff	74

Various

VenaPure®	76
Annie Black	77
Lana Black	77
Flute Black	77
Grace Seafoam	77
Protector Black	78
Performance Black/Grey	78
Classic Black	78
Winner Black	78
Brian Black	78

Terms and Conditions	79
----------------------	----

UPPER LIMB

AVAILABLE ON REQUEST



UltraSling Quadrant™

DONJOY®

INDICATIONS

Post-operative immobilisation in 15° of abduction or greater for shoulder hemiarthroplasty, total shoulder replacement, shoulder reconstruction, superior labral repair and shoulder debridement, protocols involving stabilisation following fractures to the proximal humerus, elbow and forearm, Bankart procedures, rotator cuff repair, biceps tendon surgery, shoulder dislocation and elbow ligament and tendon procedures.

CHARACTERISTICS

- TempGuard™ - Cool to touch, anti-microbial material and optional padded shoulder strap for extended duration of use
- Optimal versatility provides four angles of immobilisation in one brace; internal, external, abduction or adduction
- ProTherapy™ - Precision based system which allows the patient to flex and extend their arm, promoting therapeutic movement throughout recovery
- Elbow opening to eliminate pressure on the olecranon bursa
- Patient friendly, single handed application for patient compliance and optional shoulder strap
- Optional extension belt available separately

ORDERING INFORMATION

Part Number	Arm	Waist Circumference
11-0211	Left	61cm - 114cm
11-0210	Right	61cm - 114cm
25-0212	Extension Belt	114cm - 152cm

Measurement taken at waist circumference 1" (2.5cm) above belly button

COMING SOON



UltraSling® Pro

DONJOY®

INDICATIONS

Universally designed to fit the right or the left arm, the UltraSling PRO features an envelope that can be adjusted and sized to accommodate a wide range of patients. Plus, the revolutionary shoulder strap system allows for fast, easy application.

CHARACTERISTICS

- Post-operative rotator cuff repairs
- Post-operative Bankhart procedures
- Capsular shifts
- Glenohumeral Dislocation/Subluxation
- Soft tissue strains/repairs

ORDERING INFORMATION

Part Number	Measurement
11-0447-9	UltraSling PRO
11-0471-9	UltraSling PRO ER, 15 Degrees
11-0472-9	UltraSling PRO ER, 30 Degrees
11-0473-9	UltraSling PRO AB



UltraSling® III

DONJOY®

INDICATIONS

Immobilisation for rotator cuff repairs, capsular shifts, Bankart repairs, glenohumeral dislocations/subluxation and soft tissue repairs/strains.

CHARACTERISTICS

- Easy open front panel allowing forearm exercise and exercise ball to stimulate circulation
- Moisture wicking mesh fabric for greater comfort to reduce risk of secondary infections, particularly in longer term use
- Allows arm to be positioned in a variety of positions dependent on post-operative requirements (either more posterior or anterior)
- Quick release buckles make it easy to take on and off for exercise/therapy
- Fits left or right arm

ORDERING INFORMATION

Part Number	Measurement	Size
11-0449-2	<28cm	S
11-0449-3	28cm - 33cm	M
11-0449-4	33cm - 38cm	L
11-0449-5	>38cm	XL

Measurement taken from elbow crease to base of index finger



UltraSling® III ER

DONJOY®

INDICATIONS

Maintaining the shoulder and arm in an externally rotated position, providing immobilisation for rotator cuff repairs, capsular shifts, Bankart repairs, glenohumeral dislocation/subluxation and soft tissue strains/repairs.

CHARACTERISTICS

- Available with 15° or 30° foam pillow to suit your protocol
- De-rotational shoulder strap to keep product in place and quick release buckles aid application

ORDERING INFORMATION

Part Number		Rotation	Measurement	Size
Left Arm	Right Arm			
11-1344-2	11-1343-2	15°	<28cm	S
11-1344-3	11-1343-3	15°	28cm - 33cm	M
11-1344-4	11-1343-4	15°	33cm - 38cm	L
11-1346-2	11-1345-2	30°	<28cm	S
11-1346-3	11-1345-3	30°	28cm - 33cm	M
11-1346-4	11-1345-4	30°	33cm - 38cm	L

Measurement taken from elbow crease to base of index finger



UltraSling® III AB

DONJOY®

INDICATIONS

Immobilisation in neutral to an internally rotated position for rotator cuff repairs, capsular shifts, Bankart repairs, glenohumeral dislocations/subluxation and soft tissue repairs/strains.

CHARACTERISTICS

- Large comfortable pillow gives 45° or 60° abduction where required
- Easy open front panel allowing forearm exercise and exercise ball to stimulate circulation
- Promotes auxiliary air exchange to reduce risk of secondary infections and moisture wicking mesh fabric for greater comfort, particularly in longer term use
- Thumb straps allowing functional hand positioning
- Fits left or right arm

ORDERING INFORMATION

Part Number	Measurement	Size
11-0450-2	<28cm	S
11-0450-3	28cm - 33cm	M
11-0450-4	33cm - 38cm	L
11-0450-5	>38cm	XL

Measurement taken from elbow crease to base of index finger



Quick-Fit Shoulder Immobilizer

AIRCAST®

INDICATIONS

Support for post-operative rotator cuff repairs, Bankart procedures, capsular shifts, glenohumeral, dislocation/subluxation, soft tissue strains/repairs. The pillow would be considered when the tissues across the top of the glenohumeral joint need to be in a more relaxed position during healing

CHARACTERISTICS


- Constructed of breathable mesh for a cool, comfortable fit
- Release clip design on the waist allows for easy removal without changing the fit adjustment
- Custom fit by folding the sling onto itself, to adjust to the arm length of the patient
- 15° abduction pillow for proper arm placement in internal rotation
- Thumb loop inside the sling to hold the hand in place

ORDERING INFORMATION

Part Number	Min Length	Max Length	Depth
06AB	26.5 cm	43 cm	18 cm
06A (Sling Only)	26.5 cm	43 cm	18 cm

Measurement taken from elbow crease to base of index finger

Healthcare Bear Arm Sling



INDICATIONS

Injuries to the arm, wrist or hand of paediatric patients.

CHARACTERISTICS

- Made from a cotton/poly envelope with a cotton web strap with quick release buckle for easy application
- Paediatric friendly, bear print design on the sling
- Suitable for use in conjunction with a cast
- Fits left or right arm

ORDERING INFORMATION

Part Number	Depth x Length	Size
79-99132	11.5cm x 28cm	XS / Paediatric
79-99133	16.5cm x 33cm	S / Child

Clinical Shoulder Immobiliser



INDICATIONS

Preventing shoulder rotation, post trauma after shoulder (sub)luxations and for over use shoulder injuries.

CHARACTERISTICS

- High quality, durable poly cotton sling
- Contact closure adjustment on the shoulder and body straps
- Shoulder pad to provide additional comfort
- Fits left or right arm

ORDERING INFORMATION

Part Number	Depth x Length	Size
79-84012	15cm x 28cm	XS
79-84013	18cm x 33cm	S
79-84015	20.5cm x 38cm	M
79-84017	23cm x 46cm	L
79-84018	25.5cm x 51cm	XL

X-Act ROM Elbow



INDICATIONS

Fixed or controlled range of motion post injury or post operatively, stable fractures, biceps repair, mild to severe soft tissue injuries.

CHARACTERISTICS


- User-friendly range of motion control with flexion settings from -10° to 120° in 10° increments and extension settings from -10° to 90° in 10° increments
- Quick-lock function from 0° to 90° in 15° increments
- Cuffs telescope independently allowing a precise fit on any patient, adjustability from 36cm to 46cm
- Durable lateral aluminium hinge maximises protection and minimises brace weight
- Unique semi-rigid thermoformed and mouldable liners with anti-migration properties

ORDERING INFORMATION

Part Number	Arm	Size
11-9122	Left	U
11-9121	Right	U

Part Number	Accessory	Size
25-6755	Brace Lock Kit (Pack of 10)	U

Sully® Shoulder Stabiliser



INDICATIONS

Shoulder immobilisation for anterior, multi-directional, inferior, and posterior instabilities, rotator cuff deceleration, shoulder separations, and muscle strains.

CHARACTERISTICS

- Stabilises, assists and restricts movement of the shoulder post-operatively
- Can pull the humerus posteriorly as it restricts external rotation and abduction
- Elastic straps attach with Velcro at any point, in any direction for improved patient fit
- Smooth, controlled assistance or restraint follows the natural movement of the muscles and joints
- Fits left or right arm

ORDERING INFORMATION

Part Number	Chest Circumference	Bicep Circumference	Size
11-0525-1	86.5cm - 101.5cm	22.5cm - 26.5cm	XS
11-0525-2	91.5cm - 101.5cm	28cm - 33cm	S
11-0525-3	96.5cm - 106.5cm	33cm - 38cm	M
11-0525-4	106.5cm - 116.5cm	35.5cm - 43cm	L
11-0525-5	116.5cm - 129.5cm	38cm - 43cm	XL

Measurement taken at bicep in non-flex position

AVAILABLE ON REQUEST

DONJOY®

Humeral Stabilising System (HSS)



INDICATIONS

Immobilisation of stable midshaft humeral fractures, stable supracondylar fractures and biceps/triceps repairs.

CHARACTERISTICS

- Lightweight design offers three point stabilisation
- 900g weight for axial traction
- Full circumference shells and padding with the ability to break down the upper arm cuff to suit all patients
- Alternative to serial casting

ORDERING INFORMATION

Part Number	Arm	Size
11-0113-9-13066	Left	U
11-0112-9-13066	Right	U

AVAILABLE ON REQUEST

DONJOY®

Spider™ Pad Elbow



INDICATIONS

The unique honeycomb design contours for comfortable movement while providing cushioning protection. The Spider Pad is ideal for football and volleyball.

CHARACTERISTICS

- High impact resistant foam for protection during contact sports and job related activities
- Neoprene for warmth and compression
- Fits left or right elbow

ORDERING INFORMATION

Part Number	Elbow Circumference	Size
11-0392-1-06060	18cm - 20cm	XS
11-0392-2-06060	20cm - 22.5cm	S
11-0392-3-06060	22.5cm - 25cm	M
11-0392-4-06060	25cm - 28cm	L
11-0392-5-06060	28cm - 30.5cm	XL
11-0392-6-06060	30.5cm - 33cm	XXL
11-0392-7-06060	33cm - 35.5cm	XXXL

Measurement taken at elbow joint

AVAILABLE ON REQUEST

PROCARE®

Humeral Fracture Brace/Over the Shoulder



INDICATIONS

Indicated for management of mid-shaft and distal third humeral fractures.

CHARACTERISTICS

- Pre-formed orthopaedic grade thermoplastic cast-like construction for the humerus
- Hook and loop straps for easy adjustment
- Deltoid cap extension, providing a greater alignment control and prevents slippage
- Supplied with two double thickness stockinettes
- Fits left or right arm

ORDERING INFORMATION

Part Number	Bicep Circumference	Length	Size
79-97953	20.5cm - 28cm	25.5cm	S
79-97955	25.5cm - 33cm	26cm	M
79-97957	30.5cm - 38cm	28.5cm	L
79-97958	35.5cm - 43cm	32.5cm	XL

AIRCAST®

Pneumatic Armband



INDICATIONS

Lateral epicondylitis (tennis elbow), medial epicondylitis (golfer's elbow).

CHARACTERISTICS

- Single pre-inflated aircell concentrates compression on the extensor muscle, providing more support, less constriction
- Latex free
- Fits left or right elbow

ORDERING INFORMATION

Part Number	Forearm Circumference	Size
05A-B	20cm - 35.5cm	U

Part Number	Description
0501	Armband Aircell

WRIST & HAND

A Great Alternative to Conventional Casting!

Exos is the only bracing system offering a removable, adjustable, reformable and waterproof solution for the treatment of fractures and other injuries requiring stabilisation.



REMOVABLE

ALTHOUGH THE EXOS BRACE IS CUSTOM FITTED TO YOUR PATIENT, IT CAN BE EASILY REMOVED, ALLOWING:

- Examination of the fracture site
- Examination of a wound or incision
- Assessment of muscle atrophy
- The ability to initiate rehabilitation therapy earlier in the healing process
- Skin and brace cleaning
- Sawless removal

REFORMABLE

THE EXOS BRACE IS COMPLETELY REFORMABLE, ALLOWING:

- Modification to accommodate anatomical changes as healing progresses
- Serial changes of position
- One brace for duration of treatment
- Less material waste


ADJUSTABLE

EXOS USES THE EXCLUSIVE BOA® REEL AND LACE DYNAMIC CLOSURE SYSTEM, ALLOWING:

- Expansion to allow for initial swelling and to reduce the risk of Compartment Syndrome
 - Tension modification for custom fit despite atrophy or changes in soft tissue
 - Improved patient comfort and fit
-
- Exos braces are waterproof, allowing patients to continue their regular daily activities such as showering, swimming
 - Anti-microbial treated material to reduce odour
 - Radiolucent, allowing x-rays to be taken without removal of the brace
 - Latex free

EXOS SIZING CHART

Hand Circumference	Size
14.5cm - 16.5cm	XXXS
16.5cm - 18.5cm	XXS
18.5cm - 20cm	XS
20cm - 23cm	S
23cm - 25.5cm	M
25.5cm - 28cm	L
>28cm	XL



Wrap tape around hand to measure size



1. Simply heat brace in our oven for 3-5 minutes.



2. Open and slide onto patient's limb.



3. Gently mould and adjust for a comfortable fit.



4. Use the reel and lace to tighten as needed.



Short Arm Fracture Brace with Boa® - Open Thumb



INDICATIONS

Fractures to the distal radius or distal ulna, and can be employed for acute injuries as well as post-operatively. Support and protection for major injuries to the wrist.

CHARACTERISTICS

- The open thumb design allows for enhanced ease of fit and patient comfort

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
312-31-1111	312-32-1111	18.5cm - 20cm	XS
312-41-1111	312-42-1111	20cm - 23cm	S
312-51-1111	312-52-1111	23cm - 25.5cm	M
312-61-1111	312-62-1111	25.5cm - 28cm	L
312-71-1111	312-72-1111	>28cm	XL



Wrist Brace with Boa®



INDICATIONS

Carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapholunate dislocations, triangular fibrocartilage complex tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures.

CHARACTERISTICS

- The openings around the base of the thumb and across the palmar crease enable range of motion and provide greater comfort without loss of stability
- Detachable nylon strap across the web space provides comfortable, adjustable control distally

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
221-21-1111	221-22-1111	16.5cm - 18.5cm	XXS
221-31-1111	221-32-1111	18.5cm - 20cm	XS
221-41-1111	221-42-1111	20cm - 23cm	S
221-51-1111	221-52-1111	23cm - 25.5cm	M
221-61-1111	221-62-1111	25.5cm - 28cm	L
221-71-1111	221-72-1111	>28cm	XL



Boxer's Fracture Brace with Boa®



INDICATIONS

Stabilises fractures to the 4th and/or 5th metacarpal bones, provides stabilisation to the 4th and 5th interphalangeal and wrist joints. May be used for head/neck fracture of the 4th and/or 5th metacarpal.

CHARACTERISTICS

- Dorsal hook-and-loop secure tab secures the overlapping (top) portion of brace for better stability
- Hook-and-loop fastener straps across the web space and around the 4th and 5th digits can be attached anywhere to the brace surface for better stability
- Trimmable finger section accommodates desired range of motion

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
325-31-1111	325-32-1111	18.5cm - 20cm	XS
325-41-1111	325-42-1111	20cm - 23cm	S
325-51-1111	325-52-1111	23cm - 25.5cm	M
325-61-1111	325-62-1111	25.5cm - 28cm	L
325-71-1111	325-72-1111	>28cm	XL



Radial Gutter Fracture Brace with Boa®



INDICATIONS

Stabilises fractures to the 2nd and/or 3rd metacarpal bones and the respective metacarpophalangeal joints of the hand. Radial gutter portion allows the caregiver to achieve greater than 70° of flexion at the metacarpal joints of the 2nd and 3rd digits.

CHARACTERISTICS

- Full hook-and-loop release ribbon allows for easier application
- Dorsal hook-and-loop secure tab secures the over-lapping (top) portion of brace for better stability
- Hook-and-loop straps across the web space and around the 2nd and 3rd digits can be attached anywhere to the brace surface for better stability
- Trimmable finger section accommodates desired range of motion

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
326-31-1111	326-32-1111	18.5cm - 20cm	XS
326-41-1111	326-42-1111	20cm - 23cm	S
326-51-1111	326-52-1111	23cm - 25.5cm	M
326-61-1111	326-62-1111	25.5cm - 28cm	L
326-71-1111	326-72-1111	>28cm	XL



Wrist Brace

INDICATIONS

Carpal bone injuries such as lunate, pisiform, or triquetral fractures, scapho-lunate dislocations, triangular fibrocartilage complex tears, radio-carpal ligament injuries, or minimally displaced or stable distal radius or distal styloid fractures. May be used to support stable wrist fractures or control wrist motion for injuries that require stabilisation.



CHARACTERISTICS

- Openings around the base of the thumb and across palmar crease enable range of motion without loss of stability
- Detachable nylon strap across the web space provides comfortable, adjustable control distally
- Ulnar styloid relief holes help reduce pressure on bony prominence

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
220-31-1111	220-32-1111	18.5cm - 20cm	XS
220-41-1111	220-42-1111	20cm - 23cm	S
220-51-1111	220-52-1111	23cm - 25.5cm	M
220-61-1111	220-62-1111	25.5cm - 28cm	L
220-71-1111	220-72-1111	>28cm	XL



Thumb Spica Fracture Brace with Boa®

INDICATIONS

Stabilisation to the thumb, first metacarpophalangeal joint, the first carpometacarpal joint, wrist joint, distal radius and distal ulna. Useful for all injuries to the thumb, such as Bennett's or complex Rolando's fractures, suspected or confirmed navicular (scaphoid) fractures, and for other carpal bone injuries.



CHARACTERISTICS

- Bivalve thumb opening accommodates swelling/dressings
- Extended, trimmable length of the thumb for maximum protection
- Full hook-and-loop release ribbon for easier application

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
315-31-1111	315-32-1111	18.5cm - 20cm	XS
315-41-1111	315-42-1111	20cm - 23cm	S
315-51-1111	315-52-1111	23cm - 25.5cm	M
315-61-1111	315-62-1111	25.5cm - 28cm	L
315-71-1111	315-72-1111	>28cm	XL



Long Thumb Spica II with Boa®

INDICATIONS

Stabilisation to the thumb in adduction, metacarpophalangeal joint, carpometacarpal joint, distal radius, distal ulna and conditions where immobilisation of the wrist in the sagittal and frontal plane is required. May be used pre-operative, post-operative, post-traumatic as well as for carpal tunnel syndrome.



CHARACTERISTICS

- Increased forearm length for better immobilisation and a Velcro release tab allows the brace to open more fully
- Adjustable thumb strap provides better immobilisation

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
232-31-1111	232-32-1111	18.5cm - 20cm	XS
232-41-1111	232-42-1111	20cm - 23cm	S
232-51-1111	232-52-1111	23cm - 25.5cm	M
232-61-1111	232-62-1111	25.5cm - 28cm	L
232-71-1111	232-72-1111	>28cm	XL



Short Thumb Spica II

INDICATIONS

Stabilisation to the first metacarpophalangeal joint and may be used pre-operative, post-operative, post-traumatic as well as for degenerative conditions.



CHARACTERISTICS

- Adjustable thumb circumference accommodates a variety of sizes
- Thumb strap provides better immobilisation
- Lock/loop contact closure

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
212-31-1111	212-32-1111	18.5cm - 20cm	XS
212-41-1111	212-42-1111	20cm - 23cm	S
212-51-1111	212-52-1111	23cm - 25.5cm	M
212-61-1111	212-62-1111	25.5cm - 28cm	L
212-71-1111	212-72-1111	>28cm	XL



Hand Based Radial Gutter Brace

INDICATIONS

Stabilise fractures to the 2nd and/or 3rd metacarpals or phalanges. Radial gutter portion allows the caregiver to achieve greater than 70° of flexion at the metacarpal joints of the 2nd and 3rd digits.



CHARACTERISTICS

- Web straps with hook-and-loop contact closure
- Maintains free range of motion at the wrist
- Entire brace can be trimmed for desired range of motion

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
126-31-1111	126-32-1111	18.5cm - 20cm	XS
126-41-1111	126-42-1111	20cm - 23cm	S
126-51-1111	126-52-1111	23cm - 25.5cm	M
126-61-1111	126-62-1111	25.5cm - 28cm	L
126-71-1111	126-72-1111	>28cm	XL



Hand Based Ulnar Gutter Brace

INDICATIONS

Stabilises fractures to the 4th and/or 5th metacarpals or phalanges. The ulnar gutter portion of the brace allows the caregiver to achieve greater than 70° of flexion at the metacarpal joints of the 4th and 5th digits.



CHARACTERISTICS

- Web straps with hook-and-loop contact closure
- Maintains free range of motion at the wrist
- Entire brace can be trimmed for desired range of motion

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
125-31-1111	125-32-1111	18.5cm - 20cm	XS
125-41-1111	125-42-1111	20cm - 23cm	S
125-51-1111	125-52-1111	23cm - 25.5cm	M
125-61-1111	125-62-1111	25.5cm - 28cm	L
125-71-1111	125-72-1111	>28cm	XL



Paediatric Short Arm Fracture Brace with Boa® - Open Thumb



INDICATIONS

Support and protection for major injuries to the wrist. Primarily for fractures to the distal radius or distal ulna and can be employed for acute injuries as well as post-operatively.

CHARACTERISTICS

- Open thumb design allows for enhanced ease of fit
- Lycra top cover allows for faster heating and cooling, speeding the application process

ORDERING INFORMATION

Part Number		Hand Circumference	Colour	Size
Left	Right			
311-11-2293	311-12-2293	14.5cm - 16.5cm	Spider Print	XXXS
311-21-2293	311-22-2293	16.5cm - 18.5cm	Spider Print	XXS
311-11-3285	311-12-3285	14.5cm - 16.5cm	Polka Dot	XXXS
311-21-3285	311-22-3285	16.5cm - 18.5cm	Polka Dot	XXS
311-11-1111	311-12-1111	14.5cm - 16.5cm	Black	XXXS
311-21-1111	311-22-1111	16.5cm - 18.5cm	Black	XXS





Extended Short Thumb Spica

INDICATIONS

Immobilises the base of the thumb and the interphalangeal joint to stabilise injuries to the ulnar collateral ligament (UCL) (Gamekeeper's thumb), fractures or avulsion injuries to the proximal first phalanx or distal first metacarpal, severe sprains of the thumb, or for arthritic conditions involving the metacarpophalangeal joint, or the interphalangeal joint. May also be used to immobilise the thumb joints following surgical repairs and soft tissue injury.



CHARACTERISTICS

- Bivalve thumb opening accommodates swelling or dressings
- Hook-and-loop contact closure

ORDERING INFORMATION

Part Number		Hand Circumference	Size
Left	Right		
211-31-1111	211-32-1111	18.5cm - 20cm	XS
211-41-1111	211-42-1111	20cm - 23cm	S
211-51-1111	211-52-1111	23cm - 25.5cm	M
211-61-1111	211-62-1111	25.5cm - 28cm	L
211-71-1111	211-72-1111	>28cm	XL



Sheet Materials

INDICATIONS

Multiple layered thermoformable material in a single construct ideal for custom fabrication for musculoskeletal stabilisation.



CHARACTERISTICS

- 45.5cm x 61cm sheet size
- Length of UBL edging tape is available as an accessory (Size 2.5cm x 60cm) - pack of 10
- Latex free

ORDERING INFORMATION

Part Number	Sheet Thickness
800-40-1111	1.0mm
800-50-1111	1.3mm
800-60-1111	1.5mm

Part Number	Description
800-1111	UBL Edging Tape (Pack of 10)



Humeral Fracture Brace with Boa®

INDICATIONS

Stabilise fractures to the mid-shaft of the humerus and some injuries to the more proximal humerus.



CHARACTERISTICS

- Brace can be easily trimmed to the appropriate circumference
- Edge tape provided in the package can be easily applied to the trimmed edge of the brace

ORDERING INFORMATION

Part Number	Width	Arm Circumference	Size
330-40-1111	10.5cm	15.5cm - 40.5cm	S
330-50-1111	12.5cm	23cm - 45.5cm	M
330-60-1111	15cm	28cm - 48cm	L



Boa Adjustable Tensioning Kit

INDICATIONS

Universal tensioning kit for use in combination with Exos Sheet Material for creating custom made splints.



CHARACTERISTICS

- Boa reel and lace dynamic closure system allows for expansion to accommodate any initial swelling

ORDERING INFORMATION

Part Number	Size
800-ATK	U



Straps 12.5cm

INDICATIONS

Universal straps for use in combination with Exos Sheet Material for creating custom made splints.

CHARACTERISTICS

- Web material straps with hook/loop contact spots at both ends
- Length 12.5cm
- Supplied in a pack of 10



ORDERING INFORMATION

Part Number	Length
800-50-STRAP	12.5cm



Top Loading Oven

INDICATIONS

Conveniently sized oven used to heat Exos lower and upper extremity products, to gently mould the brace to the patient's limb.

CHARACTERISTICS

- 230v, UK three pin plug
- Electronic thermostat controls temperature range
- Temperature settings (low 82°C, medium 93°C and high 104°C)
- 0 - 30 minutes timer control



ORDERING INFORMATION

Part Number	Description
800-OVEN-EURO	OVEN, EXOS, 230V, EURO



Oven Extender Ring Kit

INDICATIONS

The thermal oven Extension Ring Kit is intended to be used for heating larger braces or brace components, so the heating element does not come in direct contact with any portion of the product while being heated.

CHARACTERISTICS

- 8.5cm extender ring
- Includes one thermal wrap for improved heat retention within the oven



ORDERING INFORMATION

Part Number	Size
800-ERK	U



Edging Tape

INDICATIONS

Use Edge Tape to cover cut edges of Exos products. It can also be used on areas of irritation to soften hard edges.

CHARACTERISTICS

- Provides a finished look and increases patient comfort.
- Available in Lycra or UBL.



ORDERING INFORMATION

Part Number	Description
800-02-1111	EXOS™ Edging tape Lycra-Black 1,9 x 30,48 cm
800-1111	EXOS™ Edging tape UBL-Black 2,54 x 60,96 cm



RespiForm® Wrist

DONJOY®

INDICATIONS

Wrist sprains and trauma. Tendinopathy. Pre- and post-operative immobilization

CHARACTERISTICS

- Dedicated for either right or left hand use
- Ergonomic, flexible lightweight design
- Anatomically designed, provides tailored individual fit
- Easy and simple application via adjustable straps and "Stop Strap" system
- Incorporates two semi-rigid dorsal stays and an aluminum palmar stay for overall wrist immobilisation
- Advanced Ozone™ material, promotes air flow
- Breathable in general use, providing good skin care and around-the-clock comfort in different heat and humidity conditions

ORDERING INFORMATION

Part Number		Wrist Circumference	Size
Left	Right		
82-96160	82-96150	11 - 14 cm	Pediatric/XS
82-96162	82-96152	14 - 16 cm	S
82-96163	82-96153	16 - 19 cm	M
82-96164	82-96154	19 - 21 cm	L
82-96165	82-96155	21 - 24 cm	XL



RespiForm+® Wrist

DONJOY®

INDICATIONS

Wrist & thumb sprains and trauma. Tendinopathy, rhizarthrosis. Pre- and post-operative immobilization

CHARACTERISTICS

- Ergonomic, flexible lightweight design
- Anatomically designed, provides tailored individual fit
- Easy and simple application via adjustable straps and "Stop Strap" system
- Incorporates two semi-rigid dorsal stays and an aluminum palmar stay for overall wrist immobilisation
- Advanced Ozone™ material, promotes air flow
- Breathable in general use, providing good skin care and around-the-clock comfort in different heat and humidity conditions
- Modular designed malleable aluminum thumb stay for optimum positioning
- Hook and loop system for easy and secure distal thumb fitting

ORDERING INFORMATION

Part Number		Wrist Circumference	Size
Left	Right		
82-96180	82-96170	11 - 14 cm	Pediatric/XS
82-96182	82-96172	14 - 16 cm	S
82-96183	82-96173	16 - 19 cm	M
82-96184	82-96174	19 - 21 cm	L
82-96185	82-96175	21 - 24 cm	XL



DigiForm® Wrist

DONJOY®

INDICATIONS

All indications requiring immobilization of three or four fingers. Sprains, subluxations, dislocations, fractures. Pre- / post-surgery wrist immobilization.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface
- Contact closure straps

ORDERING INFORMATION

Part Number	Wrist Circumference	Size
82-96200	11-15 cm	Pediatric/XS
82-96202	15-19 cm	S/M
82-96203	19-23 cm	L/XL



DigiForm+® Wrist

DONJOY®

INDICATIONS

All indications requiring immobilization of three or four fingers. Sprains, subluxations, dislocations, fractures. Pre- / post-surgery wrist immobilization

CHARACTERISTICS

- Easy application via four adjustable straps and "Stop Strap" system
- Fits either right or left hand
- Lightweight, low profile design
- Prepositional aluminum palmar stay provides tailored support
- Removable modular components for ease of washing
- Finger tips protective cover
- Splint adjustable for combinations of three fingers, or for all four fingers
- Breathable 3D material
- Available in five sizes

ORDERING INFORMATION

Part Number	Wrist Circumference	Size
82-96210	11 - 14 cm	Pediatric/XS
82-96212	14 - 16 cm	S
82-96213	16 - 19 cm	M
82-96214	19 - 21 cm	L
82-96215	21 - 24 cm	XL

PolyForm® Wrist

INDICATIONS

Immobilization of the wrist/thumb/fingers as part of pain management
Pain relief associated with rheumatoid arthritis



CHARACTERISTICS

- Supportive and comfortable thumb-immobilization piece
- Easy application via five adjustable straps and "Stop Strap" system
- Removable straps for thumb and finger support
- Adjustable aluminum dorsal stay for improved comfort
- Conformable hand stay for optimal adjustment
- Lightweight (180g), low profile design
- Soft, breathable, latex-free material
- Removable modular components for ease of washing
- Fits either right or left hand
- Available in five sizes (XS to XL)

ORDERING INFORMATION

Part Number	Wrist Circumference		Size
	CM	Inches	
82-96300	11 - 14cm	4 ¼" - 5 ½"	XS
82-96302	14 - 16cm	5 ½" - 6 ½"	S
82-96303	16 - 19cm	6 ½" - 7 ½"	M
82-96304	19 - 21cm	7 ½" - 8 ½"	L
82-96305	21 - 24cm	8 ½" - 9 ½"	XL

CarpaForm® Wrist

INDICATIONS

Carpal tunnel syndrome.
Tendinopathy, Teno-synovitis



CHARACTERISTICS

- Ergonomic, ultra-lightweight design
- Supportive exoskeleton design for decompression of carpal nerve
- Malleable aluminum structure provides tailored support
- Supports the wrist in neutral position, can be used at night
- Dedicated for either right or left hand use
- Easy and simple to adjust and fit
- Available in three sizes

ORDERING INFORMATION

Part Number		Wrist Circumference	Brace Length	Size
Left	Right			
82-96232	82-96222	13 - 16.5 cm	19cm	S
82-96233	82-96223	16.5 - 19.5 cm	19cm	M
82-96234	82-96224	19.5 - 25 cm	20cm	L

AVAILABLE ON REQUEST

ErgoForm® Wrist

INDICATIONS

All indications requiring immobilization of the thumb joints. Indicated for sprains of the thumb, subluxations, dislocations, rhizarthrosis



CHARACTERISTICS

- Ergonomic design
- Universal model, fits either right or left hand
- Easy application via adjustable strap and "Stop Strap" system
- Malleable aluminum thumb stay for optimum positioning
- Distal strap to secure thumb
- Advanced Ozone™ material, promotes air flow
- Breathable in general use, providing good skin care and around-the-clock comfort in different heat and humidity conditions
- Anatomically designed, providing tailored individual fit
- Lightweight, low profile design
- Available in two sizes


ORDERING INFORMATION

Part Number	Wrist Circumference	Size
82-96191	14 - 18 cm	T1
82-96192	18 - 22 cm	T2

RhizoForm® Wrist

INDICATIONS

Indicated for rhizarthrosis of the thumb



CHARACTERISTICS

- Ergonomic universal design, fits both right and left hands
- Ultra-lightweight supportive exoskeleton structure
- Malleable structure provides tailored support
- onSteam® breathable microfiber lining promotes comfort
- Easy and simple to both adjust and fit
- Available in two sizes

ORDERING INFORMATION

Part Number		Wrist Circumference	Size
Left	Right		
82-96292	82-96242	13.5-16.5 cm	S/M
82-96293	82-96243	16.5-20 cm	L/XL

ComfortFORM™ Wrist

INDICATIONS

Treatment of sprains, strains, and control of carpal tunnel syndrome symptoms.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface
- Contact closure straps



ORDERING INFORMATION

Part Number		Wrist Circumference	Brace Length	Size
Left	Right			
81-97892	81-97882	11.5cm - 14cm	19cm	XS
81-97893	81-97883	14cm - 16.5cm	19cm	S
81-97895	81-97885	16.5cm - 19cm	20cm	M
81-97897	81-97887	19cm - 21.5cm	21.5cm	L
81-97898	81-97888	>21.5 cm	21.5cm	XL

ComfortFORM™ Wrist & Thumb

INDICATIONS

Treatment of sprains, strains, and control of carpal tunnel syndrome symptoms.

CHARACTERISTICS

- Durable, lightweight Lycra lined foam laminate for breathability and patient comfort
- Dual mouldable aluminium stays support the palmar surface
- Contact closure straps



ORDERING INFORMATION

Part Number		Wrist Circumference	Wrist Length	Size
Left	Right			
81-87312	81-87302	< 14 cm	17.5 cm	XS
81-87313	81-87303	14-16 cm	17.5 cm	S
81-87315	81-87305	16-19 cm	17.5 cm	M
81-87317	81-87307	19-21 cm	17.5 cm	L
81-87318	81-87308	>21 cm	17.5 cm	XL

Quick Fit Wrist II

INDICATIONS

The Quick-Fit Wrist II is a universal, functional wrist orthosis for the treatment of wrist and hand injuries and trauma.

CHARACTERISTICS

- Simple, single pull lace design for ease of application
- Soft nylon/ foam laminate materials offer maximum comfort
- Adjustable, malleable, removable dorsal stay for efficient immobilization
- Unique contoured palmar stay provides proper anatomical fit and can be contoured as needed to address a variety of indications
- Universal size, available in right or left



ORDERING INFORMATION

Part Number		Wrist Circumference	Length	Size
Left	Right			
79-87560	79-87570	< 25,4 cm	20 cm	U
79-87561	79-87571	> 25,4 cm	20 cm	U/XL

Quick Fit WTO

INDICATIONS

Ideal for addressing de Quervain's syndrome, Gamekeeper's thumb, Scaphoid injuries, sprains and strains.

CHARACTERISTICS

- Comfortable, lightweight nylon/foam laminate construction
- Adjustable dorsal stay to accommodate various sized patients
- Pre-contoured aluminum palmar and thumb stays may be adjusted for proper angle to address wrist and thumb pathologies
- Elastic straps with D-ring closure provide compressive support and comfort



ORDERING INFORMATION

Part Number		Wrist Circumference	Size
Left	Right		
79-87480	79-87490	< 33 cm	U
79-87591	79-87581	> 30 cm	U/XL



Frog Splints / Baseball Splints

PROCARE®

INDICATIONS

Helps maintain phalangeal joints in position without bulky dressing.

CHARACTERISTICS

- Padded aluminum splint.
- Easily contours for proper fit.
- Foam Padded.
- 12 pack.

ORDERING INFORMATION

Frog Splints		
Part Number	Length	Size
79-71963	7 x 6 cm	S
79-71965	6.5 x 6.5 cm	M
79-71967	9 x 7.5 cm	L



BaseBall Splints		
Part Number	Length	Size
79-72093	10 cm	S
79-72095	11.4 cm	M
79-72097	12.7 cm	L

PROCARE®

10" Lace up Wrist

INDICATIONS

Ideal for wrist sprains, strains, and Carpal Tunnel Syndrome.

CHARACTERISTICS

- Perforated suede/flannel construction.
- Single pull laced closure for easy application.
- Length: 10" length.



ORDERING INFORMATION

Part Number		Wrist Circumference	Size
Right	Left		
79-87222	79-87232	≤14 cm	XS
79-87223	79-87233	14 - 16,5 cm	S
79-87225	79-87235	16,5 - 19 cm	M
79-87227	79-87237	19 - 21,5 cm	L
79-87228	79-87238	21,5 cm +	XL

PROCARE®

Finger Splints

INDICATIONS

Repair of ruptured tendons, Fingertip and nail bed injuries and mallet Finger.

CHARACTERISTICS

- Lightweight plastic splints
- Available in clear or skin tone
- A pack of padded finger straps is available separately for securing the splints
- Starter set of clear finger splints available (various sizes)



ORDERING INFORMATION

Part Number	Length	Diameter	Quantity	Size
79-72251	5.3cm	1.5cm	Pack of 12	1
79-72252	5.2cm	1.7cm	Pack of 12	2
79-72253	5.8cm	1.9cm	Pack of 12	3
79-72254	6.3cm	1.9cm	Pack of 12	4
79-72255	5.9cm	1.9cm	Pack of 12	5
79-72256	6.1cm	2.3cm	Pack of 12	5.5
79-72257	6.4cm	2.3cm	Pack of 12	6
79-72258	6.3cm	2.5cm	Pack of 12	7


Part Number	Description	Quantity
79-72260	Clear Splint Starter Set	Pack of 30
79-72250	Beige Starter Set	Pack of 30



**BACK & NECK
SUPPORTS**

 **DJO**® POWERING MOTION™

Product Advisor Chart





Exos FORM™ II



Exos Form II Spine Braces


Exos Form II is designed with patented thermoformable VST-X™ material, Comfort CORE™ Foam that contours to the patient's body, Diamond Grid construct for a more durable and longer-wearing brace and Boa® Technology's advanced closure system.

The end result is a spine brace that delivers superior decompression and support, a quick and easily customizable fit, and provides patient comfort and satisfaction by reshaping form, fit and function.

FEATURES AND BENEFITS


With a simple, lightweight, low profile design, the Exos Form™ II range is designed with patient comfort in mind without sacrificing functionality Adjustable (angle) belt wings accommodate varying patient body structures for optimal fit

The modular step-up/step-down design enables you to adjust and configure your brace as you progress through rehabilitation and recovery A semi-universal design minimizes inventory management requirements



COMMON INDICATIONS

	SACRUM	L1 - L5	T9 - S1	T9 - S1
	621	626	631	637
Sacroiliac Pain	●			
Acute and Chronic Lower Back Pain		●	●	●
Lower Back Sprains / Strains		●	●	●
Lumbar Disc Displacement		●	●	●
Osteoporosis		●	●	●
Disc Herniation and Degeneration		●	●	●
Post-Operative Laminectomy		●	●	●
Post-Operative Disectomy		●	●	●
Spondylolisthesis		●	●	●
Post-Operative Fusion			●	●
Spinal Stenosis			●	●
Stable / Non-Displaced Fractures			●	●



AVAILABLE ON REQUEST

Exos FORM™ II 637

Provides a higher degree of spinal support and relief for a wide range of indications such as: acute and chronic back pain, post-operative support, activity related back pain, disc herniation, stable, non-displaced fractures.

Removable rigid anterior and posterior panels combined with lateral support panels, provide superior posterior support from T9-S1.



ORDERING INFORMATION

Part Number	Description	Size
300621-40	Exos™ Form™ II 621	S/M
300621-60	Exos™ Form™ II 621	L/XL
300631-40	Exos™ Form™ II 631	S/M
300631-60	Exos™ Form™ II 631	L/XL
300626-40	Exos™ Form™ II 626	S/M
300626-60	Exos™ Form™ II 626	L/XL
300637-40	Exos™ Form™ II 637	S/M
300637-60	Exos™ Form™ II 637	L/XL

Exos FORM™ II 631

The Exos Form II 631 provides a higher degree of spinal support and relief from a wide range of indications from acute lower back pain to post-operative support. The Exos Form II 631 provides superior sagittal control and support from T9-S1.



Accessories		
Part Number	Description	Size
900-631-40	Posterior T9 Panel	S/M
900-631-60	Posterior T9 Panel	L/XL
900-637-40	Lateral Extensions	S/M
900-637-60	Lateral Extensions	L/XL
900-AP-40	Anterior Panel	S/M
900-AP-60	Anterior Panel	L/XL
900-1-WINGS-40	Belt Wing Set	S/M
900-1-WINGS-60	Belt Wing Set	L/XL
900-1-BWING-40	Bridge Wing	S/M
900-1-BWING-60	Bridge Wing	L/XL
900-1-INSERT-40	Lordotic Inserts	S/M
900-1-INSERT-60	Lordotic Inserts	L/XL
900-JITT	SI Belt Pad Set	Universal
900-POSTPAD-40	Posterior T9 Panel Padding	S/M
900-POSTPAD-60	Posterior T9 Panel Padding	L/XL
900-ANTPAD-40	Anterior Panel Padding	S/M
900-ANTPAD-60	Anterior Panel Padding	L/XL
900-HOTCOLD-40	Hot/Cold Pack Accessory	S/M
900-HOTCOLD-60	Hot/Cold Pack Accessory	L/XL

Exos FORM™ II 626

Provides relief for chronic back pain, post-operative support, activity related back pain, disc herniation. Superior sagittal control and support from L-1 to L-5.

Features removable customisable inserts for additional support.



Exos FORM™ II 621

Provides stabilisation of the sacroiliac (SI) joint, in case of sacroiliac pain and dysfunction.

Superior stabilization of the sacroiliac (SI) joint through direct compression. Silicon grip strips prevent the brace from migrating.



Sizing Information - Waist Circumference

Description	Exos FORM™ II 637	Exos FORM™ II 626	Exos FORM™ II 621
S/M	63cm - 125cm	63cm - 125cm	< 89cm
L/XL	125cm - 152cm	125cm - 152cm	91cm - 140 cm

NEW



DonJoy Back Brace II TLSO

INDICATIONS

Control and relieve acute pain often associated with anomalies of the spine such as osteoporosis, compression fractures, spinal stenosis, strain, and excessive kyphosis.

CHARACTERISTICS

- Padded shoulder straps with quick release buckles enhance patient comfort and compliance
- Mechanical advantage pulley system allows for easy thoracic adjustment
- Midline cutout avoids contact with surgical incision
- Rigid overlapping side panels for lateral stability
- Rigid anterior plate helps to resist flexion

ORDERING INFORMATION

Part Number	Waist Circumference	Back Panel Height	Size
11-1585-2	76cm - 89cm (30" - 35")	43cm	S
11-1585-3	89cm - 101cm (35" - 40")	43cm	M
11-1585-4	101cm - 114cm (40" - 45")	48cm	L
11-1585-5	114cm - 127cm (45" - 50")	48cm	XL
11-1585-6	127cm - 140cm (50" - 55")	48cm	XXL



Immostrap™

INDICATIONS

Acute or subacute lower back pain, lumbar arthrosis, lumbar sciatica, post disc herniation.

CHARACTERISTICS

- Quadruple bilateral straps to adjust compression level
- Independent superior and inferior compression adjustment for optimal anatomical fit
- Ergonomic front panel and a strong compressive elastic material enhances posture control
- Four semi rigid dorsal stays for stabilisation of the sacro-lumbar area
- Rolling buckles for effortless tightening
- Brace height 26cm (10")

ORDERING INFORMATION

Part Number	Waist Circumference	Size
82-0229-1	56 cm - 68 cm	XS
82-0229-2	68 cm - 82 cm	S
82-0229-3	82 cm - 98 cm	M
82-0229-4	98 cm - 116 cm	L
82-0229-5	116 cm - 136 cm	XL



DONJOY®
Bodystrap




INDICATIONS

Acute, sub-acute or chronic lower back pain. Return to daily living activities, lightweight tasks.

CHARACTERISTICS

- Strong compressive elastic material
- 4 Flexible dorsal stays
- Bilateral elastic straps
- Finger loop for ease of application
- Ergonomic front panel provides optimal comfort while sitting or standing
- Low profile design



ORDERING INFORMATION

Part Number	Waist Circumference	Size
82-0260-1	56 cm - 68 cm	XS
82-0260-2	68 cm - 82 cm	S
82-0260-3	82 cm - 98 cm	M
82-0260-4	98 cm - 116 cm	L
82-0260-5	116 cm - 136 cm	XL



Actistrap

DONJOY®

INDICATIONS

Acute, subacute or chronic lower back pain, lumbar sciatica.

CHARACTERISTICS

- Strong compressive elastic material to enhance posture control with double bilateral straps to adjust compression level
- Four rigid dorsal stays aid stabilisation of the sacro-lumbar area and lordosis correction
- Anatomical design for improved comfort and fit
- Two flexible abdominal stays
- Hand loop for simple application with finger pockets for ease of adjustment
- Brace height 26cm (10")

ORDERING INFORMATION

Part Number	Waist Circumference	Size
82-0231-1	56 cm - 68 cm	XS
82-0231-2	68 cm - 82 cm	S
82-0231-3	82 cm - 98 cm	M
82-0231-4	98 cm - 116 cm	L
82-0231-5	116 cm - 136 cm	XL
82-0231-6	136 cm - 156 cm	XXL



Abdostrap II

DONJOY®

INDICATIONS

Abdominal support for temporary abdominal weakness, post-partum, post-op use following abdominal surgery.

CHARACTERISTICS

- Moderate compressive elastic material
- 2 Flexible abdominal stays for increased abdominal support
- Anterior panel: trimmable textile for probe or catheter insertion
- 2 Flexible dorsal stays
- Removable additional abdominal strap included for enhanced support

ORDERING INFORMATION

Part Number			Waist Circumference	Size
18 cm*	25 cm	32 cm*		
82-0256-1	82-0257-1	82-0258-1	56 cm - 68 cm	XS
82-0256-2	82-0257-2	82-0258-2	68 cm - 82 cm	S
82-0256-3	82-0257-3	82-0258-3	82 cm - 98 cm	M
82-0256-4	82-0257-4	82-0258-4	98 cm - 116 cm	L
82-0256-5	82-0257-5	82-0258-5	116 cm - 136 cm	XL

***AVAILABLE ON REQUEST**

MyBabystrap™ Fusion

NEW



INDICATIONS

Pain and discomfort associated with pregnancy

CHARACTERISTICS

- Modular design
- Adjustable back pad
- Four posterior rigid stays
- Bilateral straps
- Anatomical shape
- Stylish maternity belt
- Thin, soft and breathable fabric
- Ergonomic front panel
- Easy-to-use with simple, hook-loop closure
- Removable elastic strap

ORDERING INFORMATION

Part Number	Waist Circumference
82-0276	80 - 150 cm

ComfortFORM™ Back Support



INDICATIONS

Support for lower back injuries. Prevents spinal overload during daily activity.

CHARACTERISTICS

- Durable mesh elastic support with moulded back compression pad
- Low profile design for comfort and fit
- Double pull straps to promote spinal alignment

ORDERING INFORMATION

Part Number	Waist Circumference	Size
81-89352	51 cm - 64 cm	XS
81-89353	64 cm - 76 cm	S
81-89355	76 cm - 86 cm	M
81-89357	86 cm - 97 cm	L
81-89358	97 cm - 107 cm	XL
81-89359	107 cm - 142 cm	XXL
81-89359-10	142 cm - 158 cm	XXXL



Premium Panel Elastic Binder

PROCARE®

INDICATIONS

Compression and support for post natal or post abdominal surgery, abdominal strains and weakness.

CHARACTERISTICS

- Constructed with panel elastic and soft cotton type lining for patient comfort
- Universal style support with contact closure
- Multi panel design helps prevent panel rolling
- Available in 22.5cm (9") 3 panel and 30cm (12") 4 panel heights

ORDERING INFORMATION

Part Number		Size	Waist Circumference
9" 3-Panel Binder	12" 4-Panel Binder		
79-89070	79-89090	S	76-114 cm (30-45")
79-89071	79-89091	M	114-157 cm (45-62")
79-89210	79-89220	L	157-188 cm (62-74")
79-89270	79-89280	XL	188-213 cm (74-82")
79-89250	79-89260	XXL	213-254 cm (82-100")



Hernia Belt

DONJOY®

INDICATIONS

Provides adjustable, gradual compressive support to help reduce inguinal hernias and alleviate discomfort.

CHARACTERISTICS

- Easy to adjust, soft waist band and perineal straps provide comfortable compression
- Removable foam compression pads for gentle, gradual compression
- Treats right, left or double hernias
- Comfortable, lightweight nylon material
- Latex free

ORDERING INFORMATION

Part Number	Size	Waist Circumference
81-97903	S	70 cm - 80 cm
81-97905	M	80 cm - 95 cm
81-97907	L	95 cm - 110 cm
79-89298	XL	110 cm - 120 cm

AVAILABLE ON REQUEST



Form Fit Collar

PROCARE®

INDICATIONS

This collar provides semi-rigid support for use in soft tissue damage.

CHARACTERISTICS

- Contoured for anatomical fitting
- Comfortable and secure fit

ORDERING INFORMATION

Part Number		Height x Length	Neck Circumference	Size
Firm-Density	Med-Density			
79-83003	79-83013	7.5 x 47 cm	11 - 16 cm	S
79-83005	79-83015	10 x 51 cm	13 - 18 cm	M
79-83007	79-83017	11.5 x 57 cm	15 - 20 cm	L

AVAILABLE ON REQUEST



California Collar

PROCARE®

INDICATIONS

Ideal for cervical spine immobilization and extraction.

CHARACTERISTICS

- Pre-formed, lightweight, two-piece polyethylene foam with rigid plastic supports front and rear
- Collar maintains head in neutral position
- May be used to provide support in combination with a cervical or scoliosis cast

ORDERING INFORMATION

Part Number	Height	Neck Circumference	Size
79-83123**	5,7 cm	25,5 - 33 cm	S
79-83125**	5,7 cm	33 - 40,5 cm	M
79-83127**	5,7 cm	40,5 - 48 cm	L
79-83133**	8,3 cm	25,5 - 33 cm	S
79-83135**	8,3 cm	33 - 40,5 cm	M
79-83137**	8,3 cm	40,5 - 48 cm	L
79-83143**	10,8 cm	25,5 - 33 cm	S
79-83145**	10,8 cm	33 - 40,5 cm	M
79-83147**	10,8 cm	40,5 - 48 cm	L
79-83153**	13,3 cm	25,5 - 33 cm	S
79-83155**	13,3 cm	33 - 40,5 cm	M
79-83157**	13,3 cm	40,5 - 48 cm	L

Height: sternal notch to chin **Items in Stock ***Available on request

AVAILABLE ON REQUEST

A male athlete is captured in a dynamic running pose on a green field. He is wearing a grey t-shirt, dark shorts, white socks, and white sneakers with blue accents. A grey and white knee brace is visible on his right knee. The background is a blurred stadium setting. The text 'KNEE LIGAMENT BRACING' is overlaid on the left side of the image.

KNEE LIGAMENT BRACING

4-Points-of-Leverage™ Bracing System

The 4-Points-of-Leverage is a dynamic bracing system made up of a rigid cuff and strap configuration, whereby a posterior force is applied to the tibia, preventing anterior movement. This reduces the strain on the anterior cruciate ligament (ACL).

The 4-Points-of-Leverage brace design:

- Significantly reduces ACL strain*
- Provides a constant dynamic load



Point of Leverage #1

The tibia is anchored by a cuff placed at the bottom of the posterior calf

Point of Leverage #2

The femur is anchored by a cuff placed at the top of the anterior thigh

Point of Leverage #3

A strap across the back of the lower thigh pushes the femur anteriorly

Point of Leverage #4

A strap placed on the anterior tibia applies an active constant load to prevent anterior tibia translation



**Fleming et al. (2000), The influence of functional Knee Bracing on the ACL Strain Biomechanics in Weightbearing and Nonweightbearing Knees, The American Journal of Sports Medicine, Vol., Nr. 6, p. 815-824*

Measuring Guide

For DonJoy knee braces throughout the catalogue when not otherwise specified on the product page



Sizes Available:

XS-XXXL Measurements taken 15cm above mid-patella

Part number code "x"	Size	Circumference (cm)	Circumference (Inches)
1	XS	33cm - 39cm	13" - 15½"
2	S	39cm - 47cm	15½" - 18½"
3	M	47cm - 53cm	18½" - 21"
4	L	53cm - 60cm	21" - 23½"
5	XL	60cm - 67cm	23½" - 26½"
6	XXL	67cm - 75cm	26½" - 29½"
7	XXXL	75cm - 81cm	29½" - 32"

How the DonJoy® 4-Point Brace Protects the ACL Deficient Knee


Patients with an ACL deficient knee have an anterolateral instability that causes the knee to give way.

It is very difficult to control instability in an ACL deficient knee in a direct way using a knee brace.

Indirect stabilisation of the knee is achieved with the biomechanical 4 point design of the DonJoy brace.

The 4-point system is clinically proven to be the most effective way to prevent anterior translation of the tibia.

A DonJoy brace will provide a high level of stability for the ACL deficient knee by effectively reducing or preventing the anterior tibial translation through the use of the 4-points-of-leverage.



FourcePoint™ Hinge Technology

Aggressive eccentric contractions of the quadriceps produce significant anterior tibial translation and can injure the ACL or permanently stretch out an ACL graft (DeMorat et al., AJSM, 2004). Additional studies show quadricep-intensive activities should be avoided at joint angles, particularly the last 30 degrees of extension, when strains are at a maximum. DJO has developed a hinge to address maximising ACL and ACL graft protection in these "at risk" positions.

Small knee flexion angles at initial contact during landing tasks may contribute to ACL injuries. The FourcePoint Hinge can be used to control joint range-of-motion by increasing the knee flexion angle during landing, resulting in a decrease of anterior shear force on the ACL by 9-13% (Garret and Yu, AJSM, 2004).

Flexion



Extension



The DonJoy FourcePoint utilises a spring mechanism in the hinge to apply gradually increasing resistance during knee extension. The resistance engages in the last 25 degrees of extension (relative to the extension stop).

The hinge has three levels of resistance with the option to turn off resistance if appropriate.

The resistance serves three critical roles:

- First, it reduces the time spent near full extension or in the "at risk" position
- Second, it increases the posterior load applied by the brace on the tibia to prevent anterior tibial translation
- Third, it eliminates the extension shock felt when a patient extends into a 10 degree standard rigid stop

FourcePoint Hinge Technology works to enhance DonJoy's 4-Points-of-Leverage design by damping knee joint extension, which improves the mechanical performance of the brace and reduces shear forces at the knee. Addressing the "at-risk" position, the hinge resistance engages in the last 25 degrees of knee extension.

The FourcePoint hinge combined with the 4-Points-of-Leverage cuff and strapping design provides a more comfortable brace that reduces anterior shear forces at the knee. This is beneficial for patients as it reduces strain on the healing ACL graft, providing stability to an unstable, ACL deficient knee or for prophylactic use.



DONJOY®
Defiance™



FourcePoint hinge technology

DonJoy Defiance sets the industry standard for custom made knee braces. A host of options allows customisation of the fit and functions, to ensure maximum protection and confidence. It has become the market leader in technology, materials and design. Defiance incorporates DonJoy's patented 4-Points-of-Leverage for effective support, and offers a choice of patented hinge designs.

INDICATIONS


Moderate to severe ACL, PCL, CI, MCL and LCL instabilities. The OA Defiance is suitable for the active patient with osteoarthritis.

CHARACTERISTICS

- Custom fabricated, strong, lightweight, low profile carbon composite frame - up to 7 degrees of correction can be built into the frame for osteoarthritis correction
- Anti-migration band and supra-condyle pad for enhanced brace suspension
- Flexion and extension control
- Minimum 5 years warranty on frame and hinges, 6 months on soft goods

OPTIONS

- Hinge options; FourcePoint, Standard Polycentric, DropLock or Adjustable OA hinges - adjustable OA hinge allows for 0 - 3° of "fine tuning" load correction (see page 10)
- Strap configurations; ACL, PCL, or Combined Instability (CI)
- Frame options; swooping thigh cuff, short thigh cuff, short calf cuff
- Over forty frame colours and graphic options (see page 10)
- Additional layer of reinforced carbon composite in the frame for maximum strength



- Swivelling strap tabs for comfort or PowerCaps (metal strap tab caps recommended for contact sports)
- High strength stainless steel gears (recommended for contact sports)
- Chamois or suede liners and strap pads
- The Lifetime Warranty on Frame and Hinges now applies only to Braces with FourcePoint Hinges & Steel High Strength Gears only.
- 5 Year Warranty on Aluminum Geared Hinges (Standard & OA)
- OA and Standard hinges can be upgraded to Steel High Strength Gears with a \$20 up-charge from aluminum gears to steel = Buying a lifetime warranty

ORDERING INFORMATION





Size	Thigh Circumference	Calf Circumference
Minimum	31.8 cm	26 cm
Maximum	97 cm	66.5 cm

Each Defiance is custom made. Measurements taken with CCMI Mark III Individual Custom Order. 2-3 weeks waiting period.

Defiance™ Colour Options



Defiance™ Hinge Options

<p>STANDARD POLYAXIAL</p>  <p>ACL/PCL/CI</p> <p>DESCRIPTION Polycentric hinge motion stabilises the knee while matching natural knee kinematics. Polycentric hinges also maximise the anterior/posterior and varus/valgus stabilising forces of the brace.</p> <p>INDICATION Patient type: ACL, PCL, Combined, MCL and LCL</p>	<p>DROP LOCK LOCKING</p>  <p>M. QUADRICEPS DEFICIENCY</p> <p>DESCRIPTION The DropLock polycentric hinge is a spring-loaded system which automatically engages the locking device when the patient fully extends their leg. The DropLock hinge has an easy push button release to unlock, allowing free ROM for assisted gait. There are no additional Flexion or Extension stops available with DropLock hinge design.</p> <p>INDICATION Patient type: Symptoms and/or physical findings of m. quadriceps deficiencies, standing assist, ost-polio syndrome</p>	<p>FOURCEPOINT DAMPENING</p>  <p>ACL ONLY</p> <p>DESCRIPTION This patented hinge design incorporates a spring leaf and dampening mechanism that provides increasing resistance in the last 25° of extension to reduce extension shock and increase posterior tibial loads. The FourcePoint hinge is shipped as standard with a 10° extension stop to engage resistance at 25°.</p> <p>INDICATION • The only hinge clinically proven to protect the ACL • Reduces time spent in "at risk" position • Not indicated for PCL instabilities</p>	<p>OA UNLOADING</p>  <p>OA</p> <p>DESCRIPTION This polyaxial hinge comes with a telescoping condyle pad that allows patient-controlled additional unloading on the base load of 1 to 7° built into the frame. A condyle pad is included for increased patient comfort and compliance.</p> <p>INDICATION • Patient type: Unicompartmental OA with or without combined ACL, PCL, Combined, MCL and LCL instabilities and activities of daily living</p>
---	--	---	--



Armor™ CI

DONJOY®

INDICATIONS

Moderate to severe ACL, MCL and LCL instabilities. For active lifestyles involving extreme sports e.g. skiing, motocross, contact sports.

CHARACTERISTICS

- Enhanced 4-Points-of-Leverage with ACL strapping
- Available in 7 sizes, XS - XXXL
- Low profile FourcePoint hinge with steel outer plate with optional flexion and extension control
- Mouldable, medial swooping frame for custom fit
- Short calf for boot clearance
- Strong aircraft quality aluminium frame with internally mounted strap tabs
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

Part Number		Thigh Circumference	Size
Left	Right		
11-1030-1	11-1029-1	33cm - 39cm	XS
11-1030-2	11-1029-2	39cm - 47cm	S
11-1030-3	11-1029-3	47cm - 53cm	M
11-1030-4	11-1029-4	53cm - 60cm	L
11-1030-5	11-1029-5	60cm - 67cm	XL
11-1030-6	11-1029-6	67cm - 75cm	XXL
11-1030-7	11-1029-7	75cm - 81cm	XXXL
13-8925		FourcePoint Allen Key	

Measurement taken 15cm above mid-patella



FullForce™ with Force Point

DONJOY®

INDICATIONS

Moderate to severe ACL, MCL & LCL instabilities. For active lifestyles involving mid-range sports e.g. tennis, netball, cricket.

CHARACTERISTICS

- Enhanced 4-Points-of-Leverage with ACL strapping
- Available in 7 sizes, XS - XXXL
- Low profile FourcePoint hinge
- Strong, mouldable aircraft quality aluminium frame with carbon fibre look finish
- Internally mounted swivelling strap tabs to help accommodate leg movement
- Swooping thigh cuff for inner thigh clearance & bilateral use
- Optional short or standard calf length
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

Part Number		Description
Left	Right	
11-0265-x	11-0264-x	ACL, Standard Calf
13-8925		FourcePoint Allen Key

Measurement taken 15cm above mid-patella. Sizing information on page 6

*Accessories available on request



Playmaker Xpert

DONJOY®

INDICATIONS

For mild medial/lateral support of the knee. For use in warmer environments.

CHARACTERISTICS

- Advanced high tenacity nylon/spandex material promotes a cooling effect while delivering a lighter and slimmer profile
- Removable polycentric hinges and aluminium uprights
- Popliteal cut-out for extra comfort in flexion
- Fits left or right knee
- Available in 2 versions: Pull-On Sleeve and Wraparound model
- 6 month warranty

ORDERING INFORMATION

Part Number		Thigh Circumference	Size
82-0022-1	82-0023-1	33cm - 39cm	XS
82-0022-2	82-0023-2	39cm - 47cm	S
82-0022-3	82-0023-3	47cm - 53cm	M
82-0022-4	82-0023-4	53cm - 60cm	L
82-0022-5	82-0023-5	60cm - 67cm	XL
82-0022-6	82-0023-6	67cm - 75cm	XXL
82-0022-7	82-0023-7	75cm - 81cm	XXXL

Measurement taken 15cm above mid-patella



Quickfit Hinged Knee

DONJOY®

INDICATIONS

Protection, compression and support for knee sprains, mild to moderate ligament lesions, ligament instabilities, meniscus injuries, patella subluxation, laxity or misalignment, mild arthritic pain, return to mobility or sports activities.

CHARACTERISTICS

- Patella donut maintain patella in correct position.
- Donut can be trimmed to provide tailored patella misalignment correction.
- Removable aluminium polycentric hinges offer medial / lateral stability.
- Encased hinges for protection from lateral impacts (ideal for sports use).
- Soft, anti-pinch posterior Lycra panel enhances comfort.
- Breathable, form-fitting 5 mm spacer mesh material promotes airflow.
- Fit Left or Right knee.

ORDERING INFORMATION

Part Number	Thigh circumference	Size
82-0010-1	26,5 cm - 33 cm	XXS
82-0010-2	33 cm - 39 cm	XS
82-0010-3	39 cm - 47 cm	S
82-0010-4	47 cm - 53 cm	M
82-0010-5	53 cm - 60 cm	L
82-0010-6	60 cm - 67 cm	XL
82-0010-7	67 cm - 75 cm	XXL
82-0010-8	75 cm - 81 cm	XXXL

Measurement taken 15cm above mid-patella



Reddie Brace

DONJOY®

INDICATIONS

3 mm neoprene wraparound configuration and strap tabs provide ease of application, especially for the older patient. Posterior strap adjustment allows for one-time adjustment and proper positioning of hinges. Removable dual axis polycentric hinges provide medial/lateral support. Available with popliteal cut-out only.

CHARACTERISTICS

- Wraparound configuration and strap tabs for easy application
- Posterior strap adjustment for proper hinge positioning
- Removable Hinges
- Available in 7 sizes

ORDERING INFORMATION

Part Number	Measurement	Size
79-82392	34 cm - 39 cm	XS
79-82393	39 cm - 45,5 cm	S
79-82395	45,5 cm - 52 cm	M
79-82397	52 cm - 58,5 cm	L
79-82398	58,5 cm - 64,5 cm	XL
79-82399	64,5 cm - 71 cm	XXL
79-82399-10	71 cm - 78,5 cm	XXXL

Measurement taken 15cm above mid-patella



Stabilized Knee Support

PROCARE®

INDICATIONS

Ideal for mild patella tendonitis, chondromalacia or subluxations.

CHARACTERISTICS

- 1/8" neoprene with medial/lateral spiral stays for added support
- Superior and inferior compression straps help prevent migration and provide additional support
- Universal neoprene buttress may be adjusted to allow for proper patella management and stabilization

ORDERING INFORMATION

Part Number	Measurement	Size
79-82722	34 cm - 39 cm	XS
79-82723	39 cm - 45,5 cm	S
79-82725	45,5 cm - 52 cm	M
79-82727	52 cm - 58,5 cm	L
79-82728	58,5 cm - 64,5 cm	XL
79-82729	64,5 cm - 71 cm	XXL
79-82729-10	71 cm - 78,5 cm	XXXL

Measurement taken 15cm above mid-patella

KNEE OA





Measuring Guide

For DonJoy knee braces throughout the catalogue when not otherwise specified on the product page




Sizes Available:
XS-XXXL Measurements
taken 15cm above mid-patella

Part number code "x"	Size	Circumference (cm)	Circumference (Inches)
1	XS	33cm - 39cm	13" - 15½"
2	S	39cm - 47cm	15½" - 18½"
3	M	47cm - 53cm	18½" - 21"
4	L	53cm - 60cm	21" - 23½"
5	XL	60cm - 67cm	23½" - 26½"
6	XXL	67cm - 75cm	26½" - 29½"
7	XXXL	75cm - 81cm	29½" - 32"

Osteoarthritis (OA)

We are at a unique stage in our healthcare system where patients are being diagnosed with knee OA at a younger age, and there is increasing demand for conservative treatment options. OA is the fourth leading disability in the world with 8.75 million people in the UK seeking treatment for OA*. This combined with a growing active aging population sees the healthcare system not being able to continue to offer total knee replacements as a primary intervention. OA knee bracing is a conservative management alternative.

*Reference - Arthritis Research UK



Adjustable OA Defiance™ **DONJOY®**

INDICATIONS

The Adjustable OA Defiance™ is a custom-made brace designed for the active OA patient who does not want to compromise their high level of activity. The Adjustable OA Defiance is an extremely lightweight brace that provides durable support for mild to severe levels of osteoarthritis. The bilateral hinge bars offer support for ligament instabilities, as well as a patented telescoping condyle pad that puts the load control in the hands of the patient.

CHARACTERISTICS

- Custom-made for a perfect fit utilising the DonJoy CCMI measurement system
- Base level of corrective load built into the frame
- Patient controlled adjustable condyle applies 0° to 3° of load correction
- Lightweight, low-profile, carbon composite frame moulded to accommodate unique varus or valgus angulations
- Micro-sphere condyle pad for increased patient compliance
- Patented 4-Points-of-Leverage system (ACL, PCL, CI strapping options available)
- Extension stops: 0°, 10°, 20°, 30° and 40° (10° installed)
- Flexion stops: 45°, 60°, 75° and 90°
- Suede frame and strap pads standard on each brace
- Available in standard or short lengths.
- Arthritis-friendly features

ORDERING INFORMATION

Size	Thigh Circumference	Calf Circumference
Minimum	31.8 cm	26 cm
Maximum	97 cm	66.5 cm

Each Defiance is custom made. Measurements taken with CCMI Mark III Individual Custom Order. 2-3 weeks waiting period.



OA Nano™

DONJOY®

Combining Nano MAG™, a new hi-tech material, with DonJoy's proven adjuster technology, OA Nano is the world's lightest functional OA knee brace.

INDICATIONS

For mild to moderate unicompartmental OA, the OA Nano is ideal for patients who want to stay active and move freely.

CHARACTERISTICS

- The flexibility of the magnesium frame provides an intimate fit and increased comfort, while Nano MAG (on the loading side) provides strength for off-loading
- DonJoy's clinically proven three point off-loading technology shifts the compressive loads from the degenerative area to the healthy compartment
- 4 Points of Leverage and contoured frame design enhances brace suspension for all levels of activity
- Available in 7 sizes, XS - XXXL
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

Part Number		Affected Compartment
Left	Right	
11-1215-x	11-1214-x	Medial Compartment
11-1217-x	11-1216-x	Lateral Compartment

Note: Off-loading hinge is on opposite side of affected compartment. Sizing information on page 42

Part Number	Description
25-1376-0-00000	Stop kit, extension 20°, 30°, 45° flexion 60°, 75° and full lockout

DONJOY®

OA Reaction™ Web

INDICATIONS

Suitable for early to mild unicompartmental osteoarthritis of the knee & patella. OA Reaction Web is also a good option for patients wishing to use it as a lighter alternative to their usual rigid brace.

CHARACTERISTICS

- Four way stretch wrap-around design
- Silicone web technology
- Low profile, lightweight, discreet
- Single aluminium upright
- Silicone condyle pad
- Fitted undersleeve
- Available in 7 sizes, XS - XXXL
- 6 month warranty

ORDERING INFORMATION

Part Number	Affected Compartment
11-7426-x	Right Medial/Left Lateral
11-7427-x	Left Medial/Right Lateral

Note: Off-loading hinge is on opposite side of affected compartment. Sizing information on page 42



DONJOY®

OA Adjuster™ 3

INDICATIONS

For patients with moderate to severe unicompartmental OA, needing significant off-load and ligament stability in a tough durable frame.

CHARACTERISTICS

- Nano MAG™ technology on the off-loading side, combined with the durability of the aluminum frame, provides lightweight function with optimal strength
- 4 Points of Leverage system - provides clinically proven ligament stability combined with low-profile hinges and double upright design allow a contoured, personalised fit for all levels of activity
- Available in 7 sizes, XS - XXXL
- 1 year warranty on frame and hinge, 6 month warranty on soft goods

ORDERING INFORMATION

Part Number		Affected Compartment
Left	Right	
11-1591-x	11-1590-x	Medial Compartment
11-1593-x	11-1592-x	Lateral Compartment

Note: Off-loading hinge is on opposite side of affected compartment. Sizing information on page 42

Part Number	Description
25-1376-0-00000	Stop kit, extension 20°, 30°, 45° flexion 60°, 75° and full lockout





POST-OP

 **DJO**® POWERING MOTION™

NEW



X-ROM

INDICATIONS

Easy-to-use post-op ROM protection with fast application and comfortable stability.

CHARACTERISTICS

- Updated hinge design: Precisely control ROM with settings in 10° increments from -10° to 90° extension and -10° to 120° flexion
- Quick-lock design: Easily lock at 10° increments from -10° to 30°
- Easy-to-use design: Strap, tab and buckle design make post-op application and adjustment fast and simple
- Comfortably fits a wide range of patients: Easily accommodate patients from 4'11" to 6'5"
- Push button slider settings: Adjustments are easy and precise
- Ready to fit: Quickly fit patients right out of the box
- Optimal fit & stability: Bendable hinge bars help enhance varus/valgus contouring
- Secure contact points: Enhanced ankle pads and strategically placed contact points between the leg and brace help reduce slippage
- Lightweight strength: Strong yet lightweight design helps reduce brace migration
- Anti-microbial pads: Helps neutralize odor while increasing comfort and durability

ORDERING INFORMATION

Part Number	Description	Size
11-2151-9	X-ROM Post-op Knee Brace	U
25-7530	X-ROM Replacement Straps Kit	
25-7531	X-ROM Replacement Pad Kit	



X-Act™ ROM Lite

INDICATIONS

Designed for immobilisation and protected range of motion associated with post-traumatic treatment, ACL, PCL, LCL and MCL surgeries, patella realignment and meniscal or articular repairs.

CHARACTERISTICS

- Economy version of the X-Act ROM knee brace with easy to use strap system
- Adjustable uprights can be contoured to the patient's leg to provide varus and valgus support
- Extension settings from -10° to 90°
- Flexion settings from -10° to 120°
- Easy quick-lock at 0°, 15° and 30°
- Full circumferential foam
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

Part Number	Description	Size
11-2161-9	X-Act ROM Lite	U*

*Length adjustability: 46cm to 65cm *Maximum thigh circumference: 76cm

**INDICATIONS**

For warmer environments, the Cool IROM provides immobilisation and precise range of motion control for ACL, PCL, MCL and LCL surgeries.

CHARACTERISTICS

- 10° incremental range of motion control
- Mouldable hinge bars allow for varus/valgus positioning
- Hinge bar breakdown capability
- 4 straps and abbreviated foam at key control points
- Integrated strap/hinge bar design
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

Part Number	Description	Leg Inseam
11-0175-2-06000	Cool IROM Short	<80cm
11-0175-4-06000	Cool IROM Long	>80cm

DONJOY®

TROM™ Advance

INDICATIONS

Offering range-of-motion control and locking capabilities following ACL, PCL, MCL, and LCL surgeries, meniscal repairs, patella realignment, regenerative chondroplasty, stable femora fractures, total knee replacements and high tibial osteotomies.

CHARACTERISTICS

- Streamlined hinge with ROM control and quick lock from -10° to 90° to accommodate all post-operative protocols
- Telescopic sliding system to allow customised strap position above and below knee
- Adjustability of length from 45.5 cm to 58.5 cm to adapt to patient's leg length
- Quick-release buckle for fast and easy application
- Wide strap covers for additional support and comfort
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

Part Number	Description	Thigh Circumference	Size
11-9114-9	TROM Advance, Cool	<66cm	U
11-9116-9	TROM Advance, XL Cool	66cm - 81cm	XL

Length adjustability: 45cm to 60cm, measurements taken 15cm above mid-patella



DONJOY®

Immo AT4 3V

INDICATIONS

Pre- or post-operative total immobilization (knee ligament injuries, patella or meniscus injuries, knee surgery, emergency immobilization or delayed diagnosis).

CHARACTERISTICS

- Universal 3-panel design, adapts to patient thigh circumference.
- Perforated foam for enhanced breathability.
- Malleable & removable aluminium bilateral stays.
- Aluminium posterior stay allows for strict immobilisation.
- Anti-migration strap avoids brace slipping.
- Bilateral model, fits Left or Right leg.
- Available in 4 heights: 30, 40, 50, 60 cm.

ORDERING INFORMATION

Part Number	Brace Length	Size
82-43-30	30 cm	Ped.
82-43-01	40 cm	S
82-43-02	50 cm	M
82-43-03	60 cm	L
82-43-04	65 cm	XL

For correct height: measure the inseam of the leg and deduct 20 cm
3-Panel design adapts to thigh circumference





VersaROM™ Hip

DONJOY®

INDICATIONS

DonJoy® VersaROM™ Hip is indicated for use with post-op patients who have undergone a surgical procedure and need an orthosis until the component is integrated and soft tissue heals.

CHARACTERISTICS

- Integrated belt compression with cinching straps
- Malleable hip and thigh frames
- Precise ROM Adjustments
- Range of Motion allows controlled hip flexion and extension adjustments in 15° increments
- Telescoping Height Adjustment
- Laterally placed abduction paddle with three-point pressure
- Customized Fit
- Lightweight, slim-profile design

ORDERING INFORMATION

Part number	Description
11-3335	DonJoy® VersaROM™ Hip



Tibial Fracture Brace

DONJOY®

INDICATIONS

Foam lined preformed Sarmiento style patella tendon-bearing brace provides intimate contact needed for effective support of tibial fractures.

CHARACTERISTICS

- A bilateral upright foot component to add medial/lateral support as well as to prevent distal slippage
- Each brace comes complete with two (2) double-thick stockinettes
- Measure length from superior border of patella to mid-level malleoli

ORDERING INFORMATION

Part Number		Calf Circ.	Length	Size
Right	Left			
79-97703	79-97693	28 cm - 33 cm	33 cm - 42 cm	S
79-97704	79-97694	33 cm - 38 cm	38 cm - 51 cm	M/Long
79-97705	79-97695	33 cm - 38 cm	42 cm - 47 cm	M
79-97707	79-97697	38 cm - 46 cm	38 cm - 51 cm	L
79-97706	79-97696	38 cm - 46 cm	49.5 cm - 54.5 cm	L/Long

Measure length from superior border of patella to mid-level malleoli.


AVAILABLE ON REQUEST

**PATELLA
FEMORAL**



DONJOY®

Reaction™ Web



The New Reaction Knee Brace is an active, responsive, webbed approach to general anterior knee pain, and a distinct alternative to basic knee sleeves.

INDICATIONS

- General anterior knee pain
- Chondromalacia Patella
- Osgood-Schlatter Disease
- Quadriceps and Patella Tendonitis/Tendinosis

CHARACTERISTICS

- The elastomeric web design helps reduce pain by dispersing energy across the knee
- The flexible membranes absorb shock and shift peak loads away from the painful area
- The contoured fit, compression and stretch of the membranes provide increased proprioception
- The web design dynamically stabilises the patella on all sides
- Dual axis hinges for additional support
- Fitted under-sleeve for comfort
- Open framework provides intimate fit around the patella and captures the anatomic features of the knee
- The mesh backing creates a comfortable fit and a breathable solution
- The adjustable straps help the brace remain in place during daily or sports activities



ORDERING INFORMATION

Part Number	Thigh Circumference	Size
11-0215-2	33cm - 46cm	XS/S
11-0215-3	47cm - 59cm	M/L
11-0215-4	60cm - 75cm	L/XL
11-0215-5	>75cm	XXL/XXXL

Measurement taken 15cm above mid-patella

Accessories	Description	Size
11-0541-x	Replacement Undersleeve	S to XXXL



DONJOY®
Tru-Pull® Lite



The new Tru-Pull Lite™ is DonJoy's most comfortable, low profile and shortest patellofemoral brace, yet strong enough to relieve pain.

INDICATIONS

Recommended for acute mild to moderate anterior knee pain, support for patella instability, lateral subluxation and maltracking.

CHARACTERISTICS

- Patented dual durometer (silicone) buttress system provides consistent corrective force on the patella
- Bifurcated strap dynamically pulls the patella into proper alignment
- Semi-flexible removable hinges included
- 21 cm long support giving maximum comfort in a low profile design
- Fabric is hypoallergenic, breathable, and lightweight
- 6 month warranty



ORDERING INFORMATION

Part Number		Thigh Circumference	Size
Left	Right		
11-0261-1	11-0260-1	33cm - 39cm	XS
11-0261-2	11-0260-2	39cm - 47cm	S
11-0261-3	11-0260-3	47cm - 53cm	M
11-0261-4	11-0260-4	53cm - 60cm	L
11-0261-5	11-0260-5	60cm - 67cm	XL
11-0261-6	11-0260-6	67cm - 75cm	XXL
11-0261-7	11-0260-7	75cm - 81cm	XXXL

Measurement taken 15cm above mid-patella



Lateral "J"

DONJOY®

INDICATIONS

The tubular Lateral "J" buttress is ideal for supporting lateral patella subluxation, dislocations or lateral patellofemoral malalignment.

CHARACTERISTICS

- Lateral to medial strapping for additional stabilisation
- Drytex material for patients with allergies or sensitivity to neoprene
- Available in left or right configuration
- 6 month warranty

ORDERING INFORMATION

Part Number		Thigh Circumference	Size
Left	Right		
81-06602	81-06592	33cm - 39cm	XS
81-06603	81-06593	39cm - 47cm	S
81-06605	81-06595	47cm - 53cm	M
81-06607	81-06597	53cm - 60cm	L
81-06608	81-06598	60cm - 67cm	XL
81-06609	81-06599	67cm - 75cm	XXL

Measurement taken 15cm above mid-patella

AVAILABLE ON REQUEST



Infrapatellar Band

AIRCRAFT®

INDICATIONS

Osgood Schlatter disease; patella tendonitis.

CHARACTERISTICS

- An aircell focuses compression on the patella tendon, helping to decrease stress at the tibial tubercle
- Fits left or right knee
- Latex free
- 6 month warranty

ORDERING INFORMATION

Part Number	Circumference	Size
08A-B	25.5cm - 43cm	U

Measurement taken just below the kneecap

A person is walking barefoot on a beach, splashing water. The person is wearing a light-colored short-sleeved shirt and dark shorts. The background shows a beach with waves and mountains in the distance. The scene is captured in a warm, golden light, suggesting a sunset or sunrise. The water splashes are frozen in time, creating a dynamic and energetic feel. A yellow square is visible on the left side of the image.

WALKERS

AirSelect™ Range

Designed for Comfort, Engineered for Healing

Our lightest Aircast walking boot ever, AirSelect brings clinically tested innovations that improve comfort, compliance, and clinical outcomes.

Integrated Inflation System for Ease of Use


PUTS THE PATIENT IN CONTROL FOR AN INDIVIDUALIZED FIT.



Ultra Light Technology for Compliance


THE LIGHTEST FULL-SHELL WALKING BOOT AVAILABLE.

Semi-rigid shell integrates vents that maximize airflow and minimize mass, without compromising support.



Softstrike Technology & Optimized Rocker Sole Design for Comfort


A patented grid design provides cushioning that reduces heel loading and helps prevent impact-related stresses. Shape and height of the rocker sole provide a more natural gait and less hip differential.



Heel strike is absorbed and dissipated by the patented SoftStrike technology.

Mid-stance section is designed to minimize hip differential and knee flexion movement.

Toe-off section continues the rocking motion, enhancing mobility.




Duplex Technology for Faster Healing

Multiple aircells are layered within the semi-rigid shell to provide pulsating compression with each step.

This design is clinically proven to reduce edema 3x faster, relieve pain, and speed healing^{1, 2}


During ambulation, the overlapping aircells fluctuate – resulting in a massaging motion that resorbs the edema.



¹ Conservative Therapy for Acute Lateral Ligament Lesions – Single Chamber vs. Two-Chamber Orthosis Systems; Schmidt, Mainers, Reintges, Lipke, Benesch, Gerngross- Surgery Dept of the Federal Arm Hospital, Ulm, Germany, 1999
² A New Concept in Fracture Immobilization ; Dale PA, Bronk JT, O'Sullivan ME, et al – Clinical Orthopedics and Related Research, 1993



AIRCAST®
AirSelect™ Elite



AirSelect puts comfort, control and healing in the hands of the patient with the most advanced pneumatic walking boot available. Beneath its sleek appearance are patented technological advances that work together to enhance treatment and improve outcomes. The integrated inflation system lets patients inflate each aircell individually, to maximise comfort and minimise oedema.

INDICATIONS

Specifically designed for stable fractures of the lower leg, foot, and ankle; severe ankle sprains; and post-operative use. Clinically proven effective in healing fractures and managing oedema.

CHARACTERISTICS

- The most advanced walking boot, engineered to provide the ultimate in protection, comfort, and oedema control
- Overlapping posterior and distal duplex aircells provide intermittent compression, reducing oedema 3 times faster
- Pre-inflated front panel aircell and custom inflation system provide a full-contact fit
- Enclosed toe cover provides protection from hazards and the elements
- Latex free
- Fits left or right leg
- 6 month warranty



ORDERING INFORMATION

Part Number	Shoe Size		Size
	EU	UK	
01EP-XS	<35.5	<3	XS
01EP-S	35.5 - 39.5	3 - 6	S
01EP-M	39.5 - 43	6 - 9	M
01EP-L	43 - 47	9 - 12	L
01EP-XL	>47	>12	XL



AIRCAST®
AirSelect™ Standard



The Standard walking boot offers a moderate level of support and Aircast's quality, effectiveness and comfort in an economical model.

INDICATIONS

Stable fracture of foot/ankle/lower leg, severe ankle sprain, post-operative immobilisation.

CHARACTERISTICS

- Semi-rigid shell that supports the limb while providing protection
- Two adjustable distal aircells that provide compression and support
- Integrated inflation system for ease of use
- Lab-tested rocker sole encourages natural gait
- Latex free
- Fits left or right leg



ORDERING INFORMATION

Part Number	Shoe Size		Size
	EU	UK	
01EF-XS	<35.5	<3	XS
01EF-S	35.5 - 39.5	3 - 6	S
01EF-M	39.5 - 43	6 - 9	M
01EF-L	43 - 47	9 - 12	L
01EF-XL	>47	>12	XL



AirSelect™ Short

AIRCAST®

INDICATIONS

Metatarsal fractures, stress and stable fracture of the foot, bunionectomy, hallux valgus, soft tissue injury, severe ankle sprain.

CHARACTERISTICS

- The Short walking boot has a durable, semi-rigid shell that helps support the limb while providing protection
- The pre-inflated front panel and two customisable distal aircells provide compression and support
- Integrated inflation system for ease of use
- SoftStrike technology absorbs and dissipates shock, while the lightweight and lab-tested rocker sole encourages natural gait
- Latex free
- Fits left or right leg

ORDERING INFORMATION

Part Number	Shoe Size		Size
	EU	UK	
01ES-XS	<35.5	<3	XS
01ES-S	35.5 - 39.5	3 - 6	S
01ES-M	39.5 - 43	6 - 9	M
01ES-L	43 - 47	9 - 12	L
01ES-XL	>47	>12	XL



XP Diabetic Walker™

AIRCAST®

INDICATIONS

Stable foot and/or ankle fracture; severe ankle sprain; post-operative; Charcot foot; ulcer in malleolus, hind-foot and heel.

CHARACTERISTICS

- Specifically designed to meet the needs of the higher risk diabetic patient with a high rocker sole for maximum off-loading of the diabetic foot
- Dual density insole to help eliminate pressure points
- Graduated compression to help reduce oedema
- Regulated aircell inflation with the specialised hand bulb with pressure gauge
- Latex free
- Fits left or right leg
- 6 month warranty

ORDERING INFORMATION

Part Number	Shoe Size		Size
	EU	UK	
01PD-S	35 - 38	4 - 5.5	S
01PD-M	39 - 42	6 - 8	M
01PD-L	>42	>8	L
01PD-XL	>45	>10.5	XL



Heel Wedges

AIRCAST®

INDICATIONS

Designed to be used with the Aircast walkers. Indicated for post Achilles tendon surgery, the heel supports can be placed on top of each other to achieve the desired heel rise.

CHARACTERISTICS

- Supplied as a pack of five heel wedges of varying size to suit most walker boots
- Up to three wedges can be combined in a single walker to achieve the following levels of plantarflexion:
 - 3 heel wedges: 22 degrees
 - 2 heel wedges: 16 degrees
 - 1 heel wedge: 10 degrees
- Latex free

ORDERING INFORMATION

Left	Right	Walking Brace Size
01K-L	01K-R	S/M/L
01K-XLL	01K-XLR	XL



Nextep™ Contour Walker

PROCARE®

INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg or forefoot.

CHARACTERISTICS

- Secure wraparound shell housing provides immobilisation and support for ligament and tendon repair
- Eliminates need for serial casting
- Brushed nylon/foam liner with detachable toe cover
- Rocker bottom with non-skid sole
- Fits left or right leg

ORDERING INFORMATION

Part Number	UK Shoe Size		Size
	Mens	Womens	
79-95273	3 - 6	2 - 5	S
79-95275	6.5 - 9	5.5 - 8	M
79-95277	9.5 - 12.5	8.5 - 13.5	L

AVAILABLE ON REQUEST

COMING SOON



XcelTrax™ Tall

PROCARE®

INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

CHARACTERISTICS

- Lightweight design
- Unique pre-shaped ergonomic frame helps to modify the level of pressure on the injury at both the malleoli and lower leg
- Fits left or right leg

ORDERING INFORMATION

Part Number		UK Shoe Size	Size
Pack of 1	Pack of 16		
79-95492-1		<3.5	XS
79-95493-1	79-95493-16	4 - 6.5	S
79-95495-1	79-95495-16	7 - 10	M
79-95497-1	79-95497-16	10 - 12	L
79-95498-1		>12	XL

COMING SOON



XcelTrax™ Ankle

PROCARE®

INDICATIONS

Ideal following metatarsal fractures, bunionectomies, acute ankle sprains, stress and stable fractures of the lower leg and forefoot.

CHARACTERISTICS

- Short version
- Lightweight design
- Unique pre-shaped ergonomic frame helps to modify the level of pressure on the injury at both the malleoli and lower leg
- Fits left or right leg

ORDERING INFORMATION

Part Number		UK Shoe Size	Size
Pack of 1	Pack of 16		
79-95502-1		<3.5	XS
79-95503-1	79-95503-16	4 - 6.5	S
79-95505-1	79-95505-16	7 - 10	M
79-95507-1	79-95507-16	10 - 12	L
79-95508-1		>12	XL



PROCARE®
MiniTrax™



SMALL

MEDIUM

LARGE

The Short walking boot offers a solution to a number of indications related to forefoot, mid-foot and hind foot injuries, as well as some ankle injuries.

INDICATIONS

Ideal support and comfort for your smallest patients with mild to severe foot and ankle injuries. MiniTrax™ is the only pediatric walking boot to feature kid-friendly design elements that make wearing a boot less intimidating. Beneath the fun exterior though, the MiniTrax is a serious boot with big benefits.

CHARACTERISTICS

- Lightweight design features anatomically shaped uprights to facilitate fit
- Low-profile rocker sole of the MiniTrax features a durable, non-slip tread tough enough to withstand the test of time - and of a kid!
- Easy to apply and adjust with a simple strapping and liner system



ORDERING INFORMATION

Part Number	Age	Foot Length	Leg Length	Child Shoe Size	Size
79-95510	1 - 2 ½	<15cm	<23cm	<7.5	S
79-95520	2 ½ - 5	15cm - 17.5cm	<28cm	7.5 - 10	M
79-95530	5 - 8	17.5cm - 20cm	<36cm	10.5 - 1 adult	L



ANKLE

AIRCAST®

AirSport+

AVAILABLE ON REQUEST



The AirSport+ Ankle Brace is a novel design, intended to help injured or weakened ankles regain their full potential, by providing support, compression and stability in a comfortable and easy to apply brace. Ideal for active individuals or athletes recovering from ankle injuries and keen to return to sports activities.

INDICATIONS

- Mild to moderate sprains (Grade I & II)
- Ankle ligament distortions, tears
- Chronic ankle instabilities
- Prophylactic use, return to sports activities
- Post-op rehab use

CHARACTERISTICS

- Innovative ATOP/MOZ Lacing System provides adjustable, uniform compression without creating pressure points
- Bilateral anatomically-shaped shells minimize the risk of ankle rollover (inversion/eversion) and provide maximum medial & lateral support for unstable ankles
- Posterior high ankle strap provides additional stability to the upper ankle
- Aircast clinically-proven semi-rigid shells and aircells provide compression and support
- Foam-filled aircells cushion and protect the ankle
- ATF (anterior talofibular) cross strap and forefoot wrap provide additional compression and stability to the joint
- Step-in design makes applying the brace easy
- Elastic fabric at heel ensures a proper fit
- Breathable material for coolness & comfort during daily or sports activities



ORDERING INFORMATION

Part Number		UK Shoe Size	Size
Left	Right		
82-0266-1	82-0265-1	<3	XS
82-0266-2	82-0265-2	3 - 5	S
82-0266-3	82-0265-3	6 - 8	M
82-0266-4	82-0265-4	9 - 11	L
82-0266-5	82-0265-5	>11	XL



Stabilising Speed Pro

DONJOY®

INDICATIONS

Chronic ankle sprains and prophylactic use.

CHARACTERISTICS

- Improved low profile design reduces the risk of ankle injury by up to 61%
- Wide (7.5cm) circumferential elastic strap for compression and stability
- Removable medial/lateral stays for additional support and stability
- Speed lacing system for fast application
- Fits left or right ankle
- 6 month warranty

ORDERING INFORMATION

Part Number	Heel/Ankle Circumference	Size
11-3234-0	23cm - 25cm	XXS
11-3234-1	25cm - 28cm	XS
11-3234-2	28cm - 30cm	S
11-3234-3	30cm - 33cm	M
11-3234-4	33cm - 36cm	L
11-3234-5	36cm - 38cm	XL
11-3234-6	38cm - 41cm	XXL
11-3234-7	41cm - 43cm	XXXL

Measurement taken from base of the heel around the ankle whilst wearing a sports sock



Stabilising Pro

DONJOY®

INDICATIONS

Chronic ankle sprains and prophylactic use.

CHARACTERISTICS

- Improved low profile design reduces the risk of ankle injury by up to 61%
- Wide (7.5cm) circumferential elastic strap for compression and stability
- Removable medial/lateral stays for additional support and stability
- Speed lacing system for fast application
- Fits left or right ankle
- 6 month warranty

ORDERING INFORMATION

Part Number	Heel/Ankle Circumference	Size
11-3235-0	23cm - 25cm	XXS
11-3235-1	25cm - 28cm	XS
11-3235-2	28cm - 30cm	S
11-3235-3	30cm - 33cm	M
11-3235-4	33cm - 36cm	L
11-3235-5	36cm - 38cm	XL
11-3235-6	38cm - 41cm	XXL
11-3235-7	41cm - 43cm	XXXL

Measurement taken from base of the heel around the ankle whilst wearing a sports sock



A60™ Ankle Support

AIRCAST®

INDICATIONS

Chronic instability, prophylaxis.

CHARACTERISTICS

- A 60° stabiliser guards against rollover and helps prevent sprains
- Breath-O-Prene® fabric keeps wearer cool and dry
- Sleek, anatomic design for a comfortable bulk-free fit
- Simple single strap application and adjustment
- Latex free
- 6 month warranty

ORDERING INFORMATION

Part Number BLACK		UK Shoe Size	Size
Left	Right		
02TSL	02TSR	<6	S
02TML	02TMR	6 - 8	M
02TLL	02TLR	>8	L

AVAILABLE ON REQUEST

AIRCAST®

Air-Stirrup®

INDICATIONS

Acute ankle distortion grade I, II & III, post-operative application, chronic ankle instability.

CHARACTERISTICS

- Semi-rigid, anatomically designed shells for protection, comfort and the prevention of inversion/eversion
- Patented Duplex aircell system to enhance circulation and reduce swelling
- Pre-inflated aircells for easy application
- Streamlined to fit in shoes for early protected weight-bearing
- Latex free
- 6 month warranty



ORDERING INFORMATION

Part Number		Patient Height	Brace Length	Size
Left	Right			
02JL	02JR	Small children	15cm	XS (Paediatric)
02CL	02CR	<157cm (<5'2")	22cm	S
02BL	02BR	157cm - 167cm (5'2" - 5'6")	23cm	M
02AL	02AR	>162cm (>5'4")	27cm	L

AIRCAST®

Air-Stirrup® Universe™

INDICATIONS

Ankle sprains and chronic ankle instability.

CHARACTERISTICS

- Patented Duplex aircell system to reduce swelling and pain by enhanced circulation
- Universal semi-rigid shell design fits left or right ankle
- Universal size
- Latex free
- 6 month warranty



ORDERING INFORMATION

Part Number	Size
02E	U



**DIABETIC,
OFFLOADING
& FOOT
SUPPORTS**



Squared Toe Post-Op Shoe

PROCARE®

INDICATIONS

Post-operative treatment of hallux valgus, toe deformities and protocols where maximum immobilisation and protection of the forefoot are required.

CHARACTERISTICS

- The squared design allows additional room and comfort for the patient whilst keeping the dressing clean and toes dry
- High ankle strapping keeps the foot in the proper location, for ease of ambulation and pressure reduction
- Quality construction provides support and minimal movement of the foot while ambulating
- Rigid rocker sole design allows for ease of ambulation and the skid resistant sole provides traction
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size		Size
	Mens	Womens	
79-81232	3 - 4.5	2 - 3.5	XS
79-81233	5 - 6.5	4 - 5.5	S
79-81235	7 - 8.5	6 - 7.5	M
79-81237	9 - 11.5	8 - 10.5	L
79-81238	>11.5	>10.5	XL



Off-Loading Diabetic Shoe

PROCARE®

INDICATIONS

Pressure relief during the treatment of diabetic ulcers.

CHARACTERISTICS

- Tri-laminate IMPAX grid insole designed to redistribute load away from ulcers with easy to remove IMPAX foam blocks
- IMPAX grid allows for a continuous plantar surface
- Lightweight breathable material enhances patient's comfort and promotes compliance
- Squared toe box and wide foot bed to leave ample room for wound dressings
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size Mens	UK Shoe Size Womens	Size
79-81513	3.5 - 5.5	4 - 6	S
79-81515	6 - 8	6.5 - 8.5	M
79-81517	8.5 - 10.5	9 - 10	L
79-81518	11 - 13.5	>10.5	XL

REPLACEMENT INSOLE

Part Number	UK Shoe Size Mens	UK Shoe Size Womens	Size
79-81523	3.5 - 5.5	4 - 6	S
79-81525	6 - 8	6.5 - 8.5	M
79-81527	8.5 - 10.5	9 - 10	L
79-81528	11 - 13.5	>10.5	XL



Podalux II

DONJOY®

INDICATIONS

Post-op, post-trauma foot and toe applications and for patients requiring extra stability, more specifically after hallux valgus and ensuring the foot, toes and forefoot are held in a proper position during rehabilitation.

CHARACTERISTICS

- Wide, solid skid resistant sole for extra stability
- Removable rigid foot plate for additional firmness and arch support
- Lined with breathable micro-mesh textile for patient comfort
- Two contact closure straps for easy application and adjustment
- Toe piece included to cover up the toes if necessary
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size Mens	UK Shoe Size Womens	Size
82-0042-1	3 - 4	2.5 - 3.5	S
82-0042-2	4.5 - 5.5	4 - 5	M
82-0042-3	6 - 7	6 - 6.5	L
82-0042-4	7.5 - 8	7 - 7.5	XL
82-0042-5	8.5 - 10.5	8 - 10.5	XXL
82-0042-6	11 - 13	11 - 13	XXXL



Remedy Pro™

PROCARE®

INDICATIONS

Reducing forefoot pressure for forefoot trauma, post-surgical healing and diabetic ulcers or wounds.

CHARACTERISTICS

- 15° wedge design reduces forefoot pressure by shifting weight to midfoot and heel
- Easy forefoot closure
- Durable outer sole, padded liner and adjustable ankle strap
- Seamless, padded heel
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size		Size
	Mens	Womens	
79-81722	<4.5	<4	XS
79-81723	4.5 - 6.5	4 - 6	S
79-81725	6.5 - 8.5	6 - 8	M
79-81727	8.5 - 10.5	8 - 10	L
79-81728	>10.5	>10	XL



Medical/Surgical Shoe

PROCARE®

INDICATIONS

Post-op, post-trauma applications where some mobility of the foot is allowed.

CHARACTERISTICS

- Loop/lock contact closure straps
- Comfortable contoured heel with foam mesh upper and semi-flexible skid resistant sole
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size Mens	Size
79-81132	<4	XS
79-81133	4 - 8	S
79-81135	8 - 11	M
79-81137	11 - 14	L
79-81138	>14	XL

Part Number	UK Shoe Size Womens	Size
79-81142	<1	XS
79-81143	1 - 3.5	S
79-81145	3.5 - 5.5	M
79-81147	5.5 - 7.5	L



Super-Lite A.F.O.

PROCARE®

INDICATIONS

Correction of drop foot.

CHARACTERISTICS

- Low profile
- Trim to fit ankle orthosis
- Female (narrower) and male (wider) sizes
- 6 month warranty

ORDERING INFORMATION

Part Number	Description	Size
79-81150	Right, Male	U
79-81151	Left, Male	U
79-81160	Right, Female	U
79-81161	Left, Female	U

AVAILABLE ON REQUEST

PROCARE®

Podous Boot

INDICATIONS

Prevent lower extremity disorders associated with trauma or immobility including: pressure necrosis, ankle contractures and post-operative hip abduction or rotation.



CHARACTERISTICS

- Lightweight, quality plastic shell with breathable poly/pile liner
- Adjustable toe extension helps protect toes from bedding
- Anti-rotation bar may be repositioned to assist with internal or external rotation
- Latex free
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size		Size
	Mens	Womens	
79-90550	<9.5	<9	U
79-90551	>9.5	>9	XL

Replacement Liners	Pack Size	Size
79-90560	Pack of 3	U
79-90561	Pack of 3	XL

PROCARE®

Plantar Fasciitis Night Splint

INDICATIONS

Controlled application of a prolonged stretch on the plantar fascia and Achilles tendon.



CHARACTERISTICS

- Lightweight, sturdy one-piece flexible shell with a terry liner
- Adjustable bilateral dorsiflexion straps
- Padded calf and ankle straps for additional comfort
- Includes optional toe wedge for added stretch of the plantar fascia
- Fits left or right foot

ORDERING INFORMATION

Part Number	UK Shoe Size		Size
	Mens	Womens	
79-97753	<5	<6	S
79-97755	5 - 8	6 - 8.5	M
79-97757	8 - 10.5	8.5 - 11	L
79-97758	>10.5	11 - 13.5	XL



Silicone Heel Cups

DONJOY®

INDICATIONS

Plantar fasciitis, heel spurs and Achilles tendonitis.

CHARACTERISTICS

- Moulded medical-grade soft silicone to reduce shock and vibration loads to joints
- Softer centre in heel region to reduce strain to area sensitive to pain
- Sold in pairs

ORDERING INFORMATION

Part Number	UK Shoe Size		Size
	Mens	Womens	
79-81102	Children/youth	Children/youth	XS
79-81103	4.5 - 8.5	3 - 7	S/M
79-81107	>8.5	>7.5	L/XL



Podalib

AIRCAST®

INDICATIONS

Dynamic AFO, indicated for neurological weakness or paralysis of the foot-lifting muscles, resulting in drop foot.

CHARACTERISTICS

- Dynamic and functional approach to address drop foot.
- Helps maintain the foot at 90° angle, prevents lateral imbalance, helps improve ambulation.
- Allows the patient to regain a more supple and natural walking pattern.
- Lightweight, supple textile structure makes for a comfortable and discreet walking aid.
- Ankle cuff with adjustable elastic strap and hook helps maintain foot at 90° angle.
- System adapts to all types of shoes or can be used barefoot.
- Breathable, durable and comfortable materials.
- Bilateral model, available in 4 sizes.

ORDERING INFORMATION

Part Number	Ankle Circumference	Ankle Cuff Height	Size
82-0013-1	13-19 cm	10 cm	XS
82-0013-2	13-20 cm	13 cm	S
82-0013-3	20-27 cm	13 cm	M
82-0013-4	27-34 cm	13 cm	L

AVAILABLE ON REQUEST



Hallux Valgus

AIRCAST®

NEW

ACTYTOE™

INDICATIONS

Hallux Valgus - or bunion - is one of the most common foot deformities, especially for women. Because of its weaker connective tissue, the big toe can become misaligned, resulting in an unsightly and painful protrusion from the foot.

CHARACTERISTICS

- Low-profile design
- Memory foam cushion
- Adjustable, thin and soft straps
- Anatomically contoured splint
- Breathable and lightweight materials
- For day and night wear
- Available in 3 sizes (S, M, L)
- For left and right foot
- Latex free

ORDERING INFORMATION

Part Number			Size
ActyToe	ActyToe Night	ActyToe Lift	
-	-	82-06	Universal
82-04S	82-05S	-	S
82-04M	82-05M	-	M
82-04L	82-05L	-	L



AVAILABLE ON REQUEST

ACTYTOE™ NIGHT



AVAILABLE ON REQUEST

ACTYTOE™ LIFT

COLD THERAPY





IceMan® CLEAR³™

DONJOY®

INDICATIONS

The IceMan CLEAR³ cold therapy unit helps reduce pain and swelling, speeding up rehabilitation. It provides extended cold therapy for a variety of indications and protocols. It utilizes DonJoy's patented semi-closed loop recirculation system, which maintains more consistent and accurate temperatures than other cold therapy units, in a preset configuration.

CHARACTERISTICS

- Sure fit, stackable, locking lid with watertight seal
- Low profile design, small cooler for easy portability and storage
- Transparent cooler design for easy refill indication
- Preset inlet temperature for simple and convenient, plug and play use
- Practical, integrated hose hanger for ease of use
- Compatible with a variety of pads

ORDERING INFORMATION

Part Number	Description
11-1639	IceMan CLEAR ³ (device without cuff).
11-1640	IceMan CLEAR ³ (device with cuff).

DONJOY®

Universal Cuff

INDICATIONS

Delivers even cold therapy throughout the pad.

CHARACTERISTICS

- Designed with a hook engageable material and delivered with a wrap to secure the cold pad around the affected area.
- The universal size is ideal for the shoulder, knee and ankle.
- Designed to be used in conjunction with a moisture barrier.
- Available in Regular Hose or Extended Hose.



ORDERING INFORMATION

Part Number	Description	Dimensions
11-0679-9-00000	RH Universal	27.5 cm x 29 cm
11-0680-9-00000	EH Universal	27.5 cm x 29 cm
11-0681-9-00000	RH, XL Universal	30 cm x 29 cm
11-0682-9-00000	EH, XL Universal	30 cm x 29 cm

DONJOY®

McGuire Cuff

INDICATIONS

Delivers even cold therapy throughout the pad.

CHARACTERISTICS

- Designed with a hook engageable material and delivered with a wrap to secure the cold pad around the affected area.
- Designed to be used in conjunction with a moisture barrier.



ORDERING INFORMATION

Part Number	Description	Dimensions
11-1268-9-00000	McGuire Knee Wrap w/ open patella	27.5 cm x 29 cm



Cryo/Cuff™ IC Cooler

AIRCAST®

INDICATIONS

Integrated cold and compression clinically proven to help reduce post-operative swelling and regain range of motion. Ideal for post-operative recovery, trauma, athletic training rooms and home use.

CHARACTERISTICS

- Provides both gravity and motorised cold and compression therapy for flexibility of care
- Integrated pneumatic pump within the cooler lid provides automated cold and compression therapy, with simplicity of design and ease of use operation
- Compatible with all Cryo/Cuff cuffs
- 6 month warranty

ORDERING INFORMATION

Part Number	Description
52A	Cryo/Cuff IC Cooler
2125	Cryo/Cuff Cooler

Replacement Parts	Description
25-0238	Cryo/Cuff IC Lid
25-2740	Cryo/Cuff IC Power Supply

AIRCAST®

Knee Cuff

INDICATIONS

Minimise haemarthrosis, swelling and reduce pain.

CHARACTERISTICS

- Anatomically designed cuff provides complete knee coverage for optimum treatment
- Suitable for use with Cryo/Cuff Gravity Cooler (Item 2125) and Cryo/Cuff IC Cooler (Item 52A)
- Paediatric knee cuff also available (can be used as an elbow cuff)
- Latex free
- Fits left or right knee
- 6 month warranty



ORDERING INFORMATION

Part Number	Circumference	Size
11C01	25cm - 50cm	S
11A01	45cm - 60cm	M
11B01	50cm - 80cm	L
11P01	Under 7 years old	Paediatric

Measurement taken 15cm above mid-patella

AIRCAST®

Ankle Cuff

INDICATIONS

Ideal for acute sprains, trauma, post-operative and rehabilitation, reducing haemarthrosis, swelling and pain.

CHARACTERISTICS

- Anatomically designed to completely fit the ankle to provide maximum cryotherapy
- Suitable for use with Cryo/Cuff Gravity Cooler (Item 2125) and Cryo/Cuff IC Cooler (Item 52A)
- Paediatric version available, suitable for ages 1 - 7 years
- Latex free
- Fits left or right ankle
- 6 month warranty



ORDERING INFORMATION

Part Number	Size
10A01	U
10P01	Paediatric

*All of the other cuffs are available on request

Shoulder Cuff



INDICATIONS

Minimise haemarthrosis, swelling and reduce pain.

CHARACTERISTICS

- Anatomically designed cuff provides complete shoulder coverage for optimum treatment
- Suitable for use with Cryo/Cuff Gravity Cooler (Item 2125) and Cryo/Cuff IC Cooler (Item 52A)
- Available with an extra long strap to adapt to a larger chest circumference
- Latex free
- Fits left or right shoulder
- 6 month warranty

ORDERING INFORMATION

Part Number	Chest Circumference	Size
12A01	81cm - 122cm	M/L
12AXL01	107cm - 137cm	XL Strap

*All of the other cuffs are available on request

Elbow Cuff

INDICATIONS

Ideal for epicondylitis (tennis elbow), post-operative and throwing injuries.

CHARACTERISTICS

- Minimises haemarthrosis, swelling and pain
- Anatomically designed cuff provides complete elbow coverage for optimum treatment
- Suitable for use with Cryo/Cuff Gravity Cooler (Item 2125) and Cryo/Cuff IC Cooler (Item 52A)
- Paediatric version available, suitable for ages 1 - 7 years
- Latex free
- Fits left or right elbow
- 6 month warranty



ORDERING INFORMATION

Part Number	Size
15A01	U
11P01	Paediatric

*All of the other cuffs are available on request

*SINGLE patient use only

Hand/Wrist Cuff

INDICATIONS

Ideal use for carpal tunnel syndrome, tendonitis, post-operative, and sports injuries.

CHARACTERISTICS

- Anatomically designed to cover the hand and wrist
- Versatile design allows repositioning for desired hand coverage
- Minimises swelling and pain
- Latex free
- Fits left or right hand
- 6 month warranty



ORDERING INFORMATION

Part Number	Size
16A01	U

*All of the other cuffs are available on request

VARIOUS





VenaPure anti-embolism stockings have been designed with Venous Thromboembolism (VTE) prevention in mind. DJO Global has been able to use the years of clinical experience gained in the arena of VTE prophylaxis to combine all of these positive attributes into VenaPure.

INDICATIONS

For use in the prophylaxis of deep vein thrombosis and pulmonary embolism in hospitalised patients and for use in other conditions as directed by the clinician.

CHARACTERISTICS

- Graduated compression according to Sigel profile (18mmHg at the ankle, decreasing to 8mmHg at the thigh) to minimise venous stasis and venous dilation
- Inspection toe opening on top of the foot to minimise risk of pressure point and ensure that skin integrity is maintained
- Popliteal break is easily identified by stitch change to ensure that pressure is relieved around the popliteal vein
- Pressure relieving panel at the inner thigh to ensure reduced compression of the femoral vein
- Interrupted silicone band at upper thigh provides excellent fixation to the skin, without compromising compression
- Inlay circumferential knit maintains correct positioning of stocking every time
- Colour coded toe and top ensure that the correct size is selected according to patient measurement
- Tape measure included in the dispenser to make measuring of the patient easier
- Available in two styles with lengths and sizes to fit all patients



ORDERING INFORMATION - THIGH LENGTH (PACK OF 6 PAIRS)

Thigh Circumference	Calf Circumference	Leg Length	Code	Size
<63.5cm	<30.5cm	<74cm	401S-2	S
		74cm - 84cm	401R-2	
		>84cm	401L-2	
	30.5cm - 38cm	<74cm	401S-3	M
		74cm - 84cm	401R-3	
		>84cm	401L-3	
38cm - 44.5cm	<74cm	401S-4	L	
	74cm - 84cm	401R-4		
	>84cm	401L-4		
63.5cm - 81.5cm	44.5cm - 55cm	<74cm	401S-5	XL
		74cm - 84cm	401R-5	
		>84cm	401L-5	
81cm - 91cm	55cm - 66cm	<74cm	401S-6	XXL
		74cm - 84cm	401R-6	
		>84cm	401L-6	

*Thigh with belt available on request

ORDERING INFORMATION - KNEE LENGTH (PACK OF 12 PAIRS)

Calf Circumference	Leg Length	Code	Size
<30.5cm	<41cm	400R-2	S
	>41cm	400L-2	
30.5cm - 38cm	<43cm	400R-3	M
	>43cm	400L-3	
38cm - 44.5cm	<46cm	400R-4	L
	>46cm	400L-4	
44.5cm - 51cm	<46cm	400R-5	XL
	>46cm	400L-5	
51cm - 58.4cm	<46cm	400R-6	XXL
	>46cm	400L-6	
58.5cm - 66cm	<46cm	400R-7	XXXL
	>46cm	400L-7	



Dr.Comfort®

Dr Comfort proudly carries the American Podiatric Medical Association seal of acceptance on products as indicated.

At Dr Comfort we understand the importance of foot care and are committed to helping people get back on their feet. With over 100 years combined experience in foot care, footwear knowledge and pedothic training, our team of educators know a thing or two about proper foot health.

Did you know? Revolutionizing the way you walk, our Stability outsole is designed to help provide enhanced support and prevent over pronation. Athletically-inspired shoes that are bold, comfortable and performance-focused. (PDAC approved).

- Wider sole helps provide enhanced stability
- Increased width of toe box provides extra comfort and accommodation
- EVA and rubber outsole provides a layer under foot for added support and comfort
- Firm medial rearfoot helps prevent overpronation



Annie Black

ORDERING INFORMATION

Part Number	Description
4510	Annie Black



Lana Black

ORDERING INFORMATION


Part Number	Description
110410	Lana Black



Flute Black

ORDERING INFORMATION

Part Number	Description
2315	Flute Black



Grace Seafoam

ORDERING INFORMATION

Part Number	Description
10325	Grace Seafoam



Protector Black

ORDERING INFORMATION

Part Number	Description
9310	Protector Black



Performance Black/Grey

ORDERING INFORMATION

Part Number	Description
7610	Performance Black/Grey



Classic Black

ORDERING INFORMATION

Part Number	Description
8410	Classic Black



Winner Black

ORDERING INFORMATION

Part Number	Description
6710	Winner Black



Brian Black

ORDERING INFORMATION

Part Number	Description
6510	Brian Black

Terms and Conditions

1. PARTIES

- a) The "Company" means DJO Orthopaedic South Africa (Pty) Ltd.
- b) The "Customer" means the person, firm or company buying goods from the Company to whom the goods are to be invoiced (unless the context otherwise requires) and any person to whom or for the benefit of whom the goods may be delivered.
- c) The "Goods" means the goods described on our invoice.

2. ORDERS

- a) In the Republic of South Africa, credit account customer orders with a value of less than R6,500 (excluding VAT, net discounts) will incur a R150.00 delivery charge. All orders received from credit account customer in neighboring countries will incur a delivery charge based on a detailed quotation on the products ordered.
- b) Any orders submitted or variations made orally must be confirmed in writing within 5 days.
- c) All orders are to be forwarded to the address below.
- d) V.A.T. Exemption – Buyers who qualify for exemption from paying V.A.T. should include a valid V.A.T. exemption certificate (or other documentation approved by the Customs and Excise for the purpose) with each order.
- e) Orders to non credit account customers require payment to accompany the order before the Goods can be dispatched.

3. QUOTATIONS

Any estimate or quotation shall comprise an invitation to treat and will remain open for a period of 30 days provided that the Company has not previously withdrawn it.

4. SPECIFICATION

Every effort is made to ensure our catalogue descriptions and illustrations are clear and accurate. All colours and sizes are approximate. The reproduction of colours is as accurate as photographic and printing processes allow. Due to our policy of continual improvement, some products may vary from those described. The company reserves the right to amend the specifications, construction and/or design without prior notice.

5. PRICES

- a) All prices quoted are inclusive of Value Added Tax.
- b) The Company reserves the right to vary the price of the Goods to take account of any increases in the cost of raw materials, manufacture, packing, transport, wages or otherwise.
- c) Prices charged will be those effective at the date of receipt of Customer order.

6. PAYMENT

- a) Unless the Company otherwise agrees, settlement is due within 30 days of date of invoice to DJO Orthopaedic South Africa (Pty) Ltd at the address below.
- b) The Company reserves the right to charge interest on all overdue payments and outstanding balances at the rate of 3% per month.
- c) The Customer shall not be entitled to withhold payment either wholly or in part once the same has become due except where expressly agreed with the Company.
- d) Non credit account customers and sales to members of the public require either payment of a pro-forma invoice or payment to accompany the order before the Goods will be dispatched.

7. DELIVERY / CARRIAGE CHARGES

- a) Unless the Company otherwise agrees, delivery shall be made to the delivery point specified in the order.
- b) Stated delivery times are an estimate on the part of the Company and shall not be binding upon the Company which shall not be liable for any costs or damages caused by reason of any delay in delivery.
- c) The Company reserves the right to make delivery by installments and to tender a separate invoice in respect of each installment.
- d) Orders received before 15:00 Monday to Friday will only be dispatched the same day for

Republic of South Africa customers with credit accounts that do not have overdue payments. This is subject to products being in stock and does not apply to custom manufactured products and non stock items. Delivery will be on a next day basis with some exceptions. Deliveries to remote areas may take up to 72-96 hours.

e) Standard Delivery Charges:

- (i) All deliveries to Republic of South Africa residential addresses including webshop orders and drop ship addresses will incur a delivery charge.
- (ii) Timed and Saturday deliveries can be arranged at a cost. This will be detailed at the time of ordering.
- (iii) Oversized items are subject to additional delivery and handling charges, regardless of the value of the order.
- f) Oversized products are delivered to a ground floor location and are not installed in-situ. Should additional delivery services be required, they must be specified at the time of order and may be subject to additional costs.
- g) Delivery for custom manufactured products may vary from 2 weeks to 6 weeks from date of order.
- h) Spare parts, unusual sized products and special order items can take between 15 and 30 working days for delivery.

8. CANCELLATIONS

- a) No notice of cancellation of an order shall be valid unless given in writing and agreed to by an authorised representative of the Company.
- b) If written notice of cancellation is received by the Company, the Company reserves the right to levy a cancellation charge up to 50% of the amount of the order.
- c) If the contract provides for delivery by installments, each installment shall be deemed to be the subject of a separate contract, and non-delivery or delay in the delivery of any installment shall not entitle the Customer to repudiate the contract as regards the balance.

9. PASSING OF TITLE

The legal title to the Goods shall not pass to the Customer until all sums due or payable by the Customer to the Company whether in respect of the Goods or otherwise have been received by the Company.

10. PASSING OF RISK

The risk in the Goods shall pass to the Customer when the Goods are tendered for delivery at the delivery point specified in the order. In the event that the Goods are collected by the Customer from the Company, the risk of the Goods shall pass to the Customer when they are loaded onto the vehicle of the Customer or of his carrier or other agent.

11. WARRANTY

The Company warrants that the Goods shall at the time of delivery correspond to the specifications published by the Company when used for purposes for which goods of that type are ordinarily used. In the absence of any such published specifications applicable to the Goods, the Company warrants that the Goods shall at the time of delivery be of merchantable quality. If the Goods do not conform to that warranty, the Company will, at its option:

- a) Replace the Goods found not to conform to the warranty;
 - b) Bring the Goods into conformity with the published specifications of the Company or (as the case may be) merchantable quality; or
 - c) Take back the Goods found not to conform to the warranty and refund the appropriate part of the purchase price.
- Replacement of DrComfort shoes – within 3 months of purchase.
Replacement of BSV products within 6 months of purchase.

12. FORCE MAJEURE

The Company shall not be liable for any loss or damage caused by circumstances beyond the Company's control, including but not limited to, an act of God, war, civil disturbance, governmental restrictions, import or export regulations, industrial disputes and difficulties in obtaining labour or materials. Should any such event occur, the Company might cancel or suspend this contract without incurring any liability for any loss or damage caused.

13. NOTICE OF DAMAGE, NON-DELIVERY OR WRONG DELIVERY

- a) The Company gives no guarantee of the exact quantity of goods delivered, and in the event that quantity of Goods delivered does not correspond with the quantity stated on the order the customer shall only be liable to pay for the quantity delivered.
- b) If upon receipt the Goods appear damaged or incomplete, please indicate on the carrier's documentation and notify our Customer Services Department of the situation no later than 3 days after delivery. If in doubt, please sign the carrier's documentation "Unchecked". Should an expected consignment not arrive, please notify our Customer Services Department within 21 days from date of invoice, failing which the company shall be under no liability for non-delivery.
- c) Any surplus goods delivered shall remain the property of the Company, and the Customer shall take all reasonable precautions for the safe custody and protection of such surplus goods until the time of their removal by the Company.
- d) In the event of goods being wrongly delivered, under no circumstances shall any goods be returned to the Company without its prior consent, otherwise liability is not accepted. The Company will issue a Returns Authorisation Number (RAN) and this must be marked on the exterior of every package returned. Returned Goods must be suitably packaged for transportation by the carrier. Goods received without suitable packaging will not be credited. Where the Company has made the error, full credit will be given on uplifted goods in good condition.
- e) The Company is under no obligation to accept goods for return, e.g. in the case of Customer ordering error, surplus stock, etc., but where it agrees to do so the following handling charges will apply:-

05 days of delivery, no restocking fee. (Excluding Dr.Comfort & custom product orders) 30 days of delivery, a 15% restocking/packaging fee.
30 days – 365 Days of delivery, a 25% restocking/packaging fee. Over 12 months of delivery, no credit will be given.

The Customer will arrange carriage and pay all freight charges. The Customer should obtain a Returns Authorisation Number, which should be marked clearly on each package being returned. Only complete units of sale may be returned and all goods are to be in good condition. The Customer will be liable for a 25% restocking fee for any goods returned in a non saleable condition or with missing accessories, regardless of the number of days since delivery, with the balance of the goods being paid for in full by the Customer.

- f) Refused Deliveries – We reserve the right to make a carriage charge plus 15% handling charge on all consignments which are tendered for delivery and refused without justification. Goods returned in a non-saleable condition will be paid for in full by the Customer. Goods should not be returned without the prior agreement of the Company.
- g) Only unused and complete goods may be returned, and they must be undamaged, complete with all accessories and less than 12 months old. Goods not conforming to these criteria will not be accepted for credit. The Company reserves the right to dispose of any non saleable and non creditable returned goods.
- h) Custom manufactured products and non stock products are non returnable.

14. DISPUTES AND SET-OFF

Any liability of the Company under the contract shall be subject to and conditional upon the due performance and observance by the Customer of all its obligations under these conditions. Subject to these conditions, the Customer shall not be entitled to withhold or delay payment or exercise any right or set-off whatsoever and howsoever arising or arisen which might otherwise be available to it.

15. COPYRIGHT

All copyrights shall remain the property of the Company alone unless otherwise expressly agreed in writing with the Customer. All product images remain the property of the Company and may not be reproduced in print or electronic format without the express written permission of the Company.

16. RETAIL/RESELLING

Reselling of DJO branded products is forbidden without written consent from the Company.

17. PATENT RIGHTS

The Customer warrants that any instructions furnished or given by the Customer shall not be such as will cause the Company to infringe any letters patent or trademark or other intellectual property right in execution of the Customer's order and shall indemnify the Company against all actions, costs, claims and such demands arising therefrom to the extent that such instructions lead to such infringements.

18. GOVERNMENT CHARGES

The Customer will reimburse the Company for any increase in tax or governmental charge, or for any new tax or governmental charge hereafter becoming effective which the Company may be required to pay to any government upon the sale, production, or transportation of the Goods and which has the effect of increasing the cost of the Goods sold hereunder.

19. ASSIGNMENT

This contract is personal to the Customer and may only be assigned by it with the prior written consent of the Company.

20. BANKRUPTCY / LIQUIDATION

If the Customer being an individual (or, when the Customer is a firm, any partner in that firm) shall at any time become bankrupt, or having a receiving order made against them or shall make any composition or arrangement with or for the benefit of the creditors or if the Customer, being a company, an order for winding up is made or a resolution for winding up is passed or if a receiver or manager is appointed on behalf of a creditor or otherwise, the Company shall have the right to stop Goods in transit, to suspend further deliveries and to determine any subsisting contract with the Customer and (in any case where title to the Goods has to be passed to the Customer but payment is due and has not been made) the Company shall have a lien over the Goods and right of resale, provided that any such action shall be without prejudice to any right or claim the Company may have against the Customer.

21. CUSTOMER CARE

In our quest to "Never Stop Getting Better", we may from time to time monitor and record customer communications. This is for internal training purposes only and the confidentiality of our Customers and their transactions will be fully protected at all times.

22. PERSONAL CALLERS

We regret that we are unable to deal directly with personal callers to our office.

23. TERMS AND CONDITIONS

Terms and conditions quoted in this catalogue supersede all previous terms and conditions unless otherwise stated.

24. COMPANY DETAILS

DJO Orthopaedic South Africa (Pty) Ltd DJO Global

Unit 6, 2 on London, 2 London Circle, BrackenGate Business Park, Brackenfell, 7560, South Africa

T +27 (0) 21 276 2968 F +27(0) 86 609 8891

Company Registration Number: 2009/010124/07 VAT Registration Number: 4080257027

orders.zaf@djoglobal.com

Facebook: DJO South Africa

Instagram: @djo_southafrica

