



PHILIPS

Patient monitors

Affordable, reliable patient monitoring

Philips Goldway GS20
patient monitor



Key benefits of the GS20 at a glance

- **Easy to use and move around** thanks to lightweight, ergonomic design
- **Delivers accurate measurements** with an extensive range of market-proven functions and features, including Philips FAST SpO₂ and Philips ECG algorithms
- **Allows 2.5 hours of continuous monitoring** thanks to 4800 mAh, 11.1V rechargeable lithium ion battery
- **Bed-to-bed overview function** - including waveforms - for up to 32 monitors on the same local network
- **Supports a wide range of patient acuity levels** for neonatal, pediatric and adult patients
- **Facilitates a low total cost of ownership** with affordable servicing, parts and accessories

A compact, cost-effective approach to patient monitoring

A patient's condition can change rapidly, and clinicians need to respond swiftly to any changes that occur. Having a highly reliable, versatile patient monitoring solution on hand saves precious time, and provides caregivers with the insights they need to make informed clinical decisions.

Built according to our exacting quality standards the Philips GS20 preconfigured patient monitor provides advanced functionalities for various clinical environments and patient acuity levels. In addition to offering market-proven essential parameters, the GS20 system is cost-effective and easy to learn and use.

Ergonomic design, market-proven capabilities


The lightweight, ergonomic design of the GS20 monitor makes it highly practical to use. The monitor can be quickly and easily mounted on a roll stand and moved from bed to bed.

The GS20 offers essential measurements for daily use, such as the oxyCRG and Mini-trend software functions that provide basic clinical support. It can be connected to the Philips Efficia CMS200 Central Station, and includes an HL7 data output function for EMR/HIS integration. Moreover, the monitor features the Philips FAST SpO₂ and Philips ECG algorithms.

These essential market-proven capabilities facilitate accurate measurements, helping to guide clinical decision making.

A closer look at the **GS20 monitor**

Affordable and reliable patient monitoring solution with low total cost of ownership



An affordable patient monitoring solution

The GS20 monitor delivers market-proven Philips technology at an affordable price point, making it an attractive option for healthcare environments seeking to deliver outstanding patient care with a manageable level of capital investment and low total cost of ownership.

GS20 patient monitoring supports a wide range of patient acuity levels and provides customers with a cost-effective, easy-to-learn, and preconfigured bedside monitoring solution, with proven Philips technology.

