



System One

Category: CPAP

REMstar Auto



Optimizing CPAP prescription pressures at home

System One features

- Advanced event detection*
- Opti-Start feature
- A-Flex for comfort
- System One Resistance Control
- Mask fit check
- Waveform reports
- Oximetry module
- On-board memory
- Smart ramp
- Encore patient management system
- SD card and wired, wireless, and broadband modem capable

Humidification

- System One humidification control
- Pre-heat humidifier option
- Heated Tube humidifier option

REMstar Auto advantages

- CPAP, Auto, CPAP-Check, Auto-Trial modes
- A-Flex for comfort
- Opti-Start technology
- Advanced event detection
- Modem and SD card capable
- Optional, integrated oximetry available
- Heated Tube option

*Philips Respironics devices with advanced event detection detect, respond to (some models), and log these advanced events: Apnea Hypopnea Index (AHI), respiratory-related arousal (RERA), snore, flow limitation (FL), leak, vibratory snore (VS), obstructive airway apnea (OA), large leak (LL), clear airway apnea (CA), periodic breathing (PB), and hypopnea (H).

Home Healthcare
Respiratory Solutions

★ Company of the Year ★
2012-2014

awarded by
FROST & SULLIVAN

PHILIPS


RESPIRONICS

REMstar Auto specifications

Part number	Description
N561S	REMstar Auto A-Flex with SD card
IN561HS	REMstar Auto A-Flex with humidifier and SC card
IN561TS	REMstar Auto A-Flex with Heated Tube humidifier and SD card

Specifications	
Pressure range	4-20 cm H ₂ O
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	4 to CPAP; patient adjustable (CPAP mode) 4 to min pressure (Auto mode)
Dimensions	18 cm L x 14 cm W x 10 cm H
Weight	1.36 kg (without humidifier)
Filters	Foam and optional ultra-fine
Device set-up	LCD/control wheel/push button
Data storage capacity (minimum)	Display: compliance 30-day averages SD card > 1 year On-board: up to 6 months
Compliance meter	Breathing detection
Altitude compensation	Automatic
Electrical requirements	100-240 VAC, 50/60 Hz
Humidification	Integrates with System One 60 Series heated humidifier and System One 60 Series Heated Tube humidifier
Warranty	2 years
Miscellaneous	Lighted LEDs

Sleep therapy accessories	
Humidification	Part number
System One 60 Series humidifier	IN6H
System One 60 Series Heated Tube humidifier	IN6T
Water chamber (includes both 50 and 60 Series lids)	1063785
Heated Tube humidifier upgrade kit	6T-CONVKIT
Tubing	Part number
System One Heated Tube	SYS1HT15
System One Performance Tube (15mm)	1074400
Performance Tubing (22mm)	1032907
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1035443
Ultra-fine filter, disposable (1/pack)	1063096
Accessories	Part number
Device side panel	1099592
Carrying case	1063857
Power cord Europe	1039014
Power cord UK	1039026



Respironics, A-Flex, Auto-Trial, C-Flex, CPAP-Check, Encore, REMstar, and System One are trademarks of Koninklijke Philips N.V. and its affiliates. All rights reserved.

Please visit www.philips.com/respironics



© 2014 Koninklijke Philips N.V. All rights are reserved.

Philips Healthcare is part of Royal Philips

Philips Healthcare reserves the right to make changes in specifications and/or to discontinue any product at any time without notice or obligation and will not be liable for any consequences resulting from the use of this publication.

Broudy NB 12/19/13 MCI 4105940 PN 1113946

Philips Respironics Asia Pacific
+65 6882 5282
Philips Respironics Australia
+61 (2) 9947 0440
Philips Respironics Benelux
+33 1 47 28 30 83
Philips Respironics Deutschland
+49 8152 93 06 0
Philips Respironics France
+33 2 51 89 36 00
Philips Respironics Iberia
+34 91 566 90 00

Philips Respironics Italy
+39 039 203 1
Philips Respironics Sweden
+46 8 120 45 900
Philips Respironics Switzerland
+41 6 27 45 17 50
Philips Respironics United Kingdom
+44 800 1300 845
Philips Respironics
1010 Murry Ridge Lane
Murrysville, PA 15668 USA