








Clear, concise labeling on the package and pads guides experienced and infrequent caregivers alike through pad selection, placement, and connection.

## Technical specifications - ALS



	<b>989803166021</b>	<b>M3713A</b>	<b>M3716A</b>	<b>M3718A</b>	<b>M3501A</b>	<b>M3717A</b>	<b>M3719A</b>	<b>M3504A</b>	
<b>Description/name</b>	Adult/child preconnect	Adult/child Plus	Adult/child radiolucent	Adult/child radiotransparent	Adult/child with barrel connector	Infant Plus	Infant radiotransparent	Infant with barrel connector	
<b>Pads per case</b>	10	10	10	10	10	5	5	5	
<b>Indications</b>	Monitoring, external pacing, defibrillation, and cardioversion								
<b>Patient</b>						Infant	Infant	Infant	
Infant < 10 kg (22 lbs)									
Child 10 kg – 25 kg	Child	Child	Child	Child	Child				
Adult > 25 kg (55 lbs)	Adult	Adult	Adult	Adult	Adult				
<b>Conductor gel surface area (cm<sup>2</sup>)</b>	75	75	75	97	75	32	50	32	
<b>Approximate lead wire length (cm)</b>	120	120	120	120	60	60	60	60	
<b>Number of hours on skin</b>	24	24	24	24	24	24	24	24	
<b>Number of defibrillation shocks</b>	50	50	50	50	50	50	50	50	
<b>Duration of pacing (hours) at 100mA, 80 bpm</b>	8	8	8	3	8	8	1.5	8	
<b>Temperature range for storage</b>	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	15° – 35°C	
<b>Shelf life (in months) from date of manufacture (may be less upon receipt)</b>	24	24	24	24	24	24	18	24	
<b>Philips adapter compatibility</b>	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200		N/A	05-10000 05-10100	05-10000 05-10100	N/A

## Defibrillator compatibility

<b>Intrepid*</b> 867172; M3508A plug connector cable	●	●	●	●		●	●	
<b>Intrepid*</b> 867172; M3507A barrel connector cable					●			●
<b>DFM100*</b> 866199; M3508A plug connector cable		●	●	●		●	●	
<b>DFM100*</b> 866199; M3507A barrel connector cable					●			●
<b>XL+**</b> 861290; M3508A plug connector cable	●	●	●	●		●	●	
<b>XL+**</b> 861290; M3507A barrel connector cable					●			●
<b>MRx**</b> M3535A, M3536A; M3508A plug connector cable	○ <sup>2</sup>	○ <sup>2</sup>	○ <sup>2</sup>	○ <sup>2</sup>		○ <sup>1</sup>	○ <sup>1</sup>	
<b>MRx**</b> M3535A, M3536A; M3507A barrel connector cable					○ <sup>2</sup>			○ <sup>1</sup>
<b>FR3</b> 861388, 861389	○ <sup>3,4</sup>	○ <sup>3,4</sup>	○ <sup>3,4</sup>					
<b>FRx</b> 861304, 861305	○ <sup>5</sup>	○ <sup>5</sup>	○ <sup>5</sup>					
<b>FR2/FR2+</b> M3840A, M3841A, M3860A, M3861A, 861458, 861459	○ <sup>3,4</sup>	○ <sup>3,4</sup>	○ <sup>3,4</sup>	○ <sup>3,4</sup>				


\* Not available for sale in the USA.

\*\* MRx and XL+ have been discontinued and no longer for sale.


Not all products are available in all geographies. Please contact your Philips representative.

● Fully compatible  
○ Partially compatible


## Technical specifications - AEDs



989803158211 (1 set)  
989803158221 (5 sets)



989803139261



989803149981 (1 set)  
989803149991 (5 sets)



M3870A

Description/name	Defibrillator pads	SMART Pads II	SMART Pads III	FR2 Infant/child
<b>Pads per case</b>	1 or 5	1	1 or 5	1
<b>Indications</b>	Monitoring, external pacing, defibrillation, and cardioversion			Defibrillation
<b>Patient</b>				
Infant < 10 kg (22 lbs)	Infant	Infant	Infant	Infant
Child 10 kg – 25 kg	Child	Child	Child	Child
Adult > 25 kg (55 lbs)	Adult	Adult	Adult	
<b>Conductor gel surface area (cm<sup>2</sup>)</b>	85	80	80	85
<b>Approximate lead wire length (cm)</b>	120	120	120	120
<b>Number of hours on skin</b>	24	24	24	24
<b>Number of defibrillation shocks</b>	50	50	50	50
<b>Duration of pacing (hours) at 100mA, 80 bpm</b>	8	8	8	N/A
<b>Temperature range for storage</b>	0° – 50°C	0° – 50°C	0° – 50°C	0° – 32°C
<b>Shelf life (in months) from date of manufacture</b> (may be less upon receipt)	30	30	30	30
<b>Philips adapter compatibility</b>	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	05-10000 05-10100 05-10200	N/A

## Defibrillator compatibility

<b>Intrepid*</b> 867172; M3508A plug connector cable	●	●	●	
<b>Intrepid*</b> 867172; M3507A barrel connector cable				
<b>DFM100*</b> 866199; M3508A plug connector cable	●	●	●	
<b>DFM100*</b> 866199; M3507A barrel connector cable				
<b>XL**</b> 861290; M3508A plug connector cable	●	●	●	
<b>XL**</b> 861290; M3507A barrel connector cable				
<b>MRx**</b> M3535A, M3536A; M3508A plug connector cable	○ <sup>2</sup>	○ <sup>2</sup>	○ <sup>2</sup>	
<b>MRx**</b> M3535A, M3536A; M3507A barrel connector cable				
<b>FR3</b> 861388, 861389	○ <sup>4</sup>	●	●	
<b>FRx</b> 861304, 861305		●		
<b>FR2/FR2+</b> M3840A, M3841A, M3860A, M3861A, 861458, 861459	○ <sup>3,4</sup>	○ <sup>3,4</sup>	○ <sup>3,4</sup>	○ <sup>4</sup>

\* Not available for sale in the USA.

\*\* MRx and XL+ have been discontinued and no longer for sale.

Not all products are available in all geographies. Please contact your Philips representative.

1 Not in AED mode, as the energy automatically delivered is 150 Joules

2 Not for use with infants/children (10–25 kg / 22–55 lbs) in AED Mode

3 Infant/Child Key is required for use with infants/children (10–25 kg / 22–55 lbs) in Manual and AED Mode

4 Not preconnectable and cannot be left attached to device after use, when in Standby Mode

5 Infant/Child Key is required for use with infants/children (10–25 kg / 22–55 lbs) in AED Mode

© 2021 Koninklijke Philips N.V. All rights are reserved.

● Fully compatible  
○ Partially compatible

4522 991 60851 \* JAN 2021