OXYGEN ADMINISTRATION

Start Time:				
Stop Time:	Date:			
Candidate's Name:				
Evaluator's Name:			Points Possible	Points Awarded
Takes, or verbalizes, b	ody substance isolation precautions		1	
Assembles the regulator to the tank			1	
Opens the tank			1	
Checks for leaks			1	
Checks tank pressure			1	
Attaches non-rebreather mask to oxygen			1	
Prefills reservoir			1	
Adjusts liter flow to 12 liters per minute or greater			1	
Applies and adjusts the mask to the patient's face			1	
Note: The examiner i	must advise the candidate that the patient is	not tolerat	ing the	
non-rebreather mask	. The medical director has ordered you to ap	ply a nasa	l cannula	
to the patient.				
Attaches nasal cannula to oxygen			1	
Adjusts liter flow to 6 liters per minute or less			1	
Applies nasal cannula to the patient			1	
Note: The examiner i	must advise the candidate to discontinue oxy	gen thera	ру	
Removes the nasal cannula from the patient			1	
Shuts off the regulator			1	
Relieves the pressure within the regulator			1	
		Total:	15	
Critical Criteria				
Did not tai	ke, or verbalize, body substance isolation precau	utions		
Did not as	semble the tank and regulator without leaks			
Did not pr	efill the reservoir bag			
Did not adjust the device to the correct liter flow for the non-rebreather mask (12 liters per minute or greater)				
	ljust the device to the correct liter flow for the national minute or less)	sal cannula	ı	