



BHM-Tech Produktionsgesellschaft mbH  
 Grafenschachen 242  
 7423 Grafenschachen, Austria  
[www.bhm-tech.at](http://www.bhm-tech.at)

**CE** meets the requirements of  
 93/42 EEC | EN ISO 60601-1 | EN ISO 10993 | EN ISO 14971  
 Quality management system according to DIN EN ISO 13485:2016



*Perfection made in Austria*

*Your glasses as a hearing module*

*"Better hearing at last –  
 with my smart pair  
 of glasses!"*



Version 1.5 - 05/2021



It's always great to hear from you! Even for people with hearing loss, access to reliable hearing should be a matter of course, preferably elegant and inconspicuous. After all, the ability to communicate unhindered and be able to orient oneself is essential!

## *pan – as individual as you are*

Hear better – with glasses. **pan** is an aesthetic system that captures sound via two directional microphones and transfers it to the ear through a small loudspeaker and is easy to integrate into any pair of glasses. The experts at BHM have created an innovative product that is not discernible as a hearing module, thus setting completely new standards in terms of innovation and sophisticated technology.

Maximum choice – the utmost satisfaction. Today, a pair of glasses is more than just a visual aid: as a fashion accessory, it is, in addition to providing a great fit and excellent vision, an expression of your personality and lifestyle. So it's good to see how **pan** has made it possible to integrate a hearing aid into your glasses for the first time, in a way that provides the best



hearing quality without attracting any attention. Regardless of whether you would like to keep using your existing glasses, or prefer to opt for a new pair entirely.

## *pan can do more*

Good to see indeed. **pan** also provides further reassurance in a variety of ways. You can focus on the person you're talking to and enjoy the exchange.

**pan** automatically adapts to the acoustic level of the environment. Thanks to the adaptive directivity feature, you can even make out whispers very clearly against a loud background, without the ambient noise disappearing from your auditory sensation area. This also applies whilst on the phone: you will be able to hear the person on the other end without isolating yourself acoustically from the environment. Evenings at the opera or theatre will once again become an unforgettable event, as **pan** has an inductive signal input that allows you hear the sounds in all their beauty.

**pan** also provides relief for tinnitus. Silence is pretty much the worst for tinnitus patients. Annoying ringing in the ears can be treated with the support of **pan** in the context of a tinnitus therapy.



For single-sided deafness or severe impairment, we recommend CROS (contralateral routing of signals) or BiCROS, which ensures the safe transfer of sound from the affected ear to the healthy or less affected one.

Overall, **pan** provides a variety of operation modes, which are specially tailored to your requirements together with your hearing aid specialist. The individualised programme range allows you to adjust **pan** to your current specific hearing requirements flexibly at any time. This digital hearing module is not only customised on a visual level, but is also precisely tailored to your personal hearing requirements in sophisticated quality.

## *Advantages at a glance:*

- See and hear safely and comfortably – maximum comfort
- Easy to switch – fits flexibly to your glasses
- Perfect hearing experience even with ambient noise
- Sound in all its glory – also provides support in the treatment of tinnitus
- pan-CROS or BiCROS for safe transfer of sound – particularly for unilateral deafness