


BAO[®] **HEALTH**

BODY ARMOR OUTLET

PROTECTING THE FRONT LINE NO MATTER THE THREAT



Due to the market volatility pricing is subject to change without notice. Please contact your Government Account Manager for up to date pricing and availability.

BODY ARMOR OUTLET • 888-884-2427 • SALES@BODYARMOROUTLET.COM



SANITIZING WIPES

Hand Wipe Pack

- Hand Sanitizing Wipe Pack
- 70% Alcohol / Benzalkonium Chloride
- EMS Grade

ITEM CODE

EPHW-01

CERTIFICATES

FDA

MOQ

10,000

CARTON QUANTITY

80 PACKS/CTN

CARTON / PALLET

90



SANITIZING WIPES

Hand Wipe Canister

- Hand Sanitizing Wipe Pack
- 70% Alcohol / Benzalkonium Chloride
- EMS Grade

ITEM CODE

EPHW-02

CERTIFICATES

FDA

MOQ

10,000

CARTON QUANTITY

24 CANISTERS / CTN

CARTON / PALLET

60



SANITIZING WIPES

Surface Wipe Pack

- Hand Sanitizing Wipe Pack
- 70% Alcohol / Benzalkonium Chloride
- EMS Grade

ITEM CODE

EPSW-01

CERTIFICATES

FDA

MOQ

10,000

CARTON QUANTITY

40 PACKS/CTN

CARTON / PALLET

90



SANITIZING WIPES

Surface Wipe Canister

- Hand Sanitizing Wipe Pack
- 70% Alcohol / Benzalkonium Chloride
- EMS Grade

ITEM CODE

EPSW-02

CERTIFICATES

FDA

MOQ

10,000

CARTON QUANTITY

24 CANISTERS / CTN

CARTON / PALLET

60



SANITIZING SPRAY

Wexford Labs CleanCide RTU - Spray Bottle

- Kills 99.9% of household germs
- Eliminates odor-causing bacteria
- Sanitizes in as little as one minute
- Low pH
- Virucidal, germicidal, fungicidal and tuberculocidal efficacy
- Effective against antibiotic resistant organisms
- Effective against Norovirus

ITEM CODE

WXL-CCRTU

CERTIFICATES

N-List

MOQ

10,000

CARTON QTY

25 BOTTLES / CTN

CARTON / PALLET

80



SANITIZING SPRAY

Wexford Labs Thymo-Cide RTU - Spray Bottle

- Streak-free on glass, mirrors, stainless steal
- Leaves no residue
- Free from harsh chemicals
- Virucidal, germicidal and fungicidal efficacy in the presence of 5% blood serum
- Effective against antibiotic resistant organisms including MRSA & MRSE
- Effective against Norovirus
- Sanitizes in as little as one minute
- No rinse required (even on toys)

ITEM CODE

WXL-TCRTU

CERTIFICATES

EPA

MOQ

10,000

CARTON QTY

25 BOTTLES / CTN

CARTON / PALLET

80



HAND SANITIZER

2 fl oz

2 fl oz ScrewCap Pouch
Alcohol Antiseptic 70% Gel
Alcohol Antiseptic 80% Liquid

QUANTITY

20 /BOX

CASE SIZE

12.3 + 11.3 + 17.4

PALLET

640 BOXES (12,800 UNITS)



HAND SANITIZER

8 fl oz

8 fl oz Pump Top Bottle
Alcohol Antiseptic 70% Gel
Alcohol Antiseptic 80% Liquid

QUANTITY

36 /CASE

CASE SIZE

112.1 + 12.1 + 8.4

PALLET

60 CASES (2,160 UNITS)



HAND SANITIZER

16 fl oz

16 fl oz Pump Top Bottle
Alcohol Antiseptic 70% Gel
Alcohol Antiseptic 80% Liquid

QUANTITY

16 /CASE

CASE SIZE

112.1 + 12.1 + 8.4

PALLET

64 CASES (1,024 UNITS)



HAND SANITIZER

Gallon

Gallon Jug w/Screw Cap / Pump Tool
Alcohol Antiseptic 70% Gel
Alcohol Antiseptic 80% Liquid

QUANTITY

4 /CASE

CASE SIZE

12 + 12 + 11.88

PALLET

48 CASES (192 UNITS)



FACE MASKS

Surgical Mask

3-ply SurgicalMask
Size: 175 x95mm
2Layers of Non-woven Fabric
MedicalMelt-blownFabric
Medical Surgical
1-year Validity Period
Surgical Use

FACE MASKS

LEVEL 1

ITEM CODE

3PLY-ASTM L1

DESCRIPTION

3PLY LEVEL 1 FACE MASK
W/EAR LOOP

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic
blood minimum pressure in mm Hg
for pass result . 80mm Hg

Differential pressuremmH20/cm2
(Breathability):<4.0

Bacterial filtration efficiency: 2 ≥ 95%

Sub-micron particulates filtration
efficient at 0.1 micron: ≥ 95%

CASE QTY

50 pcs/box
2000 pcs/carton

PACKAGE

Box: 7"x4"x4"
Carton: 21"x 16"x16"

FACE MASKS

LEVEL 2

ITEM CODE

3PLY-ASTM L2

DESCRIPTION

3PLY LEVEL 2 FACE MASK
W/EAR LOOP

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic
blood minimum pressure in mm Hg
for pass result . 120mm Hg

Differential pressuremmH20/cm2
(Breathability):<5.0

Bacterial filtration efficiency: 2 ≥ 98%

Sub-micron particulates filtration
efficient at 0.1 micron: ≥ 98%

CASE QTY

50 pcs/box
2000 pcs/carton

PACKAGE

Box: 7"x4"x4"
Carton: 21"x 16"x16"

FACE MASKS

LEVEL 3

ITEM CODE

3PLY-ASTM L3

DESCRIPTION

3PLY LEVEL 3 FACE MASK
W/EAR LOOP

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic
blood minimum pressure in mm Hg
for pass result . 160mm Hg

Differential pressuremmH20/cm2
(Breathability):<5.0

Bacterial filtration efficiency: 2 ≥ 98%

Sub-micron particulates filtration
efficient at 0.1 micron: ≥ 98%

CASE QTY

50 pcs/box
2000 pcs/carton

PACKAGE

Box: 8.2"x4.5"x4.5"
Carton: 23"x20"x20"

FACE MASKS

LEVEL 3

ITEM CODE

3PLY-ASTM L3T

DESCRIPTION

3PLY LEVEL 3 FACE MASK
W/HEAD TIES

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic
blood minimum pressure in mm Hg
for pass result . 160mm Hg

Differential pressuremmH20/cm2
(Breathability):<5.0

Bacterial filtration efficiency: 2 ≥ 98%

Sub-micron particulates filtration
efficient at 0.1 micron: ≥ 98%

CASE QTY

50 pcs/box
2000 pcs/carton

PACKAGE

Box: 8.2"x4.5"x4.5"
Carton: 23"x20"x20"



FACE MASKS

KN95Mask

- Adult/Child
- PP Non-woven Fabric Material
- 5 Layers of Protection
- High-efficiency Filtering
- Adjustable Nose Bridge

ITEM CODE

KN-95

STANDARD

GB2626-2006

MOQ

100,000 pcs

CASE QTY

**50 pcs/box
30 box/ctn**

PACKAGE

**Carton Size: 27"x19"x16"
Gross Weight: 27lbs**



FACE MASKS

N95 NIOSH Mask

- Adult
- PP Non-woven Fabric Material
- 5 Layers of Protection
- High-efficiency Filtering
- Adjustable Nose Bridge

ITEM CODE

N95 L-188

STANDARD

NIOSH

MOQ

100,000 pcs

CASE QTY

**20 pcs/box
20 box/ctn**

PACKAGE

**Carton Size: 16"x16"x14
Gross Weight: 11lbs**



FACE MASKS

N95 Hard Shell Molded Mask

- Adult
- PP Non-woven Fabric Material
- 5 Layers of Protection
- High-efficiency Filtering
- Adjustable Nose Bridge

ITEM CODE

N95 L-288

STANDARD

NIOSH

MOQ

100,000 pcs

CASE QTY

**10 pcs/box
40 box/ctn**

PACKAGE

**Carton Size: 26"x11"x16.5"
Gross Weight: 14lbs**

Brands Offered

- Makrite
- Great Hawk
- Kimtech
- Indiana Face Mask





GOWN

LEVEL 1

ITEM CODE

GOWN-L1

DESCRIPTION

**PEVA 40 GSM,
ISOLATION GOWN**

CASE QTY

**1 pcs/polybag
200 bags/carton**

NW / GW

42lbs / 45lbs

PACKAGE

**Bag: 12"x12"x16"
Carton: 22"x12"x19"**



GOWN

LEVEL 2

ITEM CODE

GOWN-L2

DESCRIPTION

**PP+PE 45 GSM,
ISOLATION GOWN**

CASE QTY

**1 pcs/polybag
110 bags/carton**

NW / GW

42lbs / 45lbs

PACKAGE

**Bag: 12"x12"x16"
Carton: 22"x12"x19"**



GOWN

LEVEL 3

ITEM CODE

GOWN-L3

DESCRIPTION

**45 GSM SMS
ISOLATION GOWN**

CASE QTY

200 bags/carton

NW / GW

40lbs / 45lbs

PACKAGE

Carton:24"x14"x20"



GOWN PERFORMANCE

LEVEL 2 - STERILIZED

ITEM CODE

FQS172

DESCRIPTION

45GSM.PP+PE GOWN

CASE QTY

100 pcs/carton

NW / CW

32lbs / 35lbs

PACKAGE

Carton: 24"x15"x20"



GOWN PERFORMANCE

LEVEL 3- STERILIZED

ITEM CODE

FQS173

DESCRIPTION

45GSM.PP+PE GOWN

CASE QTY

100 pcs/carton

NW / CW

32lbs / 35lbs

PACKAGE

Carton: 24"x15"x20"



GOWN PERFORMANCE

LEVEL 2 - NON STERILIZED

ITEM CODE

FQS174

DESCRIPTION

45GSM.PP+PE GOWN

CASE QTY

100 pcs/carton

NW / CW

32lbs / 35lbs

PACKAGE

Carton: 24"x15"x20"



GOWN PERFORMANCE

LEVEL 3 - NON STERILIZED

ITEM CODE

FQS175

DESCRIPTION

45GSM.PP+PE GOWN

CASE QTY

100 pcs/carton

NW / CW

32lbs / 35lbs

PACKAGE

Carton: 24"x15"x20"



PROTECTIVE COVERING

Disposable Chemo Gown

ITEM CODE

FQS181

DECRPTION

65-75 GSM, CHEMO GOWN
SMS COATING

CASE QTY

100 pcs /carton

NW / GW

40lbs /45lbs

PACKAGE

Carton: 24"x15"x20"



PROTECTIVE COVERING

Reusable Surgical Gown

- Sizes:S/M /L/XL/XXL/XXXL
- Color: Blue/ Yellow
- Base Weight: 22-65 gsm
- Wrist: Thumb Loop / Elastic
- Medical Use

ITEM CODE

FQS199

STANDARD

LEVEL 2



PROTECTIVE COVERING

Disposable Coverall

- Sizes: 170/M, 175/L, 180/XL, 185/XXL
- Radiation Sterilization: Type 5/6, Weight: 62-65g
- Waterproof against blood, bacteria, virus, antistatic, breathable, whole body protection
- Medical Use

COVERALL

LEVEL 3

ITEM CODE

NDDG-20465

STANDARD

LEVEL 3

MOQ

50,000 pcs

CASE QTY

Vacuum Packed, 40 pcs /ctn

PACKAGE

Carton Size: 25"x15"x20"
Gross Weight: 33lbs



MEDICAL GLOVES

Nitrile Gloves

- Specifications: Powder-Free Non-Sterile
- Disposable, Universal for Right and Left Hand
- Sizes: S/M/L - 3 Year Shelf Life
- Water Impervious: G-1 / AQL 1.5
- Weight: 0.16oz (Size M) / Length: 9.5in
- Medical Use

ITEM CODE

GLOVES-NSUP

STANDARD

ASTM D6319

MOQ

50,000 pairs



MEDICAL GLOVES

Latex Gloves

- Specifications: Powder-Free Non-Sterile
- Disposable, Universal for Right and Left Hand
- Sizes: S/M/L - 3 Year Shelf Life
- Water Impervious: G-1 / AQL 1.5
- Weight: 0.16oz (Size M) / Length: 9.5in
- Medical Use

ITEM CODE

GLOVES-LX

STANDARD

ASTM D3577

MOQ

50,000 pairs



MEDICAL GLOVES

Vinyl Gloves

- Specifications: Powder-Free Non-Sterile
- Disposable, Universal for Right and Left Hand
- Sizes: S/M/L - 3 Year Shelf Life
- Water Impervious: G-1 / AQL 1.5
- Weight: 0.16oz (Size M) / Length: 9.5in
- Medical Use

ITEM CODE

GLOVES-VL

STANDARD

ASTM D3577

MOQ

50,000 pairs



SANITIZER DISPENSER

33 fl oz

- Automatic
- Size: 5 x 11 x 10 in
- Weight: 2.1 lbs

SPECIFICATION

- Suitable for 3 different liquids - Foam, Soap, Sanitizer
- Refillable bottle (100ml) or disposable bag (800ml)
- Lockable with window design
- With UV disinfection can kill virus on hands
- Can use DC power or battery
- Battery supports up to 30,000 pushes
- Light indicator for low battery and leak control



SANITIZER DISPENSER

40 fl oz

- Automatic
- Available with or without Water Tray
- Size: 6 x 11 x 11 in / 6 x 11 x 18 in
- Weight: 2.1 lbs

SPECIFICATION

- Suitable for 3 different liquids - Foam, Soap, Sanitizer
- Refillable bottle (1200ml) or disposable bag (1000ml)
- Lockable with window design
- With UV disinfection can kill virus on hands
- Can use DC power or battery
- Battery supports up to 30,000 pushes
- Light indicator for low battery and leak control



SANITIZER DISPENSER

Adjustable

- Size: 12 x 12 x 64 in
- Weight: 10.8 lbs

SPECIFICATION

- Suitable for 3 different liquids - Foam, Soap, Sanitizer
- Refillable bottle (1000ml / 1200ml) or disposable bag (800ml / 1000ml)
- Lockable with window design
- With UV disinfection can kill virus on hands
- Can use DC power or battery
- Battery supports up to 30,000 pushes
- Light indicator for low battery and leak control
- Adjustable height (0.7m to 1.2m from floor to sensor)



SANITIZER DISPENSER

Non-Adjustable

- Size: 11 x 13 x 48 in
- Weight: 10 lbs

SPECIFICATION

- Suitable for 3 different liquids - Foam, Soap, Sanitizer
- Refillable bottle (1200ml) or disposable bag (1000ml)
- Lockable with window design
- With UV disinfection can kill virus on hands
- Can use DC power or battery
- Battery supports up to 30,000 pushes
- Light indicator for low battery and leak control



EYE + FACE PROTECTION

Safety Face Shield

ITEM CODE

S-01

DESCRIPTION

MEDICAL FACE SHIELD

STANDARD

GB14866-2006

CASE QTY

10 pcs /box
300 pcs /carton

PACKAGE

Bag: 22"x15"
Carton: 27"x20"x12"



EYE + FACE PROTECTION

Safety Face Shield

ITEM CODE

S-04

DESCRIPTION

NEW FACE SHIELD

STANDARD

GB14866-2006

CASE QTY

10 pcs /box
250 pcs /carton

PACKAGE

Bag: 22"x15"
Carton: 27"x20"x12"



EYE PROTECTION

Anti Fog Goggles

ITEM CODE

ZTAF-1106

DESCRIPTION

MEDICAL GOGGLES

STANDARD

CB32166.1-2016

CASE QTY

1 pcs /box
1000 pcs /carton

PACKAGE

Bag: 7"x13"x13"
Carton: 20"x40"x45"



SHOE COVER

Disposable Shoe Cover

ITEM CODE
SHOECOV-B
DESCRIPTION
PP+PE 45GSM SHOE COVER
CASE QTY
25 pairs/ inner box 750 pairs/carton
NW / GW
39lbs /42lbs
PACKAGE
Box: 8.5"x4.5"x8.3" Carton: 25.5"x21"x18"



SHOE COVER

Disposable Shoe Cover

ITEM CODE
SQS-XT07
DESCRIPTION
35GSM SMS SHOE COVER WITH ANTI SKID
CASE QTY
500 pairs/carton
NW / GW
39lbs /42lbs
PACKAGE
Carton: 24"x15"x20"



SHOE COVER

Insulated Shoe Cover

ITEM CODE
SQS-XT07-1
DESCRIPTION
45GSM, PP+PE SHOE COVER WITH ANTI SKID
STANDARD
LEVEL 1
CASE QTY
500 pairs/carton
NW / GW
39lbs /42lbs
PACKAGE
Carton: 24"x15"x20"



BOOT COVER

Disposable Boot Cover

ITEM CODE
SQS-XT05
DESCRIPTION
35GSM, SMS SHOE COVER WITH ANTI SKID
STANDARD
LEVEL 3
CASE QTY
400 pairs/carton
NW / GW
39lbs /42lbs
PACKAGE
Carton: 24"x15"x20"



BOOT COVER

Disposable Boot Cover

ITEM CODE
SQS-XT05-1
DESCRIPTION
35GSM, SMS SHOE COVER WITH ANTI SKID
STANDARD
LEVEL 2
CASE QTY
400 pairs/carton
NW / GW
39lbs /42lbs
PACKAGE
Carton: 24"x15"x20"



HEAD COVER

Disposable Surgical Cap

ITEM CODE

SQS-Mo6

DESCRIPTION

35GSM, SMS
SURGICAL CAP

CASE QTY

600 pcs/carton

NW / GW

38 lbs/ 42 lbs

PACKAGE

Carton: 20"x14"x20"



HEAD COVER

Disposable Head Cover

ITEM CODE

BUFCAP-B

DESCRIPTION

PP+PE 45 GSM
BOUFFANT CAP

CASE QTY

35 pcs / inner box
1050 pcs / carton

NW / GW

30 lbs / 34 lbs

PACKAGE

Box: 8.5"x5"x8.5"
Carton: 25"x22"x18"



HEAD COVER

Disposable Head Cover

ITEM CODE

BUFCAP-C

DESCRIPTION

25 GSM, NON-WOVEN
BOUFFANT CAP

CASE QTY

100 pcs /carton

NW / GW

27 lbs / 31 lbs

PACKAGE

Box: 13"x8.5"x13"
Carton:25"x22"x17"



Lab Coat

ITEM CODE

FQS182

DESCRIPTION

45 GSM, PP+PE LAB COAT

STANDARD

LEVEL 1

CASE QTY

200 pcs/carton

NW / GW

40 lbs / 43 lbs

PACKAGE

Carton: 24"x15"x20"



Lab Coat

ITEM CODE

FQS183

DESCRIPTION

**45 GSM, NON-
WOVEN LAB COAT**

CASE QTY

200 pcs/carton

NW / GW

40 lbs / 43 lbs

PACKAGE

Carton: 24"x15"x20"



Seat Cover

ITEM CODE

JN-YT01

DESCRIPTION

**45 GSM, NON-WOVEN
CHAIR CUSHION**

CASE QTY

100 pcs/carton

NW / GW

40 lbs / 45 lbs

PACKAGE

Carton: 24"x15"x20"



FACE MASKS

KN95 Mask

- Child
- PP Non-woven Fabric Material
- 5 Layers of Protection
- High-efficiency Filtering
- Adjustable Nose Bridge

ITEM CODE

KN-95

STANDARD

GB2626-2006

MOQ

100,000 pcs

CASE QTY

**50 pcs/box
30 box/ctn**

PACKAGE

**Carton Size:28"x19"x17"
GrossWeight:27 lbs**



BOYS + GIRLS FACE MASKS

3 Ply

ITEM CODE

BAO-MASK3KID-50

DESCRIPTION

3 ply masks

CASE QTY

50 pcs/box



THERMOMETER

Infrared Thermometer

- Model: FTD-802, 9v Power Supply, Simple One-handed Operation, Laser Targeting, LCD with Back Light
- Measuring Range: Body Mode 32°C - 42°C
- Measuring Accuracy: Infrared Body 0.1 - 0.3°C, Surface $\pm 1.5\%$
- Measuring Distance: Infrared 3 - 5cm, Surface Distance Coefficient Ratio 12:1
- Response Time: 0.8 seconds
- Auto power off: ≤ 15 Seconds
- Transportation Storage: Ambient Temp, 20-35°C, Relative Humidity $\leq 93\%$
- Working Environment: Ambient Temp, 16-35°C, Relative Humidity $\leq 85\%$
- Medical Use

ITEM CODE

SK-T008

STANDARD

MEDICAL

MOQ

10,000 pcs

CASE QTY

40 pcs/ctn

PACKAGE

Box Size: 22"x13"x12"
Weight: 22 lbs



TESTING KIT

Universal Transport Media Nasal Swab Kit (UTM)3ML

- Virus Transport/Sampling Kit
- Includes protein for stabilization
- Antibiotics to minimize bacterial / fungal contamination
- Sterile
- Vials designed for transport of clinical sample
- Contains swabs for testing use
- Nasopharyngeal swab collection
- Oropharyngeal swab collecton

ITEM CODE	CERTIFICATES
UTM-NSK3M	MSDS
	MOQ
	10,000
	PACKQTY
	1 per pack
	CARTON / BOX
	500



SYRINGE

Safety Retractable Syringe

- Vaccine Safety Syringe
- Syringe Volume: 0.5ML/1ML or 3ML/5ML
- Needle Gauge: 23G, 25G (16mm-25mm)
- Luer connector: Integrated or Luer lock

ITEM CODE	CERTIFICATES
LH-SRS23G	FDA / CFDA
	MOQ
	5,000
	CARTONQTY
	50 SYRINGES / CTN
	CARTON / BOX
	500