



THERMOMETER

## Infrared Thermometer

- Model: FTD-802, 9v Power Supply, Simple One-handed Operation, Laser Targeting, LCD with Back Light
- Measuring Range: Body Mode 32°C - 42°C
- Measuring Accuracy: Infrared Body 0.1 - 0.3°C, Surface ± 1.5%
- Measuring Distance: Infrared 3 - 5cm, Surface Distance Coefficient Ratio 12:1
- Response Time: 0.8 seconds
- Auto power off: ≤ 15 Seconds
- Transportation Storage: Ambient Temp, 20-35°C, Relative Humidity ≤93%
- Working Environment: Ambient Temp, 16-35°C, Relative Humidity ≤85%
- Medical Use

ITEM CODE

**SK-Too8**

STANDARD

**MEDICAL**

MOQ

**10,000 pcs**

CASE QTY

**40 pcs/ctn**

PACKAGE

**Box Size: 22"x13"x12"**  
**Weight: 22 lbs**