



FACE MASKS

Surgical Mask

3-ply Surgical Mask
Size: 175 x 95mm
2 Layers of Non-woven Fabric
Medical Melt-blown Fabric
Medical Surgical
1-year Validity Period
Surgical Use

FACE MASKS

LEVEL 1

ITEM CODE

3PLY-ASTM L1

DESCRIPTION

**3PLY LEVEL 1 FACE MASK
W/EAR LOOP**

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic blood minimum pressure in mm Hg for pass result - 80mm Hg

Differential pressure mmH2O/cm2 (Breathability): <4.0

Bacterial filtration efficiency: $2 \geq 95\%$

Sub-micron particulates filtration efficient at 0.1 micron: $\geq 95\%$

CASE QTY

**50 pcs/box
2000 pcs/carton**

PACKAGE

**Box: 7"x4"x4"
Carton: 21"x 16"x16"**

FACE MASKS

LEVEL 2

ITEM CODE

3PLY-ASTM L2

DESCRIPTION

**3PLY LEVEL 2 FACE MASK
W/EAR LOOP**

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic blood minimum pressure in mm Hg for pass result - 120mm Hg

Differential pressure mmH2O/cm2 (Breathability): <5.0

Bacterial filtration efficiency: $2 \geq 98\%$

Sub-micron particulates filtration efficient at 0.1 micron: $\geq 98\%$

CASE QTY

**50 pcs/box
2000 pcs/carton**

PACKAGE

**Box: 7"x4"x4"
Carton: 21"x 16"x16"**

FACE MASKS

LEVEL 3

ITEM CODE

3PLY-ASTM L3

DESCRIPTION

**3PLY LEVEL 3 FACE MASK
W/EAR LOOP**

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic blood minimum pressure in mm Hg for pass result - 160mm Hg

Differential pressure mmH2O/cm2 (Breathability): <5.0

Bacterial filtration efficiency: $2 \geq 98\%$

Sub-micron particulates filtration efficient at 0.1 micron: $\geq 98\%$

CASE QTY

**50 pcs/box
2000 pcs/carton**

PACKAGE

**Box: 8.2"x4.5"x4.5"
Carton: 23"x20"x20"**

FACE MASKS

LEVEL 3

ITEM CODE

3PLY-ASTM L3T

DESCRIPTION

**3PLY LEVEL 3 FACE MASK
W/HEAD TIES**

STANDARD

ASTM F2100-2019

Resistance to penetration by synthetic blood minimum pressure in mm Hg for pass result - 160mm Hg

Differential pressure mmH2O/cm2 (Breathability): <5.0

Bacterial filtration efficiency: $2 \geq 98\%$

Sub-micron particulates filtration efficient at 0.1 micron: $\geq 98\%$

CASE QTY

**50 pcs/box
2000 pcs/carton**

PACKAGE

**Box: 8.2"x4.5"x4.5"
Carton: 23"x20"x20"**