



PROTECTIVE COVERING

Disposable Chemo Gown

ITEM CODE

FQS181

DESCRIPTION

**65-75 GSM, CHEMO GOWN
SMS COATING**

CASE QTY

100 pcs /carton

NW / GW

40lbs /45lbs

PACKAGE

Carton: 24"x15"x20"



PROTECTIVE COVERING

Reusable Surgical Gown

- Sizes: S/M /L/XL/XXL/XXXL
- Color: Blue/ Yellow
- Base Weight: 22-65 gsm
- Wrist: Thumb Loop / Elastic
- Medical Use

ITEM CODE

FQS199

STANDARD

LEVEL 2



PROTECTIVE COVERING

Disposable Coverall

- Sizes: 170/M, 175/L, 180/XL, 185/XXL
- Radiation Sterilization: Type 5/6, Weight: 62-65g
- Waterproof against blood, bacteria, virus, antistatic, breathable, whole body protection
- Medical Use

COVERALL

LEVEL 3

ITEM CODE

NDDG-20465

STANDARD

LEVEL 3

MOQ

50,000 pcs

CASE QTY

Vacuum Packed, 40 pcs /ctn

PACKAGE

**Carton Size: 25"x15"x20"
Gross Weight: 33lbs**