



Pipette Calibration

Calibration Levels

	BASIC	BRONZE	SILVER	GOLD
Exterior inspection & cleaning	Yes	Yes	Yes	Yes
Leak test	Yes	Yes	Yes	Yes
Volumes from .5uL, and up	Yes	Yes	Yes	Yes
Test volumes	Low/High	Low/High	Low/Mid/High	Low/Mid/High
Samples per volume	4	4	4	10
Fixed volume	Yes	Yes	Yes	Yes
Variable volume	Yes	Yes	Yes	Yes
Calibration sticker	Yes	Yes	Yes	Yes
Calibration report	No	Yes	Yes	Yes
Single channel price*	\$19	\$30	\$50	\$80
Multichannel (8, 12 or 16) price*	\$75	\$120	\$180	\$300

*Price per pipette for depot service. Price does not include return shipping and/or parts. Prices subject to change without notice

AccuCal Services 436-1641 Lonsdale Avenue, North Vancouver, BC V7M 2J5
Tel: (604) 355-3158 **Email:** accucalservices@gmail.com **Website:** www.accucalservices.com

Name: _____

Address: _____

Phone number: _____

E-mail: _____

Calibration Level: ☐ BASIC ☐ BRONZE ☐ SILVER ☐ GOLD

Make	Model & Serial Number

☐ Purchase Order: _____

☐ Credit Card (do not send credit card number. An email will be sent for credit card payment)

☐ All equipment has been decontaminated prior to shipping/pickup

AccuCal Services 436-1641 Lonsdale Avenue, North Vancouver, BC V7M 2J5
Tel: (604) 355-3158 **Email:** accucalservices@gmail.com **Website:** www.accucalservices.com