

8. Accessories

The *Flowtron* ACS900 pump should only be used with the following garments:

| CALF GARMENTS | | | |
|-----------------------|---|--------------------|------------|
| Order Code | Type | Calf Circumference | Therapy |
| DVT10 | DVT10 Standard Calf Garment | Up to 43 cm (17") | Uniform |
| DVT10S ^(a) | DVT10S Standard Calf Garment (Sterile) | Up to 43 cm (17") | Uniform |
| L501-M | L501-M Standard Calf Garment | Up to 43 cm (17") | Uniform |
| DVT20 | DVT20 Large Calf Garment | Up to 58 cm (23") | Uniform |
| DVT60L | DVT60L Bariatric Calf Garment | Up to 81 cm (32") | Uniform |
| TRP10 | Tri Pulse TRP10 Regular Calf Garment | Up to 43 cm (17") | Sequential |
| TRP20 | Tri Pulse TRP20 Large Calf Garment | Up to 58 cm (23") | Sequential |
| TRP60L | Tri Pulse TRP60L Bariatric Calf Garment | Up to 81cm (32") | Sequential |

| CALF & THIGH GARMENTS | | | |
|-----------------------|---|---------------------|------------|
| Order Code | Type | Thigh Circumference | Therapy |
| DVT30 | DVT30 Standard Thigh Garment | Up to 71cm (28") | Uniform |
| DVT30S ^(a) | DVT30S Standard Thigh Garment (Sterile) | Up to 71cm (28") | Uniform |
| L503-M | L503-M Standard Thigh Garment | Up to 71cm (28") | Uniform |
| DVT40 | DVT40 Large Thigh Garment | Up to 89cm (35") | Uniform |
| TRP30 | Tri Pulse TRP30 Regular Thigh Garment | Up to 71cm (28") | Sequential |
| TRP40 | Tri Pulse TRP40 Large Thigh Garment | Up to 89cm (35") | Sequential |

| FOOT GARMENTS | | | |
|-----------------------|----------------------------------|--|---------|
| Order Code | Type | Shoe Size | Therapy |
| FG100 | Foot Garment - Regular | UK Men/Women up to size 7 US Women up to size 9 US Men up to size 7 EURO up to size 40 | Uniform |
| FG100S ^(a) | Foot Garment - Regular (Sterile) | UK Men/Women up to size 7 US Women up to size 9 US Men up to size 7 EURO up to size 40 | Uniform |
| FG200 | Foot Garment - Large | UK Men/Women size 7½ or above US Women size 9½ or above US Men size 7½ or above EURO size 41 or above | Uniform |
| FG200S ^(a) | Foot Garment - Large (Sterile) | UK Men/Women size 7½ or above US Women size 9½ or above US Men size 7½ or above EURO size 41 or above | Uniform |

a. Please check with your local ArjoHuntleigh representative for availability of sterile garments.