

Amtryke Adaptive Tricycle Request *Form*



Please Select Rider Type: ☐ Adult ☐ Child (Under 18 must be completed by parent/guardian)

To apply for an Amtryke please complete all questions.

Recipient Name		Date of Request:	
Street Address			
City	State	Zip	County
DOB	Height	Weight	☐ Male ☐ Female
Diagnosis			

Complete if recipient is Under 18

Parent/Guardian			
Street Address			
City:	State:	Zip:	Phone:
Email:			
Secondary Contact:		Phone:	
Treating Therapist:		Clinic Location:	
Phone:		Email:	

Are you a first-time Amtryke/Adaptive Bike recipient? ☐ Yes ☐ No

Do you have a place to store the Amtryke/Adaptive Bike out of the weather? ☐ Yes ☐ No

Do you have a vehicle to transport the Amtryke/Adaptive Bike? ☐ Yes ☐ No

*If NO, do you have a safe place to ride the Amtryke/Adaptive Bike? ☐ Yes ☐ No

Tell us about the recipient:

Iron City Trykes/Birmingham AMBUCS has a suggested \$100 donation, but this is not required in order to receive your bike. We will provide each accepted recipient with a fundraising link that can be shared with family, friends, and your community to assist us cover the costs of your bike and others. Our bikes average between \$900-\$1900, so any amount contributed through your efforts is greatly appreciated and will help us impact more riders.

**Our bikes are distributed based on available funds, current inventory, and need. Individual placements of the bikes are at the discretion of the local chapter.*

Individuals will not be considered for placement until this application is complete and submitted including: request form, waiver, assessment & Bike Selection Forms (completed by PT, OT, or other Qualified Professional).

EMAIL COMPLETED APPLICATION TO: birminghamambucs3@gmail.com

AmTryke Adaptive Tricycle Waiver Form

(Must be filled out completely by adult rider or parent/guardian)

AMBUCS members nationwide are dedicated to creating opportunities for mobility and independence by providing Amtryke adaptive tricycles, offering educational scholarships to therapy students and performing various forms of community service.

Purpose: The Amtryke adaptive tricycle creates a feeling of freedom, builds self-esteem, strengthens muscles and improves motor coordination and range of motion—all while making exercise fun.

Steering: Initially, the rider may have difficulty turning or changing directions. Encourage the rider to go straight ahead, back up and slowly turn around. On many models there are three steering options for the Amtryke. On the front column of the tricycle, you will find two holes for the steering pin. The top hole is straight steering, the bottom hole allows a 20-degree turning radius. Leaving the pin out gives the rider free steering.

Safety Cautions

- Fast speeds and sharp turns can cause the Amtryke adaptive tricycle to tip or turn over.
- Always wear a helmet when riding an Amtryke or any bike. Use of other protective gear is highly recommended.
- Adult supervision required if used by younger or developmentally delayed riders.
- Use caution near vehicles, swimming pools and other bodies of water, hills, alleys and sloped driveways.
- Always wear closed toe shoes.
- Never allow more than one rider.
- Use of the steering pin is recommended to prevent over-steering or possible tip-overs.

The information contained in this service is not intended nor implied by National AMBUCS™, Inc. to be professional medical advice by National AMBUCS, Inc. Always seek the advice of your physician, therapist or other qualified healthcare provider prior to starting any treatment or with any question you may have regarding a medical condition. Nothing contained in this document is intended by National AMBUCS, Inc. to be for medical diagnosis or treatment by National AMBUCS, Inc. or on behalf of National AMBUCS, Inc.

In no event shall National AMBUCS, Inc. or Iron City Trykes/Birmingham AMBUCS be liable for any direct, indirect, incidental, consequential, special, exemplary, punitive, or any other monetary or other damages, personal injury or property damages, fees, fines, costs, attorney fees, or liabilities of any kind arising out of or relating in any way to this service or use of the Amtryke® adaptive tricycle or any other bike provided, and/or content or information provided herein.

☐ I agree that the rider listed on the front of this document may be photographed. I also agree that the photos and recipient's name may be used in promotional efforts for National AMBUCS, Amtryke or the local AMBUCS chapter. I further grant AMBUCS the ability to use the photos and name for advertising/publicity purposes without additional compensation, except where prohibited by law. If anybody in my party does not want to be photographed under these same terms, I will let the photographer know as soon as possible.

*Please email a photo of the recipient to ironcitytrykes@gmail.com (By sending a photo, you are giving consent for Iron City Trykes/Birmingham AMBUCS to use the photo online, in print, digitally to help us promote our work.

By signing below, I acknowledge that I have read and understood this liability waiver.

Signature

If Recipient is Under Age 18 Parent /Guardian MUST sign below

Signature

Individuals will not be considered for placement until all three forms are returned to AMBUCS: this form filled out by the recipient/guardian, as well as the Assessment Form and Tryke Selection Form filled out by the therapist.

This form, Request Form & Liability Waiver Form and Assessment Form must be received before Amtryke placement is considered.

(205) 823-7800
www.ironcitytrykes.org

EMAIL COMPLETED APPLICATION TO: birminghamambucs3@gmail.com

AmTryke Assessment Form


(Must be filled out completely by therapist or physician)

Recipient's Name	Height	Weight
Diagnosis		

*This private information is only used to help appropriately fit the rider.

Is the recipient ambulatory? ☐ Yes ☐ No Right Handed ☐ Left Handed ☐

RIDER'S MEASUREMENTS				
Arm Measurements (inches)				Total Length
Left	A to B:		B to C:	
Right	A to B:		B to C:	
Trunk	A to D:			
Leg Measurements (inches)				Total Length
Left	D to E:		E to F:	
Right	D to E:		E to F:	



A	Center of Shoulder
B	Center of Elbow
C	Center of Digit Crease
D	Center of Hip
E	Center of Knee
F	Bottom of Foot

Arm Length & Leg Length Measurements are critical to correct Amtryke Selection Sizing Chart is available online: www.ambucs.org/riders/wish-list/sizing-chart/

Notes on Provided Measurements (if any):

--

Helmet Sizing	
Size	Measurement (head circumference)
Toddler (XS)	17.7" – 19.3"
Child (S)	20.5" – 21.7"
Youth (L)	20.9" – 22.4"
Adult (XL)	22.4" – 23.6"

Therapist Name	Is this the treating therapist? ☐ Yes ☐ No
Credentials	Facility
Phone	Email

	Yes	No	Notes
Does the recipient have good head control?	☐	☐	
Does the recipient sit w/o support?	☐	☐	
Does the recipient transfer independently?	☐	☐	
Does the recipient stand independently?	☐	☐	
Does the recipient walk independently?	☐	☐	
Does the recipient have muscle tone issues?	☐	☐	
Does the recipient have orthopedic concerns?	☐	☐	
Does the recipient have sensory needs?	☐	☐	
Has the recipient used an adaptive Tryke? Yes _____ No _____ If yes, which model and how much assistance was needed? _____			
Any additional info? _____			

By signing below, you are signifying that in your professional opinion this rider would benefit from an Amtryke. You assume no liability.

Therapist Signature: _____ Date: _____

(205) 823-7800
ironcitytrykess.org

EMAIL COMPLETED APPLICATION TO: birminghamambucs3@gmail.com

Recipient's Name:

GENERIC ACCESSORIES(not model specific)

Leg & Foot Items	☐ Foot Cups (pair): ☐ Small ☐ Medium	☐ Knee Adductor Strap	☐ Pedal Block (1 = ¾") _____ qty
Hand Items	☐ Variable Range of Motion Kit (only for Hand & Foot Cycles)	☐ Wrist Wraps (Includes right & left) ☐ X-Small ☐ Medium ☐ Small ☐ Large	☐ Wrist Brace Mitt: ☐ Right ☐ Left ☐ Small ☐ Medium

FOOT TRYKES

☐ **1410 FOOT TRYKE** – (10" wheels, fixed drive) 50-FC-0 100 Arm Length 12-20", Leg Length 15-21"

Standard with Rear Steering

☐ **Option #1 – Standard Seating System: Blue Bucket Seat**

☐ **Option #2 – Alternate Seating System: Snappy Seat System**

Snappy Seat Accessories: ☐ Laterals ☐ Head Rest

1410 Accessories: ☐ Separator Cube ☐ H-Harness - 11.5"

☐ **1412 ProSeries FOOT TRYKE** (12" wheels, fixed drive) – 50-FC-1412 Arm Length 16-24", Leg Length 23-28"

Standard with Rear Steering

☐ **Option #1 – Standard Seating System: Medium Pommel Saddle Seat/1600 Simple Seatback**

Seat Bottom Alternates: ☐ Bench Seat ☐ Large Pommel Saddle Seat ☐ SaddleSeat

☐ **Option #2 – Alternate Seating System: 1400 ProSeries Seat Back System** (pushgrip, backpad, med. pommel & 2 laterals)

SeatBottomAlternates: ☐ Bench Seat ☐ Large Pommel Saddle Seat ☐ SaddleSeat

Trunk Support Options: ☐ ProSeries Full Padded Back ☐ ProSeries Head Rest ☐ ProSeries Lumbar Pad
☐ Recumbent 10 Degree Seat Post

☐ **Option #3 – Alternate Seating System: Gray Bucket Seat**

1412 Accessories: H-Harness - 11.5" ☐ Knee Separator: ☐ 2" ☐ 5" ☐ 9" Extender Tube
☐ ½" Exercise Pedals ☐ Steering) Hi-Rise Handlebars: ☐ 5" (BMX) ☐ 8"
☐ ½" Expanding Pedals

☐ **1416 ProSeries FOOT TRYKE** (16" wheels, fixed or freewheel drive) 50-F C-141 6 Arm Length 16-24", Leg Length 26-32"

☐ **Option #1 – Standard Seating System: 1400 ProSeries Seat Back System** (push grip, backpad, med. pommel & 2 laterals)

SeatBottomAlternates: ☐ BenchSeat ☐ TractorSeat with bracket ☐ LargePommel Saddle Seat ☐ Saddle Seat

Trunk Support Options: ☐ ProSeries Full Padded Back ☐ ProSeries Head Rest ☐ ProSeries Lumbar Pad
☐ Recumbent 10Degree Seat Post

☐ **Option #2 – Alternate Seating System: Medium Pommel Saddle Seat/1600 Simple Seat Back**

Seat Bottom Alternates: ☐ Bench Seat ☐ TractorSeatwithbracket ☐ LargePommel Seat ☐ Saddle Seat

☐ **Option #3 – Alternate Seating System: Gray Bucket Seat**

1416 Accessories: ☐ RearSteeringKit ☐ Calf & Leg Supports ☐ 9"extender tube
☐ H-Harness - 11.5" ☐ Hi-Rise Handlebars: ☐ 8" ☐ 12"
☐ ½" Exercise Pedals ☐ Knee Separator: ☐ 2" ☐ 5"

(205) 823-7800

ironcitytrykes.org

EMAIL COMPLETED APPLICATION TO: birminghamambucs3@gmail.com

Recipient's Name: _____

Foot Trykes (continued)

1420 ProSeries FOOT TRYKE (20" wheels, fixed or freewheel drive) –5 0-F C-142 0 Arm Length 20-30", Leg Length 29-35"
☐ Standard w/large pommel seat

☐ **Option #1 – Standard Seating System: 1400 ProSeries Seat Back System** (push grip, backpad, large pommel & 2 laterals)
SeatBottomAlternates: ☐ BenchSeat ☐ TractorSeatw/bracket ☐ Medium Pommel Seat ☐ Saddle Seat

Trunk Support Options: ☐ ProSeries Head Rest ☐ ProSeries Full Padded Back ☐ ProSeries Lumbar Pad
☐ Recumbent 10 Degree Seat Post

☐ **Option #2 – Alternate Seating System: Large Pommel Saddle Seat/1600 Simple Seatback**
Seat Bottom Alternates: ☐ Bench Seat ☐ TractorSeatw/bracket ☐ Medium Pommel Seat ☐ Saddle Seat

☐ **Option #3 – Alternate Seating System: Gray Bucket Seat**

1420 Accessories: ☐ RearSteeringKit ☐ Calf & LegSupports
☐ H-Harness - 11.5" ☐ Hi-Rise Handlebars: ☐ 8" ☐ 12"
☐ ½" Exercise Pedals ☐ Knee Separator: ☐ 2" ☐ 5" ☐ 9"extender tube

1420XL ProSeries FOOT TRYKE (20" wheels, fixed/freewheel drive) – 50 -FC -14 20-X L Arm Lgth 22-34", Leg Lgth 30-45"
☐ Standard w/large pommel seat

☐ **Option #1 – Standard Seating System: 1400 ProSeries Seat Back System** (push grip, backpad, large pommel & 2 laterals)
SeatBottomAlternates: ☐ BenchSeat ☐ TractorSeatw/bracket ☐ Medium Pommel Seat ☐ Saddle Seat
Trunk Support Options: ☐ ProSeries Head Rest ☐ ProSeries Full Padded Back ☐ ProSeries Lumbar Pad
☐ Recumbent 10 Degree Seat Post

☐ **Option #2 – Alternate Seating System: Large Pommel Saddle Seat/1600 Simple Seat Back**
Seat Bottom Alternates: ☐ BenchSeat ☐ TractorSeatw/bracket ☐ Medium Pommel Seat ☐ Saddle Seat

☐ **Option #3 – Alternate Seating System: Large Wheelchair Seat**
Wheelchair Seat Accessories: ☐ SwingAwayArms ☐ Wheelchair SeatBar Ends

1420 XL Accessories: ☐ RearSteeringKit ☐ Calf & Leg Supports
☐ H-Harness - 11.5" ☐ Hi-Rise Handlebars: ☐ 8" ☐ 12"
☐ ½" Exercise Pedals ☐ Knee Separator: ☐ 2" ☐ 5" ☐ 9" extender tube

☐ **High Functioning Set Up** (Both Options Allow Coasting) Can be used with 1416, 1420, 1420XL
(Rider must have the ability to steer, apply brake and make a complete pedal rotation on their own.)
Comes Standard with Pommel Seat (1416 & 1420 – L), 1600 Simple Seat Back, Hi-Rise Handlebars, Exercise pedals and safety flag.
MUST choose Option 1 or 2.

Choose Tryke: ☐ 1416 ☐ 1420 ☐ 1420XL

☐ **Option #1:**(switch to existing freewheel hub and add Disk Brake Kit - necessary for safety) Adds the ability to coast
☐ Dual HandBrake(Option1only)

☐ **Option #2:** (3-Speed Kit and coaster brake) Adds the ability to coast and switch gears for hilly terrain —Replaces the 1600 Series

Standard SeatingSystem (for either option): Pommel Seat/1600 Simple Seat Back

Seat Bottom Alternates: ☐ Bench Seat ☐ TractorSeat with bracket ☐ 1400Seat Back Set (w/ laterals)

Adaptive Accessories: ☐ ½" Expanding Pedals ☐ Telescoping Loop Handlebars

Recipient's Name:	
-------------------	--

HAND TRYKES & HAND-FOOT TRYKES

☐ AM-10 HAND TRYKE (10" wheels & 2.5" crank arms) 50-HFC-0105 Arm length 13-17" standard with Rear Steering (converts to push bar by inserting pin) and safety flag.			
☐ AM-10 HAND-FOOT TRYKE (10" wheels & 2.5" crank arms) 50-HFC-0105 Arm length 13-17", leg length 15-21" standard with Rear Steering.			
☐ Option #1 – Standard Seating System: Blue Bucket Seat			
☐ Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: ☐ Head Rest ☐ Laterals			
AM-10 Accessories: ☐ H-Harness - 11.5" ☐ Vertical Hand Grips ☐ Separator Cube			

☐ AM-12Small HAND TRYKE (12" wheels & 3" crank arms) 50-HFC-0110 Arm length 15-20" standard with Saddle Seat, Rear Steering (converts to push bar by inserting pin) and safety flag.			
☐ AM-12Small HAND-FOOT TRYKE (12" wheels & 3" crank arms) 50-HFC-0110 Arm length 15-20", leg length 19-24" standard with Saddle Seat, Rear Steering (converts to push bar by inserting pin) .			
☐ Option #1 – Standard Seating System: Saddle Seat & 1600 Simple Seat Back			
☐ Option #2 – Upgraded Seating System: Snappy Seat System Snappy Seat Accessories: ☐ Head Rest ☐ Laterals			
☐ Option #3 – Upgraded Seating System: Blue Bucket Seat			
AM-12S Accessories: ☐ H-Harness - 11.5" ☐ Separator Cube ☐ Vertical Hand Grips ☐ 9/16" Exercise pedals (hand + foot only)			

☐ AM-12 HAND TRYKE (12" wheels & 4" crank arms) 50-HFC-0210 Arm length 14-23": Standard w/rear steering			
☐ AM-12 HAND-FOOT TRYKE (12" wheels & 4" crank arms) 50-HFC-0210 Arm length 14-23", leg length 21-29" standard with Rear Steering.			
☐ Option #1 – Standard Seating System: Saddle Seat/1600 Simple Seat Back Seat Bottom Alternates: ☐ Bench Seat ☐ Pommel Seat ☐ Medium (11.5x10") ☐ Large (14x13")			
☐ Option #2 – Upgraded Seating System: ProSeries Seat Back System (push grip, backpad, saddle seat & 2 laterals) Seat Bottom Alternates: ☐ Bench Seat ☐ Pommel Seat ☐ Medium (11.5x10") ☐ Large (14x13") Trunk Support Upgrades: ☐ ProSeries Full Padded Back ☐ Head Rest ☐ ProSeries Lumbar Pad			
☐ Option #3 – Upgraded Seating System: Gray Bucket Seat			
AM-12 Accessories: ☐ H-Harness-11.5" ☐ 9" Extender Tube ☐ 9/16" Exercise ☐ Knee Separator: ☐ 2" ☐ 5" ☐ Vertical Hand Grips pedals (hand + foot only)			

AmTryke Sizing Chart

TRYKE TYPE (if?) (How will the tryke be used?)	RIDER LEG LENGTH (to bottom of shoe) of hip (greater trochanter)	RIDER ARM LENGTH (to center of digit crease) of shoulder to center of middle of hip	MODEL	RIDER WEIGHT (Pounds)	RIDER MAX HEIGHT (Inches)	TRYKE WEIGHT (Pounds)	WHEEL SIZE (Inches)	TTRYKE HEIGHT (inches)	TRYKE LENGTH (inches)	TRYKE WIDTH (Inches)
Hand & Foot	15-21	13-17	AM-10	55	40	45	10	24	38	21
	19-24	15-20	AM-12S	150	40	45	12	27	38	24
	21-29	14-23	AM-12	150	47	45	12	36	60	32
	24-36	18-27	AM-16	175	66	55	16	36	68	33
Foot	15-21	12-20	1410	55	40	45	10	24	38	21
	23-28	16-24	1412	125	42	72	12	13	43	27
	26-32	16-24	1416	175	60	74	16	49	58	30
	29-35	20-30	1420	250	68	74	20	50	64	30
	30-45	22-34	1420XL	275	76	89	20	43	72	29
	30-41	20-28	JT-2000/ JT2300USS	250	74	80	20	48	72	32
	36-45	17-29	TP-3000	300	75	47	20	33	63-70	33
Hand	up to 37	19-30	1020	250	67	74	20	41	69	30
	up to 41	22-26	1024	250	72	85	24	45	75	32
	All trykes in the Hand & Foot section can be converted to Hand trykes.									

A Center of Shoulder
B Center of Elbow
C Center of Digit Crease

D Center of Hip (greater trochanter)
E Center of Knee
F Bottom of Foot



RIDER'S MEASUREMENTS

Arm Measurements (inches) Total Length

Left A to B + B to C = _____

Right A to B + B to C = _____

Trunk A to D = _____

Leg Measurements (inches) Total Length

Left D to E + E to F = _____

Right D to E + E to F = _____

HELMET SIZING

Sizes Head Circumference Inches

Toddler (XS) 17.7" – 19.3"

Child (S) 20.5" – 21.7"

Youth (L) 20.9" – 22.4"

Adult (XL) 22.4" – 23.6"

EMAIL COMPLETED APPLICATION TO: birminghamambucs3@gmail.com