


Cannulas



CASTROVIEJO


Cyclodialysis Cannula
Gently curved Spatulated
Angled Tip, Top opening

- WSC-640** 10 mm (20 gauge, 0.5 mm opening)
- WSC-640 A** 15 mm (20 gauge, 0.5 mm opening)
- WSC-640 B** 10 mm (22 gauge, 0.3 mm opening)
- WSC-640 C** 15 mm (22 gauge, 0.3 mm opening)




RANDOLPH

Cyclodialysis Cannula
13 mm Spatulated Gently Curved
Angled Tip, 19 gauge
WSC-641




WESTERN

Lasik Cannula
13 mm Angled Spatulated
Gently Curved Tip, 26 Gauge
WSC-642





BANAJI

Lasik Cannula
25 Gauge Gently Curved
Angled Tip with Six Irrigating Ports
WSC-643




BISHOP-HARMON

Anterior Chamber Wash Cannula
10 mm Angled Spatulated Tip
Sizes: 18G, 19G, 20G
WSC-644




WESTERN

Viscoelastic Cannula
10 mm Angled 22 Gauge, 25 Gauge
WSC-645



WESTERN

Irrigating Silicone Bulb
with Luer Lock Adaptor
WSC-646
WSC-646 A Spare Luer Lock Adaptor





Cannulas





DALJIT

Blast Cannula Used for Washing The Cortex
WSC-647




WESTERN

Lachrymal Cannula in
23 Gauge, 24 Gauge, 25 Gauge Long Hub
WSC-648 Straight
WSC-648 A Half Curved
WSC-648 B Full Curved




RYCROFT

Air Injection Cannula
7 mm Angled Blunt Tip
WSC-649 25 Gauge, 26 Gauge,
27 Gauge, 30 Gauge




RAININ

Air Injection Cannula 7mm
Angled Blunt Spatulated Ttip
WSC-650 26 Gauge, 27 Gauge, 30 Gauge



TROUTMAN


Alpha-Chymotrypsin Cannula
Delicate Olive Tip 25 Gauge with
1 mm Smooth Bulbous Tip
WSC-651




Cannulas





ATKINSON
Retrolubar Needle Round Bevelled
Tip 35 mm Long
WSC-652 23 Gauge, 25 Gauge, 27 Gauge




WESTERN
Needle 30 Gauge for Aspiring
Intra-Vitreous Fluid
WSC-653 12 mm
25 mm



PIERCE
Co-axial Irrigating and Aspiring
Cannula, 19 Gauge Outer Cannula,
23 Gauge Inner Cannula with 0.3 mm
Aspiring Port
WSC-654 Spare Inner Needle





GILLS-WELSH
Olive Tip Cannula 10 mm Angled
Shaft with 1 mm Smooth Bulbous Tip
WSC-655




Cannulas





SIMCOE
Irrigating And Aspiring Cannula
23 Gauge Curved Shaft with 0.3 mm
Top Aspiring Port with Silicone Tubing
WSC-656 Irrigation Through Main Hub and
Aspiration Through Side Tubing




WESTERN
Irrigating and Aspiring Cannula Same As
Nut with Front Aspiring Port Beveled Anterior
For Proper Visualization During Aspiration
WSC-656 A Irrigation Through Main Hub and
Aspiration Through Side Tubing




SIMCOE (REVERSE)
Irrigating and Aspiring Cannula 23 Gauge
Curved Shaft with 0.3mm Top Aspiring Port
with Silicone Tubing
WSC-657 Irrigation Through Side Tubing and
Aspiration Through Main Hub



WESTERN (REVERSE)
Irrigating and Aspiring Cannula Same as
But with Front Aspiring Port
Beveled Anterior for Proper Visualization
During Aspiration
WSC-657 A Irrigation Through Side Tubing and
Aspiration Through Main Hub



WESTERN
'J' shaped 12 O' Clock Irrigating and Aspiring
Cannula 23 Gauge Side Curved Shaft with Aspiring
Port Facing 12 O' Clock for Easy Aspiration with
Silicone Tubing
WSC-658 Irrigation Through Main Hub
Aspiration Through Side Tubing
WSC-658 A Irrigation Through Side Tubing and
Aspiration Through Main Hub (Reverse)





Cannulas




GILL-WELSH

Olive Tip Irrigating and Aspiring Cannula
23 Gauge Irrigation and 25 Gauge Aspiration
Angled Shaft with Silicone Tubing
WSC-659 Right
WSC-659 A Left





DALJIT

Irrigating and Aspiring Cannula Curved
Co-axial Shaft with Outer Tubing 18 Gauge
And Inner Tubing 24 Gauge Aspiring Port
Bent Upwards for Proper Visualization During
Cortex Aspiration without Rupturing the
Posterior Capsule, with Silicone Tubing
WSC-660





O'GAWA

Irrigating and Aspiring Cannula
Straight 20 Gauge Irrigating and
Aspiring Tube with Silicone Tubing
WSC-661



GESS

Irrigating and Aspiring Cannula
23 Gauge U-Shape Tip with 0.3 mm
Aspiring Port for Aspiring
Between 11 and 1 O' Clock Through
Silicone Tubing
WSC-662 Right
WSC-662 A Left





Cannulas




SILICONE TUBING

Spare Silicone Tubing for Irrigating
And Aspiring Cannula
WSC-663




WESTERN

Stromal Hydration Cannula
27 Gauge 1mm Angled Tip
WSC-664





GILLS

Aspiration Cannula
25 Gauge Angled Blunt Tip
WSC-665 5 mm, 7 mm, 8 mm, 10 mm, 12 mm Tip





GILLS-WELSH

Flat Olive Tip Aspiration Cannula
25 Gauge Angled Shaft
WSC-666 Straight
WSC-666 A Right
WSC-666 B Left




BINKHORST

Irrigating Cannula
25 Gauge Angled Shaft
WSC-667 Straight
WSC-667 A Right
WSC-667 B Left





KELMAN

Irrigating Cystotome Cannula
22 gauge
WSC-668 Straight
WSC-668 A Angled



BLUMENTHAL

Irrigating Cystotome Cannula
27 Gauge, Tip Angled 90°, Then
Rotated 90°, Angled 9 mm
Gently Curved Shaft
WSC-669




Cannulas




WESTERN

Irrigating Cystotome Cannula
25 Gauge, Tip Angled 90°,
Gently Curved 19 mm Shaft
WSC-670




KRAFT

Posterior Capsule Polishing Curette
Irrigating, 1.75 mm Diameter Cup
25 Gauge, Gently Curved Shaft
WSC-671




KNOLLE

Posterior capsule polishing irrigating
2mm angled sandblast tip 23 gauge
WSC-672





JENSON

Posterior and anterior capsule
polisher, irrigating with carbide
impregnated olive tip 1 by 2mm
long 23 gauge gently curved shaft
WSC-673




WESTERN

Capsulorhexis and Posterior Capsule
Polisher, Irrigating Vacuum Action
Left Hand Aspirating Angled Shaft
23 Gauge with 0.3mm Reverse Aperture
WSC-674




WESTERN

Silicone Oil Aspirating Cannula During
Vitreotomy
Sizes: 18 Gauge by 35 mm Oval Aspirating Ports
WSC-675




WESTERN

Hydro-Dissection Cannula Angled
11mm Shaft
Sizes: 25, 26, 27 Gauge
WSC-676



WESTERN

Hydro-Dissection Cannula J-Shaped
10 mm Angled Shaft for 12 O' Clock
Hydro-Dissection 25 Gauge
WSC-677




Cannulas




BLUMENTHAL

Hydro-Dissection and Hydro-Delineator
Cannula Tapered Tip 4 mm Tip Angled Shaft
WSC-678




WESTERN

Hydro-Dissection and Hydro-Delineator
Cannula Beveled Sharp Blocked Tip.
10 mm Angled Shaft with Two Irrigating Ports
WSC-679




CORYDON

Hydro-Dissection and Hydro-Expression
Cannula Angled 30°, 1.8 mm Tip Gently
Curved Shaft
WSC-680 25 Gauge for Visco-Expression
WSC-680 A 27 Gauge for Hydro-Expression




KNOLLE

Irrigating Vectis 23 Gauge 5 mm Wide by
9 mm Long Loop with Three Irrigating Ports
WSC-681 3.2 mm Wide by 5 mm Long Loop




DREWS

Irrigating Vectis 23 Gauge 5 mm wide by
9 mm Long Loop with Three Irrigating Ports
One at The Front and Two at The Back
WSC-682




ANIS

Irrigating Vectis 25 Gauge 5 mm wide by
7mm Long Loop with a Dimple at The Top
For Retraction of Upper Capsular Flap
With Two Irrigation Ports
WSC-683



WESTERN

Irrigating Vectis for Small Incision
non Phaco Serrated Tip one Irrigating Port
WSC-684




Cannulas




SIMCOE

Cortex Extractor Cannula
23 Gauge Gently Curved Shaft
0.3 mm Aspiration Port
WSC-685 13mm
WSC-685 A 21mm





BLUMENTHAL

Cortex Extractor Cannula
21 Gauge, 2 mm Sand Blasted Tip
Gently Curved Shaft 0.3 mm
Aspiration Port
WSC-686





ANIS

Cortex Extractor Cannula
23 Gauge, 2 mm Carbide Impregnated
Angled Tip Straight Shaft 0.3 mm
Aspiration Port
WSC-687





KRATZ

Lens Manipulating Cannula
0.2mm tip
WSC-688




McINTYRE

Irrigating Iris Hook 23 Gauge
Angled Shaft
WSC-689




YOUENS

Irrigating Lens Manipulating Cannula
30 Gauge Angled Shaft 1.5 mm Bent Tip
WSC-690



OSHER

Irrigating Lens Manipulating Cannula
24 Gauge C-Shaped Hook at The Tip
Straight Shaft
WSC-691




Cannulas





GREATHER

Irrigating Iris Retractor and Posterior
Capsule Polisher 23 Gauge
WSC-692 Straight
WSC-692 A Angled




SHAPIRO

Irrigating Posterior Capsule Polisher
2.5 mm Long Carbide Impregnated Tip
0.3 side Port 23 Gauge
WSC-693





LEWICKY

Chamber Maintainer 20 Gauge Threaded
Tip for Fixation Very Useful for Manual
Small Incision and Phaco
WSC-694





BLUMENTHAL

Chamber Maintainer 20 Gauge Smooth
Oval Tip for Fixation very useful for
Manual Small Incision and Phaco
WSC-695



WESTERN

Infusion Cannula 20 Gauge Beveled Tip with
Luer Lock adapter for Vitreoretinal Surgery
WSC-696 2.5mm
WSC-696 A 4mm
WSC-696 B 6mm



WESTERN

Suction Cannula
WSC-697

