

Category	Test Item	Specimen Type	Specimen Volume	Reaction Time	Measuring Range	Clinical Reference	Qualification
COVID-19	SARS-CoV-2 Ab ^{New}	WB/Serum/Plasma	10 µL	10 min	Qualitative	Negative	CE
	SARS-CoV-2 IgM ^{New}	WB/Serum/Plasma	10 µL	10 min	Qualitative	Negative	CE
Vitamin	Vitamin D	Serum/Plasma	75 µL	15 min	5-100 ng/mL	30-100 ng/mL	CE
Diabetic & Renal Injury Markers	HbA1c	WB	10 µL	5 min	4.0-14.5%	4-6.5%	CE IFCC NGSP
	MAU	Urine	75 µL	3 min	5-300 mg/L	0-20 mg/L	CE
	CysC	WB/Serum/Plasma	10 µL	5 min	0.2-10 mg/L	0.5-1.1 mg/L	CE
	NGAL	Urine	75 µL	15 min	10-1500 ng/mL	≤131.7 ng/mL	CE
	β2-MG	WB/Serum/Plasma	10 µL	10 min	0.3-20 mg/L	1.0-2.7 mg/L	CE
Thyroid Function	T3	WB/Serum/Plasma	75 µL	15 min	0.61-9.22 nmol/L	1.23-3.07 nmol/L	CE
	T4	WB/Serum/Plasma	75 µL	15 min	12.87-300 nmol/L	66-181 nmol/L	CE
	TSH	WB/Serum/Plasma	75 µL	15 min	0.1-100 mIU/L	0.3-4.2 mIU/L	CE
	fT3	WB/Serum/Plasma	75 µL	15 min	0.40-50.00 pmol/L or 0.26-32.55 pg/mL	2.8-7.1 pmol/L (1.82-4.61 pg/mL)	CE
	fT4	WB/Serum/Plasma	75 µL	15 min	1.00-100 pmol/L or 0.078-7.77 ng/dL	12-22 pmol/L (0.94-1.72 ng/dL)	CE
Fertility	β-HCG	WB/Serum/Plasma	20 µL	15 min	2-200,000 mIU/mL	5 mIU/mL	CE
	LH	WB/Serum/Plasma	75 µL	15 min	1-100 mIU/mL		CE
	FSH	WB/Serum/Plasma	75 µL	15 min	1-100 mIU/mL		CE
	Prog	WB/Serum/Plasma	75 µL	15 min	1.5-60 ng/mL		CE
	PRL	WB/Serum/Plasma	75 µL	15 min	1-200 ng/mL		
	E2	WB/Serum/Plasma	75 µL	15 min	9-3000 pg/mL		
	Testosterone	WB/Serum/Plasma	75 µL	15 min	0.2-15 ng/mL		
	AMH	WB/Serum/Plasma	75 µL	15 min	0.1-16 ng/mL		
Cardiac Markers	NT-proBNP	WB/Serum/Plasma	75 µL	15 min	18-35000 pg/mL	<75 years old 0-300 pg/mL ≥75 years old 0-450 pg/mL	CE
	cTn I	WB/Serum/Plasma	75 µL	15 min	0.1-50 ng/mL	0-0.3 ng/mL	CE
	Myo	WB/Serum/Plasma	75 µL	15 min	2.0-400 ng/mL	0-58 ng/mL	CE
	CK-MB	WB/Serum/Plasma	75 µL	15 min	0.3-100 ng/mL	0-5 ng/mL	CE
	H-FABP	WB/Serum/Plasma	75 µL	15 min	1-120 ng/mL	0-7 ng/mL	CE
	3 in 1 (cTn I/Myo/CK-MB)	WB/Serum/Plasma	75 µL	15 min	Same with single items	Same with single items	CE
	BNP	WB/ Plasma	75 µL	15 min	5-5000 pg/mL	0-100 pg/mL	CE
	2 in 1 (cTn I/ NT-proBNP)	WB/Serum/Plasma	75 µL	15 min	Same with single items	Same with single items	CE
	hsCRP	WB/Serum/Plasma	10 µL	3 min	0.1-10 mg/L	≤1.0 mg/L	CE
	TnT ^{New}	WB/Serum/Plasma	100 µL for whole blood 75 µL for plasma/serum	15 min	0.03-10 ng/mL	0-0.1 ng/mL	
Coagulation	D-Dimer	WB/Plasma	15 µL for whole blood 10 µL for plasma	5 min	0.1-10 mg/L	0-0.5 mg/L	CE
	D-Dimer one step	WB/Plasma	100 µL for whole blood 75 µL for plasma	15 min	0.1-10 mg/L	0-0.5 mg/L	CE
Inflammation Markers	PCT	WB/Serum/Plasma	75 µL for whole blood 50 µL for plasma/serum	15 min	0.1-100 ng/mL	0-0.5 ng/mL	CE
	CRP (hsCRP+CRP)	WB/Serum/Plasma	8.5 µL for whole blood 5 µL for plasma/serum	3 min	0.5-200 mg/L	CRP:0-10 mg/L hsCRP:0-1 mg/L	CE
	2 in 1 (CRP+PCT)	WB/Serum/Plasma	75 µL	15 min	CRP:0.5-150 mg/L PCT:0.1-100 ng/mL	CRP:<10 mg/L PCT:<0.5 ng/mL	CE
	SAA	WB/Serum/Plasma	10 µL	15 min	1.0-300 mg/L	0-10 mg/L	CE
	IL-6 ^{New}	WB/Serum/Plasma	75 µL	15 min	3 ~ 4000 pg/mL	≤10 pg/mL	CE
Tumor Markers	AFP	WB/Serum/Plasma	75 µL	15 min	5-400 ng/mL	0-20 ng/mL	CE
	PSA	WB/Serum/Plasma	75 µL	15 min	2-100 ng/mL	0-4 ng/mL	CE 0123
	CEA	WB/Serum/Plasma	75 µL	15 min	1-500 ng/mL	0-5 ng/mL	CE
	fPSA	WB/Serum/Plasma	75 µL	15 min	0.2-30 ng/mL	< 1 ng/mL	CE 0123
Craniocerebral Injury Markers	S100β ^{New}	WB/Serum/Plasma	75 µL	15 min	0.05-10 ng/mL	< 0.2 ng/mL	
Others	Cortisol	WB/Serum/Plasma	75 µL	15 min	50-1000 nmol/L	7:00-10:00 a.m.: 201.31-536.54 nmol/L 4:00-8:00 p.m.: 65.78-330.85 nmol/L	

Ongoing: Ferritin, FOB, H.Pylori antigen, Vit-B12, ST2



Wave Biotech Inc

No. 54, 7th Cross, HMT Layout, R.T. Nagar, Bengaluru - 560 032.
Phone : 9980996930, 9844141515, 8553135135
E-mail : wavebiotechinc@gmail.com, www.wavebiotechinc.in



Wave Biotech Inc



Fincare™

FIA Meter III Plus (FS-205)

Fluorescence Immunoassay Analyzer

CE NGSP ifcc

Finecare™ FIA METER III PLUS

Multi-Channel

20 channels, detect more items one time

Reliable Accuracy

Good correlation with CLIA method

Temperature control

Provide better reaction condition

Rapid Test

High speed of detection, enable 200 tests per hour

FS-205

- Dimensions: 374*396*480.5 mm
- Weight: Approx. 19.5Kg

ID chip inlet

- Automatically reading of batch number and calibration information in the ID chip

Test Channel

- 20 channels
- Random access for samples
- Automatic internal incubation

Built-in Thermal Printer

- Real-time auto-printing
- Support for external printers

Designed Sampler

- Specially-designed sampler for small volume of specimen

Large Display

- 10" LCD Touch Screen
- Intuitive user interface

Connection

- Direct connection to LIS/HIS


Cartridge Tray

- Cartridges will be automatically ejected into the tray when test is finished


Test Reagent

- Room temperature storage
- Suitable for whole blood, serum, plasma and urine
- Get results in 3-15 minutes


TEST PROCEDURE




Sampling



Mixing



Dropping



Reading



Laboratory



ER



ICU



Clinics



GP Office

APPLICATION