

CARDIAC AND CANCER SURGERY SOLUTIONS

Turnkey OT Setup - Devices & equipments Instruments - Consumables & disposables

PRODUCT CATALOG CARDIO VASCULAR INSTRUMENTS



Heart Lung Pack • Perfusion Cannulas • Oxygenator Heart Valves • IABP Machine • Heart Lung Machine Aortic Punch • Surgical Instruments • Chest Drain













WYLIE Abdominal Aorta Clamps

320-291 9 1/2" (24cm), angled shanks, curved jaws, jaw length 2 3/4" (7cm) **320-290** 10 1/2" (26.7cm), angled shanks, curved jaws, jaw length 3 1/8" (8cm)



WYLIE Hypogastric Clamp

320-289 9 1/2" (24.1cm), angled shanks, curved jaws, jaw length 2 3/4" (7cm)



CHERRY Supraceliac Clamp, 60° Handle, Double Curve

320-287 12 1/2" (31.7cm), angled shanks, curved jaws, jaw length 5 3/4" (14.5cm)



LEMOLE-STRONG Aorta Clamps

320-257 8 1/4" (21cm), straight shanks, jaw length 2" (5cm), depth 7/8" (2.2cm) 320-258 8 1/4" (21cm), angled shanks, jaw length 2" (5cm), depth 7/8" (2.2cm)





BAHNSON Aortic Clamp

320-263 9 1/8" (23.2cm), angled shanks, curved jaws, jaw length 2 3/4" (7cm)



LAMBERT Abdominal Aortic Aneurysm Clamp

320-266 9 1/4" (23.5cm), angled shanks, curved jaws, jaw length 2" (5cm)



GRANT Abdominal Aortic Aneurysm Clamps ("Texas Twister") 2x3 DeBAKEY Jaws Serrations

320-269 10" (25.5cm), jaw length 1 3/4" (5.5cm), 2" (5cm) beyond curve **320-264** 10 1/2" (26.5cm), jaw length 2 5/8" (7cm), 2 1/2" (6.5cm) beyond curve





DeBAKEY Aorta Clamp 2x3 DeBAKEY Jaw Serrations



JARIT Aortic Aneurysm Clamp, DeBAKEY Jaw Serrations

320-357 10 3/4" (27.3 cm), jaw length 2" (5cm), 2 1/4" (5.7cm) beyond curve



MORRIS Aorta Clamps 2x3 DeBAKEY Jaw Serrations

320-267 8 1/4" (21cm), curved shanks, angled jaws, jaw length 2 3/4" (6.7cm) wide jaws 8 1/2" (21.5cm), curved shanks, angled jaws, jaw length 3 3/4" (9.5cm) delicate jaws 11" (27.9cm), curved shanks, angled jaws, jaw length 3" (7.6cm), std jaws 320-259 320-268



COOLEY Iliac Clamps

8 5/8" (22cm), jaw length 2 1/2" (6.5cm) 320-494 320-495





COOLEY Aorta Clamps

320-243 10 7/8" (27.6cm), jaw length 2 1/4" (5.5cm), depth 5/8" (1.5cm) **320-244** 10 1/4" (26cm), jaw length 3" (7.5cm), depth 3/4" (1.7cm)



OCHSNER Aorta Clamp

320-242 9 3/8" (24cm), angled shanks, jaw length 2 1/8" (5.4cm), depth 3/4" (1.9cm)



DIETHRICH Aorta Clamp

320-246 8 1/4" (21 cm), angled shanks, jaw length 2" (5 cm), depth 5/8" (1.6 cm)



KAY Aorta Clamp

320-247 7 3/4" (20cm), angled shanks, jaw length 2 1/4" (5.5cm), depth 1 1/8" (3cm)





LAMBERT-KAY Aorta Clamp, Titanium

320-728TI 8" (20.3cm), jaw length | 1/4" (3cm), depth 3/4" (2cm)



LAMBERT-KAY Aorta Clamp

320-248 7 3/4" (20cm), angled shanks, jaw length 2" (4.9cm), depth 3/4" (2cm)



WYLIE "J" Clamp

320-249 8 1/2" (21.5cm), curved shanks, jaw length | 3/4" (4.5cm), depth | 3/8" (3.5cm)





DeBAKEY Tangential Occlusion Clamps

320-320 7 7/8" (20cm), jaw length | 5/8" (4cm), depth | 1/2" (1.2cm)
320-322 9 7/8" (25cm), jaw length 2 1/8" (5.5cm), depth 3/4" (1.8cm)
320-324 | 10 1/4" (26cm), jaw length 2 1/2" (6.5cm), depth 3/4" (1.8cm)
320-326 | 10 3/8" (26.5cm), jaw length 2 5/8" (6.7cm), depth 3/4" (1.8cm)
320-328 | 10 5/8" (27cm), jaw length 3" (8cm), depth 3/4" (2cm)



COOLEY Tangential Occlusion Clamps

320-447 7 3/4" (19.7cm), jaw length 3/4" (1.9cm), depth 1/4" (6mm) 7 3/4" (19.7cm), jaw length 1" (2.5cm), depth 1/4" (7mm) 7 3/4" (19.7cm), jaw length 1 3/8" (3.5cm), depth 3/8" (10mm)





COOLEY Partial Occlusion Clamps

 320-335
 8" (20.3cm), jaw length 1" (2.7cm), depth 1/4" (7mm)

 320-336
 8 1/8" (20.6cm), jaw length 1 1/2" (4cm), depth 3/8" (8mm)

 320-337
 9" (22.9cm), jaw length 2" (5.2cm), depth 1/2" (1.2cm)

 320-338
 9 1/4" (23.5cm), jaw length 2 3/8" (6cm), depth 1/2" (1.2cm)

 320-339
 11" (27.9cm), jaw length 2 7/8" (7.2cm), depth 1/2" (1.2cm)







BECK Aorta Clamps

 320-296
 8" (20.3cm), straight shanks, jaw length | 1/4" (3cm), depth | 1/4" (8mm)

 320-297
 8 3/8" (21.3cm), straight shanks, jaw length | 3/4" (4.5 cm), depth 3/8" (1cm)

 320-298
 8 1/2" (21.6cm), straight shanks, jaw length | 2 5/8" (6.7cm), depth 3/8" (1cm)


SATINSKY Vena Cava Clamps

320-390 9 3/4" (24.5cm), jaw length 2" (5cm), depth 3/8" (1cm), large **320-391** 9 3/4" (24cm), jaw length 1 5/8" (4cm), depth 3/8" (1cm), medium **320-392** 9 3/4" (23.5cm), jaw length 1 3/8" (3.5cm), depth 1/4" (8mm), small



SATINSKY Vena Cava Clamps

310-204 6 1/4" (15.9cm), cross-serrated, pediatric size **310-200** 9 1/4" (23.5cm), cross-serrated **310-202** 10 1/4" (26cm), cross-serrated, heavy





BECK Aorta Clamps

320-296 8" (20.3cm), straight shanks, jaw length | 1/4" (3cm), depth 1/4" (8mm) **320-297** 8 3/8" (21.3cm), straight shanks, jaw length | 3/4" (4.5 cm), depth 3/8" (1cm) **320-298** 8 1/2" (21.6cm), straight shanks, jaw length 2 5/8" (6.7cm), depth 3/8" (1cm)


SATINSKY Vena Cava Clamps

320-390 9 3/4" (24.5cm), jaw length 2" (5cm), depth 3/8" (1cm), large **320-391** 9 3/4" (24cm), jaw length | 5/8" (4cm), depth 3/8" (1cm), medium **320-392** 9 3/4" (23.5cm), jaw length | 3/8" (3.5cm), depth 1/4" (8mm), small



SATINSKY Vena Cava Clamps

310-204 6 1/4" (15.9cm), cross-serrated, pediatric size **310-200** 9 1/4" (23.5cm), cross-serrated **310-202** 10 1/4" (26cm), cross-serrated, heavy





GLOVER Auricular Appendage Clamp

320-375 10 7/8" (27.5cm), jaw length 1 7/8" (4.7cm), depth 1/4" (7mm)



REYNOLDS Clamps

320-483	50°	7" (18cm), angled shanks, jaw length 1 1/8" (3cm), depth 1/8" (5mm), small
320-484	50°	7" (18cm), angled shanks, jaw length 1 3/8" (3.5cm), depth 1/4" (7mm), large
320-485	85°	6 1/2" (16.5cm), angled shanks, jaw length 1 1/8" (3cm), depth 1/8" (5mm), small
320-486	85°	6 1/2" (16.5cm), angled shanks, jaw length 1 3/8" (3.5mm), depth 1/4" (7mm), large





GLOVER Auricular Appendage Clamp

320-375 10 7/8" (27.5cm), jaw length 1 7/8" (4.7cm), depth 1/4" (7mm)



REYNOLDS Clamps

320-483	50°	7" (18cm), angled shanks, jaw length 1 1/8" (3cm), depth 1/8" (5mm), small
320-484	50°	7" (18cm), angled shanks, jaw length 1 3/8" (3.5cm), depth 1/4" (7mm), large
320-485	85°	6 1/2" (16.5cm), angled shanks, jaw length 1 1/8" (3cm), depth 1/8" (5mm), small
320-486	85°	6 1/2" (16.5cm), angled shanks, jaw length 1 3/8" (3.5mm), depth 1/4" (7mm), large





HARKEN Auricle Clamps



GLOVER Clamps


320-370 320-372 320-360 320-361 320-362 320-363	8 1/4" (21cm), straight shanks, spoon-shaped jaws, jaw length 1 1/2" (4cm) 10 1/4" (26cm), straight shanks, spoon-shaped jaws, jaw length 2 1/8" (5.5cm) 8 1/4" (21cm), straight shanks, standard curve, jaw length 1 3/4" (4.5cm) 9 1/2" (24cm), straight shanks, extra long curve, jaw length 2 1/4" (5.7cm) 7 1/4" (18.5cm), straight shanks, acutely curved, jaw length 1 5/8" (4cm), 8 3/4" (22cm), straight shanks, extra long acute curve, jaw length 2 1/8" (5.3cm)	
		00



DeBAKEY-BAHNSON Clamps

320-275 320-276 320-277

10 7/8" (27.6cm), jaw length 3" (7.6cm), light curve 11 3/8" (28.9cm), jaw length 2 1/2" (6.4cm), medium curve 10 5/8" (27cm), jaw length 2 3/8" (6cm), acute curve



DeBAKEY Vascular Clamps

320-270 9 5/8" (24.5cm), curved shanks, curved jaws, jaw length 2 3/4" (7.2cm) 320-272 10 3/8" (26.5cm), curved shanks, fully curved jaws, jaw length 4 3/8" (11cm)



COOLEY Multipurpose Clamps

8" (20.3cm), jaw length | 1/4" (3.2cm), depth 5/8" (1.6cm) 8" (20.3cm), jaw length | 1/2" (3.8cm), depth 3/4" (1.9cm) 320-472 320-473





DeBAKEY Multi-Purpose Clamps, 60° Angle

 320-302
 8" (20.5cm), jaw length 2 1/8" (5.5cm)

 320-304
 9 3/8" (24cm), jaw length 2 3/4" (7cm)


 320-300
 9 7/8" (25cm), jaw length 1 5/8" (4cm)

 320-306
 12" (30.5cm), jaw length 4" (10cm)



LELAND JONES Peripheral Vascular Clamps

320-215 7 5/8" (19.5cm), straight jaws, jaw length 2 3/4" (7cm) **320-216** 7 5/8" (19.5cm), 15° angle jaws, jaw length 2 3/4" (7cm) **320-217** 7 1/2" (19cm), 30° angle jaws, jaw length 2 3/4" (7cm)





COOLEY Peripheral Vascular Clamps

320-209	6 1/4" (15.9cm), straight shanks, straight jaws, jaw length 1 3/4" (4.5cm)
320-210	6 1/2" (16.5cm), curved shanks, 55° angle jaws, jaw length 1 1/4" (3.2cm)
320-211	6 3/4" (17.1 cm), curved shanks, 55° angle jaws, jaw length 1 1/2" (3.8 cm)
320-212	7 1/4" (18.4cm), curved shanks, 55° angle jaws, jaw length 2 3/8" (6cm)
320-213	7" (17.8cm), curved shanks, 45° angle jaws, jaw length 3/4" (4.5cm)
320-214	6 3/8" (16.2cm), straight shanks, curved jaws, jaw length 1/4" (3cm)







DeBAKEY Peripheral Vascular Clamps

320-250 7 1/2" (19cm), curved shanks, angled jaws, jaw length 3 1/8" (8cm) 320-255 7 7/8" (20cm), curved shanks, curved jaws, jaw length 3 1/2" (9cm), S-shape



DeBAKEY Aortic Exclusion Clamp

320-260 7 1/2" (19cm), curved shanks, curved jaws, jaw length 3 3/4" (9.4cm), S-shape



DeBAKEY Peripheral Vascular Clamp

320-265 6 1/4" (16cm), straight shanks, curved jaws, jaw length 2 1/4" (5.6cm)







BAINBRIDGE Clamps, Longitudinal Serrations

 105-300
 6" (15.2cm), straight jaws, jaw length 1 3/4" (4.5cm)

 105-302
 6" (15.2cm), curved jaws, jaw length 1 3/4" (4.5cm)


 105-304
 7 1/4" (18.4cm), straight jaws, jaw length 2 1/2" (6.5cm)

 105-306
 7 1/4" (18.4cm), curved jaws, jaw length 2 1/2" (6.5cm)



DeBAKEY-BAINBRIDGE Clamps, DeBakey Serrations

320-228 7 7/8" (20cm), straight jaws, jaw length 2 3/8" (6cm), delicate **320-229** 7 7/8" (20cm), curved jaws, jaw length 2 3/8" (6cm), delicate





DeBAKEY Bulldog Clamps, Ring Handle

320-180	5 1/8" (13cm), straight shanks, straight jaws, jaw length 1 3/4" (4.4cm)
320-182	5" (12.7cm), straight shanks, curved jaws, jaw length 1 3/4" (4.4cm)
320-185	4 3/4" (12.1cm), curved shanks, 45° angle jaws, jaw length 1 1/2" (3.8cm)
320-187	4" (10.2cm), straight shanks, 90° angle jaws, jaw length 1 1/2" (3.8cm)
320.100	5" (12.7cm) curved shaple curved jour jour length 3/4" (4.4cm) S shape












DeBAKEY Coarctation Clamps

320-230 8 5/8" (22cm), straight shanks, straight jaws, jaw length 2 1/8" (5.5cm) 320-231 8 5/8" (22cm), angled shanks, straight jaws, jaw length 2 1/8" (5.5cm)



GLOVER Coarctation Clamps

320-238	8 5/8" (22cm), straight shanks, straight jaws, jaw length 2 5/8" (6.8cm)
320-239	8 5/8" (22cm), angled shanks, straight jaws, jaw length 2 5/8" (6.8cm)
320-240	10 1/8" (25.7cm), straight shanks, straight jaws, jaw length 3" (7.5cm)
320-241	10 1/8" (25.7cm), angled shanks, straight jaws, jaw length 3" (7.5cm)





COOLEY Coarctation Clamps

320-397 8 3/4" (22.2cm), straight shanks, straight jaws, jaw length 2 1/4" (5.5cm)
320-398 8 3/4" (22.2cm), angled shanks, straight jaws, jaw length 2 1/4" (5.5cm)
320-399 9 1/4" (23.5cm), angled shanks, straight jaws, jaw length 2 3/4" (7cm)

CRAFOORD Coarctation Clamp

310-214 9 3/4" (24.8cm), longitudinal serrations

COOLEY Patent Ductus Clamps

320-395 8 3/4" (22.2cm), straight shanks, straight jaws, jaw length 2 1/4" (6cm) **320-396** 8 3/4" (22.2cm), angled shanks, straight jaws, jaw length 2 1/4" (6cm)




DeBAKEY Patent Ductus Clamps

320-233 7 5/8" (19.5cm), straight shanks, straight jaws, jaw length | 5/8" (4.2cm) **320-234** 7 5/8" (19.5cm), angled shanks, straight jaws, jaw length | 1 5/8" (4.2cm)



GLOVER Patent Ductus Clamps

320-235 8" (20.5cm), straight shanks, straight jaws, jaw length | 1/4" (3cm) **320-236** 8" (20.5cm), angled shanks, straight jaws, jaw length | 1/4" (3cm)





GREGORY Profunda Clamps

5" (12.8cm), jaw length 1 1/2" (4cm), delicate 6" (15.2cm), jaw length 2" (5cm) 7 3/8" (18.7cm), jaw length 2 3/8" (6cm) 320-218 320-220 320-221

7 3/4" (19.7cm), jaw length 2 3/4" (7cm) 320-222



DALE Femoral-Popliteal Clamps


7" (17.8cm), depth 1/2" (1.3cm), small 7 1/8" (18.1cm), depth 3/4" (1.9cm), medium 320-225





RENAL Artery Clamps

320-480 7 1/4" (18.5cm), jaw length 2 1/4" (5.8cm) **320-481** 7 3/4" (19.5cm), jaw length 2 5/8" (6.5cm) **320-482** 10 3/4" (27.3cm), jaw length 3 1/2" (9cm)



HENLY Vascular Clamp

320-262 8" (20.3cm), double angled jaws, jaw length 3 1/2" (8.9cm)





DEMOS Tibial Artery Clamp

320-223 8 1/2" (21.6cm), jaw length 1 1/2" (3.8cm), perpendicular ring handles



DeBAKEY Multi-Purpose Clamps, 90° angle


320-312 7" (18cm), jaw length 2 1/8" (5.5cm)
320-314 8 1/2" (21.5cm), jaw length 2 3/4" (7cm)
320-310 9" (23cm), jaw length 1 5/8" (4cm)
320-316 10 7/8" (27.5cm), jaw length 4" (10cm)






LEE Right Angle and Bronchus Clamps

320-200 9" (23cm), right angle clamp, jaw length 1 3/8" (3.3cm) beyond angle **320-205** 9 3/4" (25cm), bronchus clamp, jaw length 1 3/8" (3.3cm) beyond angle



SAROT Bronchus Clamps

305-450 9" (22.9cm), jaw length 1 1/2" (3.8cm), right **305-451** 9" (22.9cm), jaw length 1 1/2" (3.8cm), left





HARKEN Auricle Clamps

310-190 9 1/2" (24.1 cm), cross-serrated, No. 1 curve **310-191** 9 1/2" (24.1 cm), cross-serrated, No. 2 curve **310-192** 9 1/2" (24.1 cm), cross-serrated, No. 3 curve



STILLE Vessel Clamp

475-500 9" (22.9 cm), horizontal serrations



MAYO (GUYON) Vessel Clamp

475-505 9 1/4" (23.5cm), horizontal serrations





HERRICK Kidney Clamp

475-510 9 1/2" (24.1 cm), longitudinal serrations



WERTHEIM Pedicle Clamp

475-520 10" (25.4cm), horizontal serrations



WERTHEIM-CULLEN Pedicle Clamp

475-525 8 1/4" (21cm), longitudinal serrations





OCHSNER Forceps

105-235 4 1/2" (11.4cm), horizontal serrations, 90° angle, very fine tips





BAINBRIDGE Forceps

105-300 6" (15.2cm), longitudinal serrations, straight jaws 105-302 6" (15.2cm), longitudinal serrations, curved jaws



MIXTER Pediatric Forceps

5 1/4" (13.3cm), curved 140-216 140-217 5 1/8" (13cm), right angle





MIXTER Forceps, Petit-Point, Right Angle





GEMINI Clamp

140-329 5 1/2" (14cm)



JARIT Vascular Forceps, Finely Serrated Jaws


320-630	5 1/2" (14cm), straight	
320-631	5 1/2" (14cm), curved	
320-632	5 1/3" (13.5cm), full curve	
320-633	5 1/2" (14cm), straight, 1x2 teeth	
320-634	5 1/2" (14cm), curved, 1x2 teeth	





BAILEY Forceps, Very Fine Tips

320-650 7" (17.8cm), 30° angle **320-651** 7" (17.8cm), 60° angle **320-652** 7" (17.8cm), 90° angle



SELMAN Peripheral Blood Vessel Forceps, Cooley Serrations

 320-195
 4 3/4" (12.1cm), straight jaws, jaw length 1/2" (1.3cm)

 320-196
 5 1/4" (13.3cm), straight jaws, jaw length 1" (2.5cm)

 320-197
 4 3/4" (12.1cm), angled jaws, jaw length 1/2" (1.3cm)

 320-198
 5 1/4" (13.3cm), angled jaws, jaw length 1" (2.5cm)





DeBAKEY Multi-Purpose Clamps, Jaw Length I 3/8" (3.5cm)



COOLEY Multi-Purpose Clamps, Jaw Length I 3/8" (3.5cm)

	COOLET Pluid-Purpose Clamp	s, Jaw Length 1 3/6 (3.3cm)	
320-450 320-451 320-452	6 1/2" (16.5cm), 30° angle 6" (15.2cm), 60° angle 5 1/2" (14cm), 90° angle		
\			

CASTANEDA Multi-Purpose Clamps, Jaw Length I 3/8" (3.5cm)

	CASTAIREDA Fluid-Furpose Cial	inps, Jaw Length 1 3/6 (3.3cm)	
320-600 320-601 320-602 320-603	4 7/8" (12.4cm), curved shanks 5" (12.7cm), 30° angle 4 7/8" (12.4cm), 50° angle 4 1/4" (10.8cm), 90° angle		
			3
-			



COOLEY Pediatric Multi-Purpose Clamps, 90° Angle

320-416 4 1/2" (11.4cm), jaw length 3/4" (1.9cm) **320-417** 4 1/2" (11.4cm), jaw length 1" (2.5cm) **320-418** 4 1/2" (11.4cm), jaw length 1 3/8" (3.5cm)





DARDIK Multi-Purpose Clamps - "D" CLAMP

320-800 5 7/8" (15 cm), jaw length 1" (2.7cm), straight 320-801 5 7/8" (15 cm), jaw length 1 1/4" (3cm), 15° angle 320-802 5 7/8" (15 cm), jaw length 1 1/4" (3cm), 30° angle 5 7/8" (15 cm), jaw length 1 1/4" (3cm), 45° angle 320-803 5 1/2" (14 cm), jaw length | 1/4" (3cm), 90° angle 6 3/4" (17 cm), jaw length | 1/2" (4cm), straight 6 3/4" (17 cm), jaw length | 1/2" (4cm), 15° angle 320-804 320-810 320-811 6 3/4" (17 cm), jaw length 1 1/2" (4cm), 30° angle 320-812 6 3/4" (17 cm), jaw length 1 1/2" (4cm), 45° angle 320-813 320-814 6 1/3" (16 cm), jaw length 1 1/2" (4cm), 90° angle









MICRO CLAMPS, Pediatric, DeBakey Serrations, Titanium

320-700TI 4 3/4" (12.1cm), straight shanks, jaw length 1 1/4" (3.2cm)

320-701TI 4 3/4" (12.1 cm), curved shanks, jaw length 1 1/4" (3.2 cm) **320-702TI** 4 1/4" (10.8 cm), angled shanks, jaw length 7/8" (2.2 cm)

320-708TI 4 5/8" (11.7cm), acutely angled shanks, jaw length 1 1/4" (3.2cm)









MICRO CLAMPS, Pediatric, DeBakey Serrations, Titanium

320-703TI 320-709TI 320-704TI 320-710TI 4 1/4" (10.8cm), straight shanks, ,45° angled jaw, jaw length 1" (2.5cm) 4 1/2" (11.4cm), straight shanks, 45° angled jaw, jaw length 1 1/4" (3.2cm) 4 1/8" (10.5cm), straight shanks, ,60° angled jaw, jaw length 7/8" (2.2cm) 4 1/4" (10.8cm), straight shanks, 60° angled jaw, jaw length 1 1/4" (3.2cm)











MICRO CLAMPS, Pediatric, DeBakey Serrations, Titanium

 320-706TI
 4 1/8" (10.5cm), straight shanks, jaw length 1/2" (1.3cm), depth 1/4" (6mm)

 320-705TI
 4 1/8" (10.5cm), slightly angled shanks, jaw length 1/2" (1.3cm), depth 1/4" (6mm)

 320-711TI
 4 3/8" (11.1cm), curved shanks, jaw length 7/8" (2.2cm), depth 1/4" (6mm)





.



.



MICRO CLAMPS, Pediatric, DeBakey Serrations, Titanium

4 1/8" (10.5cm), slightly angled shanks, jaw length 5/8" (1.6cm), depth 3/8" (1cm) 4 3/8" (11.1cm), slightly angled shanks, jaw length 3/4" (2cm), depth 3/8" (1cm) 4 5/8" (11.7cm), straight shanks, spoon-shaped jaw length 1" (2.5cm), depth 1/4" (6mm) 320-707TI

320-712TI

320-713TI









DeBAKEY Patent Ductus Clamps, Delicate

320-199 6" (15.2cm), straight shanks, jaw length 1" (2.5cm), 2.5mm tips 320-453 6 3/4" (17.1cm), straight shanks, jaw length 2" (5.2cm), 2mm tips 320-454 6 3/8" (16.2cm), angled shanks, jaw length 1 7/8" (4.7cm), 2mm tips

COOLEY Patent Ductus Clamps

320-400 6 3/4" (17cm), straight shanks, straight jaws, jaw length 1 1/8" (3cm) 320-401 6 1/4" (16cm), angled shanks, straight jaws, jaw length 1 1/8" (3cm)

DeBAKEY Coarctation Clamps, Delicate, 2mm Tips

320-455 7" (17.8cm), straight shanks, straight jaws, jaw length | 3/8" (6cm) 7" (17.8cm), angled shanks, straight jaws, jaw length 1 3/8" (6cm) 320-456 320-457

7" (17.8cm), straight shanks, spoon-shaped jaws, jaw length 1 3/4" (4.4cm)

COOLEY Coarctation Clamps

320-405 6 1/4" (16.cm), angled shanks, straight jaws, jaw length 1 1/8" (3cm) 320-403 6 7/8" (17.5cm), straight shanks, straight jaws, jaw length 1 3/4" (4.5cm) 320-406 6 3/4" (17.cm), straight shanks, curved jaws, jaw length 1 3/4" (4.4cm) 7" (18cm), straight shanks, spoon-shaped jaws, jaw length 2 3/8" (6cm) 320-410




COOLEY Pediatric Clamps

320-411 5 1/2" (14cm), straight shanks, straight jaws, jaw length 1 1/8" (3cm) 320-412 5 1/4" (13.3cm), straight shanks, angled jaws, jaw length 1 1/8" (3cm)

5 1/4" (13.3.cm), angled shanks, straight jaws, jaw length 1 1/8" (3cm) 5 1/2" (14cm), straight shanks, spoon-shape jaws, jaw length 1 1/8" (3cm) 320-413

320-414



DeBAKEY Peripheral Vascular Clamp

320-265 6 1/4" (16cm), curved jaws, jaw length 1 3/4" (4.4cm)



BECK Vessel Clamps

320-470 6" (15.2cm), straight shanks, jaw length 1" (2.5cm), depth 3/8" (1cm)





MULTI-PURPOSE CLAMPS, Cooley Serrations, Angled Shanks


320-475 6 1/2" (16.5cm), jaw length 1 1/2" (3.8cm), depth 5/8" (1.5cm) **320-476** 6 1/2" (16.5cm), jaw length, 2" (5cm), depth 3/4" (2cm)



KITZMILLER Clamps, Cooley Serrations, Jaw Length I 1/8" (3cm), Stainless Steel

320-540 5 1/4" (13.3cm), left angled **320-542** 5 1/4" (13.3cm), right angled





KITZMILLER Micro Clamps, DeBakey Serrations, Jaw Length 7/8" (2.3cm), Titanium

320-540TI 4 1/4" (10.8cm), left angled 4 1/4" (10.8cm), right angled








COOLEY Neonatal Vascular Clamps

5 1/8" (13cm), straight shanks, straight jaws, jaw length 1 1/2" (3.8cm) 320-500 5 1/8" (13cm), angled shanks, straight jaws, jaw length 1 5/8" (4.2cm) 5 1/8" (13cm), curved shanks, angled jaws, jaw length 1 5/8" (4.2cm) 320-501 320-502 320-503 5 1/8" (13cm), straight shanks, curved jaws, jaw length 1 5/8" (4.2cm) 5 1/8" (13cm), curved shanks, jaw length 3/4" (2cm), depth 1/4" (5mm) 5 1/8" (13cm), jaw length 1 1/2" (3.8cm), carotid clamp 320-504

320-505







CALMAN Carotid Closure Clamp

320-520 4 7/8" (12.4cm), jaw length 1 1/4" (3.2cm), depth 1/4" (7mm)



THOMPSON Carotid Closure Clamp

320-525 6" (15.2cm), jaw length 1" (2.5cm), depth 1/8" (4mm)



EXAMPLE Clamp

320-530 6 1/4" (15.9cm), jaw length 13/8" (3.5cm), depth 3/8" (9mm)



KARTCHNER Internal Carotid Clamp

320-535 6 1/2" (16.5cm), jaw length 2" (5cm)





GREGORY Profunda Clamp

320-218 5" (12.8cm), jaw length 1 1/2" (4cm), delicate



BECK Miniature Aorta Clamp

320-487 5 1/2" (14cm), jaw length 1" (3cm), depth 1/4" (6mm)



SUBRAMANIAN Miniature Aorta Clamp

320-490 5 1/2" (14cm), jaw length 1 1/8" (3cm) beyond curve





COOLEY-SATINSKY Clamp

320-415 5 1/2" (14cm), jaw length 3/4" (2cm), depth 1/4" (5mm)



COOLEY Pediatric Vascular Clamps

320-506 5 1/4" (13.3cm), curved shanks, jaw length 5/8" (1.6cm), depth 1/4" (5mm) 5 1/4" (13.3cm), curved shanks, jaw length 1/2" (1.2cm), depth 1/4" (5mm) 5 1/4" (13.3cm), curved shanks, jaw length 1/3" (8mm), depth 1/4" (5mm)



COOLEY Vascular Clamps, 90° Angle

320-440 6 1/2" (16.5cm), jaw length | 1/8" (2.7cm), depth 3/8" (1cm) **320-441** 6 1/2" (16.5cm), jaw length | 1/4" (3.4cm), depth | 1/4" (7mm)






REYNOLDS Clamps, 85° Angle

320-485 6 1/4" (15.9cm), jaw length 1" (2.6cm), depth 1/4" (5mm) 320-486 6 1/4" (15.9cm), jaw length 1 1/8" (3cm), depth 1/4" (6mm)



ANASTOMOSIS Clamps

320-420 6 3/4" (17cm), angled shanks, jaw length 7/8" (2.3cm), depth 3/8" (1cm), DeBakey serrations **320-422** 6 1/2" (16.5cm), angled shanks, jaws angled left, jaw length 3/4" (2cm), depth 3/8" (1cm), Cooley serrations





JARIT Anastomosis Clamp

320-425 6 3/4" (17cm), jaw length 1 5/8" (2.5cm), depth 3/8" (1cm), Cooley serrations




CASTANEDA Anastomosis Clamps

320-610 5 3/4" (14.6cm), jaw length 5/8" (1.5cm), depth 1/4" (5mm) 6" (15.2cm), jaw length 6/8" (1.8cm), depth 1/4" (6mm) 6 1/8" (15.6cm), jaw length 7/8" (2.1cm), depth 1/3" (7mm)

DERRA Anastomosis Clamps

320-430 6 3/4" (17cm), jaw length 3/4" (1.9cm), depth 3/8" (1cm), small 320-431 6 3/4" (17cm), jaw length 1" (2.5cm), depth 1/2" (1.3cm), medium 320-432 6 3/4" (17cm), jaw length 1 1/4" (3.3cm), depth 5/8" (1.5cm), large





COOLEY-DERRA Anastomosis Clamps

6 1/2" (16.5cm), jaw length 5/8" (1.7 cm), depth 1/4" (8mm), small 6 1/2" (16.5cm), jaw length 7/8" (2.2cm), depth 3/8" (1 cm), medium 320-435 320-436 320-437 6 1/2" (16.5cm), jaw length 1 1/8" (3cm), depth 1/2" (1.3cm), large



DERRA Partial Occlusion Clamps

4 1/2" (11.4cm), angled shanks, jaw length 3/8" (1cm), short, small curve 4 3/8" (11.1cm), angled shanks, jaw length 5/8" (1.6cm), short, medium curve 320-433 320-463 4 I/2" (11.4cm), angled shanks, jaw length 1 I/8" (2.8cm), short, large curve 5 I/2" (14cm), curved shanks, jaw length 3/4" (1.7cm), small curve 5 I/2" (14cm), curved shanks, jaw length 7/8" (2.2cm), medium curve 320-464 320-465 320-466



CASTANEDA Partial Occlusion Clamp

320-604











LILLEHEI-POTTS Scissors, 7" (17.8cm)

100-240	Straight blades, blunt tips			
100-241	Curved blades, blunt tips			
100-242	Straight blades, sharp tips			
100-243	Curved blades, sharp tips			
100-244	Straight blades, tenotomy tips			
100-245	Curved blades, tenotomy tips			





LILLEHEI-POTTS Scissors, 9" (22.9cm), SUPERCUT

102-305	Curved blades, sharp tips
102-306	Curved blades, blunt tips
102-307	Curved blades, tenotomy tips









REYNOLDS Scissors, Narrow Tip, Curved

 100-295
 5 1/2" (14cm), plain edge, very narrow tips

 100-296
 5 1/2" (14cm), serrated edge, very narrow tips

 100-297
 7" (17.8cm), plain edge, narrow tips

 100-298
 7" (17.8cm), serrated edge, narrow tips



JAMISON-METZENBAUM Scissors, Narrow Tip, Curved

 100-300
 6" (15.2cm)

 100-302
 7" (17.8cm)

 100-304
 8" (20.3cm)

 100-306
 9" (22.9cm)



METZENBAUM (LAHEY) Scissors, CARB-EDGE

101-251 5 3/4" (14.6cm), curved



METZENBAUM Scissors, CARB-EDGE

 101-258
 7" (17.8cm), curved, blunt/blunt tips

 101-259
 7" (17.8cm), curved, delicate, sharp/sharp tips

 101-260
 7" (17.8cm), straight, delicate, blunt/blunt tips

 101-261
 7" (17.8cm), curved, delicate, blunt/blunt tips





DISSECTING Scissors, "Flatback" SUPERCUT

102-120	5 3/4" (14.6cm), curved
102-122	7" (17.8cm), curved
102-310	7" (17.8cm), curved, delicate
102-311	8" (20cm), curved, delicate
102-125	9" (22.9cm), curved
102-126	10" (25.4cm), curved
102-128	11" (27.9cm), curved
102-129	12" (30.5cm), curved
102-127	14" (35.6cm), curved
102-245	5" (12.7cm), curved







METZENBAUM (NELSON) Scissors

305-100	9" (22.9cm), straight
305-101	9" (22.9cm), curved
305-105	I I'' (27.9cm), straight
305-106	11" (27.9cm), curved
305-110	14 1/2" (36.8cm), curved



METZENBAUM (NELSON) Scissors, CARB-EDGE

101-266	9'' (22.9cm), straight, blunt tips
101-267	9" (22.9cm), curved, blunt tips
101-275	9" (22.9cm), straight, delicate tip
101-276	9" (22.9cm), curved, delicate tip
101-278	9" (22.9cm), straight, sharp tips
101-279	9" (22.9cm), curved, sharp tips
101-269	10" (25.4cm), curved, blunt tips
101-270	I I'' (27.9cm), curved, blunt tips
101-272	12" (30.5cm), curved, blunt tips
101-256	14" (35.6cm), curved, blunt tips





DEMARTEL Vascular Scissors, Curved

310-302 6 7/8" (17.5cm), blunt tips **310-303** 6 7/8" (17.5cm), sharp tips



JARIT Endarterectomy Scissors, SUPERCUT

102-315 6 3/4" (17.1 cm), curved



JARIT Endarterectomy Scissors

310-330 6 3/4" (17.1 cm), curved



DeBAKEY Endarterectomy Scissors

310-349 7" (17.8cm), full curve





REYNOLDS, Petit-Point Scissors, SUPERCUT

102-245 5" (12.7cm), curved, tenotomy tips



REYNOLDS, Petit-Point Scissors, CARB-EDGE

101-255 5" (12.7cm), curved, one blade serrated



JAMISON Scissors, SUPERCUT

102-300 5 1/2" (14cm), curved **102-302** 7" (17.8cm), curved





JARIT DISSECTING Scissors, Tenotomy Tips, SUPERCUT

102-313 7" (17.8cm), curved



LINCOLN Scissors

310-369	6 3/4" (17.1cm), straight
310-370	6 3/4" (17.1cm), curved
310-367	10" (25.4cm), straight
310-368	10" (25.4cm), curved



MAYO Scissors, SUPERCUT

102-115	9'' (22.9cm), straight
102-116	9" (22.9cm), curved
102-113	10" (25.4cm), straight
102-118	11" (27.9cm), curved
102-114	12" (30.5cm), straight



MAYO-HARRINGTON Scissors, CARB-EDGE

101-225	9" (22.9cm), straight
101-226	9" (22.9cm), curved
101-230	14" (35.6cm), straight
101-232	14" (35.6cm), curved







MAYO-HARRINGTON Scissors

 305-140
 9" (22.9cm), straight

 305-141
 9" (22.9cm), curved

 305-145
 11" (27.9cm), straight

 305-146
 11" (27.9cm), curved





JORGENSON Scissors

305-120 8 3/4" (22.2cm), curved



SATINSKY Scissors

310-300 9 3/4" (24.8cm), curved shanks, acutely curved tips





COOLEY Scissors

305-125 7 1/2" (19.1cm), curved



WILLAUER Scissors

 305-130
 10" (25.4cm), straight

 305-131
 10" (25.4cm), curved

 305-135
 12" (30.5cm), straight

 305-136
 12" (30.5cm), curved



POTTS-SMITH Dissecting Scissors, Curved

310-290 5 1/2" (14cm) **310-291** 7 1/2" (19.1cm) **310-292** 9" (22.9cm)



CHURCH Scissors, 7" (17.8cm)

310-375 Straight **310-376** Curved





MATTOX-POTTS Scissors

310-304 5 1/2" (14cm), 45° angle, 22mm blades





JARIT Peripheral Vascular Scissors, 5" (12.7cm), 16.5mm Blades

25° angle
45° angle
60° angle
90° angle












DeBAKEY-POTTS Scissors, 6 3/4" (17.1cm), 23mm Blades

310-332	25° angle
310-333	45° angle
310-334	60° angle





POTTS-SMITH Scissors, SUPERCUT, 7" (17.8cm), 21mm Blades



POTTS-SMITH Scissors, 7" (17.8cm), 22mm Blades



DeBAKEY Scissors, 23mm Blades

	DEDARLI SCISSOIS,	Lonnin Diades		
310-340 310-341 310-342 310-343 310-344 310-345	7" (17.8cm), 25° angle 7" (17.8cm), 45° angle 7" (17.8cm), 60° angle 9" (22.9cm), 25° angle 9" (22.9cm), 45° angle 9" (22.9cm), 60° angle			
				3
		-		











DIETRICH Circumflex Artery Scissors, SUPERCUT, 7" (17.8cm), 13.5mm Blades

130° angle, with stabilizing bar





310-360 I 30° angle, with stabilizing bar



FAVALORO Scissors, SUPERCUT, 6 3/4" (17.1cm), 13.5mm Blades

102-374 | 125° angle, curved shanks



FAVALARO Coronary Scissors, 6 3/4" (17.1cm), 13.5mm Blades

310-374 I 30° angle, curved shanks





KARMODY In-Situ Valve Scissors

310-535 6 1/4" (15.9cm) shaft, serrated blades, blunt tips



VALVE Leaflet Excision Scissors

310-540 4 1/4" (10.8cm), 45° angle up










WEBSTER Needle Holders, CARB-BITE Jaws

121-118 4 1/2" (11.4cm), smooth

121-119 4 1/2" (11.4cm), serrated, 3,600 teeth per square inch



CRILE-WOOD Needle Holders, 4,000 Teeth Per Square Inch, Straight, CARB-BITE Jaws

	CRILE-WOOD Needle	Holders, 4,000 feeth Fer Square Inch, Straight	, CARD-DITE Jaws
121-130 121-131 121-132 121-133 121-134 121-129 121-128	6" (15.2cm) 7" (17.8cm) 8" (20.3cm) 9" (22.8cm) 10 1/4" (26cm) 12" (30.5cm) 14" (35.6cm)		
		And the set of	

MAYO-HEGAR Needle Holders, 2,500 Teeth Per Square Inch, CARB-BITE Jaws

		, ,		,	
121-135 121-140 121-145 121-146 121-142 121-143 121-147 121-144	6" (15.2cm) 7" (17.8cm) 8" (20.3cm) 9" (22.8cm) 10" (25.4cm) 11" (27.9cm) 12" (30.5cm) 14" (35.6cm)				
			INTERNATION OF THE PARTY OF THE		



FELL Needle Holders, 15mm CARB-BITE Jaws, 2 mm Tips





JULIAN Needle Holders, 4,000 Teeth Per Square Inch, CARB-BITE Jaws

121-206 121-207 8 1/4" (2 l cm), straight 8 1/4" (2 l cm), curved

DeBAKEY Vascular Needle Holders, 7,000 Teeth Per Square Inch, Extra Light, CARB-BITE Jaws

| 121-184 6" (15.2cm) | 121-185 7" (17.8cm) | 121-186 8" (20.3cm) | 121-187 9" (22.8cm) | 121-188 10 1/2" (26.7cm) | 121-189 12" (30.5cm)







JARIT Vascular Needle Holders, 4,000 Teeth Per Square Inch, CARB-BITE, Slender Jaws



SAROT Needle Holders, CARB-BITE Jaws

 121-197
 6" (15.2cm), 4000 teeth per square inch, delicate

 121-198
 6 5/8" (16.8cm), 4000 teeth per square inch, delicate

 121-200
 7 1/8" (18.1cm), 2,500 teeth per square inch

 121-205
 10 1/8" (25.7cm), 2,500 teeth per square inch





MICROVASCULAR Needle Holders, Titanium, STARDUST, Straight

124-100TI 7 1/8" (18.1cm) **124-101TI** 8 5/8" (21.9cm) **124-102TI** 10" (25.4cm)







MASSON Needle Holder, CARB-BITE Jaws

121-180 10 1/2" (26.7cm), 2,500 teeth per square inch

STRATTE Valve Needle Holders, Double Curved, 2,500 Teeth Per Square Inch, CARB-BITE Jaws

121-195

9'' (22.8cm) 10'' (25.4cm), delicate tip 121-196








FINOCHETTO Needle Holder, Angled CARB-BITE Jaws

121-192 10 1/2" (26.7cm), 2,500 teeth per square inch



DeBAKEY Needle Holder, Extra Light CARB-BITE Jaws

121-193 9" (22.8cm), curved, 4,000 teeth per square inch



MICRO Needle Holders, STARDUST, Curved, 1.2mm Tips

 277-401
 5" (12.7cm)


 277-403
 5 3/4" (14.6 cm)

 277-406
 6 1/4" (15.9cm)

 277-409
 7 1/4" (18.4cm)

 277-413
 8 3/4" (22.2cm)

 277-415
 10" (25.4cm)





"LOW-RYDER" Needle Holders, Straight, STARDUST, I.2mm Tips

277-405 6 1/4" (15.9cm) **277-408** 7 1/4" (18.4cm) **277-412** 8 3/4" (22.2cm)



RICHARDS Micro Needle Holders, Titanium, STARDUST, Extra Fine, Straight, ImmTips

124-105TI 6 1/4" (15.9cm) **124-106TI** 7" (17.8cm)





PETIT-POINT RYDER Needle Holders, 12 mm CARB-BITE Jaws, Imm Tips

 121-171
 5" (12.7), smooth

 121-173
 6" (15.2cm), smooth

 121-175
 7" (17.8cm), smooth

 121-179
 8" (20.3cm), smooth

 121-178
 9" (22.8cm), smooth





PETIT-POINT RYDER Needle Holders, 4,000 Teeth Per Square Inch, I2mm CARB-BITE Jaws, Imm Tips

 121-170
 5" (12.7), serrated

 121-172
 6" (15.2cm), serrated

 121-174
 7" (17.8cm), serrated

 121-176
 8" (20.3cm), serrated

 121-177
 9" (22.8cm), serrated



RYDER Needle Holders, Titanium, Straight, STARDUST, I.65mm Tips

124-108TI 6" (15.2cm) 124-109TI 6 3/4" (17.1cm) 124-110TI 7 1/2" (19.1cm)

6" (15.2cm)

121-161TI



RYDER CARB-BITE Needle Holders, Titanium, Straight, 1.65mm Tips



121-166 9" (22.8cm), 4,000 teeth per square inch





RYDER Needle Holders, 4,000 Teeth Per Square Inch, 15mm CARB-BITE Jaws, 1.65mm Tips


121-160 5" (12.7)
121-161 6" (15.2cm)
121-162 7" (17.8cm)
121-163 8" (20.3cm)
121-164 9" (22.8cm)
121-165 10" (25.4cm)


JARIT Intracardiac Needle Holders, 4,000 Teeth Per Square Inch, 15 mm CARB-BITE Jaws, 2mm Tips



MICRO Needle Holders, STARDUST, Straight, 1.2mm Tips

277-400	5" (12.7cm)
277-402	5 3/4" (14.6 cm)
277-404	6 1/4'' (15.9cm)
277-407	7 1/4'' (18.4cm)
277-410	7 1/4'' (18.4cm), left handed
277-411	8 3/4'' (22.2cm)
277-414	10'' (25.4cm)







DIETRICH Microvascular Needle Holders, 7,000 Teeth Per Square Inch, CARB-BITE, I.2mm Tips

7" (17.8cm) 121-181 121-182 8" (20.3cm) 121-183 9" (22.8cm)

(IIIIIIIIII)



JARIT Wire Twister, CARB-BITE Jaws

121-257 6" (15.2cm), short delicate jaws 121-258

6" (15.2cm), short blunt jaws







JARIT Sternal Needle Holder and Wire Twister, CARB-BITE, 6mm Tips

121-260 7 3/4" (19.7cm), heavy duty





COOLEY-BAUMGARTEN Sternal Needle Holder and Wire Twister, CARB-BITE

8" (20.3cm), heavy duty, blunt jaws, textured tip





ADDITIONAL INFORMATION

We offer many different jaw surfaces which provide maximum gripping power for all suture sizes: Smooth, 1800, 2500, 4000, and 7000 teeth per square inch, **CARB-BITE** and **STARDUST**.

A unique metallurgical process bonds tungsten carbide particles to the jaws of our **STARDUST** needle holders. The jaw surfaces are lightly textured for secure handling of delicate and micro suture needles. Extra hardness assures durability and firm grip without damaging needles. These hand-crafted needle holders are light weight and superbly balanced. **STARDUST** jaws are warranted for a period of I year under normal use. The **STARDUST** surface will be restored during warranty period without charge.

All JARIT **CARB-BITE** and **CARB-EDGE** instruments are warranted for a period of 3 years under normal use. **CARB-BITE** needle holder jaw inserts will be replaced and **CARB-EDGE** scissors will be sharpened during the warranty period without charge.



DeBAKEY-DIETRICH Vascular Tissue Forceps, Straight, Micro Fine Serrations, Imm Tips

320-050 6" (15.2cm) **320-051** 7 3/4" (19.7cm) **320-052** 9 1/2" (24.1cm)



320-113 6" (15.2cm) **320-114** 7 3/4" (19.7cm) **320-115** 9 1/2" (24.1cm) **320-099** 12" (30.5cm)

DeBAKEY Forceps, Titanium, Straight, I.8mm Tips

320-111TI 7 1/8" (18.1cm)



DeBAKEY Forceps, 2mm Tips

320-110	6" (15.2cm), straight
320-111	7 3/4" (19.7cm), straight
320-112	9 1/2" (24.1cm), straight
320-109	12" (30.5cm), straight
320-108	14" (35.6cm), straight
320-105	6" (15.2cm), angled
320-106	7 3/4" (19.7cm), angled
320-107	9 1/2" (24.1cm), angled





COOLEY Forceps, Straight, 2mm Tips

320-116	6" (15.2cm)
320-117	8" (20.3cm)
320-119	9 1/2" (24 Lcm)



DeBAKEY Forceps, Titanium, Straight, 2mm Tips

320-084TI	7 3/4" (19.7 cm)
320-085TI	9 1/2" (24.1 cm)



DeBAKEY Forceps, Titanium, Straight, 2.5mm Tips

320-083TI 6" (15.2cm)





DeBAKEY Forceps, Straight, 2.7mm Tips

320-100 6" (15.2cm) **320-101** 7 3/4" (19.7cm) **320-102** 9 1/2" (24.1cm) **320-103** 12" (30.5cm) **320-104** 14" (35.6cm)



DeBAKEY Forceps, Straight, 3.2mm Tips

320-070 6" (15.2cm) **320-071** 7 3/4" (19.7cm) **320-072** 9 1/2" (24.1cm) **320-073** 12" (30.5cm) **320-074** 14" (35.6cm)



ADDITIONAL INFORMATION

Tip width is measured 1mm from the tip of the instrument. Our forceps come in 1mm, 1.2mm, 1.5mm, 1.8mm, 2mm, 2.7mm, and 3.2mm tip widths.

The most prevalent jaw types are as follows:

- STARDUST tungsten carbide particle surface
- CARB-BITE tungsten carbide inserts
- **DeBakey** jaws have a 1 x 2 rows of microscopic teeth
- Cooley jaws have two opposing rows of semi rounded microscopic teeth
- Platform tips for suturing
- The I.5mm **Ring Tip** forceps are designed for larger vessels such as the saphenous vein, femoral arteries and other I-3 mm sized vessels. The sizing pertains to the outer diameter of the ring.



GERALD Forceps, Straight, Narrow Tips, 1 x 2 Teeth

 285-125
 7" (17.8cm)

 285-127
 7 1/2" (19.1cm)

 285-129
 9 1/4" (23.5cm)



GERALD Forceps, CARB-BITE, Straight, I x 2 Teeth

122-570 7" (17.8cm), Delicate Tips



POTTS-SMITH Dressing Forceps, Serrated Tips

 130-350
 7" (17.8cm)

 130-355
 8 1/4" (21cm)

 130-360
 9 1/2" (24.1cm)

 130-367
 10 1/4" (26cm)

 130-365
 12" (30.5cm)







POTTS-SMITH Tissue Forceps, I x 2 Teeth, Serrated Tips

130-352 7" (17.8cm) 8 1/4" (21cm) 9 1/2" (24.1cm) 10 1/4" (26cm) 12" (30.5cm) 130-357 130-362 130-368 130-366



POTTS-SMITH Needle Forceps, CARB-BITE Tips

130-353 7" (17.8cm) 8 1/4" (21cm) 9 1/2" (24.1cm) 130-358 130-363 130-369 10 1/4" (26cm)



CUSHING Forceps, Semi-Sharp Dissecting End

8" (20.3cm), serrated 130-445 130-446 8" (20.3cm), I × 2 teeth, cross-serrated



CUSHING Tissue Forceps Titanium, Sharp Dissecting End, Gutsch Handle

130-425TI 7" (17.8cm), delicate 1.2mm serrated tips



CUSHING Forceps, Sharp Dissecting End, Gutsch Handle

130-425 7" (17.8cm), serrated **130-430** 7" (17.8cm), 1 × 2 teeth **130-435** 7" (17.8cm), 1 × 2 teeth, serrated

CUSHING Forceps, Semi-Sharp Dissecting End, Gutsch Handle, CARB-BITE Tips

 122-530
 7" (17.8cm), standard tips

 122-535
 7" (17.8cm), delicate tips

 130-417
 7" (17.8cm), angled









CUSHING Tissue Forceps Titanium, Sharp Dissecting End, Gutsch Handle

130-425T1 7" (17.8cm), delicate 1.2mm serrated tips



CUSHING Forceps, Sharp Dissecting End, Gutsch Handle

130-425 7" (17.8cm), serrated 7" (17.8cm), 1 × 2 teeth **130-435** 7" (17.8cm), 1 × 2 teeth, serrated



CUSHING Forceps, Semi-Sharp Dissecting End, Gutsch Handle, **CARB-BITE** Tips

 122-530
 7" (17.8cm), standard tips

 122-535
 7" (17.8cm), delicate tips

 130-417
 7" (17.8cm), angled







CUSHING Dressing Forceps, Square Dissecting End, Bayonet Handle

 130-410
 7 1/4" (18.4cm), serrated tips

 130-420
 7 3/4" (19.7cm), serrated tips

 130-415
 7 1/4" (18.4cm), 1 × 2 teeth

 130-421
 7 3/4" (19.7cm), 1 × 2 teeth

 130-418
 7 3/4" (19.7cm), CARB-BITE tips



CUSHING Forceps, Square Dissecting End

130-400 7" (17.8cm), serrated 7" (17.8cm), 1 × 2 teeth



ARREST SEEL STATE

BROWN Tissue Forceps, 8 x 8 Teeth

 130-265
 6" (15.2cm)

 130-266
 8" (20.3cm)

 130-267
 10" (25.4cm)





WHEELER-Type Plaque and Vessel Forceps, CARDIO-GRIP Tips

130-296 5 1/2" (14cm) 130-297 7 1/4" (18.4cm) 130-298 8 3/4" (22.2cm)



DeBAKEY Needle Forceps, CARB-BITE, 2 x 4 Teeth, 1.5mm Tips

 122-541
 6" (15.2cm)


 122-543
 7 3/4" (19.7cm)

 122-545
 9 1/2 (24.1cm)



SKINNER Needle Forceps, 9" (22.9cm), CARB-BITE, I x 2 Teeth

| 122-600 | Standard tips | 122-601 | Delicate tips | Micro tips |





DeBAKEY Needle Forceps, CARB-BITE, 2 x 4 Teeth, 2.0mm Tips

122-540 6" (15.2cm) 122-542 7 3/4" (19.7cm) 122-544 9 1/2 (24.1cm)





SEMKEN Needle Forceps, 6" (15.2cm), CARB-BITE, I x 2Teeth

122-590 Standard I.5mm tips 122-591 Delicate I mm tips





BONNEY Forceps, CARB-BITE Tips

122-520 6 3/4" (17.1 cm), 1 × 2 teeth





FERRIS-SMITH Forceps, CARB-BITE Tips

122-550 6 3/4" (17.1 cm), 1 × 2 teeth







NELSON Tissue Forceps, 6 x 7 Teeth

130-370	6" (15.2cm)
130-375	7" (17.8cm)
130-380	9" (22.9cm)





SELMAN Tissue Forceps, Round Cup With Teeth

135-252	6" (15.2cm)
135-253	8" (20.3cm)
135-254	10" (25.4cm)





RUSSIAN Tissue Forceps

130-390	6" (15.2cm)
130-392	8" (20.3cm)
130-395	10" (25.4cm)
130-394	12" (30.5cm)
130-396	14" (35.6cm)







SINGLEY Tissue Forceps

135-250 9" (22.9cm) **135-248** 14" (35.6cm)





WANGSTEEN Tissue Forceps, CARB-BITE Tips

 135-256
 6" (15.2cm)

 135-257
 9" (22.9cm)

 135-258
 10" (25.4cm)





REUL Coronary Tissue Forceps, Delicate, CARB-BITE, I.8mm Tips

130-437 8" (20.3cm)



EVANS Tissue Forceps, CARB-BITE Tips

130-383 10" (25.4cm)












DENNIS Micro Scissors, Round Handle, Double-Action, 7 1/4" (18.4cm), 8mm Blades

277-058 Straight blades277-059 Curved blades



JACOBSON Scissors, Bayonet Handle, 7 1/4" (18.4cm)

277-135 Straight, 8mm micro blade 277-136 Curved, 8mm micro blade





JACOBSON Micro Scissors, Round Handle, Micro Fine 9mm Blades

277-335 6 3/4" (17.1 cm), 25° angle **277-336** 6 3/4" (17.1 cm), 45° angle **277-337** 6 3/4" (17.1 cm), 60° angle **277-338** 6 1/2" (16.5 cm), 90° angle **277-339** 6 1/2" (16.5 cm), 125° angle





JACOBSON Micro Scissors, Flat Handle, Micro Fine 9mm Blades

 277-330
 6 1/2" (16.5cm), 25° angle

 277-331
 6 1/2" (16.5cm), 45° angle

 277-332
 6 1/2" (16.5cm), 60° angle

 277-333
 6 1/4" (15.9cm), 90° angle

 277-334
 6 1/4" (15.9cm), 125° angle





YOU-POTTS Scissors, Fine 10mm Blades

310-305 7 1/4" (18.4cm), 45° angle 310-306 7" (17.8cm), 90° angle 7" (17.8cm), 120° angle



JACOBSON Micro Scissors, Flat Handle, 6 7/8" (17.5cm), 12mm Blades

277-123 Straight blades, blunt tips277-124 Curved blades, blunt tips



JACOBSON Scissors, Bayonet Handle, 7 1/4" (18.4cm)

 277-130
 Straight, regular I 2mm blade

 277-131
 Curved, regular I 2mm blade

 277-132
 Angled on edge, 8.3mm blade





JARIT Microsurgical Scissors, Flat Handle, 6 7/8" (17.5cm), 14.5mm Blades

277-120 Straight blades
 277-121 Curved blades
 277-122 Angled on edge



VALVE Leaflet Excision Scissors, 45° Angled on Flat, 4 1/4" (10.8cm)

310-540



POTTS Finger Arteriotomy Scissors, 45° Angled to Side, 5 1/2" (14cm), 13mm Blades

310-542





JARIT Microsurgical Needle Holders, Extra Wide Flat Handle, Without Catch, Smooth Jaws, .5mm Micro Tip, 5 3/4" (14.6cm)

277-075 Straight jaws **277-076** Curved jaws



CASTROVIEJO Micro STARDUST Needle Holders, Round Handle, .7mm Micro Tip, 4 3/4" (12.1cm)

278-220 Straight jaws, with catch
278-221 Curved jaws, with catch
278-222 Straight jaws, without catch
278-223 Curved jaws, without catch





JARIT Microsurgical Needle Holders, Round Handle, Smooth Jaws, .7mm Micro Tip, 7" (17.8cm)

277-091Straight jaws, with catch277-090Straight jaws, without catch277-095Curved jaws, with catch277-094Curved jaws, without catch







JACOBSON Needle Holders, Flat Serrated Handle, Smooth Jaws, .8mm Micro Tip, 6 7/8" (17.5cm)

277-141 Straight jaws, with catch277-140 Straight jaws, without catch



JARIT Micro STARDUST Needle Holder, Titanium, Low Profile Box Lock, .8mm Micro Tip

278-080TI 7" (17.8cm), flat handles, straight jaws, with catch **278-082TI** 8 1/4" (21cm), flat handles, straight jaws, with catch



JARIT Micro STARDUST Needle Holder, Titanium, Low Profile Box Lock, .8mm Micro Tip

 278-141TI
 4 3/4" (12.1 cm), straight jaws, with catch

 278-093TI
 7 1/8" (18.1 cm), straight jaws, with catch

 278-095TI
 7 7/8" (20cm), straight jaws, with catch

 278-097TI
 9 7/8" (25.1 cm), straight jaws, with catch

 278-140TI
 4 3/4" (12.1 cm), straight jaws, without catch

 278-142TI
 4 3/4" (12.1 cm), curved jaws, without catch



JARIT Micro STARDUST Needle Holder, Titanium, Round Handle, Low Profile Box Lock, .8mm Micro Tip

278-113T1 5 3/4" (14.6cm), straight jaws, with catch 7" (18cm), straight jaws, with catch **278-113T1** 7 7/8" (20cm) straight jaws with catch





JARIT Microsurgical Needle Holders, Round Handle, Smooth Jaws, .9mm Micro Tip, 7" (17.8cm)

277-093 Straight jaws, with catch
277-092 Straight jaws, without catch
277-097 Curved jaws, with catch
Curved jaws, without catch



CASTROVIEJO Needle Holders, Flat Serrated Handle, Smooth Jaws, Delicate 1.5mm Tip, 5 1/2" (14cm)

355-405 Straight jaws, with catch
355-404 Curved jaws, with catch
355-404NL Straight jaws, without catch
Curved jaws, without catch



CASTROVIEJO Needle Holders, Flat Serrated Handle, Smooth Jaws, Imm Micro Tip, 5 1/2" (14cm)

355-406straight jaws, with catch355-407curved jaws, with catch355-406NLstraight jaws, without catch355-407NLcurved jaws, without catch



JACOBSON Needle Holders, Flat Serrated Handle, Delicate Imm STARDUST Tip

278-125 6 3/4" (17.1cm), straight jaws

278-125NL 6 3/4" (17.1cm), straight jaws, without catch **278-126** 6 3/4" (17.1cm), curved jaws





JARIT Microsurgical Needle Holders, Round Handle, Smooth Jaws, Imm Micro Tip, 5 3/4" (14cm)


277-072 Straight jaws, with catch 277-070 Straight jaws, without catch 277-073 Curved jaws, with catch 277-07I Curved jaws, without catch





JARIT Needle Holder, Round Handle, Imm Micro STARDUST Tip

278-110 7" (17.8cm), straight jaws, without catch 278-111 7" (17.8cm), straight jaws, with catch 278-116 7" (17.8cm), curved jaws, without catch 278-117 7" (17.8cm), curved jaws, with catch 278-112 8" (20.3cm), straight jaws, without catch 8" (20.3cm), straight jaws, with catch 278-113 278-118 8" (20.3cm), curved jaws, without catch 278-119 8" (20.3cm), curved jaws, with catch 278-114 9" (22.9cm), straight jaws, without catch 278-115 9" (22.9cm), straight jaws, with catch 9" (22.9cm), curved jaws, without catch 278-120 278-121









CASTROVIEJO Micro STARDUST Needle Holders, With Catch, Straight Jaws, Delicate 1.3mm Tip

278-040 5 3/4" Round handle 5 1/2" Flat handle 278-050







CASTROVIEJO Needle Holders, Flat Serrated Handle Straight Jaws, Delicate 1.3mm CARB-BITETip, 5 1/2" (14cm)

355-412 Smooth jaws, with catch 355-413 Serrated jaws, with catch 355-412NL Smooth jaws, without catch 355-413NL Serrated jaws, without catch







JARIT Microsurgical Needle Holders, Flat Serrated Handle, Straight Jaws, Delicate 1.3mm CARB-BITE Tip, 7" (17.8cm)

121-228 Smooth jaws, with catch 121-228NL Smooth jaws, without catch 121-229 Serrated jaws, with catch 121-229NL Serrated jaws, without catch



JARIT Microsurgical Needle Holders, Flat Serrated Handle, Straight Jaws, Delicate 1.3mm CARB-BITE Tip, 8 1/4" (21cm)

121-230 Smooth jaws, with catch 121-230NL Smooth jaws, without catch 121-231 Serrated jaws, with catch 121-23 INL Serrated jaws, without catch





JACOBSON Needle Holder, Flat Serrated Handle, Delicate 1.3mm STARDUST Tip

278-180 7 1/4" (18.4cm), straight jaws 7 1/4" (18.4cm), straight jaws 278-182 8 1/4" (21cm), curved jaws 278-181 8 1/4" (21cm), straight jaws 8 1/4" (21cm), straight jaws 278-183 278-184 8 1/4" (21cm), curved jaws 278-185 9 1/2" (24.1 cm), straight jaws 278-186







JARIT Needle Holder, Round Handle, Straight Jaws, Delicate I.3mm STARDUST Tip

278-095 8" (20.3cm), with catch



JARIT Needle Holder, Round Handle, Curved Jaws, Delicate 1.3mm STARDUST Tip

278-101 7" (17.8cm), with catch 7" (17.8cm), without catch 278-100 278-103 8 I/4" (21cm), with catch 278-102 8 1/4" (21cm), without catch 9" (22.9cm), curved jaws, with catch 278-105

278-104 9" (22.9cm), curved jaws, without catch





JARIT Micro STARDUST Needle Holder, Titanium, Round Handle, With Catch, Low Profile Box Lock, Delicate 1.3mm Tip

278-143TI 4 3/4" (12.1cm), curved jaws 6" (15.2cm), straight jaws 278-145TI 278-147TI 6" (15.2cm), curved jaws 278-069TI 278-07 ITI 278-078TI 8" (20.3cm), curved jaws





JARIT Micro STARDUST Needle Holder, Titanium, Round Handle, Without Catch, Low Profile Box Lock, Delicate 1.3mm Tip

 278-144TI
 6" (15.2cm), straight jaws

 278-146TI
 6" (15.2cm), curved jaws

 278-070TI
 7 1/4" (18.4cm), straight jaws

 278-072TI
 7 1/4" (18.4cm), curved jaws

 278-077TI
 8" (20.3cm), straight jaws

 278-079TI
 8" (20.3cm), curved jaws



JACOBSON STARDUST Needle Holder, Round Handle, With Catch, Low Profile Box Lock, Heavy Style, I.2mm Tip

 278-060
 6 1/8" (15.6cm), straight jaws

 278-061
 6 1/8" (15.6cm), straight jaws

 278-062
 6 1/8" (15.6cm), curved jaws

 278-063
 6 1/8" (15.6cm), curved jaws



JACOBSON STARDUST Needle Holder, Round Handle, With Catch, Low Profile Box Lock, Heavy Style, 2mm Tip

 278-130
 7 I/2" (19.1 cm), straight jaws

 278-131
 7 I/2" (19.1 cm), curved jaws

 278-132
 9" (22.9 cm), straight jaws

 278-133
 9" (22.9 cm), curved jaws

 278-134
 10" (25.4 cm), straight jaws

 278-135
 10" (25.4 cm), curved jaws





JACOBSON STARDUST Needle Holder, Round Handle, Without Catch, Low Profile Box Lock, Heavy Style, 2mm Tip

 278-131NL
 7 1/2" (19.1cm), curved jaws

 278-130NL
 7 1/2" (19.1cm), straight jaws

 278-133NL
 9" (22.9cm), curved jaws

 278-132NL
 9" (22.9cm), straight jaws

 278-135NL
 10" (25.4cm), curved jaws

 278-134NL
 10" (25.4cm), straight jaws



JACOBSON STARDUST Needle Holder, Titanium, Round Handle, Low Profile Box Lock, Heavy Style, 2mm Tip

 278-171TI
 7 1/8" (18.1 cm), straight jaws, with catch

 278-172TI
 7 1/8" (18.1 cm), curved jaws, with catch

 278-173TI
 7 1/8" (18.1 cm), straight jaws, without catch

 278-174TI
 7 1/8" (18.1 cm), curved jaws, without catch





AARET SESSON

COOLEY Vascular Tissue Forceps, Delicate Imm Cross-Serrated Tips

320-119 7 7/8" (20cm)

DeBAKEY-DIETRICH Vascular Tissue Forceps, Straight, Delicate Imm Tips

320-050 6" (15.2cm) **320-051** 7 3/4" (19.7cm) **320-052** 9 1/2" (24.1cm)

DeBAKEY-DIETRICH Coronary Artery Forceps, Titanium, Delicate I.3mm Tips

320-113TI 6 1/8" (15.6cm) **320-114TI** 7 1/8" (18.1cm) **320-115TI** 8 1/4" (21cm)

DeBAKEY Forceps, Titanium, 1.5mm Tips

320-080TI 6" (15.2cm) **320-081TI** 7 3/4" (19.7cm) **320-082TI** 9 1/2" (24.1cm)





DeBAKEY-DIETRICH Coronary Artery Forceps, 1.5mm Tips

320-113 6" (15.2cm) **320-114** 7 3/4" (19.7cm) **320-115** 9 1/2" (24.1cm) **320-099** 11 7/8" (30.2cm)



DeBakey Forceps, Titanium, Round Handle, I.8mm Tips

320-111TI 7 1/8" (18.1cm)



RHOTON-Type Suture Tying Forceps, Round Handle, Delicate, Imm Platform Tips

277-049 7 1/8" (18cm), straight



JARIT Micro Forceps, Titanium, Round Handle, STARDUST Platform Tips

 279-093TI
 7 1/8" (18.1cm), straight

 279-094TI
 7 1/8" (18.1cm), curved

 279-097TI
 8 1/4" (21cm), straight

 279-098TI
 10" (25.4cm), straight





JARIT Microsurgical Forceps, Round Handle, With Platform Tips, 7 1/8" (18.4cm)

 277-040
 7 I/4" (18.4cm) straight, 0.5mm micro tip

 277-044
 7 I/8" (18cm) curved, 0.5mm micro tip

 277-042
 7 I/8" (18cm) straight, Imm delicate tip

 277-046
 7 I/8" (18cm) curved, Imm delicate tip


 277-046
 7 I/8" (18cm) curved, 1mm delicate tip

 277-041
 7 I/8" (18cm) straight, 0.5mm micro tip, 1 × 2 teeth

 277-045
 7" (17.8cm) curved, 0.5mm micro tip, 1 × 2 teeth


 277-043
 7 I/8" (18cm) straight, 1mm delicate tip, 1 × 2 teeth

 277-047
 7 I/8" (18cm) curved, 1mm delicate tip, 1 × 2 teeth



JARIT Micro Suture Tying Forceps, Round Handle, With .5mm Micro Platform Tips

277-080 6" (15.2cm), straight 6" (15.2cm), curved





JARIT Micro Forceps, Round Handle, .8mm Micro DeBakey Tips

279-149 7 1/8" (18.1cm)



JARIT Micro Forceps, Round Handle, Delicate Imm STARDUST Platform Tips

 279-100
 7 1/4" (18.4cm), straight

 279-105
 7 1/8" (18cm), curved

 279-150
 7" (17.8cm), 1×2 teeth

 279-101
 8 1/4" (21cm), straight

 279-106
 8 1/4" (21cm), curved

 279-102
 9" (22.9cm), straight

 279-107
 9" (22.9cm), curved



DENNIS Micro Forceps, With Counter Balance, Round Handle, Delicate Imm STARDUST Platform Tips

 279-158
 7 I/8" (18.1 cm), straight

 279-159
 7 I/8" (18.1 cm), curved

 279-125
 8 3/8" (21.3cm), straight





MILLS Coronary Mammary Tissue Forceps

310-528 7 1/8" (18.1cm), 1mm ring tip

JARIT Micro Ring Forceps, Delicate Imm STARDUST Ring Tips

 279-110
 71/4" (18.4cm), straight

 279-130
 7 1/8" (18.1cm), curved

 279-111
 8 1/4" (21cm), straight

 279-131
 8 1/4" (21cm), curved

 279-112
 9" (22.9cm), straight



JARIT Micro Ring Forceps, Titanium, STARDUST Tips

279-083TI 6" (15.2cm), 1mm ring tips **279-088TI** 6" (15.2cm), 2mm ring tips **279-084TI** 7 1/8" (18.1cm), 1mm ring tips **279-089TI** 7 1/8" (18.1cm), 2mm ring tips







DENNIS Micro Ring Forceps, With Counter Balance, Round Handle, Imm STARDUST Tips

279-113 6" (15.2cm), straight 279-114 6" (15.2cm), curved 7 3/8" (18.7cm), straight 279-115 7 3/8" (18.7cm), curved 279-116 279-117 8 1/4" (21cm), straight 8 1/4" (21cm), curved 279-118 9 1/8" (23cm), straight 279-119 9 1/8'' (23cm), curved 9 7/8'' (25.1cm), straight 279-120 279-121 279-122 9 7/8" (25.1 cm), curved



DENNIS Micro Ring Forceps, Titanium, With Counter Balance, Round Handle, Straight, Imm STARDUST Tips

279-110TI 7 1/8" (18.1cm) **279-111TI** 8 1/8" (20.1cm) **279-112TI** 9 1/8" (23cm)



BROLI-ADSON Forceps, I.2 mm Tips

130-246	5 1/8" (13cm), straight, serrated tips
130-247	5 1/8" (13cm), curved, serrated tips
130-248	5 1/8" (13cm), straight, STARDUST tips
130-249	5 1/8" (13cm), straight, 1×2 teeth





GERALD STARDUST Forceps, Titanium, Delicate Imm Platform Tips

279-150TI 7" (17.8cm) **279-151TI** 8 1/4" (21cm)



GERALD Forceps, Titanium, Delicate Imm DeBakey Tips

285-131TI 7" (17.8cm) **285-132TI** 8 1/4" (21cm)



GERALD Forceps, Delicate 1.2mm DeBakey Tips

285-126 7" (17.8cm) **285-132** 7 7/8" (20cm) **285-130** 9 1/2" (24.1cm)



GERALD Forceps, Delicate I.3mmTips

285-120 6 7/8" (17.5cm), serrated tips **285-128** 9 1/4" (23.5cm), serrated tips **285-135** 9 7/8" (25cm), cross-serrated tips





CUSHING Tissue Forceps, Titanium, Semi-Sharp Dissecting End, Gutsch Handle

130-425TI 7" (17.8cm), delicate 1.2mm serrated tips







Certificate of Registration

This is to certify that the

Quality Management System for Medical Devices

of

BISI INTERNATIONAL PRIVATE LIMITED

at

1270, Acharyon Ka Rasta Kishanpole Bazar, Jaipur – 302 003, Rajasthan, India

has been independently assessed and is compliant with the requirements of:

ISO 13485:2016

For the following scope of activities:

Manufacturing and Trading of Non-Woven OT Gowns, Surgical Instruments, Surgical Drapes, PPE (Personal Protective Equipment) Kits, Face Mask 3 Ply, 2 Ply & N95 and COVID-19 Test Kit

Certificate Number: UQ-14385

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Certification 1st Surveillance Audit Due 2nd Surveillance Audit Due

Certificate Expiry (subject to the company maintaining its system to the required standard)

Damel ..

Authorised Signatory

20th March 2020 19th March 2021

19th March 2022

19th March 2023



This certificate is the property of UK Certification & Inspection Limited and shall be returned immediately on request.
71-75 Shelton Street, Covent Garden, London, WC2H 9JQ, United Kingdom
Website:- www.ukcertifications.org.uk, email:- info@ukcertifications.org.uk
Company No. 11847851





www.**bisiworld**.com



Support@bisiworld.com



1270, Acharyon Ka Rasta Kishanpole Bazar Jaipur 302003, Rajasthan, India





Certificate of Compliance

CE

We hereby declare that the technical file of product confirms with the requirement of Protective Equipment 89/686/EEC as Amended Regulation (EU) 2016/425 & Medical Devices 93/42/EEC as Amended 2007/47/EC.

Certificate No.: CE-4503

Manufacturer

Name : BISI INTERNATIONAL PRIVATE LIMITED

Address : 1270, Acharyon Ka Rasta Kishanpole Bazar,

Jaipur - 302 003, Rajasthan, India

Product Name : Non-Woven OT Gowns, Surgical Instruments, Surgical Drapes,
PPE (Personal Protective Equipment) Kits, Face Mask 3 Ply, 2

Ply & N95 and COVID-19 Test Kit

The manufacturer's technical documentation as required for Medical Council according to the directive 93/42/EEC has been reviewed and found to comply with the of Protective Equipment 89/686/EEC as Amended Regulation (EU) 2016/425 & Medical Devices 93/42/EEC as Amended 2007/47/EC. Any significant changes in the design or construction of the product, not agreed upon by us, this declaration will lose its validity.

This certificate is issued under the following conditions:

- It applies only to the quality system maintained in the manufacture of above referenced models and it does not substitute the design or type examination procedures, if requested.
- The certificate remains valid until the manufacturing conditions or the quality systems are changed.
- 3. The certificate validity is conditioned by positive results or surveillance audits.

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Certification 20th March 2020 1st Surveillance Audit Due 19th March 2021

2nd Surveillance Audit Due 19th March 2022

Certificate Expiry (subject to the company maintaining its system to the required standard)

Danel ..

Authorised Signatory



This certificate is the property of UK Certification & Inspection Limited and shall be returned immediately on request. 71-75 Shelton Street, Covent Garden, London, WC2H 9IQ, United Kingdom Website:- www.ukcertifications.org.uk, email:- info@ukcertifications.org.uk Company No. 11847851





www.bisiworld.com



Support@bisiworld.com

19th March 2023



1270, Acharyon Ka Rasta Kishanpole Bazar Jaipur 302003, Rajasthan, India





Certificate of Registration

This is to certify that the
Clothing for Protection against Contact with Blood and Body Fluids

BISI INTERNATIONAL PRIVATE LIMITED

at

1270, Acharyon Ka Rasta Kishanpole Bazar, Jaipur – 302 003, Rajasthan, India

has been independently assessed and is compliant with the requirements of:

ISO 16603:2004

For the following scope of activities:

Manufacturing and Trading of Non-Woven OT Gowns, Surgical Instruments, Surgical Drapes, PPE (Personal Protective Equipment) Kits, Face Mask 3 Ply, 2 Ply & N95 and COVID-19 Test Kit

Certificate Number: UQ-14384

Validity of this certificate can be verified at www.ukcertifications.org.uk/verify

Date of Certification 1st Surveillance Audit Due 2nd Surveillance Audit Due

Certificate Expiry (subject to the company maintaining its system to the required standard)

Damel ..

Authorised Signatory

20th March 2020 19th March 2021

19th March 2022

19th March 2023



This certificate is the property of UK Certification & Inspection Limited and shall be returned immediately on request. 71-75 Shelton Street, Covent Garden, London, WC2H 9JQ, United Kingdom Website:- www.ukcertifications.org.uk, email:- info@ukcertifications.org.uk Company No. 11847851





www.bisiworld.com



Support@bisiworld.com



1270, Acharyon Ka Rasta Kishanpole Bazar Jaipur 302003, Rajasthan, India



BISI INTERNATIONAL PVT LTD.. 249,BASEMENT, GIRIRAJ NAGAR, MUHANA MANDI ROAD, MANSAROVAR, JAIPUR, INDIA, RAJ-302020 PHONE NO.: +91-9460458440/9509710745 EMAIL: SUPPORT@BISIWORLD.COM BISIJPR@GMAIL.COM