


# EQUINE Blood Collection & Shipping Guidelines

**START HERE**  
**EQUINE**



Blood collection should be as **STRESS-FREE** as possible

ANIMAL ID and TUBE # REQUIRED on each Red Top Tube


IF NOT SENDING RIGHT AWAY  
1-2 Days Max in Fridge



2ml/cc of Blood


1ml/cc of Serum

Let Blood Samples Sit for at least 1 HOUR Before Putting In Fridge



WHEN READY TO SHIP

Bundle Samples in Max of 10 per Set



Cover Samples with absorbent paper



Paper Towels, Newspaper or etc...

IF SENDING RIGHT AWAY

## SHIPPING

- DO NOT send samples in a padded envelope. Samples must be sent in a box.
- Monday is the ideal time to ship samples.
- 2-day shipping required.
- Call or Email lab to provide tracking #.
- **BE AWARE OF TEST DAYS** to avoid delays.


## PAYMENT

- We accept Visa/Credit Cards, CHECKS, MONEY ORDERS and CASH. We do not advise sending cash by mail. Checks & money order should be made to 'UBRL'.


2-Day shipping required.  
GSO is recommended for California residents.



## Frozen Ice Packs & Insulation Required



TYPES OF INSULATION:  
Styrofoam Box  
Cooler Bags (Good Quality)  
Bubble Wrap  
Etc....



Check Schedule to know when tests are run



Place All Documents and Payment inside its own zip lock bag to avoid being damaged

**U** UNIVERSAL  
**B R** BIOMEDICAL  
**L** RESEARCH  
LABORATORY  
Livestock Diagnostics

**UBRL**  
1300 N. Fresno St.  
Fresno, CA 93703  
(559) 498-0820  
UBRL.ORG  
lab@ubrl.org



No PayPal Account Needed