


Mesh drain bags used for showering after surgery

Each bag holds two drains.

Cut mesh material 7 1/2 inches x 13 inches and hold in half lengthwise.

Using a 1/2 inch seam, sew 1 inch down from the top. Skip an area about 1 inch and then sew the rest of the way down the seam.



Fold so that the seam is going down the middle of the bag and sew across the bottom.

Fold over the top of the bag about 1 inch to create the channel where the ribbon will go and sew all the way around.



Cut ribbon 54 inches long, insert a safety pin to one end and thread it through the channel.

Tie the two ends of the ribbon together in a knot.




Mesh drain bags used for showering after surgery

Each bag holds two drains

Fold the mesh fabric in half lengthwise and sew the short side with a half inch se

Sew one inch. Skip one inch. Sew the rest of the seam

Turn so the seam is going down the middle of the bag



Fold down the top one inch and sew around the bag. You have created a casing:

... and then the 54 in long ribbon and tie in a knot