



# High Performance Portable Suction

## V7 G



5-Wheel Rolling Stand  
**Z260-310**

The V7 G is a high performance, ultra quiet, low-flow Continuous / Intermittent portable suction pump. Due to its ease of use and whisper quiet operation the V7 G is perfect for Tracheal and Naso-gastric applications.

- Continuous and **3 Intermittent modes:**
  - **Mode 1: 10s On / 5s Off**
  - **Mode 2: 20s On / 10s Off**
  - **Mode 3: 40s On / 20s Off**
- Digital Manometer , 0 - 150 mmHg
- Lightweight - Weighs only 4 pounds
- Power Options: 120V AC / 12V DC / Battery
- Battery Life = 120 Minutes
- Water resistant On/Off switch with flip-up protector



The V7 MAX is a premium high vacuum , high flow portable suction unit 30 L/min.

Due to its versatility it can be used in Hospitals, Ambulances, Home-Care and Emergency departments. Suitable for surgery, and out-patient surgery applications.

- Variable suction regulator , 0 - 630 mmHg
- Power Options: 120V AC / 12V DC / Battery
- LED Battery Life indicator
- Battery Life = 65 Minutes
- Water resistant On/Off switch with flip-up protector

## V7 MAX



# V7 High Performance Portable Suction



Specifications	V7 G	V7 MAX
With CR collection accessories	Z260-110	Z260-265
Power supply	120 V <sub>AC</sub> / 12 V <sub>DC</sub> / Battery	
Free air suction flow	5 L/min Continuous / Intermittent	30 L/min
Maximum vacuum	20 kPa 150 mmHg	84 kPa 630 mmHg
Power consumption	6 W (DC) 9 VA (AC)	100 VA (AC) 84 W (DC)
Internal battery	12 V – 1.3 Ah Ac-Pb (sealed lead acid)	12 V – 4.5 Ah Ac-Pb (sealed lead acid)
Battery backup	120 min	65 min
Weight	5.5 Lb	10 Lb
Bacterial filter	✓	✓
Dimensions	15" x 10" x 6" (HxWxD)	

Accessories and Parts	
Z260-302	10G Wall Bracket
Z260-310	Mobile Trolley Cart
Z260-320	Battery Charger AC/DC
Z260-319	Charging Cable 12/24V
Z260-342	Bacteria Filter
Z260-318	PC1 (1L) Reusable Collection Jar
Z260-111	Translucent Hose

For additional canister options, please contact your Vitaal Sales Associate.

Limited 2 year warranty. Excludes batteries, replacement tubing and filters.

