

11200 KEGELS IN 28 MINS

SAY NO
TO INCONTINENCE

EMSKEGEL CHIAR

- Non-invasive
- No recovery period
- Fully clothed throughout treatment




WHAT IS THE EMSKEGEL CHAIR

The EMSKEGEL chair is a non-invasive and painless device for the treatment of pelvic floor muscle weakness, urinary incontinence and vaginal muscle laxity, utilizing High Intensity Focused Electro-Magnetic to stimulate deep pelvic floor muscles and restore neuromuscular stimulation!



Key effectiveness is based on focused electromagnetic energy, in-depth penetration and stimulation of the entire pelvic floor area.

A single EMSKEGEL chair session brings more than 10,000 supramaximal pelvic floor muscle contractions compared to common Kegel exercises, which are extremely important in muscle re-education of incontinence patients.



FOCUSED ON INCONTINENCE

Repair and Tighten Pelvic Floor / Reduce Urinary Incontinence / Reduce Fecal Incontinence



- Improve Erectile Dysfunction
- Improve Men Sexual Health
- Improve Vaginal Relaxation
- Improve Women's Sexual Satisfaction
- Improve Women's Ability to Achieve Orgasm

HOW DOES EMSKEGEL CHAIR WORK?

FOR WOMEN



BEFORE

TREATMENT


AFTER

Pelvic floor muscles insufficiently support pelvic organs and affect bladder control.

Using High Intensity Focused Electro-Magnetic Technology effectively stimulates pelvic floor muscles with thousands of supramaximal contractions.

Stimulation leads to regained control over pelvic floor muscles and bladder.

FOR MEN



BEFORE

TREATMENT

AFTER

Pelvic floor muscles insufficiently support pelvic organs and affect bladder control.

EMSKEGEL effectively stimulates the pelvic floor muscles using High Intensity Focused Electro-Magnetic technology

Stimulation leads to regained control over pelvic floor muscles and bladder

ADVANTAGES



Non-invasive and no recovery time



Fda Cleared Treatment For Incontinence



75% Pad Reduction



Fully clothed throughout treatment



Walk in and walk out




Comfortable & Painless

Improvement can be observed after a single session with results improving over the treatment period.

LEATHER BACKREST

High-end and comfortable



SUPER LARGE MAGNETIC COIL

30cm large magnetic coil, maximum output magnetic force can reach 30uh

7-INCH REMOVABLE TOUCH SCREEN

The screen can be held directly in your hand, which is very convenient to use





3000W STRONG POWER SUPPLY

The unique 3000w ultra-high power supply ensures high-intensity output of the machine and the best treatment effect

STRONG MEDICAL WHEELS

The maximum load-bearing capacity can reach 200kg





TWO WORKING MODES

STRENGTH

RECOVERY

Strength mode: Strengthen the pelvic floor muscles and increase muscle strength.

Recovery mode: Relax tense muscles and reduce muscle stress.

USER FRIENDLY INTERFACE

The new intelligent operating system is easy to use and smooth to operate.

TWO MODELS, GIVING YOU MORE CHOICES



CLY-1S
Work Independently



CLY-1
Connect to the Machine to Work



FAQ

1.What is Urinary Incontinence?

Urinary incontinence is the involuntary leakage of urine, caused by weak pelvic floor muscles leading to loss of bladder control.

2.Who is the treatment suitable for?

Women or men of any age who wish to address urinary incontinence, bladder leakage and improve quality of life or as a preventive.

3.Who should not use EMSKEGEL?

Pregnant women, people with metal implants or with tumors or heart disease.

4.What does an EMSKEGEL treatment feel like?

The treatment takes about 30 minutes and you can feel the muscles contracting, like doing kegel exercises.

5.How quickly will I see results?

Usually, two or three treatment sessions.

6.How long do the results of the EMSKEGEL treatment last?

Efficacy varies from person to person, most can reach one year.

SPECIFICATIONS

Power	3000W
Output Method	Pulse Output
Control Method	Touch Control
Screen Size	7-inch Capacitive Screen
Output Frequency	1-150HZ
Net Weight	32KG
Dimension	560x620x860MM
Package Size	710x630x670MM
Package Weight	45KG
Voltage	220V/110V