





















Symbol	Explanation	Symbol	Explanation
	MD = Medical Device		Indicates that there are temperature limits for this device
	! = Consult instructions before use		Keep away from magnetic resonance equipment
	Contains a rechargeable battery		Indicates a medical device that has not been sterilized
	Keep dry and avoid excessive moisture		Do not dispose of this device in unsorted municipal waste
	CE = Signifies European technical conformity		Bluetooth® wireless or enabled technology
	FCC = Complies with United States Regulations for Radio Frequency Devices		Date of manufacture
	Serial Number		-Humidity limitation
	Keep away from heat		-Temperature limitation
	Keep dry		Class II equipment
	Follow instructions for use		Type B applied part