



## Repose 護腿套

Repose Foot Protectors Plus

減低腳步受壓

Designed specifically to minimise  
the risk of pressure damage to heels.



香港及澳門獨家代理  
Hong Kong & Macau Sole Distributor:

 健珀有限公司  
Health Park Trading Company Limited

repose®

Repose Foot Protector Range

## Repose®護腿套 Foot Protector Plus



### 靈敏充氣閥設計 Smart Valve Technology

#### 確保護腿套保持在適當充氣的狀態 Ensure the leg brace is inflated to the correct level

The 'Smart valve' ensures Repose products can never be over-inflated following the inflation instructions.

### 獨特足部設計 Offloading

#### 獨特設計減低卸載所造成壓力的風險，確保足部壓力得到完全緩解。 Unique design reduces the risk of pressure from off-loading, ensuring complete pressure relief.

Repose Foot Protector and Foot Protector Plus are designed to completely off-load pressure from the heel. A heel void area incorporated into the Repose Foot Protector provides complete heel offloading.



### XS Repose 護腿套-加細碼 (附氣泵) Repose Foot Protector Plus Extra Small

產品 Product	護腿套-加細碼 x1 & 氣泵 x1 1 x Foot Protector Plus Extra Small (Pair) & 1 x Pump
尺寸 Dimensions	L 290mm x W 200mm x H 180mm

### S Repose 護腿套-細碼 (附氣泵) Repose Foot Protector Plus Small

產品 Product	護腿套-細碼 x1 & 氣泵 x1 1 x Foot Protector Plus Small (Pair) and 1 x Pump
尺寸 Dimensions	L 330mm x W 200mm x H 205mm

## Repose®護腿套 Foot Protector Plus

! 適合不同程度有壓瘡患者  
Repose products protect patients at all levels of risk

- ✓ **VERY HIGH RISK**
- ✓ **HIGH RISK**
- ✓ **AT RISK**

### 腳跟留空設計

#### Ventilation & Visual Inspection Panel

腳跟留空設計，可增加空氣流通及容易觀察腳部傷口情況。

The open heel void of Repose Foot Protector Plus increases ventilation of the foot and so contributes to patient comfort. The open heel void also allows for easier visual inspection of the heel for wound observation and monitoring.



### 腳跟 留空設計 Ventilation & Visual Inspection Panel

### 可調校磁石貼帶 Magnaffix Adjustable Strap

專利設計配合可調校的固定磁石貼帶，磁石裝嵌在氣囊內，可以減少和病人皮膚造成的磨擦，又可穩定地保護腳跟部位。

Magnaffix is a unique, patented multi-adjustable air filled strap incorporating magnets. This retains the product quickly and easily but is also easy to clean and remove.

The magnets have been sealed within the air cell pocket barrier so that they do not come into contact with the patient's skin to minimise the risk of any pressure or skin tear damage.



### STD Repose 護腿套-中碼 (附氣泵) Repose Foot Protector Plus Standard

產品 Product	護腿套-中碼 x1 & 氣泵 x1 1 x Foot Protector Plus Standard (Pair) & 1 x Pump
尺寸 Dimensions	L 380mm x W 200mm x H 230mm.

# 尺碼指南 Sizing Guide

在選擇合適尺寸的Repose 護腿套時，您應確保在患者坐著時從腳後跟底部到膝蓋後部的中點進行測量。


When selecting an appropriate size Repose Foot Protector, you should ensure to measure from the base of the heel to the mid-point of the back of the knee, whilst the patient is sitting.

產品尺寸 Product Size	腿部測量 Leg Measurement
XS 加細碼 Extra Small	460mm - 410mm
S 細碼 Small	500mm - 450mm
STD 中碼 Standard	550mm - 490mm

最大小腿外圍尺寸 : 330mm  
Maximum calf circumference 330mm

如下圖所示。將您的測量結果與腿部尺寸表相對應，以選擇適用的產品。

This is shown in the diagram below. Correspond your measurements to the leg size chart to select an applicable product.



Description	Size (L/W/H)	Order Code
Repose Foot Protector Plus & Pump - Extra Small (Pair)	L 290 x W 200 x H 180mm	6551100
Repose Foot Protector Plus & Pump - Small (Pair)	L 330 x W 200 x H 205mm	6541100
Repose Foot Protector Plus & Pump - Standard (Pair)	L 380 x W 200 x H 230mm	6511100

如需查詢更多資訊 For more information:

香港及澳門獨家代理 Hong Kong & Macau Sole Distributor:

 健珀有限公司  
Health Park Trading Company Limited

九龍新蒲崗大有街1號勤達中心502室  
Rm 502, Midas Plaza, 1 Tai Yau Street, San Po Kong, Kowloon  
電話 Tel: (852) 2724 6880  
傳真Fax: (852) 2724 6566  
Email: [info@hptltd.com](mailto:info@hptltd.com)  
Website: [www.frontier-group.co.uk](http://www.frontier-group.co.uk)  
Repose(1)/202305/[4]V1

 <sup>®</sup>

Frontier Medical Group, Block G Newbridge Road Industrial Estate,  
Blackwood, Caerphilly NP12 2YN United Kingdom.

 © Copyright Frontier Therapeutics Limited. 2018  
All rights reserved. GP0076\_A73\_06