



VISCENTOM BIO LIMITED

1/F Hung Wai Industrial Building, No. 3 Hi Yip Street,
Yuen Long, N.T., Hong Kong

+852 3905-2260 <https://viscentom.com>



Genephyrs

珍妮芙絲

HA Filler with Lidocaine

Beauty Miracle
Easy Way to Stay Youthful

Genephyrs

珍妮芙絲



Korea Biotechnology · Beauty Miracle · Easy Way to Stay Youthful

1

Certificated
Hyaluronic Acid
Ingredients

HA ingredients are
certificated by
FDA, KFDA,
DMF, KDMF

2

High Purity of
Hyaluronic Acid

Endotoxin test results:
<0.0025IU/mg

3

High
Stability


Experienced
in product
registration and
declaration

SAFETY ZERO Technology

The BDDE (1,4-Butanediol diglycidyl ether) residue removal process is the most important part to consider because BDDE residue left after cross-linking HA, can cause serious side effects. Genephyrs uses ZERO technology in order to remove all residual BDDE which use **8-time** process can minimizes possible side effects and endures safety.



Korea High Technology · Perfectly Perfect




PRODUCTS Information

Genephyrs offers 3 kinds of products to suit different users. Customers can choose according to your needs.




01 STREAM LINE

For Superficial wrinkles, Crows' feet



02 MIDDLE LEVEL

For Deep Skin Folds, Facial Contouring



03 PROFOUND


For Facial Contouring, Chin and Cheek Augmentation



Genephyrs STREAM LINE

Genephyrs MIDDLE LEVEL

Genephyrs PROFOUND



Superficial to mix dermis

Mid to deep dermis

Deep to subcutaneous dermis

Genephyrs product has 0.3% lidocaine and different sizes of needles which helps to improve the comfort of the injection.

Genephyrs	STREAM LINE	MIDDLE LEVEL	PROFOUND
HA Concentration	24mg/mL	24mg/mL	24mg/mL
Lidocaine	0.3%	0.3%	0.3%
Needle Size	30G	27G	25G
Packing	1.1ml x 1syr.	1.1ml x 1syr.	1.1ml x 1syr.
Duration	8-12 months	12-18 months	12-18months
Storage	2-25°C	2-25°C	2-25°C
Shelf Life	24months	24months	24months

RECOMMENDED Parts for Injection


Customers can choose according to your needs.



DOCTOR'S CHOICE TSK Needles

High viscous dermal filler gels in combination with small sized needles cause high pressure inside the syringe. TSK needle allow for an improved injection performance with less risk of injection malfunction.

CONVENTIONAL NEEDLE



Genephyrs
ULTRA LIGHT



TSK has developed a special skin filling needle which has thin wall. TSK increased the inside diameter of the 27g needle to 24% and 30g needle to 17%.

- 01 No bending
- 02 No leaking
- 03 Reduce needle ejection
- 04 Compatible with all syringes
- 05 Injection comfort
- 06 Easy control



Besides, when applying high pressure or force, TSK Needle can help to prevent the problem of needle bending, leaking of filler and fall of needle. TSK Needle is the nest needle for doctor during surgery.



Korea Biotechnology · Made in Hong Kong
Keep Your Health and Natural Beauty

Viscentom Bio Ltd is a global enterprise that dedicated to demonstrate the high quality of life to all mankind by maintaining personal attractiveness. In the past 10 years, Viscentom has taken the lead in the international filler market. In order to enter the global market, we established Viscentom Bio Ltd in the global trade center – Hong Kong. Here, we not only demonstrate Korea's advanced technology, but also facilitate the communication with customers.



Obtained KIBO Certification of Venture Business for the patented BDDE removal technology



Obtained patent in Korea for the Vitamin A HA Filler



Registered PCT Patent for Long-lasting Hyaluronic acid gel composition



Obtained approval from MFDS



Meripple
Patent certification Vitamin A HA Filler

Viscentom's Meripple have immediate effect on users, which brings unexpected filler effects to customers.

Cellnoc
Natural anti-aging DERMA brand

CELLNOC can help to improve skin's health from the inside out, helping the long-lasting of hyaluronic acid.

