



variamed

TEARAWAY PATIENT DRAPE

TEARAWAY DRAPE Item 302411

TearAway patient drape for intraoperative application and sterile operative site maintenance during 3D imaging process.


The Variamed TearAway Patient Drape is designed to protect the intraoperative patient during 3D imaging and help reduce the risk of contamination during removal. Visualization is maintained through a transparent, sealed window that allows a clear view of the work site. The drapes are removed safely by pulling each drape apart at the center seam and away from the patient, not over the patient.

CLEAR VISUALIZATION

Clear polyethylene window of the drape allows visualization of work site and image-guided hardware.

UNIQUE DESIGN

Cover separates from embedded seam for safe and easy removal. Design helps reduce the risk of contamination during removal.



STERILE APPLICATION

The TearAway Drape can be placed in seconds, saving valuable operating room time. The drape is specifically designed to maintain the sterile field while being applied

SAFE REMOVAL

Drape is removed without bringing the section of the drape that falls below the sterile field above the sterile field.



More Information:

variamed.com/tearaway-drape

PO Box 22546
Chattanooga, TN
37422-2546
www.variamed.com
Phone 423.893.4060
Fax 423.893.4020

