


BUR BLOCKS



BUR BLOCK W/ SILICONE

24 Hole - \$20

16 Hole - \$16

10 Hole - \$12



ENDO BUR BLOCK

22 Hole - \$14

12 Hole - \$10



MAGNETIC BUR BLOCK W/ COVER

14 Hole - \$14

8 Hole - \$10



ALUMINUM BUR BLOCK \$20

Options: Gold, Purple, Blue

6 RA Holes/18 FG Holes



BUR BLOCK SILICONE - LOCKING

24 Hole - \$20


16 Hole - \$16



ENHANCE

\$65 FOR 30

Options: Cup, Disc and Cone



Roundflex discs \$65 for 100

Uses standard So-flex mandrels
Each pack comes with one size and grit only



GENERIC SOF-LEX DISCS

\$65 / PACK OF 100

MATRIX SYSTEMS AND POLISHING DISCS




CONTOURED MATRICES

\$30 FOR 30 PCS

Two Options:

Thickness: 0.030mm


Thickness: 0.038mm



PALODENT MATRIX – NON-STICK

\$40 (50 PCS)

Comes in 4 sizes: 3.5mm,
4.5mm, 5.5mm and 6.5mm



Sectional Matrix \$32 (50 pcs)

Each box comes with one size only

Bicuspid: 4.5mm

Molar: 5.2mm

Tall Molar: 6.4mm

Molar with Extension: 6.4mm

DECIDUOUS
3.5MM

PREMOLAR
4.5MM

MOLAR
5.5MM

MOLAR WITH
EXTENSIONS
6.5MM

PALODENT MATRIX \$37 (50 PCS)

Each box comes with one size
only



MATRIX FORCEPS
\$60 EACH



PROPHY CUPS – LATCH TYPE
\$18/Pack of 100



PROPHY CUPS – SCREW IN
\$18 FOR 100 PCS



PROPHY CUPS – SNAP ON
\$18 FOR 100 PCS

FENDER WEDGES



FENDER WEDGES \$16(10)/\$145(100)

Yellow – Large
Green – Medium
Orange – Small
Purple – Extra Small



WEDGE GUARD \$10(10)/\$85(100)

Purple – Large
Orange – Medium
Green – Small



PLASTIC V-WEDGES
\$15 FOR 100 PCS

Purple – Large
Yellow – Medium
Green – Small
Blue – Extra Small



RESIN WEDGES
\$20 FOR 100 PIECES


Options: Extra small,
Small, Medium and
Large

WEDGES



WOODEN WEDGES \$12 FOR 100

Purple – Large
Blue – Medium
White – Small
Orange – Extra Small




TULWAR WEDGES

\$19 for pack of 36
Dark Blue – Small
Light Blue – Medium
White – Large



DIAMOND WEDGES

\$18 for kit of 122 pieces



ADAPTIVE WEDGES

\$15 for 40 pieces



SECTIONAL MATRIX KIT \$40

Included in Kit: Clamps x 2
Bi-cuspid x 5, Molar x 10,
Tall Molar x 10, Sub-gingival x 5




Forceps - Ring Placement \$60 each



GARRISON SUBSTITUTE RINGS

\$60 each
Blue – premolar
Orange - molar



PALODENT STYLE RING CLAMPS

\$70 EACH



COTTON BUFF \$15 FOR 5



SHEEP SKIN BUFF \$15 FOR 5



POINTED BRUSHES \$15 FOR 5

Comes in Fine, Medium and
Coarse



RUBBER POLISHERS \$7 FOR 10

Options: Fine, Medium and
Coarse

SPIRAL POLISHERS



ASAP DISC MEDIUM \$12/PCS

SR204M – Regular Size
SR204MS – Small Size



ASAP DISC FINE \$12/PCS

SR204F – Regular Size
SR204FS – Small Size



ASAP DISC ULTRA FINE \$12/PCS

SR204UF – Regular Size
SR204UFS – Small Size



ONE GLOSS POLISHER

50 pcs + 2 mandrels \$55
25 pcs + 1 mandrel \$30



1:1 IMPRESSION GUN
\$34



1:1 IMPRESSION TIPS
\$25 FOR 50



COMPOSITE GUN
\$26



SNAP-ON POLISHER \$55

Included in Kit: Mandrels x 4
Coarse x 25, Medium x 25,
Fine x 25, Extra Fine x 25



GOAT HAIR FINE
\$15 FOR 5

Three sizes: Small, Medium,
Large



GOAT HAIR MEDIUM
\$15 FOR 5

Two sizes: Medium, Large



GOAT HAIR COARSE
\$15 FOR 5

Two sizes: Medium, Large



FENDER PREP
\$35 FOR 10
\$320 FOR 100



WEDGE GUARD (INTER GUARD)
\$15 FOR 20

Two Options:
#1 – Height: 4mm
#2 – Height: 5.5mm



MATRIX #1 REGULAR
\$20/Pack of 144



MATRIX #1 DEAD SOFT
\$30/Pack of 144
















MATRIX #2
\$20/Pack of 144



















SURGICAL DIAMOND BURS – 25MM LENGTH - \$30 (PK OF 5)

												
Figure No.	801	801	801	801	801	801	801	801	859	368	879K	10839
ISO Shape No.	001	001	001	001	001	001	001	001	166	257	299	
ISO Ø 1/10(mm)	012	014	016	018	023	029	033	035	010	023	016	016
Head Length (mm)									11	5	10	
Standard	OSA2	OSA2L	OSA3	OSA4	OSA5L	OSA5LL	OSA5L3	OSA5L4	OSC6L	OSC15	OSS3	OSA18R3
Coarse	OSA2C	OSA2LC	OSA3C	OSA4C	OSA5LC	OSA5LLC	OSA5L3C	OSA5L4C	OSC6LC	OSC15C	OSS3C	OSA18R3C
Extra Coarse	OSA2XC	OSA2LXC	OSA3XC	OSA4XC	OSA5LXC	OSA5LLXC	OSA5L3XC	OSA5L4XC	OSC6LXC	OSC15XC	OSS3XC	OSA18R3XC

MUELLER BURS - \$22 (PK OF 6)

													
Size	½	1	2	3	4	5	6	1-6				2	6
ISO Ø 1/10(mm)	006	008	010	012	014	016	018	ASSORTED				010	018
Length (mm)	28	28	28	28	28	28	28	28				34	34
Crosstech Code	MU-1/2	MU-1	MU-2	MU-3	MU-4	MU-5	MU-6	MU-1-6				MU-2L	MU-6L
Quantity Per Pack	6	6	6	6	6	6	6	6				6	6

GATE GLIDDEN BURS - \$22 (PK OF 6)

																																																																																																																																																																																																																																																																			
--	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	---	--