




TrollDental General

	<p>Aiming Ring Colour-coded, five different positions</p> <ul style="list-style-type: none">● Fits TrollByte Plus and TrollByte Kimera● Easy to apply and remove● Fully autoclavable
	<p>Bite Block For Anterior PA's</p> <ul style="list-style-type: none">● Fits TrollByte, TrollByte Plus and TrollByte Kimera holders● Easy to apply and remove● Fully autoclavable
	<p>Image Plate #1 TROLLBAG IMAGE PLATE #1 500pcs</p> <ul style="list-style-type: none">● Cover for size #1 Image Plate● Double seal● Box with 500 pcs

**Troll Shade**

Use Trollshade Filter to prevent premature curing of material during procedures. Easily attaches to operatory lights by a magnet.

**Trollbag #1**

TROLLBAG CORDLESS #1, 500 PCS

- Sensor cover for size #1 cordless sensor
- Double seal
- Box with 500 pcs

**Trollbag #2**

TROLLBAG CORDLESS #2, 500 PCS

- Sensor cover for size #2 cordless sensor
- Double seal
- Box with 500 pcs

**Trollbag #1-500's**

TROLLBAG #1, 500 PCS

- Sensor cover for size #1 sensor
- Box with 500 pcs



TrollBag #2-500's
TROLLBAG #2, 500 PCS

- Sensor cover for size #2 sensor
- Box with 500 pcs



Trollbag PSP #1
ULTRASOFT #1

- Covers for image plates size #1
- 500 pcs



Trollbag PSP #2
ULTRASOFT #2

- Covers for image plates size #1
- 500 pcs



Trollbag XL, 500's

- Covers both the sensor and the sensor holder
- For size #1 or #2 sensor
- Box with 500 pcs

**Trollbib**

Easy to use TrollBib provides cover for the patient's neck, shoulders and front - covers a wide area for better infection control.

Highly absorbent with a water-resistant backing, with ties that fasten around the patient's neck.

**TrollByte**

The TrollByte is a thin, sturdy one-piece holder, so there is no assembling of different pieces. The TrollByte has a level pin – you will always get the sensor in the right angle against the tooth and obtain pictures with accuracy.

The TrollByte is made to fit different sensors. We deliver a complete system with three models of TrollByte to suit your brand of sensor.

You only need three holders for a complete intraoral x-ray exam.

TrollByte is fully autoclavable.

**TrollByte Kimera - Twin Polymer**

The key feature of the TrollByte Kimera is the use of two different high grade polymers, joined together in a single holder.

The result is a one-piece holder with a stiff arm and a softer, adjustable claw. The stiff arm ensures accuracy of sensor placement in the patient's mouth. The softer claw allows you to finely adjust the tension to reduce the risk of tearing the sensor cover.



TrollByte Kimera Endo

For indicator images.

The height of the block guarantees that the sensor holder does not interfere with the reamers or with the rubber dam clamp.

Also for ortho radial pictures of implants.




TrollByte Lynx YELLOW #2

A yellow holder used for posterior PA and horizontal bitewing.



TrollByte Lynx BLUE #2

A blue holder used for anterior PA and vertical bitewing.



TrollByte Plus

TrollByte Plus is made in one piece. Three holders for an FMX. It's colour coded for easy identification of the image type. The aiming device of the holder shows the position of the sensor in bitewing images.

The aiming ring is a valuable tool, especially helpful when taking anterior and posterior periapicals.



TrollByte Plus Endo

A smart solution for indicator x-rays. Safe, easy and comfortable for the patient.

TrollByte Plus Endo.

The height of the block guarantees that the sensor holder does not interfere with the reamers or with the rubber dam clamp



TrollByte Refills

TrollByte Refill Kits (3pcs * accessories)

- HBW
- VBW/AF
- AMR
- STARTER KIT



TrollDapp

The TrollDapp is easy to work with. An optional grip makes it easy to hold, even in connection with using other instruments.

The TrollDapp has a serrated edge to make the cleaning of the instruments easier.

The TrollDapp is made of material that is resistant to all types of substances – chloroform, fluoride, polishing paste and more.



TrollFoil, 100pc

TrollFoil can be used under a wide variety of clinical situations including wet or dry teeth, limited opening, limited vestibular space, gaggers, and metal and non-metallic restorations.

The double-sided foil is only 8 microns thick, and it has no problem marking wet surfaces, dry surfaces or highly polished surfaces such as cast gold or BruxZir.



TrollFoil, 500pc

TrollFoil can be used under a wide variety of clinical situations including wet or dry teeth, limited opening, limited vestibular space, gaggers, and metal and non-metallic restorations.


The double-sided foil is only 8 microns thick, and it has no problem marking wet surfaces, dry surfaces or highly polished surfaces such as cast gold or BruxZir.



TrollMount, TR 1 B


For child's #0 which can also be mounted in TrollMounts for film size #2

- Pedo insert
- Box of 100 pcs



TrollMount, TR 2+2 (200)


- 4 horizontal size #2 film
- Outside measurement: 58 x 190mm
- Box of 200 pcs



TrollMount, TR 2


For occlusal pictures.

- 2 size #4 film
- Outside measurement: 120 x 140 mm
- Box of 100 pcs



TrollMount, ZR 2

- 2 vertical size #2 film
- Outside measurement: 70 x 80 mm
- Box of 200 pcs



TrollMount, TR 6

For bitewing or apical pictures.


- 6 vertical alternatively 6 horizontal size #2 films
- Outside measurement: 120 x 120mm
- Box of 100 pcs



TrollMount, TR 6 B

For bitewing or apical pictures.


- 6 vertical alternatively 6 horizontal size #0 films
- Outside measurement: 120 x 120mm
- Box of 100 pcs



TrollMount, TR 6 BWL


For bitewing or apical pictures

- 8 vertical alternatively 8 horizontal size #2 films
- Outside measurement: 120 x 140mm
- Box of 100 pcs




TrollMount, TR 6 BWL

- 6 horizontal size #3 film
- For long bitewing
- Outside measurement: 132 x 120mm
- Box of 100 pcs



TrollMount, TR 12


- 3-year bitewing history
- 12 horizontal size #2 film
- Outside measurement: 133 x 200mm
- Box of 100 pcs



TrollMount, TR 14

For bitewing and apical pictures

- 14 size #2 film
- Outside measurement: 133 x 200mm
- Box of 100 pcs



TrollMount, TR 16

For full mouth series with 2 horizontal bitewings


- 16 size #2 film
- Outside measurement: 133 x 200 mm
- Box of 100 pcs



TrollMount, TR 16L


Full mouth series with 4 horizontal bitewings

- 18 size #2 film
- Outside measurement: 122 x 293 mm
- Box of 100 pcs.



TrollMount, TR 1 H/V

TrollMount protects both sides of x-ray film longer than you'll ever need.



TrollMount, TR 16 BW


- 4-year bitewing history
- 16 vertical alternatively 16 horizontal size #2 film
- Outside measurement: 133 x 200 mm
- Box of 100 pcs



TrollMount, OPG 6X12

For assembly of Panoramic film


- Box of 100 pcs



TrollMount, PARO 22

Full mouth series with 4 vertical bitewings.


- 12 size #2 film plus 10 size #1 film
- Outside measurement: 5 3/8" x 8 3/8"
- Box of 100 pcs



TrollMount, TR 18

Full mouth series with 4 horizontal bitewings


- 18 size #2 film
- Outside measurement: 4 3/4" x 11 1/2"
- Box of 100 pcs



TrollMount, TR 20V

Full mouth series with 4 vertical bitewings


- 18 size #2 film
- Outside measurement: 4 3/4" x 11 1/2"
- Box of 100 pcs



TrollMount, TR 20 H

Full mouth series with 4 horizontal bitewings.


- 12 size #2 film plus 8 size #1 film
- Outside measurement: 4 3/8" x 13 1/8"
- Box of 100 pcs



TrollMount, TR 21 F


Full mouth series with 4 horizontal bitewing.

- 12 size #2 film plus 9 size #1 film
- Outside measurement: 4 3/8" x 13 1/8"
- Box of 100 pcs




TrollMount, TR 7-1

- 7 size #2 film
- Outside measurement: 2 3/4" x 9 1/2"
- 100 pcs




TrollMount, TR 7-2

- 4 size #2 film plus 3 size #1 film
- Outside measurement: 2 3/4" x 8 3/4"
- 100 pcs



TrollMount, TR 7-3


- 7 size #2 film
- Outside measurement: 2 3/4" x 11"
- 100 pcs



TrollMount, TR 8 C


Complete status for pedo.

- 2 size #2 film plus 6 size #0 film
- Outside measurement: 4 3/8" x 5 1/2"
- Box of 100 pcs



TrollMount, TR 20 BW


- 5-year bitewing history
- 20 horizontal size #2 film
- Outside measurement: 8 1/2" x 7 1/4"
- Box of 100 pcs



TrollMount, TR 27

For bitewing or apical pictures.


- 27 size #2 film
- Outside measurement: 7 5/8" x 11 5/8"
- Box of 100 pcs



TrollStrip CUV

Matrix bands for use with light cured composite.

- Curved, 300 pcs



TrollStrip ST

Matrix bands for use with light cured composite.

- Straight, 300 pcs