

# SURGICAL

## CATALOGUE



### PRODUCING QUALITY GOODS FOR 15 YEARS

We are Producers of quality surgical, dental, beauty and general stainless steel products in singleuse and multiuse type. We guarantee on our products and facilitate our customers with easy return policies.



For more information visit [www.mjwintl.net](http://www.mjwintl.net)

Manufacturers and Exporters of Surgical, Dental & Beauty Instruments

#### Head Office

BI24S 88/8 Bara Pathar  
Christian Town, Sialkot  
Pakistan (51310)  
0092 524 264829

#### Branch Office

Flat 1, Kinsale Grange  
57 Langley Park Road,  
Sutton-Surrey SM2-5GZ UK.  
Mob: 00447701084895

#### Contact Details

Email: [info@mjwintl.net](mailto:info@mjwintl.net)  
[mjwinternational2005@gmail.com](mailto:mjwinternational2005@gmail.com)  
Mobile: 0092 321 712 1510  
Mobile: 0092 333 371 3443



## FACTORY OVERVIEW



# INTRODUCTION

for fifteen years MJW Surgical has been supplying the Healthcare Industry with finest quality surgical instruments & innovative products for surgery. With our extensive knowledge and experience in the medical field, our products are specifically designed to offer both exceptional quality and value for money to Healthcare Professionals and Organisations. During the last few years there have been numerous developments within our company, with the Introduction of many new innovative products and services to our portfolio, which has allowed MJW Surgical to make a significant impact on the healthcare market. Our head office is based in Sialkot, Pakistan housing all Customer Service and Administrative functions. Additionally, we have an office in the UK covering Europe and the USA. Our range of specialised surgical instrumentation is uniquely manufactured in Sialkot to the highest standards of craftsmanship. Our goal is to continue to develop new, innovative medical products that will enhance the quality of surgical procedures.

### OUR SERVICE COMMITMENT

QUALITY of All surgical instruments and products manufactured and distributed by MJW Surgical conform to the highest European and International Standards with all products bearing the relevant CE mark for their classification.

### GUARANTEE OF SATISFACTION

All products undergo rigorous quality control assessment and are fully guaranteed. Any product proving unsatisfactory in normal use will be replaced or repaired free of charge following inspection.

### VALUE FOR MONEY

We believe our surgical instruments and products to be competitively priced representing excellent value for money. We are happy to offer discounts for contracts.

### QUOTATION SERVICE

We offer a no obligation quotation service whereby we will specify, code and price your requirements using our extensive product library which enables us to cross-reference most products.

### PRODUCT INFORMATION

Sectional catalogues are available for most surgical specialities.

### INNOVATION

Our objective is to work in close conjunction with all medical and nursing personnel to develop new and innovative products and services.



19-204 diam 8 mm-19-205 diam 10 mm  
BOZEMANN-FRITSCH  
Uterine Flushing Tube  
350 mm



19-206  
SCHULZE  
Salpingograph, for utero-salpingography  
in wooden case, consisting of:  
19-207 Manometer only  
19-208 Connection only, without manometer  
19-209 Uterine cannula with tape without cone  
19-210 Cone only, small  
19-211 Cone only, medium  
19-212 Cone only, large  
19-213 Record syringe, 10CC, luer-lock mount.



19-209



19-210  
19-211  
19-212



19-213

19-214  
SCHULZE  
Salpingograph, equipment like(19-206)  
but without manometer.



20-101  
THOLE  
for internal and external measure-  
ments, graduated in centimeters

20-102  
DE LEE  
graduated in centimeters  
and inches

20-103  
COLLIN  
graduated in  
centimeters

20-104  
MARTIN  
graduated in  
centimeters



for internal measurements  
between symphysis and  
promontory

20-105  
SCHNEIDER  
graduated in centimeters  
The additional straight arm  
can be replaced. The pelvimeter  
with straight arm is used for  
localisation of the portio vaginalis  
uteri and of vaginal metastasis.

20-106  
BREISKY  
graduated in  
centimeters



20-107 4 mm  
20-107A 6 mm  
20-107B 8 mm  
20-107C 10 mm  
20-107D 12 mm  
20-107E 14 mm  
20-107F 16 mm

Gyneco



Smellie  
20-107G  
33 cm / 13"



20-107J 12 mm  
20-107K 14 mm  
20-107L 16 mm  
20-107M 20 mm  
Blunt

Wallich  
20-107H  
42 cm / 16 1/2"



CUZZI  
30 cm / 12"



Toothed back



20-108  
14.5 mm  
Blunt



20-108A  
18 mm  
Blunt



20-108B  
21 mm  
Blunt



DANIS  
20-108C  
24 cm / 9 1/2"



BOER  
20-109  
29 cm  
Craniotomy forceps



SELLHEIM  
20-110  
30 cm  
Obstetric lever



DUBIOS  
27.5 cm  
Decapitating scissors  
(embryotomy)



BRAUN  
20-111  
32.5 cm  
Decapitating hook



20-112



20-113



20-114  
NAEGELE  
355 mm



20-115  
NAEGELE  
400 mm



KIELLAND  
20-115A  
42 cm



BARNES-NEVILLE  
20-115B  
43.5 cm  
With traction handle



PIPER  
20-115C  
44.5 cm



SIMPSON-BRAUN  
20-116  
305mm



KIELLAND-LUIKART  
20-116B  
39.5 cm



MCLEAN  
20-116C  
39.5 cm / 15 1/2"  
Solid blades,  
short shank.



BARTON  
20-116D  
35.5 cm / 14"  
One hinged blade



DEMLIN SUZOR  
20-117  
32 cm



BOERMA  
20-117A  
28 cm  
For caesarean section



SIMPSON  
20-118  
23.5 cm  
For caesarean section



WRIGLEY  
20-118A  
25 cm  
For caesarean section



DEWEY  
20-119A  
35.5 cm  
Obstetrical forceps  
with axis traction

MURLESS  
20-119  
25.5 cm  
Head extractor  
with hinged shaft and sliding collar lock



