

I-DEW GREEN



Entod Research Cell UK Ltd.

Quality, Innovation & Research

IS IT TIME TO CONSIDER LISSAMINE GREEN OPHTHALMIC STRIPS


1.5 mg Lissamine Green Strips

The Conjunctival Staining Expert

- Allows observation of **Bulbar Staining**, **Circumlimbal Staining**, Lid Wiper Epitheliopathy & Lid Parallel Conjunctival Folds
- Useful for **Testing Early Dry Eyes and Contact Lens-related Dryness**
- Can **Evaluate Contact Lens Fitting & Comfort**

What's Different?

- **Firmer tip** for a more convenient application
- Sterile & **CE certified**
- **Optimised Impregnation & Faster Spread** of the Lissamine Green dye for maximum performance in Ocular Staining and Surface Retention
- **Fully Automated Production** with "NO HUMAN TOUCH"
- **UK Produced & Manufactured under GMP**



Slit Lamp Set-Up:

**White Light,
Low Illumination**

ORDER TODAY



threesixty

positive-impact

