




LARGE VOLUME CONTAINERS

Genesis Industries, Inc. manufactures Back Packs which are designed to provide continuous delivery of fluids when used with our Vented Draw-Off Caps. Back Packs are available in 1 Liter, 2.5 Liter and 5 Liter capacity and various styles including Eurostyle, Round Pocketed, Square Pocketed or Triangle Pocketed. All Back Packs have a 38-415 Neck Finish and are made with FDA compliant HDPE resins. The 4 eyelets accept straps or hooks. Draw-Off Stems are available for using the Back Pack in an upright position.

The Vented Draw-Off Cap with Double Stepped Port is designed to attach large diameter flexible tubing for dispensing using an instrument. A compression seal on the inside diameter of the cap face prevents leakage. The Back Pack Shipper Cap has a ribbed outside diameter for a secure grip when sealing the cap onto the Back Pack. These leak-proof caps prevent spills during transportation or storage.

Back Pack Straps are made of abrasion resistant nylon webbing. This material makes them long-lasting and weather resistant. The nylon webbing can withstand wear and tear from friction and rubbing. The width of the strap (1.5") provides the user comfort while carrying the Back Pack weight. These straps also work well to hang the Back Pack on a chute or in a barn.



Genesis manufactures high density calf feeding systems specially designed to provide calves with continuous milk flow.

Calf-Maid® Bottles are available in 2, 3 and 4 Quart / Liter sizes and are used with patented nipples which promote natural suckling behaviors in calves. HANDI GRIP™ Threaded and Snap On Bottles are available in a 4 Quart size and can be used for colostrum storage, bottle feeding or esophageal feeding.

Genesis Instruments is a division of Genesis Industries, Inc.

Genesis Instruments manufactures accurate and reliable delivery instruments that are often used with Back Packs.
www.genesisinstruments.com

LARGE VOLUME CONTAINERS



1 LITER BACK PACKS

Available Style:
Eurostyle

Available Accessories:
Draw-Off Cap, Shipper Cap,
Tubing, Back Pack Straps

Standard Material:
High Density Polyethylene (HDPE)



2.5 LITER BACK PACKS

Available Styles:
Eurostyle, Round Pocketed, Square
Pocketed, Triangle Pocketed

Available Accessories:
Draw-Off Cap, Shipper Cap,
Tubing, Back Pack Straps

Standard Material:
High Density Polyethylene (HDPE)



5 LITER BACK PACKS

Available Styles:
Eurostyle, Round Pocketed, Square
Pocketed, Triangle Pocketed

Available Accessories:
Draw-Off Cap, Shipper Cap,
Tubing, Back Pack Straps

Standard Material:
High Density Polyethylene (HDPE)



BACK PACK ACCESSORIES

The **Vented Draw-Off Cap** with Double
Stepped Port is designed to attach large
diameter flexible tubing to a Back Pack.

The **Shipper Cap** is used to prevent spills
during transportation or storage.

Back Pack Straps are made of abrasion
resistant nylon webbing. In addition to wearing
the Back Pack, the Strap allows users to hang
the Back Pack onto a chute or in a barn.



CALF-MAID® BOTTLES & ACCESSORIES

Calf-Maid® Bottles are made with High Density
Polyethylene (HDPE) which makes them durable and
reliable. When used with Calf-Maid® patented Snap
On and Threaded Nipples, this innovative calf feeding
system will provide calves with continuous milk flow
with no leaks. Calf-Maid® Threaded and Snap On
Bottles are available in 2, 3 and 4 Quart / Liter.



HANDI GRIP™ BOTTLES & ACCESSORIES

HANDI GRIP™ Bottles have a unique design with flat
side walls and a large easy grip handle enabling users
to tilt the container at any angle. HANDI GRIP™ Bottles
and accessories are versatile allowing for colostrum
storage, bottle feeding with nipple, or esophageal
feeding with probe. HANDI GRIP™ Threaded and Snap
On Bottles are available in 4 Quart and made with
High Density Polyethylene (HDPE).



GALLON CONTAINER

The Gallon Container is a heavy-duty square
jug made with High Density Polyethylene (HDPE) and
a 38 mm neck finish. The convenient carry handle is
ideal for short- or long-term use. Genesis
manufactures a wide variety of Specialty Caps that are
compatible with the container including Vented Draw-
Off, Non-Vented Draw-Off, Shipper Caps and
Retaining Caps.

For over 30 years, **Genesis Instruments** has been manufacturing accurate and reliable delivery
instruments and systems used for animal health throughout the world. Genesis Instruments
specializes in products designed for Oral Delivery, Injection Delivery, Topical Delivery,
Water Dosing, Solid Dose Delivery and High-Volume Delivery. Popular brands include Power
Doser™, Hi-Flow™, Pump It™, Select Doser™, LectraVet™ and Mini-Cartridge™.

Genesis works directly with R & D managers, production managers and engineers to invent new
instruments and enhance existing products.

Genesis Instruments is a division of Genesis Industries, Inc.
For more information visit www.genesisinstruments.com.



Genesis Industries, Inc. – USA
1.800.826.3301
Sales@GenInd.com

Genesis Manufacturing Ltd. – United Kingdom
+44 (0) 1225 791 184
CustomerServices@GenesisIndustries.co.uk

Genesis Industries (Aust) Pty. Ltd.
+61 (0) 2 6372 6622
Sales@GenesisIndustries.com.au

Genesis Industries (New Zealand) Ltd.
+64 (0) 9 444 2115
NZ.Sales@GenesisIndustries.co.nz