



# CONTINUOUS THREADED BOTTLES

Genesis Industries, Inc. manufactures a variety of Continuous Threaded Bottles ranging from 8 ML to 1 Liter available in High Density Polyethylene (HDPE) or Low Density Polyethylene (LDPE).

Continuous Threaded Bottles are manufactured with threads that continuously wrap around the neck of the bottle. These bottles are available in a variety of neck finishes meaning the neck width (two-digit number in millimeters) and thread finish (three-digit number identifying the thread height).

Our gamma irradiated bottles are oriented and shrink packed into Mod Packs (modular packages). The Mod Packs are placed inside double poly bags in corrugated cartons for clean transfer into filling areas. Mod Packs can be opened and placed directly onto automatic filling lines.

Genesis Industries manufactures a variety of Specialty Caps including Vented Draw-Off, Non-Vented Draw-Off, Shipper Caps and Retaining Caps which conform to our Continuous Threaded Bottles. For more information on Specialty Caps visit [www.genesisindustriesinc.com](http://www.genesisindustriesinc.com).

All Continuous Threaded Bottles are manufactured from FDA compliant resins and industry standard neck finishes.



Neck Width (20 MM)

Thread Height (410)

30 ML Threaded Bottle with  
20-410 Neck Finish



Continuous Threaded Bottle used  
with Luer Lock Vented Draw-Off  
Cap and Luer Lock Syringe

Genesis Industries, Inc. is committed to ensuring our products are manufactured to high quality standards. We strive to ensure customer satisfaction by providing quality products in a timely manner with responsive communication.

# CONTINUOUS THREADED BOTTLES



**8 ML THREADED BOTTLES**

**Available Neck Finish:**  
15-425

**Standard Material:**  
High Density Polyethylene (HDPE)



**30 ML THREADED BOTTLES**

**Available Neck Finish:**  
20-410

**Standard Material:**  
High Density Polyethylene (HDPE)



**60 ML THREADED BOTTLES**

**Available Neck Finish:**  
20-410

**Standard Material:**  
High Density Polyethylene (HDPE)



**180 ML THREADED BOTTLES**

**Available Neck Finish:**  
38-400

**Standard Material:**  
High Density Polyethylene (HDPE)



**250 ML THREADED BOTTLES**

**Available Neck Finishes:**  
24-410, 28-410

**Standard Material:**  
High Density Polyethylene (HDPE)



**400 ML THREADED BOTTLES**

**Available Neck Finish:**  
45-400


**Standard Material:**  
Low Density Polyethylene (LDPE)  
High Density Polyethylene (HDPE)



**625 ML THREADED BOTTLES  
(Hook or Hang Tab)**

**Available Neck Finish:**  
28-410

**Standard Material:**  
High Density Polyethylene (HDPE)



**1 LITER THREADED BOTTLES  
(Hang Tab)**

**Available Neck Finish:**  
28-410

**Standard Material:**  
High Density Polyethylene (HDPE)

Genesis Industries, Inc. – USA  
1.800.826.3301  
Sales@GenInd.com

Genesis Manufacturing Ltd. – United Kingdom  
+44 (0) 1225 791 184  
CustomerServices@GenesisIndustries.co.uk

Genesis Industries (Aust) Pty. Ltd.  
+61 (0) 2 6372 6622  
Sales@GenesisIndustries.com.au

Genesis Industries (New Zealand) Ltd.  
+64 (0) 9 444 2115  
NZ.Sales@GenesisIndustries.co.nz