



March 26, 2024

Re: Preventative Maintenance

To Whom It May Concern,

The Co2mmander® ELITE (Product Code COM915) does not require any Preventative Maintenance Program.

The Instructional Manual for all devices does show the acceptable cleaning methods, which would be with a clean, damp cloth; or an alcohol pad. We also ask that the system be used in room temperature settings.

As always, AngioAdvancements strives to provide excellent Customer Service, so please feel free to contact our office should you have any questions or need additional assistance.

Thank you and have a great day.

Respectfully,

Frank L. Levy
Owner/Director
AngioAdvancements
PMDA, LLC dba Portable Medical Devices
239.731.6200