FEATURED PRODUCT

Edition: Dec. 2025

- √ High accuracy
- Convenient Operation
- Easier and safer to obtain specimen
- Specimens available with flexible choices



| Product Information | | | | | |
|---------------------|---|--------------|---------------|--|------|
| Cat. No. | Product | Specimen | Format | Cut-off | Pack |
| Xylazine(XYL) | | | | | |
| UX01-XYL | Xylazine(XYL)Test | Urine/Powder | Dipstick | 500ng/mL | 50T |
| UX02-XYL | Xylazine(XYL)Test | Urine/Powder | Device | 500ng/mL | 25T |
| UX04-XYL | Xylazine(XYL)Test | Urine/Powder | Dip Card | 500ng/mL | 25T |
| Pregabalin(PGE | 3) | | | | |
| OF02-PGB | Single-Drug Drug of Abuse Test -Pregabalin (PGB) | Oral Fluid | Device | 500ng/mL | 25T |
| OF03-PGB | Single-Drug Drug of Abuse Test -Pregabalin (PGB) | Oral Fluid | Direct Device | 500ng/mL | 25T |
| Marijuana(THC |) | | | | |
| OF02ST-THC | Single Drug Drug of Abuse Test Marijuana (THC) Home Test | Oral Fluid | Device | 3.5mg/mL | 1T |
| Ketamine(KET) | and γ-Hydroxybutyric acid(0 | GHB) | | | |
| D24-MDT(B) | DrinkSafe Check - GHB&KET | Drink | Card | KET:1.0g/250mL GHB:3.0g/250mL | 10T |
| D04-KET(B) | DrinkSafe Check - KET | Drink | Card | 1.0g/250mL | 10T |
| D04-GHB(B) | DrinkSafe Check - GHB | Drink | Card | 3.0g/250mL | 10T |
| Multi-Drug Rap | oid Test | | | | |
| QX5-MDT(X:2-5) | Multi-Drug Test | Breast Milk | Device | COC300/CAF2000/COT200 /THC50/ALC0.01% | 25T |

Featured Product Recommendation

MULTIPLE SPECIMEN OPTIONS FOR DOA DETECTION

(URINE & POWDER & BREAST MILK & ORAL FLUID &DRINK)













