



	ANTIBIOTIC	CODE
<b>UTI SensiRing™</b>	Enrofloxacin 5ug	ENF
	Cephalexin 30ug	CFX
	Augmentin (Clavamox) 30ug	AUG
	Gentamycin 10ug	GM
	Azithromycin 15ug	ATH
	Trimethoprim-Sulfa 25ug	T/S
	Cefpodoxime 10ug	CPD
	Chloramphenicol 30ug	C
	<b>Skin/Wound SensiRing™</b>	Cephalexin 30ug
Trimethoprim/Sulfa 25ug		T/S
Enrofloxacin 5ug		ENF
Azithromycin 15ug		ATH
Oxacillin 1ug		OX
Chlormphenicol 30ug		C
Amoxi/Clavulante 30ug		AUG
Cefpodoximine 10ug		CPD
<b>Gram Positive SensiRing™</b>	Cefalathin 30ug	KF
	Augmentin 30ug	AUG
	Gentamycin 10ug	GM
	Penicillin 10ug	PG
	Azithromycin 15ug	ATH
	Trmethoprim/Sulfa 25ug	T/S
	Clindamycin 2ug	CD
	Enrofloxacin 5ug	ENF
	<b>Gram Negative SensiRing™</b>	Tobramycin 10ug
Trimethoprim/Sulfa 25ug		TS
Ciprofloxacin 5ug		CIP
Augmentin 30ug		AUG
Doxycycline 30ug		DXT
Cephalexin 30ug		CFX
Chloramphenicol 30ug		C
Ticarsillin 75ug		TC
<b>Ear SensiRing™</b>		Gentamycin 10ug
	Ciprofloxacin 5ug	CIP
	Amikacin 30ug	AK
	Tobramycin 10ug	TN
	Neomycin 30ug	NE
	Polymycin 300ug	PB
	Ticarcillin 75ug	TC
	Ceftazadime 30ug	CAZ


**KACEY™**  
 DIAGNOSTICS  
 1854 A Hendersonville Rd  
 Asheville, NC 28803  
 Ph: 1(828)685-3569  
 Fax: 1(828)685-3571  
[www.kaceydiagnostics.com](http://www.kaceydiagnostics.com)