


# Veti Stx-2 Two Parameter Urine Strip w/Creatinine/Micro-Albumin Ratio



**KACEY**  
DIAGNOSTICS

HANDLE ↓



## Micro-Albumin/Creatinine Ratio Reference Chart

Microalbumin (mg/L)	Creatinine [ mg/dL(mmol/L )					
	10 (0.9)	50 (4.4)	100 (8.8)	200 (17.7)	300 (26.5)	600 (53)
10	X					
30						
80						
150						

Re-Test  $\geq 10$ mg/L

Normal  $>30$ mg/g  
 $>3.4$  mg/mmol

Abnormal  $>30-300$ mg/g  
 $>3.4-33.9$  mg/mmol

High Abnormal  $>300$ mg/g  
 $>33.9$  mg/mmol

X



Veti-Stx are a trademark of Kacey, Inc. © 2015 All Rights Reserved rev. 3.5.15

1854A Hendersonville Rd, Suite 112  
Asheville, NC 28803  
Ph 828-685-3569  
www.kaceydiagnostics.com

