



American Medical Equipment, Inc.

29450 Gratiot Ave, Roseville, MI 48066

PH: (586) 298-6161 FX: (866) 593-7841



LIFT CHAIRS



PETITE
LC-358P
3-Position, Full Recline, Chaise
Lounge

PETITE WIDE
LC-358PW
3-Position, Full Recline, Chaise
Lounge

SMALL
LC-358S
3-Position, Full Recline, Chaise
Lounge

MEDIUM
LC-358M
3-Position, Full Recline, Chaise
Lounge

LARGE
LC-358L
3-Position, Full Recline, Chaise
Lounge

X-LARGE
LC-358XL
3-Position, Full Recline, Chaise
Lounge

XX-LARGE
LC-358XXL
2-Position, Partial Recline, Chaise
Lounge

- Specifications:**
- Weight Capacity: 375 lbs.
 - User Height Range: 5'2" & Below
 - Seat to Floor: 18"
 - Seat Depth: 18"
 - Seat Width: 18"
 - Top of Back to Seat: 24"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 3-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 17.5"
 - Overall Chair Width: 34.5"
 - Chair Weight: 112 lbs.

- Specifications:**
- Weight Capacity: 375 lbs.
 - User Height Range: 5'3" & Below
 - Seat to Floor: 19"
 - Seat Depth: 19"
 - Seat Width: 24"
 - Top of Back to Seat: 23.5"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 3-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 17"
 - Overall Chair Width: 39"
 - Chair Weight: 118 lbs.

- Specifications:**
- Weight Capacity: 375 lbs.
 - User Height Range: 5'3" & Below
 - Seat to Floor: 19"
 - Seat Depth: 19"
 - Seat Width: 18"
 - Top of Back to Seat: 24"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 3-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 17"
 - Overall Chair Width: 35"
 - Chair Weight: 118 lbs.

- Specifications:**
- Weight Capacity: 375 lbs.
 - User Height Range: 5'4" - 6'0"
 - Seat to Floor: 20"
 - Seat Depth: 20"
 - Seat Width: 20"
 - Top of Back to Seat: 24.5"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 3-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 17"
 - Overall Chair Width: 35"
 - Chair Weight: 117 lbs.

- Specifications:**
- Weight Capacity: 375 lbs.
 - User Height Range: 5'9" - 6'2"
 - Seat to Floor: 21"
 - Seat Depth: 22"
 - Seat Width: 22"
 - Top of Back to Seat: 25.5"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 3-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 17.5"
 - Overall Chair Width: 37"
 - Chair Weight: 120 lbs.

- Specifications:**
- Weight Capacity: 500 lbs.
 - User Height Range: 5'6" - 6'1"
 - Seat to Floor: 20"
 - Seat Depth: 22"
 - Seat Width: 26"
 - Top of Back to Seat: 29"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 3-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 18"
 - Overall Chair Width: 39.5"
 - Chair Weight: 159 lbs.

- Specifications:**
- Weight Capacity: 600 lbs.
 - User Height Range: 5'6" - 6'1"
 - Seat to Floor: 20"
 - Seat Depth: 22"
 - Seat Width: 30"
 - Top of Back to Seat: 29"
 - Back Style: Button Back
 - Fully Padded Chaise: Yes
 - Recline Positions: 2-Position
 - Heat & Massage Option: Yes
 - Distance From Wall: 18"
 - Overall Chair Width: 43"
 - Chair Weight: 168 lbs.

Price: \$865

Price: \$935

Price: \$865

Price: \$880

Price: \$880

Price: \$1,220

Price: \$1,370

WE ACCEPT INSURANCE PAYMENTS TO BE APPLIED TOWARDS THE COST OF THE CHAIR!

Color Options:



Black Cherry Pacific Spruce Stone Walnut

**Colors on printed paper may vary slightly on chair*

**** IMPORTANT ** Please be sure to choose the appropriate chair for your needs because chairs are non-returnable & non-refundable**