

**Warning:**

It can only be used as
Mask prevent to the
liquid from overflowing

Face Shield

Direct Splash Protection

Use Instructions: Please
remove the outer plastic
protective film before using

Direct Splash Protection Face Shield

Key Features:

- the face shield is made of polymer materials and high quality PET.
- which is lightweight and comfortable to wear on.
- High transparent PET lens is anti-fog, which ensures clear visibility in working condition.
- Single use.

Application:

The face shield can be used in medical institutions to block body fluids, blood spattering or splashing during medical treatments.

Instructions for use:

- Open the package and take out the face shield.
- Place the foam on forehead and pull the webbing back to the head and fix the face shield at a proper position.
- Do not touch the transparent lens when using.
- Discard according to medical waste protocol after use.
- Store in cool dry place.
- Not for medical use

Approvals and Standards:

- EN 166:2001

Product Name: Protective Face Shield

Product model: XYM-FS001

Material: eco Anti-fog PET material

Material thickness: 0.25mm=250microns

Product Size: 33cm W x 22cm H

Date Of Production: 2020/04/14

Period Of Validity: 3 Years

Lot Number: 20200414



Ordering Information

Product Name	Model No.	Qty / opp bag	Qty/master polybag	Qty/Carton	Dimension (cm)	G.W/N.W (Ctn)
Protective Face Shield	XYM-FS001	1pc	/	200pcs	55x35x36	8.30/7.50 Kgs