

Mix 5 Presets

Handheld Applicator Settings

Pathology	Frequency (Hz)	Time	Intensity (%)	Energy (J)
Sinusitis	100	2:00	50	31.38
Trigeminal Neuralgia	100	2:00	50	31.38
Vasomotor Rhinitis	100	2:00	50	31.38
Classical Headache	400	2:00	50	35.52
*Migraine	400	2:00	50	35.52
Acne	500	2:00	50	36.91
Brachial Neuralgia	500	2:00	50	36.91
Dislocation	500	2:00	50	36.91
Hand/Foot Arthritis	500	2:00	50	36.91
Contusion	700	2:00	50	39.67
Knee Arthrosis	700	2:00	50	39.67
Lumbago	700	2:00	50	39.67
Muscular Contraction	700	2:00	50	39.67
Sprain	700	2:00	50	39.67
a Tendinitis	700	2:00	50	39.67
Cervical Arthrosis	700	3:00	50	59.51
4 Epicondylitis	700	3:00	50	59.51
Muscle Sprain Inferior Limbs	700	3:00	50	59.51
Scapulo-Humeral Periarthritis	900	2:00	50	42.44
7 Backache	900	3:00	50	63.66
Lumbosciatica	900	3:00	50	63.66
Sciatic Neuralgia	900	3:00	50	63.66
Burn	1500	2:00	50	50.73
Dorsal Arthrosis	1500	2:00	50	50.73
Haematoma	1500	2:00	50	50.73
Muscular Tear	1500	2:00	50	50.73
Varix Ulcer	1500	2:00	50	50.73
Wound	1500	2:00	50	50.73
Hip Arthrosis	1500	3:00	50	76.10
Lumbosacral Athrosis	1500	3:00	50	76.10
Trigger Point Muscle	10	0:30	50	7.5
4 TMJ	350	3:00	50	52
5 Deprogramming (Dental)	900	3:00	50	
Mike Smith's Protocols				
6 Sinusitis	200	3:00	50	
Vertigo	200	0:30	50	
3 Pulpitis	800	3:00	50% Hit Buccal and Lingual	
10 Antimycotic	1500	5:00	50	
Xerostomia	250	3:00	50	
Tinnitus	200	0:30	50	
1 Extraction (intraoral)	800	3:00	50	
Extraction (extra-oral)	900	3:00	50	
2 Implant	Same as Extraction			
Trigger point	10	0:23	25 Up to 3x per trigger	

Charlie Head Settings

Frequency (Hz)	Time	Intensity (%)	Energy (J)
100	2:00	50	90.9
100	2:00	50	90.9
100	2:00	50	90.9
400	2:30	50	117
400	2:00	50	93.59
500	2:30	50	118.1
500	4:00	50	189
500	2:30	50	118.1
500	2:30	50	118.1
700	2:30	50	120.3
700	3:30	50	168.5
700	5:30	50	264.8
700	3:30	50	168.5
700	3:30	50	168.5
700	2:30	50	120.3
700	4:00	50	192.6
700	3:00	50	144.4
700	7:30	50	361.1
900	4:00	50	196.1
900	8:00	50	392.3
900	10:00	50	490.3
900	9:00	50	441.4
1500	3:30	50	181.1
1500	8:30	50	439.8
1500	3:30	50	181.1
1500	5:00	50	258.7
1500	2:30	50	129.3
1500	2:30	50	129.3
1500	8:00	50	414
1500	9:30	50	491.6