



Caremed Alrick.

Mobilise® CH5 Layflat® care chair

An innovative early
mobilisation chair®,
made in the UK



Mobilise CH5



*patent pending

An innovative, early mobilisation chair that lays 100% flat and has a unique legrest that encourages sit to standing, this chair is manufactured in the UK, using UK materials which greatly reduces leadtimes.

Specifications

Lay flat dimensions:	1830mm x 560mm
Seat width	510 - 560mm
Height	520 - 820mm
Hi-lo range	490 - 790mm

Wheel castors: total & directional lock 125mm


Safe working load 160kg

Features

- ☒ Multi-purpose early mobilisation chair
- ☒ Lays 100% flat
- ☒ Sit to standing option
- ☒ Seat tilt with Trendelenburg
- ☒ Works with the Arjo Steady to facilitate a seamless transfer
- ☒ Retractable legrest - fully retracts beneath the seat so easier sit to standing
- ☒ Armrest width adjustable
- ☒ Headrest adjustable
- ☒ Strong, padded removable armrests

The CH5 chair has thick cushion base with different seating options – depending on patient

- ☒ Easily remove the seat and back cushions for cleaning
- ☒ Tough, soft touch Fortitude® material, with minimal stitching



Averages 10-12cm between bed & chair



Mobilise Cushion

The Mobilise CH5 Chair is built to encourage early mobilisation of patients, designed to make rehabilitation becomes easier, quicker and safer.

Some additional functions:

Features

- ☒ Optional IV pole and bracket for the Mobilise chair
- ☒ Now available with optional wings or adjustable laterals
- ☒ Neuro and mouldable headrests available



Manual Handling

Lay flat option reduces the need for hoists



Retractable legrest

Fully retracts beneath the seat so easier sit to standing



Easy to move

Allows patients to be moved outside or to other wards if necessary



Versatile

Simple and easy to use handset with icons for ease of use



Seating options

Depending on the needs of the patient, different cushioning seating options are available



No stooping

It is easier to administer care to the patient as you can get far closer when the patient is in the chair



Armrests are designed to hang catheter bags and drains on



Portable ventilators, monitors, etc., can all be hung from the back of the chair



The IV pole & bracket

It is directly attached to the patient's chair, facilitating easier patient care and monitoring



With optional lateral support

Ideal for patients needing thoracic or trunk support



Optional shoulder wings

No need for a spring; simply lift and swing them out of the way. They have 3 positions: closed, splayed, and open



<https://bit.ly/layflat-function>

Watch our video here for a simple pat slide demo using our Mobilise® CH5 chair with the Layflat® feature

You can purchase, hire or lease these chairs, making acquisition easy regardless of your funding preferences or budget situation.

Aftersales service:

We can work with your trust to provide a simple consignment stock / call off agreement for parts with pricing fixed for 12 months at a time.

This keeps your chairs in action and maximises their lifetime and availability for patients. During the lifetime of the product, we offer a full repair and maintenance service, and at the end of life, a recycle-trade-in / trade-up scheme.



Caremed Alrick.

Contact us for a demo, hire
or purchase:

updates@caremed.healthcare
0333 242 1250

Unit 4a Carlyon Road,
Carlyon Road Industrial Estate,
Atherstone
North Warwickshire, CV9 1JE

The Mobilise chair is made in the UK and therefore
has shorter leadtimes

 www.caremed-alrick.com