

#### Executive Summary

#### The UroDapter:

- A. Provides a painless and easy to use alternative to the catheter
- B. Ergonomically designed for optimal efficiency
- C. Painless and infection free alternative to patients
- D. Easy to learn for therapists for more optimal results
- E. Can treat wide range of indications
- F. More than 1 million pieces sold under Europe and USA registration and patent

+++++++++++

### A MINIMALLY INVASIVE CATHETER ALTERNATIVE

The UroDapter replaces the catheter and eliminates drawbacks

#### Issues with a Catheter

"Infection is the largest concern with catheter use, whether long term or short term. Catheter-associated urinary tract infections, abbreviated to CAUTIs, are the most commonly faced hospital-acquired infections

Cortese YJ, Wagner VE, Tierney M, Devine D, Fogarty A.

#### Catheterization...

- h Is a painful process
- + Causes microlesions
- Is a major source of infections cannot treat the urethra and the bladder simultaneously
- + Time consuming

## BASICS AND BENEFITS

Safe 

Effective 

Pain-Free


UroDapter® provides a minimally invasive method to introduce therapeutic drugs, gels and diagnostic agents into the bladder through the urethra, without the need of catheterization

Pre-curved concave, elastic isolating collar


Connection to Luer Slip or Luer Lock syringes

Radiused short tip that penetrates the urethral orifice less than 8mm

Flexible collar enables leakage-free instillation



## BENEFITS FOR PATIENTS






All patients having tried the UroDapter® have found it more comfortable than a catheter


- Pain-free treatment
- Cost and time effective treatment
- No risk of infections less complications

## BENEFITS FOR THERAPISTS

- It's possible to treat the bladder and urethra at the same time (which cannot be done with a catheter)
- No device or patient preparation a reduction in treatment cost and time
- Can be used when cath are contraindicated
- Singificantly Leas same
- 5 Using to learned



## CONDITIONS TO TREAT USING THE URODAPTER



Interstitial Cystitis / Bladder Pain Syndrome (IC/BPS)

Recurring Urinary Tract Infections (rUTIs)

Complication of Oncotherapy: Hemorrhagic Cystitis (HC)

Bladder Cancer, Female Patients

Indwelling Catheter Users

Lidocaine

Diagnostic purposes

#### A SUMMARY OF THE POSSIBLE INDICATIONS

Indication	Drugs administered	The chance of instillation needed	Can UroDapter® be used?
IC/BPS	GAG-replenishers (HA, CS etc.), DMSO, Lidocaine	High	Yes
Recurring UTIs	GAG-replenishers	Mid-high (if oral treatment and/or prophylaxis fails )	Yes
Recurrent UTIs	Antibiotics	Mid-high (if oral treatment and/or prophylaxis fails )	Probably
Chemo-cystitis, Radiation cystitis	Aminocaproic acid, alum, silver nitrate, prostaglandin, estrogens	Medium (not the only possibility)	With some drugs, further data needed
Chemo-cystitis, Radiation cystitis	HA, CS, PPS	Medium (new in guidelines)	Yes
Bladder Canc er (chemother apy)	BCG, Mitom ycin	High	Women: probably Men: unsure
Patients with an indwelling catheter	GAG-replenishers, DMSO, Lidocaine	Medium	Unsure
Diagnostic purposes	Dyes and contrast material	N/A	Probably

# BENEFITS FOR MANUFACTURERS AND DISTRIBUTORS

- 4 Class I medical device
- + CE Mark
- + FDA registered and approved
- Unique product with strong IP protection
- More than 1,000,000 sold in the European market





Packed together with IBSA's drug solution, iAluRil®, UroDapter® is available in 85 countries, under the name of iAluadapter®

